Smartline

GPC Cleanup Unit 6500 GPC-Probenreinigungssystem 6500

The **GPC Cleanup Unit 6500** is particularly developed and optimized for sample preparation in pesticide analysis:

- **AOAC Method No. 984.21,**
- **USEPA SW-846 Method 3640A,**
- **EN 12393, EN 1528 and**
- **L 00.00-34 Method in accordance with § 35 LMBG.**

Furthermore it is also the perfect choice to perform general cleanup tasks of a wide range of sample matrices including foodstuffs, tissues, plants and environmental samples. High molecular weight compounds such as lipids, polymers, copolymers, proteins, natural resins, cellular components and steroids are easily dispersed and eliminated.

The **GPC Cleanup Unit 6500** combines a highly flexible modular design with intelligent, easy-to-use software. It can be easily adapted to predetermined laboratory procedures to perform a large variety of GPC sample preparation tasks and requires only a very small amount of bench space.

The pump's versatile dual-piston design stands for pulse-free flow and outstanding flow rate precision and accuracy, achieved through highly precise mechanism and the "smart" drive control. The active back seal washing ensures for longer lifetime and less service.

For automatic injection and to avoid cross-contamination of the samples, the system is equipped with a separate sample loop for each sample (15 loops) and collects each fraction individually.

Every single sample can be dealt separately, allowing specific elution and collection protocols to be defined for each of them. It is also possible to use sample loops with various volumes on various positions (e.g. 10 loops with 1 ml and 5 loops with 5 ml).

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
A fixed wavelength detector optionally can be placed in the Smartline Assistant 6000 for calibration as for real time visualization of the elution profile.

Two combined 17-port/1-channel valves switch between up to 15 sample loops. After consecutive injection of the samples, the system will process all cleanup tasks automatically, significantly reducing labor time.

Any possible carry-over can be eliminated by programmable wash procedures for which the 16th position or sample loop is reserved. The supplied standard sample loops are by default 1.0 ml. Other volumes are possible upon request. Fractions from each sample are collected over a further 17-port/1-channel valve.

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Technical Data

System	Power supply:	90-260V, 47-63 Hz, CE-certified
	Dimensions:	380 x 460 x 390 mm (W x H x D)
	Weight:	13 kg
Pump	Delivery System:	double piston pump
	Flow Range:	up to 10 ml/min
	Flow Precision:	RSD ≤ 0.1% (1 ml/min)
	Pulsation:	< 0.5%
Detector	Type:	fixed wavelength UV detection
	Wavelength:	254 nm, Δλ 20 nm
	Lamp:	mercury vapour lamp
	Linearity:	0 - 2.5 AU
	Noise:	3 x 10 ⁻⁵ AU
	Drift:	2.5 x 10 ⁻⁴ AU
Valves	Switching time:	< 150 ms
	Materials	Stator: Stainless steel Rotor: Vespel
Sample Loops:		up to 16; 1 ml (standard), other volumes on request

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Ordering information

Order No.	Article
A5004 V20	Smartline GPC Cleanup Unit 6500 including UV detector and software
A5004 V21	GPC Cleanup Unit 6500; including software without detector
B40	BioBeads S-X3, 500 x 40 mm, glass column, (for USEPA SW-846 Method 3640A)
B41	BioBeads S-X3, 500 x 25 mm, glass column, (for USEPA SW-846 Method 3640A)
	Other column types available on request.

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Technical data are subject to change without notice.

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Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.