

Preparative Separation of Alkyl Acetates - part one: non overload, 10 mg/mL, Isocratic

Method VEV0061J
HPLC RP-Mode, preparative
Column: Eurospher II 100-10 C18, 150 x 30 mm ID
Phase: Eurospher II 100-10 C18
Conditions: Eluent: 30 % Water; 70% Methanol
 Gradient: isocratic
 Flow rate: 56 mL/min
 Temperature: ambient
 Volume: 500 µL

Order No. 15QE181E2N



Detection: UV 220 nm (+/- 2 nm), 10 mm flow cell (recommended flow cell 3 mm)

Substances: Ethyl acetate, Butyl acetate, Pentyl acetate

Keywords: Alkyl acetates

Chromatogram:

- 1 Ethyl acetate 10 mg/mL
- 2 Butyl acetate 10 mg/mL
- 3 Pentyl acetate 10 mg/mL

