

## Preparative Separation of Alkyl Acetates - part two: overload, 100 mg/mL, Isocratic

**Method** VEV0062J  
**HPLC** RP-Mode, preparative  
**Column:** Eurospher II 100-10 C18, 150 x 30 mm ID  
**Phase:** Eurospher II 100-10 C18  
**Conditions:** Eluent: 30 % Water; 70% Methanol  
 Gradient: isocratic  
 Flow rate: 56 mL/min  
 Temperature: ambient  
 Volume: 1 mL

Order No. 15QE181E2N



**Detection:** UV 210 nm (+/- 4 nm), 10 mm flow cell (recommended flow cell 3 mm)

**Substances:** Ethyl acetate, Butyl acetate, Pentyl acetate

**Keywords:** Alkyl acetates

### Chromatogram:

- 1 Ethyl acetate 100 mg/mL
- 2 Butyl acetate 100 mg/mL
- 3 Pentyl acetate 100 mg/mL

