

## Preparative Separation of Alkyl Acetates - part three: overload, 290 mg/mL, Gradient

**Method** VEV0063J  
**HPLC** RP-Mode, preparative

**Column:** Eurospher II 100-10 C18, 150 x 30 mm ID


**Order No.** 15QE181E2N

**Phase:** Eurospher II 100-10 C18

**Conditions:** Eluent: A Water; B Methanol  
Gradient:

0.0 min	45 % A	55 % B
2.0 min	40 % A	60 % B
3.0 min	40 % A	60 % B
5.0 min	35 % A	65 % B
6.5 min	35 % A	65 % B
10.0 min	20 % A	80 % B
12.0 min	20 % A	80 % B
15.0 min	45 % A	55 % B
16.0 min	45 % A	55 % B

Flow rate: 56 mL/min  
Temperature: ambient  
Volume: 1 mL



**Detection:** UV 210 nm (+/- 2 nm), 10 mm flow cell (recommended flow cell 3 mm)

**Substances:** Ethyl acetate, Butyl acetate, Pentyl acetate

**Keywords:** Alkyl acetates

**Chromatogram:**

- 1 Ethyl acetate 290 mg/mL
- 2 Butyl acetate 290 mg/mL
- 3 Pentyl acetate 290 mg/mL

