

Application Note

► Determination of sugars, ethanol and glycerol in wine



AZURA

Category	Food
Matrix	Wine
Method	HPLC
Keywords	Sugars, Alcohol, Glycerol, Wine, Eurokat Pb, AZURA Compact HPLC
Analytes ID	Saccharose, Glucose, Fructose, Ethanol, Glycerol VFD0121N_A_E


Summary

This application note presents a simple and rapid method for the simultaneous determination of Saccharose, Glucose, Fructose, Ethanol and Glycerol. The method is designed for the quality and process control of wine. All five substances are baseline separated and no hydrolysis of Saccharose was observed. An extra sample preparation is not necessary, just a dilution. This method works with water as a low cost eluent.

Introduction

Rapid identification and quantitation of Saccharose, Glucose, Fructose, Ethanol and Glycerol in wines is important. During the fermentation the Saccharose is cleaved into Glucose and Fructose. Monosaccharides are metabolized into Ethanol and Glycerol. The determination of the sugar content is necessary to obtain the endpoint of the fermentation. Therefore these five compounds need to be checked to control the quality and process of wine fermentation. Furthermore Saccharose, Glucose, Fructose, Ethanol and Glycerol influence the taste of wine. The Glycerol concentration is also a marker for premium wines.

Chemical Structures



Experimental sample preparation

The wine sample can be injected directly after micro filtration and a 3:4 dilution with 2 mmol/l phosphate buffer pH 7.0 (750 µl wine + 250 µl phosphate buffer). All standard solutions were prepared with double-distilled water.

Method parameters

Column	Eurokat Pb, 300 x 8 mm and 30 x 8 mm as precolumn
Eluent A	water
Flow rate	0.6 ml/min
Injection volume	20 µl
Column temperature	60 °C
System pressure	approx. 60 bar
Detection	Refractive Index
Run time	25 min

Results

A calibration was made for Saccharose, Glucose, Fructose, Ethanol and Glycerol. All five substances are baseline separated (figure 1). The calibration concentration was in the range between 1.5 to 75 mg/ml for each substance (figure 2). The experimental data beyond show that the linearity range is higher. The linearity r^2 for all substances was greater than 0.9999. Ethanol shows the smallest detector response and therefore the limit of detection is lower.

Figure 1
Chromatogram at
75 mg/ml





Figure 2
Calibration




Saccharose

Range: 1.5-75 mg/ml

Linearity (r^2) =
0.999960

Linearity range =
0.75-150 mg/ml

Limit of detection =
0.02 mg/ml




Glucose

Range: 1.5-75 mg/ml

Linearity (r^2) =
0.999985

Linearity range =
0.75-150 mg/ml

Limit of detection =
0.03 mg/ml




Ethanol

Range: 1.5-75 mg/ml

Linearity (r^2) =
0.999989

Linearity range =
1.0-150 mg/ml

Limit of detection =
0.06 mg/ml




Fructose

Range: 1.5-75 mg/ml

Linearity (r^2) = 0.999984

Linearity range = 0.75-150 mg/ml

Limit of detection = 0.03 mg/ml



Glycerol

Range: 1.5-75 mg/ml

Linearity (r^2) = 0.999980

Linearity range = 0.75-150 mg/ml

Limit of detection = 0.04 mg/ml

As samples Federweißer (young wine) and red wine were used. Both samples contain no Saccharose. Federweißer contains a high concentration of Fructose and Glucose and a low concentration of Ethanol and Glycerol (Figure 3). The red wine shows a different pattern. The red wine has a high concentration of Ethanol and Glycerol and a low concentration of Fructose and Glucose (Figure 4). The red wine was spiked with Saccharose to show that Saccharose can be analyzed with wine as matrix (Figure 5). The recovery rate was 99 %. The dilution with buffer is necessary, because the wine acids hydrolyse the Saccharose at 60 °C. This results in a lower recovery rate.

Figure 3
Federweißer (young wine)




Figure 4
Red Wine




Figure 5
Red Wine spiked with
saccharose (blue line)
compared to unspiked
sample (green line)




Figure 6
Results of samples

	Federweißer (young wine)	Red Wine
Saccharose	0.00 mg/ml	0.00 mg/ml
Glucose	29.9 mg/ml	1.4 mg/ml
Ethanol	21.2 mg/ml (2.7 %Vol)	93.6 mg/ml (11.9 %Vol)
Fructose	43.6 mg/ml	3.4 mg/ml
Glycerol	1.5 mg/ml	7.0 mg/ml

Method performance

Limit of detection	0.02 – 0.06 mg/ml (S/N = 3)
Linearity (r^2)	0.999960-0.999989
Linearity range	0.75-150 mg/ml
Retention time precision*	< 1 % RSD
Peak area precision*	< 2 % RSD

Conclusion

The use of this method facilitates to determine Saccharose, Glucose, Fructose, Ethanol and Glycerol in young and red wine without the loss of Saccharose. Saccharose hydrolyse very fast under acidic conditions. Therefore the sample is diluted with buffer.

If you use a sulfonated cross-linked styrene-divinylbenzene copolymer in H-form, Saccharose will hydrolyze over 25 °C. Using the H-form under this temperature leads to an incomplete baseline separation of acids and monosaccharides. Due to the high concentration range of the wine compounds, it is possible that a high concentrated substance overlaps a low concentrated one. The newly described method here circumvents this problem.

Another method to determine the sugar is the use of a HILIC (hydrophilic interaction chromatography) method. HILIC works very fine with sugar, but it is impossible to analyze Ethanol. Furthermore you need a high concentration of acetonitrile. The presented method only needs water as an eluent, which keeps the cost low.

The fact that Methanol and organic acids in wine cannot be assigned is avoidable by using KNAUER Eurokat H columns.

Authors

Dr. Stefan Weiz, Columns and Applications Department, KNAUER

Juliane Böttcher, Columns and Applications Department, KNAUER

Physical properties of recommended column

Packing materials in the Eurokat family are sulfonated cross-linked polystyrene copolymers. The lead forms are the phase of choice for the separation of sugars and sugar alcohols. The analytes are separated by ion exclusion and ligand exchange chromatography



Stationary phase	sulfonated cross-linked styrene-divinylbenzene copolymer in lead form
USP code	L34
Cross Linkage	6 %
Particle size	10 µm
Form	spherical
Dimensions	300 x 8 mm and 30 x 8 mm
Order number	30GX350EKN and 03GX350EKN

Recommended instrumentation



This application requires an isocratic HPLC system equipped with degasser, autosampler, column oven, and refractive index detector. Other configurations are also available. Please contact KNAUER to configure a system that's perfect for your needs.

Description	Order No.
AZURA Compact HPLC isocratic	AYIABACA
Autosampler KNAUER Optimas	A5007
Column Thermostat	A0585
RI Detector 2300	A5160
AZURA Eluent tray E 2.1L	AZC00

Contact information

KNAUER
Wissenschaftliche Geräte GmbH
Hegauer Weg 38
14163 Berlin, Germany

Tel: +49 030 809727-0
Fax: +49 030 8015010
E-Mail: info@knauer.net
Internet: www.knauer.net

