

## Product information - Flow splitter

The flow splitter consists of a needle valve ② and a T-piece ① which can be used in systems up to 50 bar. The T-piece ① splits the flow into two streams. The split ratio can be regulated by using the needle valve ② and capillaries with different inner diameters.

### Legend:

- ① T-piece
- ② Needle valve





Fig. 1: Flow splitter

## Mounting the flow splitter

**Prerequisites** ■ The flow splitter has been unpacked.

**Tool** ■ Capillary cutter

- Process**
1. Connect the UV detector with the T-piece ①.
  2. Connect one port with the fraction collector using 1/8" tubing.
  3. Connect the last port with the needle valve ② using 1/16" tubing. The used inner diameter depends on the flowrate and the intended split ratio. For higher flowrates or higher split ratio use smaller inner diameters.
  4. Connect the needle valve ② with the conductivity monitor using 1/16" tubing (see 3.).



**Fig. 2:** Mounting the flow splitter

## Repeat orders

Name	Order no.
Flow splitter	A5813
Adapter PEEK	A05841
Capillary PEEK, 1/16", inner diameter 0.13 mm	A2522
Capillary PEEK, 1/16", inner diameter 0.25 mm	A2524
Capillary PEEK, 1/16", inner diameter 0.5 mm	A2526