AURA

▶ **Beiblatt Ultraschallluftsensor**

Supplement Ultrasonic Air Sensor

V6876

HPLC

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Table of Contents

NOTE! For your own safety, read the manual and always observe the warnings and safety information on the device and in the manual.

Intended Use	12
Installation	13
Circuit Variants	13
Air Sensor without Valve	13
Mounting	14
Operation	16
Error Message	17
Air Bubbles	17
No Voltage Supply	17
Repeat Orders	18

Intended Use

The air sensor is a sealed ultrasonic sensor which recognizes liquid or air in a hose. It was tested for aqueous applications and a broad range of different buffers will be recognized reliably.

The intended use of the sensor is to detect air bubbles larger than 5 μL in suction hoses of AZURA and Smart-line pumps.

In case air bubbles occur, the air sensor directly causes the pump to switch off. The LEDs in the connector indicate the status of the air sensor:

- POWER LED lights green: Voltage supply is on.
- SIGNAL LED lights yellow: Error condition has been met, air bubble occurred.

Installation

On the side of the outlet a programmable electric switch can be found, which is compatible with the TTL inlets of the AZURA and Smartline product lines.

Circuit Variants

The air sensor is calibrated ex-works and is operated directly at the analogue inputs of pumps and valve drives. In pumps, which are not equipped with a low-pressure gradient, it is preferably mounted in front of or next to the pump head. In systems with several sensors, this makes it easier to determine which sensor is operated on which pump.

In principal, it is possible to achieve the program-controlled processing of the air sensor signal in HPLC, but for this a connector variant is needed, which is different from the one described in this supplement.

On most pumps, the connection is realized at the 24 V output of the device. If such a connector is part of the system, no separate power supply plug is necessary.

Air Sensor without Valve

Immediate shut-down of a pump if air bubbles are being recognized.

- The sensor is connected to the power supply. +24 V, brown on PIN5; blue (ground) on PIN6.
- The signaling cable (bw) is connected to the ERROR IN inlet of the pump: PIN6 of the pin header.

Mounting

The hose is inserted into the measuring channel of the air sensor. This prevents the air sensor from getting into contact with liquids.

Prerequisites

- The pump has been switched off.
- The air sensor has been mounted on the holding bracket


Auxiliary material

Lint-free cloth

CAUTION! Acoustic short circuit if measuring channel is impure. Avoid other buffers or moisture from entering. Wipe the measuring channel with a lint-free cloth.

Process and figure

1. Mount the air sensor ③ straight or angularly near the pump head.
2. Remove the cap of the air sensor and wipe the measuring channel with a lint-free cloth.
3. Clamp the hose ① into the measuring channel. The inlet ④ into the sensor channel is found on the bottom, the outlet ② on the top.

**Fig. 1** Hose in the air sensor**Next steps**

Bring the air sensor into service.

Operation

After commissioning the pump, the green LED POWER ② and the yellow LED SIGNAL ① light up because the hose is still filled with air. After the hose has been filled without air bubbles, the SIGNAL LED goes out and the POWER LED lights up.


Abb. 2 LEDs in the connector

Error Message

Air Bubbles

Malfunction	Solution
The yellow SIGNAL LED is alight. The pump stops to transport fluid.	The hose has to be filled entirely with fluid and the malfunction message has to be reset.

No Voltage Supply

Malfunction	Solution
The yellow SIGNAL LED and the green POWER LED are both out.	The 24 V output has to be switched off on the AZURA device.

Repeat Orders

NOTE Send the air sensor back to KNAUER for disposal. Because hoses can be contaminated, they cannot be sent back to KNAUER but have to be disposed of adequately.

A70099 - Air sensor

Air sensor for analytical pumps with pump heads < 100 ml/min.

Name	Order number
Bioline Air Sensor	A70099V1
Bioline sensor cable with angle connector M12×1	A70099V2

A70101 - Air sensor for 6 mm hose

Air sensor for preparative pumps with pump heads > 250 ml/min.

Name	Order number
Bioline Air sensor	A70099V3
Bioline sensor cable with angle connector M12×1	A70099V2

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© Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
All rights reserved.
The information in this document is
subject to change without prior notice.
Translation of the original German
edition of this supplement, version 1.0.
2013-06-21
Printed in Germany on environmen-
tally friendly paper from sustainable
forests.

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® AZURA are registered trademarks of