Handling of PEEK fittings


PEEK is an organic solvent-resistant polymer which is resistant against high temperatures and high pressures. These advantages make PEEK the material of choice for fittings. But, PEEK fittings can be shorn off more easily than stainless steel, if too much torque is applied. On the other hand, a PEEK fitting will not seal a connection correctly, if insufficient torque is used. Therefore, PEEK fittings must be handled with care in order to seal a capillary connection tightly.

We recommend you use a ¼" wrench to tighten the hex-head PEEK fittings listed below. **Do not apply more than 0.5 Nm torque**, otherwise the head of the fitting may deform and shear off. With the tubing held in place, finger-tighten the fitting completely and then wrench-tighten it another ¾ turn. If the fitting is still not swaged to the tubing or if the connection is leaky, remove the fitting and check the capillary is correctly positioned on the fitting. Do not overtighten the fitting more than one whole turn after finger-tightening.

4 different kinds of PEEK fittings can be used with your system, depending on your configuration.

Note Please refer to the table below to identify the correct fitting. Fittings can be distinguished by the marking on the head.

Figure 1

Type	Explanations
P3860 with short head and tip	P3860 and P3860V1 can be used to connect capillaries to KNAUER columns
P3860V1 with long head and tip	P3860 and P3860V1 can be used to connect capillaries to KNAUER columns
P3860V2 with long head and without tip. This screw fitting is marked by one notch ①	P3860V2 is to be used for the outlet of your PDA-1 flow cell.
P3860V3 with long head and short tip. This screw fitting is marked by two notches ②	P3860V3 is to be used for the inlet of your PDA-1 flow cell. P3860V3 is also compatible to Upchurch® scientific standard ports.

Attention! Please be aware that P3860V2 and P3860V3 **must** be used to connect the PDA-1 flow cell. Using other fittings will destroy the flow cell and/or will lead to a leakage.

Please keep in mind that all PEEK fittings are single-use consumables. If a fitting is screwed into a column (for example), the cone of the fitting will be swaged to seal this connection perfectly tight. Afterwards, the fitting can only be used with this column exclusively. It can be removed and reinstalled several times, but it can **not** be used with a different device.

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