# Knauer Interface Box IF2

## Manual/Handbuch

V7050 03/2007

Wissenschaftliche Gerätebau  
Dr. Ing. Herbert Knauer GmbH  
Hegauer Weg 38  
14163 Berlin, Germany  
Telephone: +49-(0)30-809727-0  
Telefax: +49-(0)30-8015010  
E-Mail: [info@knauer.net](mailto:info@knauer.net)  
Internet: [www.knauer.net](http://www.knauer.net)

# CONTENTS

General information	3
Contents of the accessory kit for the Interface Box IF 2	3
Technical data	3
Instrument Display and Connections	3
Front view of the instrument	3
Back panel of the instrument	4
Analog Inputs (AN.IN)	4
Analog Outputs (AN.OUT)	5
Digital Inputs (DIG.IN)	5
Digital Outputs (DIG.OUT)	5
RS-232 Interface	5
Connecting the Interface Box IF 2	6
Power supply	6
Interface Box IF 2 on the PC (RS-232)	6
Analog Inputs	6
Analog Outputs	6
Digital Connections	6
Installation of the WAGO service cable	7
Example of the Cabling	7
Warranty statement	13
Declaration of conformity	14

# INHALT

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## General information

The Interface Box IF 2 enables the recording or control of instruments that are not supported by the KNAUER software over four channels. It can be used for detector recording over the analog inputs, it can control the flow rate for HPLC pumps over the analog outputs, or it can switch valves over the digital outputs. However, the respective instrument must support this option. Please refer to the respective instrument manual regarding this option.

The Interface Box IF 2 does not have a stand-alone function; therefore, control software (KNAUER ChromGate® as of version 3.1 or KNAUER EuroChrom® 2000 for Windows as of version 3.05) is mandatory.

**The Interface Box IF 2 functions that are supported by the control software can be found in the software manual. Configuring and programming the Interface Box IF 2 control software (as a Smartline Manager 5000 Interface Module) is described in the respective manual of the software.**

### Contents of the accessory kit for the Interface Box IF 2

- 4 x WAGO strips, 9-pin
- 1 x WAGO lifting latch
- 1 x Zero modem connection cable (RS-232, socket/socket)
- 1 x Connection cable Cinch – cord
- 3 x Flat belt cable, 10-pin

### Technical data

Power Supply	90-260 V, 47 - 63Hz, 100 W
Dimensions	226 x 135 x195 mm (W x H x D)
Weight	3.9 kg

## Instrument Display and Connections

### Front view of the instrument

Figure 1 Front view of the Interface Box IF 2

The light-emitting diodes give information about the condition of the Interface Box IF 2. If the Interface Box IF 2 is connected with the power supply cable and is switched on the green “Power” LED lights

permanently. The green “Busy” LED lights shortly during calibration or when a command from the software is received. The instrument is calibrated when it’s turned on or a chromatographic run is started. If the control software regularly inquires about the status of the instrument (e.g. in second intervals), then the “Busy” LED will permanently light.

The red “Error” LED lights when a calibration error is detected in a self-test or when the software sends a command to the instrument that does not exist.

### Back panel of the instrument


Figure 2 Back panel view of the Interface Box IF 2

**Please pay attention to the maximal voltage and current of the inputs and outputs on the Interface Box IF 2 and on the instruments that will be controlled. If the maximum voltage and current are exceeded, KNAUER is not liable for defective instruments!**

### Analog Inputs (AN.IN)

To record from analog signals, up to four channels or instruments can be connected with the Interface Box IF 2. The channels are independent of each other, meaning all channels can be simultaneously operated from either one system or various systems.

Voltage Range:	-2.56 to +2.56 V (bipolar mode)
Absolute maximum ratings:	-10 V to +10 V
Input impedance:	10 MΩ
Maximum resolution:	24 bit
Minimum noise level:	7 μV (1 Hz, time constant 0.1 s)
Maximum data rate	up to 10 Hz (each channel)
Gain factor (for all channels):	1, 2, 4, 8, 16

## Analog Outputs (AN.OUT)

The outputs distribute analog voltage signals, so that instruments containing a respective input can be controlled. Depending on the number of outputs, up to four instruments can be controlled. The outputs can be independently programmed which enables the programming of high pressure gradient system pumps.

Voltage range:	0 to +10 V
Minimum voltage step: (DAC resolution)	2.5 mV (12 bit)
Maximum load resistance:	2 k $\Omega$

## Digital Inputs (DIG.IN)

These four inputs are intended as a trigger (start signal) input for the start signal that is distributed by the injection system. They work with contact closure on mass and are compatible with TTL/CMOS or the Open Collector. If the channels are used in separated systems, then every channel can be separately started.

Maximum high level input voltage:	+15 V
Minimum high level input voltage:	+2.5 V
Maximum low level input voltage:	+1 V
Minimum low level input voltage:	-15 V
Maximum input current at V (in) = 0.5 V:	10 mA

## Digital Outputs (DIG.OUT)

These four outputs are carried out in electromechanical relays with a single pole double trough (SPDT). They can be selected as NO (normally open; connection 2-3, see the back panel view) or NC (normally closed connection 1-2, see the back panel view). Pulse or continuous signals can be programmed. The length of a pulse contains approximately 1s; the length of the continuous signal can be freely programmed.

Output type:	selectable NO or NC
Maximum switching voltage:	maximum 175 V DC
Maximum switching current:	maximum 0.25 A

## RS-232 Interface

This interface is a 9-pin standard RS-232 interface (also known as a serial interface or COM port) and is connected to a PC in order to realize the software control of the interface.

### Plug Charging

Pin 2	RxD
Pin 3	TxD
Pin 5	GND

### Interface connection settings

Transfer rate:	19200 baud fix,
Data width:	8 bit
Stop bit:	1
Parity:	no

## Connecting the Interface Box IF 2

### Power supply

The Interface Box IF 2 is equipped with a default power supply. The required cable is included in the Interface Box IF 2's accessory kit.

### Interface Box IF 2 on the PC (RS-232)

In order to connect the Interface Box IF 2 with a PC, a standard zero modem cable is required and is included in the Interface Box IF 2 accessory kit. The configuration of the Interface Box IF 2 are specified under "RS-232 Interface".

### Analog Inputs

In order to connect the analog output of the respective instrument (e.g. the integrator output of a detector) with the analog input of the Interface Box IF 2, use the Cinch cable that is made by the instrument's manufacturer and delivered with the instrument. Alternatively, a suitable cable can be found in the Interface Box IF 2 accessory kit. Please pay attention to the maximal voltage. If it is exceeded, the instrument can be damaged.

**The analog input and digital output of one channel can not be used simultaneously.**

### Analog Outputs

The analog outputs must be connected with the respective control input of the instrument to be controlled. Information regarding the control input of the instrument can be found in the respective instrument's manual. Please pay attention to the maximal voltage and current of the in and outputs on the Interface Box IF 2 and on the instruments that will be controlled. If the maximum voltage and current are exceeded, KNAUER is not liable for defective instruments!

The WAGO plug strips and cable that are required for connection are included in the Interface Box IF 2 accessory kit.

**The analog input and digital output of one channel can not be used simultaneously.**

### Digital Connections

The digital inputs can be used to forward a trigger (start) from the injection unit. The trigger (start) signal must always be connected to the digital input of the same channel to which the detector for data recording is connected.

The digital outputs can be programmed over the control software.

The WAGO plug strips and cable that are required for connection are included in the Interface Box IF 2 accessory kit.


If a NO is required, pins 2 and 3 (see back panel view) must be connected. However, if a NC is required, pins 1 and 2 must be connected (see back panel view).

**If a channel's analog input is used, the digital output of this channel can not be programmed. In this case, a pulse signal is delivered over the digital output DIG.OUT when a run is started. This can be used as an autozero signal for a detector.**

## Installation of the WAGO service cable

In order to clamp the cable on the desired positions of the WAGO plug, please follow the directions below:

The cable is now firmly anchored in the plug strip.

Figure 3 Mounting Plug Strips

## Example of the Cabling

The exemplary cabling plan demonstrates how the respective connections can be switched. However, it is important that the analog input of a channel is never simultaneously operated with the digital output of the same channel, except as an autozero signal for the connected detector.


Figure 4 Example of the Wiring

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
Abb. 4 Beispiel - Verkabelungsplan

## Warranty statement

The warranty period of the Interface Box IF 2 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Defective Interface Box IF 2 should be sent to the manufacturer for repair.

Wissenschaftliche Gerätebau  
Dr. Ing. Herbert KNAUER GmbH  
Hegauer Weg 38  
D-14163 Berlin  
Tel: 030 – 809 727 – 0  
Fax: 030 – 801 50 10  
e-mail: [info@KNAUER.net](mailto:info@KNAUER.net)  
[www.KNAUER.net](http://www.KNAUER.net)

If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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## Declaration of conformity

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**Manufacturer's name and address:  
Herstellername und -adresse**

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**Interface Box IF 2  
Bestellnummer A 5330**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)  
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)  
EN 6100-3-2 (10/2006)  
EN 61326-1 (10/2006)

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