

Title	PurityChrom® 6 Release Notes			Page	1 von 4
Version:	2.0	Authors:	MaG, PS	Date:	24.10.2023

Table of Content:

1. Supported Instruments2

- Pumps
- Detectors (control and digital data acquisition)
- Interfaces
- Autosamplers
- Valves
- Fraction Collectors
- Column ovens

2. Installation3

- System Requirements
- Operating Systems
- Language Selection
- Installation under Windows
- Energy Saving Features

3. Release History4

Title	PurityChrom® 6 Release Notes			Page	2 von 4
Version:	2.0	Authors:	MaG, PS	Date	24.10.2023

1. Supported Instruments

Pumps

- AZURA Pumps P 2.1S/P 4.1S
- AZURA Pumps P 6.1L
- AZURA Pump P 2.1L
- Blueshadow Pump 80P

Detectors

- AZURA UV-Detector UVD 2.1S
- AZURA UV-Detector UVD 2.1L
- AZURA UV-Detector UV MWD 2.1L (4 channels, up to 50 Hz)
- AZURA Diode Array Detector DAD 2.1L
- AZURA Diode Array Detector DAD 6.1L
- AZURA Refractive Index Detector RID 2.1L
- AZURA Conductivity Detector CM 2.1S
- AZURA pH Monitor pH 2.1S
- BlueShadow 50D
- Mikron 81

Flowmeter

- Bronkhorst Mini Cori-Flow M12, M13, M14

Valves

- AZURA Valve Drives V2.1S, VU4.1
- VICI Valco Drives

Fraction Collectors

- ISCO Foxy R1
- ISCO Foxy R2
- LABOCOL Vario-4000, Vario-4000 Plus

Autosampler

- AZURA AS 6.1L

Heating Devices

- SCPA Eluentheater
- AZURA CT 2.1

Interfaces

- Interface Box IFU 2.1 LAN (Analog Inputs only)
- SCPA Ethernet Eventbox

Assistants

- AZURA ASM 2.1L
- AZURA ASM 2.2L

Miscellaneous

- SCPA External Pressure Sensors
- Gameport Adapter

Title	PurityChrom® 6 Release Notes			Page	3 von 4
Version:	2.0	Authors:	MaG, PS	Date	24.10.2023

2. Installation Requirements

Operating Systems	Windows 10 64 Bit Windows 11
CPU/RAM/Disc Space	Intel Core i5 6 th Generation or comparable/ 8 GB/ Disc Space 10GB
Network	Second network card is required for network integration
Graphics	Screen with minimal resolution 1024 x 768
Connectors and Slots	USB for license dongle COM, USB or LAN according to connected instruments
Input Device	Mouse and Keyboard

Language Selection

PurityChrom6 is available in English.

Installation and use under Windows:

Administrator rights are necessary during installation of PurityChrom6 and the Windows User Account Control settings should be set to “Never notify” in order to deactivate the User Account Control. It may be re-activated after the installation.

The user needs write and read permission for registry and folders C:\ProgramData\PurityChrom6 and C:\Program\PurityChrom6.

PurityChrom6 needs to be excluded from Microsoft Security Essentials. Furthermore, the Windows Defender and antivirus software should be deactivated.

Energy Saving Features

Automatic switching to energy saving mode or hibernating has to be disabled i.e. settings for “Turn Off Hardware”, “Turn off Computer” or “System Standby” must be set to “Never” in the control panels energy options.

Title	PurityChrom® 6 Release Notes			Page	4 von 4
Version:	2.0	Authors:	MaG, PS	Date	24.10.2023

3. Release History

In the latest version R2-2, the following changes and improvements have been added or bugs have been fixed from version R1-2:

- Bug Fix: incorrect calculation of Asymmetry was fixed
- Bug Fix: incorrect calculation and naming between MPa/kPa into bar
- Bug Fix: incorrect calculation of first fraction if Set Volume Limit is used in a method

- New supported devices:
 - o AZURA DAD2.1L/DAD6.1L
 - o AZURA RID2.1L
 - o AZURA CT2.1
 - o AZURA AS6.1L
 - o BlueShadow 50D

- New features:
 - o High Flow Pump Systems
 - o Stacked Injections
 - o Valve Daisy Chaining
 - o PID Controller
 - o System Documentation
 - o Read Out and Setting Pump Head Parameters for P4.1S
 - o Read Out of Pump Head Parameters for 80P/P2.1L
 - o Online Custom Annotation
 - o Valve "Next Position"
 - o "Notification/Service/Warning/Error" as method defined functions