

## Device support in KNAUER CDS / KNAUER device support in 3<sup>rd</sup> party software OpenLAB<sup>®</sup> CDS EZChrom Edition / ClarityChrom<sup>®</sup> / Chromeleon<sup>™</sup> / ChromGate<sup>®</sup> / PurityChrom<sup>®</sup>

The tables list the software package and driver version (in case multiple versions exist); the version number listed is the **first version** that supports the respective device. **Limitations are marked with superscript**, e.g. <sup>1)</sup>. Limitations may affect device functionality, control interfaces or software versions. "n/a" means, that a device is not supported, "-" that no information is available. If a device is supported via different communication ports, for most devices the desired port must be selected in the device directly. From version 7.4.1, ClarityChrom<sup>®</sup> *Prep* is discontinued; all its prep-functionality is now available in ClarityChrom<sup>®</sup> with KNAUER FRC license. For some time, ChromGate<sup>®</sup> was also available as "PLATINblue Edition" for the PLATINblue devices, while the standard version is named as "regular". From driver version 2825 all drivers are available in the regular version; the PLATINblue Edition was discontinued. ChromGate<sup>®</sup> was discontinued in 2013 at all. OpenLAB<sup>®</sup> CDS EZChrom Edition version number lists Shared Service version and EZChrom version. Drivers for OpenLAB<sup>®</sup> CDS Chemstation Edition or OpenLAB<sup>®</sup> CDS 2 are currently not available. KNAUER **OEM** devices are **not** supported by the regular KNAUER drivers. At the end of the document other software packages, e.g. KNAUER Mobile Control, are listed.

### Table of contents

Pumps.....	2	Fraction Collectors.....	16
Detectors .....	4	Mobile Control - KNAUER .....	20
Assistant.....	10	EuroChrom for Windows 3.05 .....	22
Autosamplers .....	11	Other Software.....	23
Analog/Digital Converters.....	13	LabVIEW (National Instruments) .....	23
Thermostats/Column Ovens .....	14	HyStar <sup>™</sup> .....	23
Electrical Valve Drives.....	15	Xcalibur <sup>™</sup> .....	23
Other Supported Instruments .....	15		

<b>Pumps</b>	<b>OpenLAB® CDS EZChrom Edition CDS version and Driver Version</b>	<b>ClarityChrom® CDS version and Driver Version</b>	<b>Chromeleon™ CDS version Driver version</b>	<b>ChromGate® CDS version and Driver Version</b>	<b>PurityChrom®</b>
AZURA® Pump P 6.1L	OLEE A.01.05 / A.04.05 Driver 1022 <sup>14)</sup>	CC 5.0.2 Driver 4095 <sup>1, 15)</sup>	Cm 7.2 Regular Class 3 driver 1024 <sup>1)</sup>	n/a	Version: 5.07.039 <sup>1)</sup>
AZURA® Pump P 2.1L	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3078 <sup>1)</sup>	Cm 7.2 Regular Class 3 driver 1024 <sup>1)</sup> Cm 6.8 Generic driver <sup>2)</sup>	CG 3.3.2 Driver 2825 <sup>3)</sup>	Version: 5.07.039 <sup>1)</sup>
AZURA® Pump P 2.1S/P 4.1S <sup>6)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3078	Cm 7.2 Regular Class 3 driver 1056 <sup>1)</sup> Cm 6.8 Generic driver <sup>2)</sup>	CG 3.3.2 Driver 1787	Version: 5.07.039
PLATINblue Pump P-1	OLEE A.01.05 / A.04.05 Driver 1017	CC 5.0.3 Driver 4102	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.2 PLATINblue Edition 1751; Regular 2825	n/a
Smartline Pump 1050 <sup>4)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.6.5 Driver 1100	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.2 Driver 1751	Version: 5.07.039
Smartline Pump 1000 <sup>13)</sup>	n/a	CC 2.4.4	Cm 6.8 Generic driver <sup>2)</sup>	CG 3.1	Version: 5.07.039
Smartline Pump 100 <sup>5), 6)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.6	Cm 6.8 Generic driver <sup>2)</sup>	CG 3.3.1 Driver 1640	Version: 5.07.039
Preparative Pump 1800 (also WellChrom K-1800) <sup>13)</sup>	n/a	CC 2.6	n/a	CG 3.1	Version: 5.07.039
WellChrom Pump K-1001 <sup>13)</sup>	n/a	CC 2.6	n/a	CG 2.55 <sup>7)</sup>	Version: 5.07.039
WellChrom Pump K-1000 <sup>13)</sup>	n/a	CC 2.6	n/a	CG 2.55 (not supported in version 3.1) <sup>8)</sup>	n/a
WellChrom Pump K-501 <sup>9)</sup>	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom Pump K-120	n/a	CC 2.6	n/a	CG 2.8	n/a
BlueShadow Pumps 10P/20P <sup>6)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3078	Cm 6.8 Generic driver <sup>2)</sup>	CG 3.3.2 Driver 1787	Version: 5.07.039

Pumps (continued)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
BlueShadow Pump 40P <sup>10)</sup>	OLEE A.01.05 / A.04.05 Driver 1022	CC 5.0.2 Driver 4095	n/a	n/a	Version: 5.07.039
BlueShadow Pump 80P <sup>10)</sup>	OLEE A.01.05 / A.04.05 Driver 1022	CC 5.0.2 Driver 4095	Cm 6.8 Generic driver <sup>2)</sup>	n/a	Version: 5.07.039
Kontron Pumps 320/322/325	n/a	CC 3.0.5 <sup>11)</sup>	-	CG 3.1.7 Driver 446 <sup>12)</sup>	n/a
Kontron Pumps 420/422/425	n/a	CC 3.0.5 <sup>11)</sup>	-	CG 3.1.7 Driver 446 <sup>12)</sup>	Version: 5.07.039
Kontron Pumps 520/522/525	n/a	CC 3.0.5 <sup>11)</sup>	-	CG 3.1.7 Driver 446 <sup>12)</sup>	n/a

- 1) - Constant Pressure Mode / Isobar Mode is not supported
- 2) - generic drivers do not support all of the devices features (e.g. no Isobar Mode); only control via RS-232; should also work with Chromeleon™ 7.x
- 3) - support of RS-232 port from driver version 2844
- 4) - for control via RS-232 only baud rate of 115200 is supported
- 5) - can be controlled as pump K-120 only via RS-232 without pressure sensor in former software versions
- 6) - to set a static IP address a separate tool is required
- 7) - from ChromGate® version 3.1 only with firmware version 5.xx or higher; pumps with firmware up to version 4.XX must be configured as K-1000
- 8) - from ChromGate® version 3.1 no pressure information; also K-1001 with firmware up to version 4.XX, no pressure trace
- 9) - from firmware version 1.23, model K-500 is not supported
- 10) - supported from Clarity version 4.X
- 11) - driver was developed by DataApex
- 12) - special RS-232 cable is required
- 13) - for full functionality the KNAUER RS-232 Y-cable is required (KNAUER article number A0755)
- 14) - support of switching SSVs during run from driver version 4.2.0.2161
- 15) - support of switching SSVs during run from driver version 6021

<b>Detectors</b> [for all PDA or scanning UV detectors the 3D UV license option is required to acquire and review the scan (3D) data]	<b>OpenLAB® CDS EZChrom Edition CDS version and Driver Version</b>	<b>ClarityChrom® CDS version and Driver Version</b>	<b>Chromeleon™ CDS version Driver version</b>	<b>ChromGate® CDS version and Driver Version</b>	<b>PurityChrom®</b>
AZURA® Detector UVD 2.1L <sup>1)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3078	Cm 7.2 Regular Class 3 driver 1024 Cm 6.8 Generic driver <sup>5)</sup>	CG 3.3.2 Driver 2825	Version: 5.08.034
AZURA® Detector MWD 2.1L	OLEE A.01.05 / A.04.05 driver 2068	CC 5.0.5 Driver 5011	Cm 7.2 Regular Class 3 driver 1024	n/a	Version: 5.07.067 <sup>15)</sup>
AZURA® Detector UVD 2.1S <sup>2)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3055	Cm 7.2 Regular Class 3 driver 1056	CG 3.3.2 Driver 1787	Version: 5.07.039
AZURA® Detector DAD 6.1L <sup>23)</sup>	OLEE A.01.05 / A.04.05 Driver 2047	CC 5.0.5 Driver 4110	Cm 7.2 Regular Class 3 driver 1024	n/a	Version: 5.07.102 <sup>16)</sup>
AZURA® Detector DAD 2.1L <sup>23)</sup>	OLEE A.01.05 / A.04.05 Driver 2068	CC 5.0.5 Driver 5011	Cm 7.2 Regular Class 3 driver 1024	n/a	Version: 5.07.102 <sup>16)</sup>
AZURA® Detector RID 2.1L / RID 2.1L Highflow	OLEE A.01.05 / A.04.05 Driver 2092	CC 6.0.1 Driver 5037 (compatible with CC/CCP 5.x versions)	Cm 7.2 Regular Class 3 driver 1038	n/a	Version: 5.08.024 <sup>17, 25)</sup>
AZURA® Conductivity Monitor CM 2.1S <sup>2)</sup>	OLEE A.01.05 / A.04.05 Driver 1022	CC 5.0.2 Driver 4095	n/a	Generic Driver June 2013 <sup>5)</sup>	Version: 5.07.039 <sup>18)</sup>
AZURA® ECD 2.1	n/a	Supported as Antec Decade Elite; CC 5.X, USB, LAN	Cm 7.2 – coming soon	n/a	n/a
PLATINblue Detector MW-1 <sup>3)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC6.1.0; Driver 5015	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.2 PLATINblue Edition 1751; Regular 2825	n/a
PLATINblue DAD PDA-1 <sup>23)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 5.0.2 Driver 4095 <sup>4)</sup>	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.2; PLATINblue Edition 1751; Regular 2825	n/a
Smartline UV Detector 2500 <sup>6)</sup>	n/a	CC 2.4.4	Cm 6.8 Generic driver <sup>5)</sup>	CG 3.1	Version: 5.07.039

Detectors (continued)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Smartline UV Detector 2520 <sup>3), 7)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.5 Driver 3033	n/a	CG 3.3.2 Driver 1751	n/a
Smartline UV Detector 2550 <sup>3), 7, 23)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.6.5 Driver 1079	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.1 Driver 1640	n/a
Smartline UV Detector 2600 <sup>23)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4	Cm 6.8 Generic driver <sup>5)</sup>	CG 3.1	Version: 5.07.039
Smartline UV Detector 200	n/a	CC 2.6	n/a	CG 3.1.7 Driver 1466	n/a
Smartline PDA Detector 2800 <sup>9), 23)</sup>	n/a	CC 2.4.4	n/a	CG 3.1.7 Driver 446	Version: 5.07.039
Smartline PDA Detector 2850 <sup>23)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.6.5 Driver 1100	n/a	CG 3.3.2 Driver 1751	n/a
Smartline RI Detector 2300	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4	n/a	CG 3.1	n/a
Smartline RI Detector 2400	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4	n/a	CG 3.1	n/a
Smartline Conductivity Monitor S2900	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3055	n/a	Generic Driver September 2010	n/a
WellChrom UV Detector K-2000 <sup>8)</sup>	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom UV Detector K-2001	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom UV Detector K-2500 <sup>8)</sup>	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom UV Detector K-2501 <sup>6)</sup>	n/a	CC 2.6	n/a	CG 2.8	

Detectors (continued)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
WellChrom UV Detector K-2600 <sup>22, 23)</sup>	OLEE A.01.05 / A.04.05 Driver 2026	CC 2.6	Cm 6.8 Generic driver <sup>5)</sup>	CG 2.55	n/a
WellChrom PDA Detector K-2700 <sup>9, 23)</sup>	n/a	PCI with RS-422 (PCI card driver supports up to Windows 7 32-bit) / LAN	n/a	CG 2.55, LAN from 3.1	n/a
WellChrom PDA Detector K-2800 <sup>9, 23)</sup>	n/a	PCI with RS-422 (PCI card driver supports up to Windows 7 32-bit) / LAN	n/a	CG 2.8, LAN from 3.1	n/a
WellChrom RI Detector 2300	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom RI Detector 2301	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom RI Detector 2400	n/a	CC 2.6	n/a	CG 2.8	n/a
WellChrom RI Detector 2401	n/a	CC 2.6	n/a	CG 2.8	n/a
Shimadzu fluorescence detector RF-10AxI; onboard COM or PCI/PCIe COM board required <sup>10)</sup>	n/a	CC 2.6	n/a	CG 3.1	n/a
Shimadzu fluorescence detector RF-20A/Axs serial control with special <b>firmware version 0.90</b> only from KNAUER (limited functionality); onboard COM or PCI/PCIe COM board required <sup>11)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3055	n/a	CG 3.3.2; Driver 1787	Version: 5.08.024 <sup>24)</sup>
Shimadzu fluorescence detector RF-20A/Axs <b>with CBM-20A</b>	Shimadzu Control with official Shimadzu driver; requires special license, driver not developed by KNAUER <sup>19)</sup>	CC 7.2.0, only LC-10 functionality, driver developed by DataApex	Cm 7.2 with official Shimadzu driver; requires special license, driver not developed by KNAUER <sup>20)</sup>	n/a	n/a

Detectors (continued)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Virtual Detector	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.6.5 Driver 1079	n/a	CG 3.1.6	n/a
BlueShadow Detector 10D <sup>3)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3055	n/a	CG 3.3.2 Driver 1787	Version: 5.07.039
BlueShadow Detector 40D <sup>3, 7)</sup>	OLEE A.01.05 / A.04.05 driver 1022	CC 5.0.2 Driver 4095	n/a	n/a	n/a
BlueShadow Detector 50D <sup>3, 7, 23)</sup>	OLEE A.01.05 / A.04.05 driver 1022	CC 5.0.2 Driver 4095	n/a	n/a	n/a
Alltech Conductivity Detector 650	n/a	n/a	n/a	CG 3.1.7 Driver 446	n/a
Kontron UV Detector 332/335 <sup>8), 12)</sup>	n/a	supported by DataApex; CC 5.0.2	-	CG 3.1.7 Driver 446	n/a
Kontron UV Detector 430/430A/432/435 <sup>8), 12)</sup>	n/a	supported by DataApex; CC 5.0.2	-	CG 3.1.7 Driver 446	n/a
Kontron UV Detector 530/532/535 <sup>8), 12)</sup>	n/a	supported by DataApex; CC 5.0.2	-	CG 3.1.7 Driver 446	n/a
Kontron DAD 540(/545) <sup>12, 23)</sup>	n/a	supported by DataApex; CC 5.0.2	-	CG 3.1.7 Driver 1466	n/a
Sedere Sedex ELSD 80/85LT	supported, driver was developed by MayLab Analytical <sup>13)</sup>	supported by DataApex; CC 3.0.5	Cm 7.2: 85LT supported from FW 2.7, requires Class 3 license (driver not developed by KNAUER)	supported, driver was developed by MayLab Analytical <sup>14)</sup>	n/a
Sedere Sedex ELSD 90LT	supported, driver was developed by MayLab Analytical <sup>13)</sup>	supported by DataApex; CC 5.0.2	Cm 7.2: from FW 1.0, requires Class 3 license (driver not developed by KNAUER)	supported, driver was developed by MayLab Analytical <sup>14)</sup>	n/a
Sedere Sedex ELSD 100LT	supported, driver was developed by MayLab Analytical <sup>13)</sup>	supported by DataApex; CC 5.0.2	Cm 7.2: requires Class 3 license (driver not developed by KNAUER)	n/a	n/a

Detectors (continued)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Sedere Sedex LC	supported, driver was developed by MayLab Analytical <sup>13)</sup>	supported by DataApex; CC 5.0.2	Cm 7.2: from FW A.1.5, requires Class 3 license (driver not developed by KNAUER)	n/a	Version: 5.07.039
ERC RefractoMax 520 series	n/a	supported by DataApex; CC 6.1.0, RS-232	Cm 7.2 (driver not developed by KNAUER)	n/a	n/a
Antec Leyden LC-EC Detector Decade	n/a	supported by DataApex; CC 5.0.2	Cm 6.8 (driver not developed by KNAUER)	n/a	Version: 5.07.039
Antec Decade II	n/a	supported by DataApex; CC 5.x, RS-232	n/a	n/a	n/a
Antec Decade Elite	n/a	supported by DataApex; CC 5.X, USB, LAN	n/a	n/a	n/a
Shimadzu CDD-10Avp <sup>26)</sup>	n/a	n/a	n/a	n/a	n/a
Grace/Alltech ELSD 3300	-	supported by DataApex; CC 5.x, LAN	-	Separate, license-free driver upon request or from Grace webpage, driver not developed by KNAUER	-
Elysia-raytest radioactivity detectors and monitors Gabi Nova, Gabi *, POMO, Ramona *	-	only upon request to Elysia; no driver from KNAUER available	-	n/a	n/a

- 1) - only LAN is supported
- 2) - to set a static IP address a separate tool is required
- 3) - via RS-232 the maximum data rate is limited to 10 Hz
- 4) - only in ClarityChrom® 5.0.2: Runs with data rates >10 Hz may not be stored with 3D data; recommended maximum data rate is 10 Hz, fixed in later versions from 5.0.3
- 5) - generic drivers do not support all of the devices features; only control via RS-232, data acquisition requires A/D converter; should also work with Chromeleon™ 7.x



- 6) - a trigger on Start input also starts the internal wave length program
- 7) - for control via RS-232 only baud rate 115200 is supported
- 8) - device only supports data acquisition via RS-232, no wavelength control
- 9) - LAN only via separate LAN adapter A4442; default communication adapter is PCI interface card for RS-422, PCI interface card requires a PCI slot on the computer; driver for PCI interface card only works with Windows NT/2000, Windows XP or Windows 7 32-bit. Older ISA interface card is not supported anymore.
- 10) - RS-232 only, onboard COM port or PCI/PCIe extension cards are recommended, USBtoCOM adapter and other adapters may not work; maximum data rate 0.5 Hz, A/D converter for higher data rates required; special RS-232 cable required (KNAUER part number G0743)
- 11) - RS-232 only, onboard COM port or PCI/PCIe extension cards are recommended, USBtoCOM adapter and other adapters may not work; special firmware version 0.90 required, feature set of RF-10 AXL only (1 channel only, no flow cell temperature control, no scan), maximum data rate 5 Hz via RS-232, A/D converter for higher data rates required; special RS-232 cable required (KNAUER part number G0743)
- 12) - special RS-232 cable is required
- 13) - data rate limited to 10 Hz, driver must be installed separately and is protected with a separate USB dongle, detector cannot be used for fractionation with KNAUER FRC option; KNAUER part number A2618-01
- 14) - data rate limited to 10 Hz, driver must be installed separately and is protected with a separate USB dongle, detector cannot be used for fractionation with KNAUER FRC option; KNAUER part number A1513-07
- 15) - 4 channels, up to 50 Hz data rate
- 16) - 8 channels, up to 10 Hz; PurityChrom® offers only limited 3D functionality, refer to the PurityChrom® manual
- 17) - up to 2 Hz, not in a system with CM 2.1S (only RID 2.1L or CM2.1S in a system with PurityChrom® is allowed)
- 18) - not in a system with RID 2.1L (only CM2.1S or RID 2.1L in a system with PurityChrom® is allowed)
- 19) - requires separate license, KNAUER part number A2618-02, detector cannot be used for fractionation with KNAUER FRC option
- 20) - requires separate license, KNAUER part number A1783-3
- 21) - requires separate license, KNAUER part number A1783-4
- 22) - for full functionality the KNAUER RS-232 Y-cable is required (KNAUER article number A0755)
- 23) - for acquiring 3D data (spectra) a 3D UV (PDA) license option is required.
- 24) - not in a system with IFU 2.1/IFU 2.1LAN (only RF-20/Axs or IFU 2.1/IFU 2.1LAN in a system with PurityChrom® is allowed)
- 25) - extended dynamic range for RID 2.1L is not supported in PurityChrom
- 26) - The Shimadzu CDD-10Avp requires an A/D converter for data acquisition, device control is not possible

Assistant (up to three devices in one housing) <sup>1)</sup>	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
ASM 2.1L <sup>1), 2)</sup>	OLEE A.01.05 / A.04.05 Driver 1017 <sup>3)</sup>	CC 3.0.7; Driver 3078 <sup>4) 5) 6)</sup>	Cm 7.2 Regular Class 3 driver 1056 <sup>8)</sup>	CG 3.3.2; Driver 2825 <sup>6)</sup>	Version: 5.07.039 <sup>7) 9)</sup>
ASM 2.2L <sup>1), 2)</sup>	Coming soon	CC 8.2.2; Driver 6036 <sup>4) 6)</sup>	Coming soon	n/a	Version: 5.09.036 <sup>7) 9)</sup>

- 1) - configuration depends on order; refer the allowed configuration and restrictions in our FAQ (<https://www.knauer.net/en/Support/FAQ/Assistants>)
- 2) - only one LAN port for all devices
- 3) - only one valve as fraction collector, no cascading if controlled as ASM; all included modules also can be controlled as separate devices via IP port; FRC valve cascading only is possible if a valve in the ASM is configured as KNAUER MultiValve FC and is controlled as separate device via IP port; more information about control as separate device via IP port in software manual
- 4) - ASM is only supported as complete Assistant; included modules cannot be controlled as separate devices (IP port separation is not supported)
- 5) - two pumps can be controlled either as HPG or as separate pumps from hotfix 5043, in older versions only as HPG
- 6) - only one valve as fraction collector, no cascading
- 7) - ASM is not supported as device; the included modules can be controlled as separate devices via IP port; only one valve as fraction collector
- 8) - ASM is not supported as device; the included modules can be controlled as separate devices via IP port
- 9) - ASM's Event Outs and Start In are not supported

<b>Autosamplers</b> (one autosampler in a system)	<b>OpenLAB® CDS EZChrom Edition CDS version and Driver Version</b>	<b>ClarityChrom® CDS version and Driver Version</b>	<b>Chromleon™ CDS version Driver version</b>	<b>ChromGate® CDS version and Driver Version</b>	<b>PurityChrom®</b>
AZURA® Autosampler AS 6.1L 1), 2)	OLEE A.01.05 / A.04.05 Driver 2109	CC 7.2.0.5043 <sup>6)</sup>	Cm 7.2 Regular Class 3 driver 1038 <sup>6)</sup> (requires FW 1.17 or higher)	n/a	Version: 5.08.024 <sup>6)</sup>
Smartline Autosampler 3950 <sup>1), 2)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4 / LAN from CC 2.6.5 Driver 1083 <sup>6)</sup>	Cm 7.2 Regular Class 3 driver 1024 (requires FW 1.15 or higher) Cm 6.8 SR10 Regular Class 3 driver <sup>6)</sup>	CG 3.1.7 Driver 1466 / LAN from CG 3.3.1 Driver 1640 <sup>8)</sup>	Version: 5.07.067 <sup>6)</sup>
PLATINblue Autosampler AS-1 1), 2)	OLEE A.01.05 / A.04.05 Driver 1017	CC 6.1.0, Driver 5015	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.2 PLATINblue Edition 1751 <sup>8)</sup> ; Regular 2825 <sup>8)</sup>	n/a
Autosampler KNAUER Optimas	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3055 <sup>3)</sup>	n/a <sup>9)</sup>	n/a <sup>7)</sup>	Version: 5.08.034
Smartline Autosampler 3800 (also WellChrom K-3800, Spark Basic Marathon)	n/a	CC 2.4.4 <sup>3)</sup>	Cm 6.8; Generic driver <sup>4)</sup>	CG 3.3.2 Driver 1787	Version: 5.07.039
Smartline Autosampler 3900 (also Spark Midas)	n/a	CC 2.4.4 <sup>3)</sup>	Cm 6.8; Generic driver <sup>4)</sup>	CG 2.8	Version: 5.07.039
Spark Triathlon	n/a	CC 2.4.4 <sup>3)</sup>	Cm 6.8; Generic driver <sup>4)</sup>	CG 2.55	Version: 5.07.039
Kontron Autosampler 360	n/a	CC 5.0.2 <sup>3)</sup>	-	CG 3.1.7 Driver 1466 <sup>5)</sup>	n/a
Kontron Autosampler 46x	n/a	-	-	CG 3.1.7 Driver 1466 <sup>5)</sup>	n/a
Kontron Autosampler 56x	n/a	CC 5.0.2 <sup>3)</sup>	-	CG 3.1.7 Driver 1466 <sup>5)</sup>	n/a

- 1) - some device features have been added from time to time; therefore, not all of the current features are supported in the software from beginning as well as not all autosampler firmware versions support it:

#### **AS 6.1L**

- FW 1.21 - 1000 µl syringe supported (currently only supported in ClarityChrom), User Defined Program supported (currently not supported by drivers)
- FW 1.20 - no driver-dependent changes (extra needle dry step, other needle height during first transport aspiration)
- FW 1.19 - 100 µl syringe supported: OLEE: driver 2165, CC: 8.1 driver 6034 (hotfix only), CG: AS 6.1L not supported; Cm: n/a
- FW 1.18 - 30 vial tray, 2500 µl syringe can be selected also if prep mode is not enabled: OLEE: driver 2143, CC: 7.4.2 driver 5066, CG: AS 6.1L not supported; Cm: 7.x driver 1056
- Advanced needle wash - requires FW 1.18: OLEE: driver 2158; CC: 8.1 driver 6018; CG: AS 6.1L not supported; Cm: not supported

#### **AS3950**

- LAN connection (devices with LAN connector sold from 2009): OLEE: driver 1017; CC: 2.6.5 driver 1094; Cm: 7.x driver 1024; CG: 3.3.1 driver 1640; PC: 5.07.067
- FW 1.04 - 84+3 vial tray (Aug 2011): OLEE: driver 1017; CC: 3.0.7 driver 3055; Cm: 7.x driver 1024; CG: 3.3.2 driver 1787; PC: 5.07.067
- FW 1.15 - 0.1 µl increments, tray thermostatted up to 40 °C (2014): OLEE: A.04.05 driver 1028; CC: 5.0.3 driver 4102, Cm: 7.x, driver 1024; CG: 3.3.2 driver 2846 (40 °C only); PC: 5.07.067
- TCP/IP support (is required to control more than one autosampler in a network, devices sold from Q3 2016 with TCP by default, former devices can be switched to UDP): OLEE: driver 1017; CC: 6.1.0 driver 5028 (hotfix); Cm: 7.x driver 1024; CG: n/a (not supported); PC: 5.07.067 (UDP not supported)
- FW 1.16 - 108 vial tray (Q1 2017); drivers set also 250 µl syringe and 108 vial tray as default in new method: OLEE: driver 2109; CC: 7.2 driver 5045; Cm: 7.x driver 1038; CG: n/a; PC: 05.08.034
- FW 1.18 - 2500 µl syringe selectable also if prep mode is not activate, 30 vials tray: OLEE: driver 4.2.0.2143; CC: 7.4, driver 5066; Cm: 7.x driver 1056; CG: n/a; PC: 05.09.028 - 2500 µl syringe; 5.09.069 - 30 vial tray
- 30 vial tray selectable in prep mode: OLEE: 4.2.0.2159; CC: 8.0.1.6012; Cm: 7.x driver 1056; CG: n/a; PC: 5.09.069

#### **AS-1**

- FW 1.02 - Tray thermostatted up to 40 °C: OLEE: from first version; CC: 6.1.0 driver 5015; Cm: 6.8 only; CG: 3.3.2 driver 2846; PC: AS-1 not supported

- 2) - LAN default configuration is LAN via DHCP, static IP address requires access to built-in LAN module via browser or ASM software
- 3) - driver was developed by DataApex
- 4) - upon request - driver was not developed by KNAUER
- 5) - special RS-232 cable is required
- 6) - no Stacked Injection supported
- 7) - may run with Midas/AS3900 driver
- 8) - in ChromGate TCP/IP communication is not supported, device must be set to UDP; later added features (see <sup>1)</sup>) are not supported
- 9) - generic driver for autosampler Spark Midas may work; upon request - driver was not developed by KNAUER

Analog/Digital Converters	OpenLAB® CDS EZChrom Edition CDS version and Driver Version <sup>1)</sup>	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
KNAUER Interface Box IFU 2.1 LAN <sup>2)</sup>	OLEE A.01.05 / A.04.05 Driver 2092	CC 6.0.1 Driver 5043 (compatible with CC/CCP 5.x versions)	n/a	n/a	Version: 5.08.024 <sup>6)</sup>
KNAUER Interface Box IFU 2.1 <sup>3), 4), 5)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3078, in former version as IF2 <sup>5)</sup>	n/a	CG 3.3.2 Driver 2825, in former version as IF2 <sup>5)</sup>	Version: 5.07.039 <sup>6)</sup>
KNAUER Interface Box IF2 <sup>4)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4 / 2.6 <sup>7), 8)</sup>	n/a	CG 3.1 <sup>9)</sup>	Version: 5.07.039 <sup>6)</sup>
Smartline Manager 5000 <sup>4)</sup>	n/a	CC 2.4.4/2.6 <sup>7), 8)</sup>	n/a	CG 3.1 <sup>9)</sup>	n/a
Smartline Manager 5050 <sup>4)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.6.5 Driver 1100 <sup>5)</sup>	n/a	CG 3.3.2 Driver 1751 <sup>5)</sup>	n/a
KNAUER HPLC Box Model 96 <sup>10)</sup>	n/a	CC 2.6 <sup>11)</sup>	n/a	CG 2.55 <sup>12), 13)</sup>	n/a
Kontron Interface Card <sup>14)</sup>	n/a	-	-	CG 3.1.7 Driver 1466	n/a

- 1) - A/D converter in OpenLAB® must be configured separately; menu is available via *Start - Programs - Agilent Technologies - OpenLAB® CDS EZChrom Edition*, right click on *Agilent SS420X Interface Configuration* and select *Run as Administrator*
- 2) - to set a static IP address a 3<sup>rd</sup> party tool is required
- 3) - built-in USB-to-serial converter requires separate installation of a driver in Windows operating system, can be found on the HPLC software CD or KNAUER webpage (*Download* section for IFU 2.1, Article Number AZB00); check Windows Device Manager for the COM port number of the converter, this number must be selected in the software configuration
- 4) - as the required serial number in the configuration, enter the last 5 digits of the serial number on the label; if the first digit is a 0 (zero), only enter the last four digits
- 5) - in previous software versions IFU 2.1 can be controlled as IF2 or Manager 5000 if IFU 2.1 is not listed in the module list
- 6) - analog input only
- 7) - in ClarityChrom® version 2.4.4 only data acquisition via analog input is supported (max. 10 Hz), analog out supported from ClarityChrom® version 2.6
- 8) - 50 Hz option for channel 1 supported from version ClarityChrom® version 2.6
- 9) - 50 Hz option for channel 1 supported from version ChromGate® version 3.1.4

- 10) - as the required serial number in the configuration, enter the last 5 digits of the serial number on the label, if the serial number is higher than 40000, please enter only the last four digits
- 11) - supported as "Basic Edition", ClarityChrom® will switch the box into the required mode
- 12) - supported as "Integration Box" only for data acquisition via analog input; box must be switched into the correct mode with special software
- 13) - 100 Hz option for channel 1 supported from version ChromGate® version 3.1
- 14) - requires an ISA slot on the computer; is not supported from Windows Vista (and higher)

Thermostats/Column Ovens	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Column Thermostat CT 2.1	OLEE A.01.05 / A.04.05 Driver 1027	CC 5.0.3 driver 4102	Cm 7.2 Regular Class 3 driver 1024	n/a	Version 5.09.069 <sup>4)</sup> , <sup>5)</sup>
Support of modified CT 2.1 with S/N FCC...	Coming soon	CC 8.2.3 driver 6042	Coming soon	n/a	Coming soon
PLATINblue Column Thermostat T-1	OLEE A.01.05 / A.04.05 Driver 1017	CC 5.0.2 Driver 4095	Cm 6.8 SR10 Regular Class 3 driver	CG 3.3.2 PLATINblue Edition 1751; Regular 2825	n/a
Smartline Column Oven <sup>2)</sup> / Column Thermostat JetStream <sup>1)</sup> , <sup>2)</sup> / PLATINblue T-1 Basic <sup>2)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4	Cm 6.8 SR10 Regular Class 3 driver	CG 3.1	Version: 5.07.039
Smartline Column Oven 4050 <sup>3)</sup>	n/a	CC 2.4.4	n/a	CG 3.1.7 Driver 1466	n/a
Memmert UN 55 Oven	n/a	n/a	n/a	n/a	Version 5.09.069
Eluent heating device <sup>6)</sup>	n/a	n/a	n/a	n/a	Version: 5.08.035 <sup>4)</sup>

- 1) - firmware version b1.17 or higher for JetStream required
- 2) - if connected, control panel must show "Choose Function" before the instrument will be opened
- 3) - device does not support static IP address, only LAN via DHCP from firmware version 1.11
- 4) - only CT 2.1 or eluent heater is supported in one system
- 5) - leak detector is not supported
- 6) - separate software TempLogger for Standalone operation available

Electrical Valve Drives	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
AZURA® Valve Unifier VU 4.1 <sup>1)</sup>	OLEE A.01.05 / A.04.05 Driver 2138	CC7.4.1 Driver 5066	Cm 7.2 Regular Class 3 driver 1056 <sup>2), 3)</sup>	n/a <sup>1)</sup>	Version: 5.07.039
AZURA® Valve V 2.1S <sup>1)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 3.0.7 Driver 3078 <sup>1)</sup>	Cm 7.2 Regular Class 3 driver 1056 <sup>2), 3)</sup>	n/a <sup>1)</sup>	Version: 5.07.039
WellChrom electric valve drives K-6, K-12, K-16	n/a <sup>4)</sup>	CC 2.4.4	CM 6.8 Generic driver <sup>3)</sup>	CG 2.8	Version: 5.07.039
Smartline valve drives S6, S12, S16 <sup>1)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	CC 2.4.4	CM 6.8 Generic driver <sup>3)</sup>	CG 3.1 <sup>1)</sup>	Version: 5.07.039
VICI Valco drives	OLEE A.01.05 / A.04.05 Driver 1017	CC 5.0.2 <sup>5)</sup>	Implemented (not by KNAUER)	Generic Driver <sup>3)</sup>	Version: 5.07.039

- 1) - can be controlled as KNAUER switching valve if not listed in the module list of previous software versions (limited on RS-232 and only valve types available for V 2.1S)
- 2) - generic driver for KNAUER switching valves will also work with AZURA® VU4.1 (only for valve types available for V 2.1S, RS-232 only) and V 2.1S (RS-232 only), Cm 6.8 and 7.2
- 3) - generic drivers do not support all of the devices features (e.g. RS-232 only, limited number of valve types)
- 4) - can be controlled as KNAUER switching valve
- 5) - driver was developed by DataApex

Other Supported Instruments	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Advion expressiON-S CMS MS <sup>1), 2)</sup>	n/a	ClarityChrom 5.0.3	n/a	n/a	n/a
Microsaic 4000 MiD Single Quad MS <sup>1), 2)</sup>	n/a	ClarityChrom 7.4.1	-	n/a	5.09.069
Airsensor <sup>3)</sup>	n/a	n/a	n/a	n/a	Version: 5.07.039

Other Supported Instruments (continued)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version	Chromleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Pressure control (pressure difference) <sup>4)</sup>	n/a	n/a	n/a	n/a	Version: 5.07.039
Flowmeter GJC, Model 5025000 230V	OLEE A.01.05 / A.04.05 Driver 1017	n/a	n/a	supported	n/a
Flowmeter Bronkhorst Mini Cori-Flow M12, M13, M14	OLEE A.01.05 / A.04.05 Driver 1017	CC 8.1 driver 6018	n/a	supported	Version: 5.09.028

- 1) - driver was developed by DataApex
- 2) - device control requires MS license
- 3) - communication via game port adapter
- 4) - requires KNAUER IFU 2.1/ IFU 2.1 LAN

Fraction Collectors (one fraction collector in a system, no combination of a fraction valve and a fraction collector possible)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version *)	ClarityChrom® CDS version and Driver Version **) <sup>17)</sup>	Chromleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
KNAUER MultiValve Fraction Collector (all KNAUER electric valve drives and VU 4.1 <sup>2)</sup> )	OLEE A.01.05 / A.04.05 Driver 1017 <sup>1)</sup>	CCP 2.6 <sup>1)</sup>	n/a	CG 3.1.6 Driver 423 <sup>1)</sup>	Version: 5.07.039 <sup>3)</sup>
ASM 2.2L valve	Coming soon	CC 8.2.2 6036 <sup>14)</sup>	n/a	n/a	Version: 5.09.036 <sup>16)</sup>
ASM 2.1L valve	OLEE A.01.05 / A.04.05 Driver 1017 <sup>13)</sup>	CCP 3.0.7 Driver 3078 <sup>14)</sup>	n/a	CG 3.3.2; Driver 2825 <sup>15)</sup>	Version: 5.07.039 <sup>16)</sup>
Smartline Fraction Collector 3050 <sup>4)</sup>	n/a	CCP 2.6	n/a	CG 3.3.1 Driver 1640	n/a
Buechi Fraction Collector C-684 <sup>5)</sup>	n/a	n/a	n/a	CG 3.1.6 Driver 423	n/a



<b>Fraction Collectors continued</b> (one fraction collector in a system, no combination of a fraction valve and a fraction collector possible)	<b>OpenLAB® CDS EZChrom Edition CDS version and Driver Version</b>	<b>ClarityChrom® CDS version and Driver Version <sup>17)</sup></b>	<b>Chromeleon™ CDS version Driver version</b>	<b>ChromGate® CDS version and Driver Version</b>	<b>PurityChrom®</b>
Buechi Fraction Collector C-660	OLEE A.01.05 / A.04.05 Driver 1017	CCP 2.6	n/a	CG 3.1.6 Driver 423	n/a
ISCO Fraction Collector Foxy Jr.	n/a	CCP 2.6	Implemented (not by KNAUER)	CG 3.1.6 Driver 423	Version: 5.07.039
ISCO Fraction Collector Foxy R1 <sup>6)</sup>	OLEE A.01.05 / A.04.05 Driver 1017 <sup>7)</sup>	CCP 2.6.5 Driver 1094 (RS-232 only) <sup>8)</sup> CC 7.4.2 hotfix 5075 **) also LAN	Implemented (not by KNAUER)	CG3.3.2 Driver 1751 <sup>7)</sup>	Version: 5.07.039
ISCO Fraction Collector Foxy R2 <sup>6)</sup>	OLEE A.01.05 / A.04.05 Driver 1017 <sup>7)</sup>	CCP 2.6.5 Driver 1094 (RS-232 only) <sup>8)</sup> CC 7.4.2 hotfix 5075 **) also LAN	CM 6.8 Generic driver <sup>20)</sup> (not by KNAUER)	CG3.3.2 Driver 1751 <sup>7)</sup>	Version: 5.07.039
ISCO Fraction Collector Foxy 200 <sup>9)</sup>	n/a	CCP 2.6 and CCP 2.6.5 Driver 1100 <sup>9)</sup>	Implemented	CG 3.1.7 Driver 1466 - CG3.3.2 Driver 1751 <sup>9)</sup>	n/a
Labocol Vario 2000 Fraction Collector	OLEE A.01.05 / A.04.05 Driver 1017	n/a	n/a	CG 3.1.6 Driver 423	Version: 5.07.039
Labocol Vario 4000 Fraction Collector <sup>6)</sup>	OLEE A.01.05 / A.04.05 Driver 1017 (RS-232 only) OLEE A.02.00 / A.04.06 Driver 2155 *) <sup>18)</sup> also LAN	CC 7.4.2 Driver 5075 (hotfix) <sup>19)</sup>	-	CG 3.3.1 Driver 1640 (RS-232 only)	Version: 5.07.039 (RS-232 and LAN control, free rack definition supported, mixed rack configuration supported, no USB control)

Fraction Collectors continued (one fraction collector in a system, no combination of a fraction valve and a fraction collector possible)	OpenLAB® CDS EZChrom Edition CDS version and Driver Version	ClarityChrom® CDS version and Driver Version <sup>17)</sup>	Chromeleon™ CDS version Driver version	ChromGate® CDS version and Driver Version	PurityChrom®
Labocol Vario 4000 Plus Fraction Collector <sup>6)</sup>	OLEE A.02.00 / A.04.06 Driver 2155 *) <sup>18)</sup>	CC 7.4.2 Driver 5075 (hotfix) <sup>19)</sup>	-	n/a	Version: 5.07.039 (RS-232 and LAN control, free rack definition supported, mixed rack configuration supported, no USB control)
VICI Valco MultiValve Fraction Collector	OLEE A.01.05 / A.04.05 Driver 1017 <sup>1)</sup>	n/a from KNAUER, in ASM 2.1L as FC <sup>3)</sup> , <sup>10)</sup> valve as FC is supported by DataApex <sup>11)</sup>	n/a	CG 3.3.1 Driver 1640 <sup>1)</sup>	Version: 5.07.039 <sup>3)</sup>
Virtual Fraction Collector	OLEE A.01.05 / A.04.05 Driver 1017	CCP 2.6	n/a	CG 3.1.6 Driver 423	n/a
Injection Module (combined module with pump and injection valve for injection of large injection volumes, allows for stacked injection) <sup>12)</sup>	OLEE A.01.05 / A.04.05 Driver 1017	n/a	n/a	CG 3.3.1 Driver 1640	n/a

\*) - Due to changes in OpenLAB EE by Agilent, the KNAUER FRC AddOn needed to be adjusted and is not compatible to OpenLAB EE A.01.05 / A.04.05 from driver package version 2146. All later added functionality and drivers are not available for OpenLAB EE A.01.05 / A.04.05.

\*\*) - Up to version 6.1, all KNAUER preparative functionality and fraction collector drivers are only available in ClarityChrom *Prep*. From version 7.2., the ClarityChrom *Prep* as a separate software is discontinued. The full functionality now is available in ClarityChrom with the special Knauer FRC control. If an upgrade of ClarityChrom *Prep* is required, a ClarityChrom with Knauer FRC control will be delivered (new User Code).

- 1) - up to 32 valves in a cascade
- 2) - VU 4.1 as MultiValve FC fully supported from CC 7.4.1.5066, in previous versions as V 2.1S RS-232 only
- 3) - only one valve can be used as fraction collector, no cascading is supported
- 4) - switch between RS-232 and LAN with a switch on the devices rear panel
- 5) - requires special relay module for serial communication

- 6) - touch panel is not blocked in remote, do not touch it while the device is under software control
- 7) - LAN control from firmware version 1.01.03, LAN control support in firmware versions 1.01.09 and 1.01.11 faulty
- 8) - from CCP 2.6.5 1094 only RS-232 support, from 7.4.2 hotfix 5075 also LAN control
- 9) - support is discontinued due to several incompatible firmware versions; not all devices may work with the driver
- 10) - not all available valve types can be built-in into ASM 2.1L
- 11) - driver from DataApex available; does not support KNAUER fractionation features
- 12) - several pumps were added to the module when released; however, a pump is only supported in the injection module if the mentioned pump is also available as a separate pump; in ASM2.1L pump and valve must be controlled separately, refer software manual for detailed information
- 13) - only one valve as fraction collector, no cascading if controlled as ASM; FRC valve cascading only is possible if a valve in the ASM is configured as KNAUER MultiValve FC and is controlled as separate device via IP port; more information about control as separate device via IP port in software manual
- 14) - only one valve as fraction collector in ASM, no cascading
- 15) - only one valve as fraction collector, no cascading if controlled as ASM; FRC valve cascading only is possible if a valve in the ASM is configured as KNAUER MultiValve FC and is controlled as separate device via IP port
- 16) - only one valve as fraction collector, no cascading
- 17) - DataApex supports fraction collection with an own implementation; it is available in ClarityChrom® as well as ClarityChrom® *Prep*. KNAUER FRC option is available in ClarityChrom® *Prep* and, from version 7.4.1, also in ClarityChrom® with Knauer FRC license
- 18) - from driver version 2155, additionally model 4000 Plus, LAN control and trays with 24 and 39 vials are supported, now trays with 20, 24, 39, 80 and 125 vials are available, free tray definition is not supported; all racks must have same type (no mixed configuration supported); no USB control support; former drivers only support model 4000, RS-232 control and trays with 20, 80 and 125 vials
- 19) - only trays with 20, 24, 39, 80 and 125 vials are supported; free tray definition is not supported; all racks must have same type (no mixed rack configuration supported); no USB control support
- 20) - generic driver should also work with Cm 7, generic drivers do not support all of the devices features (e.g. RS-232 only)

## Mobile Control - KNAUER

Mobile Control requires minimum device firmware versions, which may differ between the Mobile Control versions. Refer to the list of supported devices and their required firmware version below. A Mobile Control update may also require firmware updates of the controlled devices. Mobile Control download comes with a Firmware Wizard, which allows for firmware updates. For **Monitor Mode, supported in OpenLAB® EE and ClarityChrom® only**, every Mobile Control version requires a minimum KNAUER driver version for the respective CDS.

**Mobile Control, Version 5.5.36** also with data acquisition support, supports AZURA® devices:

- |  |  |  |
|--|--|--|
| - Pump P 2.1S/P 4.1S – Firmware V01.38 | - Detector DAD 2.1L* – Firmware V01.10 | - Assistant ASM 2.2L – Firmware V01.10       |
| - Pump P 6.1L – Firmware V01.06        | - Detector DAD 6.1L* – Firmware V01.23 | - Column thermostat CT 2.1 – Firmware V01.06 |
| - Pump P 2.1L – Firmware V01.09        | - Detector RID 2.1L – Firmware V01.24  | - Autosampler AS 3950 – Firmware V01.17      |
| - Detector UVD2.1S – Firmware V01.11   | - Valve V 2.1S – Firmware V05.01       | - Autosampler AS 6.1L – Firmware V01.17      |
| - Detector UVD 2.1L – Firmware V02.06  | - VU 4.1 – Firmware V05.68             | - IFU 2.1 LAN – Firmware V1.05               |
| - Detector MWD 2.1L – Firmware V01.10  | - Assistant ASM 2.1L – Firmware V01.17 |  |

\* - no 3D data acquisition, 4 channels only

Mobile Control version 5.5.36 supports Monitor Mode in OpenLAB® EE from driver version 2155 or in ClarityChrom® from version 8.1.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

**Mobile Control, Version 5.5.5** also with data acquisition support, supports AZURA® devices:

- |  |  |  |
|--|--|--|
| - Pump P 2.1S/P 4.1S – Firmware V01.37 | - Detector DAD 2.1L – Firmware V01.10  | - Assistant ASM 2.2L – Firmware V01.08       |
| - Pump P 6.1L – Firmware V01.05        | - Detector DAD 6.1L – Firmware V01.23  | - Column thermostat CT 2.1 – Firmware V01.06 |
| - Pump P 2.1L – Firmware V01.09        | - Detector RID 2.1L – Firmware V01.19  | - Autosampler AS 3950 – Firmware V01.17      |
| - Detector UVD2.1S – Firmware V01.11   | - Valve V 2.1S – Firmware V5.01        | - Autosampler AS 6.1L – Firmware V01.17      |
| - Detector UVD 2.1L – Firmware V02.06  | - VU 4.1 – Firmware V5.68              | - IFU 2.1 LAN – Firmware V1.05               |
| - Detector MWD 2.1L – Firmware V01.10  | - Assistant ASM 2.1L – Firmware V01.13 |  |

Mobile Control version 5.5.5 supports Monitor Mode in OpenLAB® EE from driver version 2155 or in ClarityChrom® from version 8.1.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

**Mobile Control, Version 5.5.0** also with data acquisition support, supports AZURA® devices:

- |  |  |  |
|--|--|--|
| - Pump P 2.1S/P 4.1S – Firmware V01.37 | - Detector DAD 2.1L – Firmware V01.10  | - Assistant ASM 2.2L – Firmware V01.08       |
| - Pump P 6.1L – Firmware V01.05        | - Detector DAD 6.1L – Firmware V01.23  | - Column thermostat CT 2.1 – Firmware V01.06 |
| - Pump P 2.1L – Firmware V01.09        | - Detector RID 2.1L – Firmware V01.19  | - Autosampler AS 3950 – Firmware V01.17      |
| - Detector UVD2.1S – Firmware V01.11   | - Valve V 2.1S – Firmware V5.01        | - Autosampler AS 6.1L – Firmware V01.17      |
| - Detector UVD 2.1L – Firmware V02.06  | - VU 4.1 – Firmware V5.68              | - IFU 2.1 LAN – Firmware V1.05               |
| - Detector MWD 2.1L – Firmware V01.10  | - Assistant ASM 2.1L – Firmware V01.13 |  |

Mobile Control version 5.5.0 supports Monitor Mode in OpenLAB® EE from driver version 2155 or in ClarityChrom® from version 8.1.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

**Mobile Control, Version 5.3.0** also with data acquisition support, supports AZURA® devices:

- |  |                                       |  |
|--|---------------------------------------|--|
| - Pump P 2.1S/P 4.1S – Firmware V01.37 | - Detector MWD 2.1L – Firmware V01.10 | - Assistant ASM 2.1L – Firmware V01.13       |
| - Pump P 6.1L – Firmware V01.05        | - Detector DAD 2.1L – Firmware V01.10 | - Column thermostat CT 2.1 – Firmware V01.06 |
| - Pump P 2.1L – Firmware V01.09        | - Detector DAD 6.1L – Firmware V01.23 | - Autosampler AS 3950 – Firmware V01.17      |
| - Detector UVD2.1S – Firmware V01.11   | - Detector RID 2.1L – Firmware V01.15 | - Autosampler AS 6.1L – Firmware V01.17      |
| - Detector UVD 2.1L – Firmware V02.06  | - Valve V 2.1S – Firmware V5.01       | - IFU 2.1 LAN – Firmware V1.05               |

Mobile Control version 5.3.0 supports Monitor Mode in OpenLAB® EE from driver version 2126 or in ClarityChrom® /Prep (5.0.5 or higher) from driver version 5053.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

**Mobile Control, Version 4.0.0** also with data acquisition support, supports AZURA® devices:

- |  |                                       |  |
|--|---------------------------------------|--|
| - Pump P 2.1S/P 4.1S – Firmware V01.37 | - Detector MWD 2.1L – Firmware V01.02 | - Assistant ASM 2.1L – Firmware V01.12       |
| - Pump P 6.1L – Firmware V01.03        | - Detector DAD 2.1L – Firmware V01.02 | - Column thermostat CT 2.1 – Firmware V01.06 |
| - Pump P 2.1L – Firmware V01.08        | - Detector DAD 6.1L – Firmware V01.23 | - Autosampler AS 3950 – Firmware V01.16      |
| - Detector UVD2.1S – Firmware V01.10   | - Detector RID 2.1L – Firmware V01.10 |  |
| - Detector UVD 2.1L – Firmware V02.05  | - Valve V 2.1S – Firmware V4.14       |  |

Mobile Control version 4.0.0 supports Monitor Mode in OpenLAB® EE from driver version 2092 or in ClarityChrom® /Prep (5.0.5 or higher) from driver version 5037.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

**Mobile Control, Version 3.0.1** also with data acquisition support, supports AZURA® devices:

- Pump 6.1L (from firmware version 1.03),
- Pump 2.1L (from firmware version 1.08),
- Pump P 2.1S/P 4.1S (from firmware version 1.32),
- UVD 2.1L (from firmware version 2.04),
- DAD 2.1L (from firmware version 1.1),
- DAD 6.1L (from firmware version 1.21),
- MWD 2.1L (from firmware version 1.1),
- UVD 2.1S (from firmware version 1.08),
- ASM 2.1L (from firmware version 1.11),
- Valve V 2.1S (from firmware version 4.14),
- CT 2.1 (from firmware version 1.06),
- AS 3950 (from firmware version 0840.133-115)

Mobile Control version 3.0.1 supports Monitor Mode in OpenLAB® EE from driver version 2087 or in ClarityChrom® /Prep (5.0.5 or higher) from driver version 5029.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

**Mobile Control, Version 2.2.0** supports AZURA® devices:

- Pump 6.1L (from firmware version 1.02),
- Pump 2.1L (from firmware version 1.07),
- Pump P 2.1S/P 4.1S (from firmware version 1.32),
- UVD 2.1L (from firmware version 2.04),
- DAD 2.1L (from firmware version 1.1),
- DAD 6.1L (from firmware version 1.21),
- MWD 2.1L (from firmware version 1.1),
- UVD 2.1S (from firmware version 1.08),
- ASM 2.1L (from firmware version 1.10),
- Valve V 2.1S (from firmware version 4.14)
- CT 2.1 (from firmware version 1.06)

Mobile Control version 2.2.0 supports Monitor Mode in OpenLAB® EE from driver version 2047 or in ClarityChrom®/Prep (5.0.5 or higher) from driver version 5015.

**Monitor Mode is not supported in any other CDS, e.g. PurityChrom® and Chromeleon™. Either CDS or Mobile Control can access the devices.**

## **EuroChrom for Windows 3.05 (KNAUER)**

Note that EuroChrom was discontinued in 2008. It will not run on Windows systems newer than Windows XP. Fraction collectors are supported exclusively in EuroChrom Prep software.

### Pumps

Smartline Pump 1000, Preparative Pump 1800, WellChrom K-1800, WellChrom K-1001, K-1000, K-501, K-120 (Smartline 100 as K-120 only via RS-232 without pressure sensor),

*Feed Delivery* (ECW Prep only): pump K-501, K-120

### Detectors

Smartline UV Detector 2500, Smartline RI Detectors 2300 / 2400, Smartline DAD 2800 (PCI/LAN), WellChrom UV Detectors K-2000/2001, K-2500/2501, K-2600, K-200, WellChrom RI Detectors K-2301/2401, WellChrom DAD K-2700 (PCI/LAN), WellChrom DAD K-2800 (PCI/LAN), Virtual Detector (WellChrom detectors 2xx0 allow only data acquisition, no control)

### Autosamplers

AS 3800 / Spark Basic Marathon, AS 3900 / Spark Midas, Spark Triathlon

### Analog/Digital Converters

KNAUER IF2 (also IFU 2.1 with built-in USBtoCOM converter), Smartline Manager 5000, KNAUER HPLC Box (Model 96)

### Thermostats/Column Ovens

JetStream (FW b1.17 or later required)

### Valves

KNAUER WellChrom K-6 / K-12 / K-16 [Smartline S 6 / S 12 / S 16 and AZURA® V2.1S and VU 4.1 (only 2, 6, 12, 16 port valves) are compatible via RS-232], Solvent Recycling: K-6 / K16

### Fraction Collectors (ECW Prep 3.05 only)

Labocol Vario 2000, ISCO Foxy Jr, KNAUER WellChrom K-6 / K-12 / K-16 (also cascade up to 8 valves), Buechi 684 via short circuits, User Defined FC (via short circuits)

## Other Software

**LabVIEW** (National Instruments) - support is discontinued, may not work with LabVIEW versions higher than 2011

Driver for AZURA® Pumps P 2.1S / P 4.1S, Smartline 100, BlueShadow Pumps 10P / 20P

### **HyStar™** (Bruker)

Compass HyStar™ 3.2 Plug-in: Drivers for PLATINblue system (Pump P-1, DAD PDA-1, UV Detector MW-1, Column Thermostat T-1, Autosampler AS-1, Manager M-1 (LPG valves only, interface module is not supported))  
Compass HyStar™ 4 Plug-in: Driver for AZURA® DAD 6.1L

### **Xcalibur™** (Thermo)

Version 2.0.7 and 2.1: Driver for PLATINblue Pump P-1, UV Detector MW-1, DAD PDA-1, Autosampler AS-1, Column Thermostat T-1 and T-1 Basic and Smartline Pump S 1050, UV Detector S 2550, Autosampler S 3950 and Column Thermostat JetStream

## TRADEMARKS

The company, product and service names used in this brochure are for identification purposes only.

All trademarks and registered trademarks are the property of their respective owners.

© KNAUER  
All rights reserved.

Version 1.17  
Last update: 2020/10/02

KNAUER Wissenschaftliche Geräte GmbH  
Hegauer Weg 38  
14163 Berlin, Germany

Tel.: +49-(0)30-809727-0  
Fax: +49-(0)30-8015010  
E-mail: [info@knauer.net](mailto:info@knauer.net)  
Internet: [www.knauer.net](http://www.knauer.net)

-----  
This document is subject to technical changes.  
Please go to [www.knauer.net](http://www.knauer.net) for latest version,  
ordering details and overview of article numbers.