


Software

▶ ChromGate® Chromatographie-Software

V7057-8



CONTENTS

Installation Guide ChromGate® 3.3.2	8
General Definitions	8
Installation	8
KNAUER FRC control option	11
Short Guide ChromGate®, KNAUER Instrument Control	13
Overview	13
Configuring a System	14
Instrument Configuration	14
System Options	14
Setting up a Method	15
The Instrument Wizard	15
The Method Window	16
Quantification of the Compound(s) of Interest	19
Overview	19
Integration	19
Generating and Using a Calibration Curve	20
Method File Commands	20
Quantification and Using Sequence Files	20
Overview	20
Single Level Calibration	21
Creating a Sample Sequence	23
Running a Sequence	25
Reporting	26
Collecting Data	28
General Data Collection Instructions with ChromGate	28
Collecting Data	28
Instrument Status of a running Control Method	29
Shutting Down the System	30
Setup and Control of Knauer HPLC Systems	31
Overview of Instrument Control	31
Configuration – device communication port	31
Configuring the Interface	31
Knauer Interface Configuration	32
Kontron Interface Configuration	34
K-2700 / 2800 Interface Card and Driver Installation	35
Configuring the KNAUER HPLC System	35
Instrument Configuration	36
Configuration – KNAUER HPLC System	36
Configuration – Knauer Pumps	38
Configuration – Kontron Pumps	44
Configuration – Knauer Detectors	45
Configuration – User Defined Detectors	53
Configuration – Kontron UV-Detectors	55
Configuration – Kontron Diode Array Detectors	56
Configuration – Virtual Detector	57
Detector Connections	57
Configuration – Assistant ASM2.1L	60
Configuration – Autosampler	63
Configuration – Autosampler 3950	64

Configuration – Autosampler Knauer Optimas.....	67
Configuration – Autosampler 3800.....	70
Configuration – Autosampler 3900.....	71
Configuration – Autosampler Triathlon / Endurance.....	72
Configuration – Kontron Autosamplers.....	76
Configuration – Miscellaneous Instruments.....	78
Configuration – Switching Valves.....	78
Configuration – Manager 5000/5050, IF2.....	79
Configuration – Column Oven 4050, Column Oven Jetstream.....	80
Configuration – Flowmeter.....	81
Creating an Instrument Control Method.....	83
Instrument Setup.....	84
Instrument Setup – Pumps.....	85
Instrument Setup – Detectors.....	91
All Detectors.....	91
Instrument Setup – RI Detectors.....	92
Instrument Setup – UV Detectors (S 200, K-200).....	92
Instrument Setup – UV Detectors (K-2000, K-2500).....	92
Instrument Setup – UV Detectors (UVD 2.1L, UVD 2.1S, S 2520, S 2500, K-2001, K-2501).....	93
Instrument Setup – Kontron Detectors.....	95
Instrument Setup – K-2600 Detector.....	95
Instrument Setup – S 2550 Detector.....	98
Instrument Setup – Diode Array Detectors (S 2600, DAD 2850, DAD 2800, K-2700 and Kontron DAD 540/545) ...	100
Instrument Setup – Fluorescence Detector RF-10Axl / RF-20A.....	103
Instrument Setup – Alltech 650 Conductivity Detector.....	104
Instrument Setup – User defined Detector.....	104
Instrument Setup – Virtual Detector.....	105
Instrument Setup – Assistant ASM2.1L.....	107
Instrument Setup – Autosamplers.....	111
Instrument Setup – Autosampler 3800.....	111
Instrument Setup – Autosampler Knauer Optimas / 3900 (Midas).....	112
Instrument Setup – Autosampler 3950 (Alias).....	116
Instrument Setup – Triathlon/Endurance Autosampler.....	121
Instrument Setup – Kontron Autosamplers.....	128
Instrument Setup – Miscellaneous Instruments.....	130
Instrument Setup – Manager 5000/5050/IF2 I/O.....	130
Instrument Setup – Switching Valves.....	131
Instrument Setup – Column Oven 4050.....	133
Instrument Setup – Column Oven Jetstream.....	133
Instrument Setup – Flowmeter.....	134
Setting up Auxiliary Traces.....	135
Setting up a Trigger.....	135
Setting up the Baseline Check.....	136
Instrument Status of a (running) Control Method.....	137
System Status.....	137
Instrument Status – Pumps.....	139
Instrument Status – Detectors.....	145
Instrument Status – RI Detectors (S 23[4]00, K-23[4]00/1).....	145
Instrument Status – UV Detectors (S 2520, S 2500, S 200, K-200, K- 2000/1, K-2500/1).....	145
Instrument Status – Kontron Detectors (3xx, 4xx, and 5xx).....	147
Instrument Status – fast scanning UV Detector (K-2600).....	148
Instrument Status – Diode Array Detectors.....	149
Instrument Status – Fluorescence Detector RF-10Axl / RF-20A.....	152
Instrument Status – Conductivity Detector Alltech 650.....	154
Instrument Status – User Defined Detector.....	154

Instrument Status – Virtual Detector	155
Instrument Status – Assistant ASM2.1L	156
Instrument Status – Autosampler	161
Instrument Status – Autosampler 3800	161
Instrument Status – Autosampler Optimas/3900	162
Instrument Status – Autosampler 3950	163
Instrument Status – Triathlon/Endurance Autosampler	165
Instrument Status – Kontron Autosamplers	166
Instrument Status – Miscellaneous Instruments	166
Instrument Status – Manager 5000/5050/IF2 I/O	166
Instrument Status – Knauer Switching Valves	167
Instrument Status – Column Oven 4050 and Jetstream	168
Instrument Status – Flowmeter	169
Knauer Instrument Control Method Options.....	170
General Settings	170
Runtime Settings	170
Download Tab / Method.....	172
Solvent Control	173
Qualification Procedures.....	175
Validation of Integration	178
Generic Drivers	183
ChromGate® System Suitability Setup.....	185
Copy & Paste	186
Suitability Calculation Selection.....	186
Running a Suitability Test.....	187
Suitability Reports	188
ChromGate® PDA Option.....	189
PDA Method Setup	189
PDA Options Library	189
PDA Options Purity	191
PDA Options Spectrum.....	192
PDA Options Multi-Chromatogram	194
PDA Options Ratio.....	195
PDA Views	196
3D View.....	196
Contour View	199
Mixed View.....	203
Chromatogram View	204
Spectrum View.....	208
Ratio View.....	210
Spectrum Similarity Table.....	211
Spectral Library Definition.....	212
How to Collect Spectra for a Library	213
How to Add Spectra to a Library.....	213
Library Search	214
Spectral Library Search	214
Custom Report.....	216
PDA Insert Graph Items.....	216
PDA Insert Report Items.....	217
Library Search Report	217
Library Definition Report	218
Spectrum Report	219
Spectral Display	221

Peak Table	221
Analysis Channel	222
PDA Analysis and Calculations	222
Chromatograms Extracted from the 3D Data.....	222
Multi-Chromatogram Channels	222
Working Chromatogram	223
Spectra Extracted from the 3D Data	223
Analysis Spectra	223
Working Spectrum.....	223
Background Correction	223
Spectrum Interpolation	224
Spectrum Smoothing	224
Spectrum Derivatives	224
Upslope and Down slope Spectra.....	225
Library Search Calculations	225
General	225
Pre-Filters.....	225
Ratio Chromatogram Calculation	226
Similarity Calculations	226
Lambda Max/Min Calculations	226
Noise Spectrum Calculations	227
Peak Purity Calculations	227
Background Correction	227
Calculating Total Purity	227
Three Point Purity	228
Spectrum Export	228
PDA Data Export.....	229
ChromGate® Preparative Option.....	230
Fraction Collector Configuration	230
Virtual Fraction Collector Configuration.....	233
Multi Valve Fraction Collector Configuration	234
Fraction Collector Setup	235
Fraction Collector Instrument Status.....	247
Fraction Annotations	249
Stacked Injection.....	250
SEC Option	257
Overview	257
SEC Calibration.....	257
Applications of SEC	258
Using ChromGate® SEC Software	258
Running an SEC Calibration Standard	259
Single Run Acquisition	260
Defining SEC Baseline.....	263
Defining SEC Result Ranges.....	263
Define SEC Peaks	263
Annotating SEC Chromatograms.....	264
SEC Setup	267
SEC Setup Narrow Standard Method	267
SEC Setup Universal Calibration	270
SEC Setup Broad Range1 Calibration.....	273
SEC Setup Broad Range 2 Calibration.....	275

SEC Custom Reports	278
ChromGate® SEC Equations	281
Typical Wiring Schemes.....	284
Index.....	291

Installation Guide ChromGate® 3.3.2

General Definitions

ChromGate is the Knauer OEM version of the original EZChrom *Elite* software from Agilent Technologies. The computers must be running under Windows 7 Professional or Ultimate 32 bit, Windows Vista Business or Ultimate 32bit or Windows XP Prof. Service Pack 2 for EZChrom Elite must be installed. Principally 64 bit or Home versions of Windows are not supported. The licenses are to be installed on the enterprise machine. The enterprise machine is also the host of all system configurations (instruments) and the user and project management data. A server is the computer that is connected to the HPLC system. A computer that controls the server in a Client/Server system is called client. The TCP/IP protocol and Microsoft's framework .Net version 3 or higher must be installed on all computers. The TCP/IP protocol can be installed using the standard procedure with the Windows installation disks, the .Net version 3 installer can be found on the ChromGate Installation Disc. For the installation of the ChromGate software, you must be logged onto the network domain in which you will be working (only for Client/Server installation) and the person performing the installation must have administrator rights. All other software, especially virus detection software, must be shut down.

ChromGate 3.3.2 is shipped on one DVD. It includes the original EZChrom Elite 3.3.2 from Agilent, the Service Pack 2 of EZChrom Elite 3.3.2 and the Knauer Add-ons as drivers and additional program options. The "XXXX" for the ChromGate version is the internal build number (e.g. 1787) and may change if a new build will be released. A new build may include patches and/or the support of additional devices. The license is stored on a USB license dongle, labeled with "ChromGate" and the serial number of the license.

Installation

First, EZChrom *Elite* 3.3.2 must be installed from the "ChromGate installation Disc". Click on the link "Install EZChrom Elite 3.3.2 software" to start the installer. Please refer to the EZChrom *Elite* installation guide to install EZChrom *Elite*.




Fig. 1

ChromGate installer start window

If EZChrom Elite is installed, after a computer restart, the Service Pack 2 of EZChrom Elite 3.3.2 must be installed. Click on the link “Install Service Pack 2 for EZChrom Elite 3.3.2” to start the installer. Finally, ChromGate as an Add-on for the EZChrom *Elite* software must be installed. If the Knauer FRC Option is purchased, the fraction collector drivers and the advanced fraction collection functionality must be installed with a separate installer (link “Install KNAUER FRC Option”). The Knauer FRC Option can only be installed, if ChromGate with the same build number has already been installed.



Do not insert the USB license dongle into the USB port before the software installation is completed, otherwise the license will not be recognized.



For all installation work login with administrator access. Before the installation, switch off all running programs, especially anti-virus software. In Windows Vista and Windows 7 you must confirm that the software should be installed. Please also login as administrator if the computer has rebooted after installation, otherwise the installation will not be completed.

To install ChromGate, click on the link “Install ChromGate software”.

The components which should be installed can be selected in the next window. As default, all components which ChromGate supports are selected.




Fig. 2

Click on “Next >” to confirm this window. In the next window click on “Install” to start the installation. To cancel the installation procedure click “Cancel”.




Fig. 3

The following window will be displayed when the installation has completed.




Fig. 4

Click on “Finish” to finalize the ChromGate installation.

KNAUER FRC control option

If you own a fraction collector supported by ChromGate you can install the KNAUER fraction collectors AddOn (KNAUER FRC control option). This option can only be installed if ChromGate v. 3.3.2 is installed. To control a fraction collector, a KNAUER FRC option license must also be subsequently installed.

Click on “Install KNAUER FRC option” to start the installation of the KNAUER FRC control option.




Fig. 5

In the next window, the fraction collectors that should be installed can be selected. As default, all supported fraction collectors will be selected.



Fig. 6

Click on “Next >” to confirm the window. In the next window click on “Install” to start the installation. To cancel the installation procedure click “Cancel”.




Fig. 7

The following window will be displayed when the installation has completed.




Fig. 8

Click on “Finish” to finalize the KNAUER fraction collectors AddOn installation.

The installer for the Generic Drivers can be found in ChromGate\GenericDrivers\Disk1. The installation is only required, if you wants to control one of the devices, supported via the Generic Driver. Generic drivers only allow for basic control of devices.

When using a DAD K-2700 or K-2800 with a PCI interface card, the driver for the interface card must be installed. This driver can be found on the ChromGate CD in ChromGate\Drivers. If the interface card has been installed, after a re-start, Windows will ask for the driver.

Start ChromGate once in Demo Mode by clicking Start – All programs – Chromatography – EZChrom Elite (default settings; refer the EZChrom Elite installation guide).

Close ChromGate. Insert the USB license dongle into a free USB port of the enterprise machine. Start ChromGate. The license will be recognized and can be used.

Short Guide ChromGate®, KNAUER Instrument Control

Overview

This section is provided to assist the operator in the day-to-day operation of Knauer HPLC systems using ChromGate® software. It focuses mainly on routine data acquisition and processing and is to be used in conjunction with the documentation provided with the ChromGate® software and the HPLC system. The full software description can be found in the online help. The online help will be installed along with the software. You can open it using the desktop link for the online help or with the help menu in a ChromGate® window. Additionally it will open in the context of the currently opened ChromGate® window, if you press the key <F1> on the computer keyboard.

This chapter assumes that the software has been successfully loaded on the computer or network and the instrumentation for which it is used has been installed. It does not include topics that are relevant to installation, interfacing of the unit to other components of the HPLC system, major maintenance and other topics that are more properly of interest to the system administrator. In addition, it is assumed that the chromatographic conditions for the separation and the overall data processing procedures are well understood.

The automated operation of the system involves the use of:

- the **Method** - The method describes instrument control, data acquisition, data processing, reporting and exporting of data for a single run. It includes all of the parameters that are used to perform and process a run.
- the **Sequence** - The sequence is a listing of the samples to be analyzed by the system on an automated basis. The user can select the method to be used for each analysis.

ChromGate® includes wizards that lead the user through the generation of the method and the sequence. In addition, the ChromGate® software contains a broad range of features to assist the user and also includes a number of security options.

This manual is provided to introduce the analyst to the ChromGate program and provides a discussion of some of the most commonly performed operations. It is not meant to present a detailed step-by-step procedure for use of the system or the data processing software.

There are a number of sources for further information, including:

- The Operator Manuals and Instruction Manuals that are provided with the modules in the system.
- The documentation provided with ChromGate® Data System, particularly the User's Guide (available on the ChromGate®-DVD in EZChrom Elite 3.3.2\English\Manual
- Release notes provided with the software (available on the ChromGate®-DVD in ChromGate\ChromGate)
- Documentation provided with Microsoft® Windows® Configuring the System

Configuring a System

Instrument Configuration

For more detailed information, please refer to the chapter Instrument Configuration on page 36. Configuring the system involves the selection of the instruments and features that are to be used (during installation, each component of the HPLC system is registered on the network). Once an instrument configuration has been made it can be easily retrieved.

There are two aspects to system configuration:

- **Indicating the modules** to be used in a given configuration. As an example, if the laboratory contains both a Model K-1001 pump and a K-1800 pump and only the Model K-1001 pump is to be used for a given separation, setting the instrument configuration includes selecting this pump and choosing various pump parameters (the type of pump head that is used, and the desired pressure units).
- **Indicating the options** that should be used with the method (e.g. if the baseline check feature should be used).

To generate a new instrument configuration:

- Select File - New - Instrument in the ChromGate® main window create a new instrument file, presented by an instrument icon on the right side of the main window.
- Right click on the instrument icon and select Configure - Instrument on the context menu that is presented to access the Instrument Configuration dialog box.
- Verify that the Instrument type field indicates Knauer HPLC system and press <Configure> to open the Knauer HPLC dialog box.
- Click on the module that you want to include in the configuration and press the green arrow to put that icon onto the right-hand side. Access the Configuration dialog box for a module by double-clicking on it. The dialog box contains information about the specific module and lets you edit the configuration information. The configuration information must match with the configuration of the device. From devices, which will be controlled via LAN connection, you can read-out the configuration.
- Add and configure all of the modules you want to control, until the appropriate configuration has been created.
- If additional system options, as SEC, SYS or PDA, are desired, access the options dialog box by pressing the <Options> button. The purchased options must be enabled to make them accessible in the instrument. Please refer the next chapter System Options for more information.
- Once the configuration and desired options have been selected, press OK to return to the Instrument Configuration dialog box, indicate the desired name of the instrument and press OK to return to the main screen.

An existing configuration can be edited by right clicking on the icon for that configuration and selecting **Configure - Instrument**. If you alter a configuration, make sure that you first make a copy of the existing instrument if you want to maintain both configurations.

System Options

The *System Options* dialog box is used to access the following:

- System Suitability - used to calculate a number of parameters pertaining to the overall efficiency of the system such as the resolution, repeatability, peak asymmetry, number of theoretical

plates, noise and drift. A detailed description of the setup of this feature you can find in the **ChromGate® System Suitability Setup** chapter below.

- SEC - used to calculate a number of parameters for SEC/GPC related chromatography. It adds special calibration methods and report options. Please find more in the **SEC Option** chapter in this manual and in the original manual *EliteSEC*, found on the installation DVD in *EZChrom Elite 3.3.2\English\Manual\Optional Software Manuals*.
- PDA - used to calculate a number of parameters that are related to the use of photodiode array detection. A detailed discussion of this feature is presented in the chapters **ChromGate® PDA Option** starting on page 185 and PDA Analysis and Calculations on page 222. The original software manual PDA Analysis can be found on the installation DVD in *EZChrom Elite 3.3.2\English\Manual\Optional Software Manuals*.
- Baseline Check - used to verify that the baseline is stable enough to provide useful chromatograms. A detailed discussion of this feature is presented in the chapter Setting up the Baseline Check on page 136.

Setting up a Method

A **method** is a set of parameters that describes the instrumental conditions for a chromatographic separation and includes information about data processing and reporting. It includes the parameters for the operation of each of the various modules in the HPLC system such as the pumps, detectors, oven and the autosampler, which are downloaded to the various modules for use when the separation run is initiated.

The Method window is used to set up a method, and is accessed via the Method wizard. The Method window includes a tab for each component of the system, and the tab presents various parameters that are available via the control panel of that module. The various tabs that are presented in the Method window are described in the chapter **The Method Window**.

Method files can be saved (*.met) and retrieved as well as imported from or exported to another instrument in the same manner as other Windows files. Please note, that for devices, which are not configured in the original instrument, the parameters will be set to default values.

The Instrument Wizard

The Instrument Wizard is used to lead you through the steps required to create a new method or modify an existing one, create a new sequence or run a method or sequence. The Instrument Wizard will always start, if you open an instrument. You can also start it from the opened instrument window with a click on the Wizard icon.



Fig. 9 *Instrument Wizard*

Select Create or modify a method to bring up the method window:



Fig. 10 *Method Wizard*

Click the appropriate button to present the Method window, which can be edited as described in the section **The Method Window**.

If you are more familiar with the menus, you can ignore the wizard. If you don't want to start the Wizard always if you open the instrument, disable the checkbox "Show at instrument startup".

The Method Window

Overview

The Method window includes a tab for each instrument in the system on which you can select the particular run parameters. The information on a given tab is essentially the same as if you were to program the component on a standalone basis. The reader should review information about parameter selection that is provided in the instruction manual for each module.

Pump Tabs

The Pump tab contains several areas where you can enter parameters related to operating mode, control pressure limits, and flow programs.

The tab for the Smartline Pump 1000 is presented below (the tabs for other pumps are very similar).




Fig. 11 Instrument setup tab (Smartline Pump 1000)

Details for the **Control Settings**, **Control Pressure Limits** as well as the creation of a **Time Table** are described in the section Instrument Setup – Pumps on page 84.

Detector Tabs

The format of the detector tab is dependent on the nature of the specific detector. A typical tab is shown below. Some of the parameters are common to essentially all detectors and are described at the beginning of the chapter Instrument Setup – Detectors on page 89, followed by details which are specific to certain detectors:

RI Detectors (S 23[4]00, K-23[4]00/1).....	page 92
UV Detectors (S 200, K-200)	page 92
UV Detectors (K-2000, K-2500)	page 92
UV Detectors (S 2500, K-2001, K-2501, UVD2.1S)	page 93
Instrument Setup – Kontron Detectors	page 95
K-2600 Detector.....	page 95
UV Detectors S 2550 / S 2520.....	page 80
Diode Array Detectors (S 2600, K-2700, DAD 2800/2850 and Kontron DAD 540)	page 100
Fluorescence Detector RF-10Axl / RF-20A/Axs	page 103
Instrument Setup – Alltech 650 Conductivity Detector.....	page 104
User defined Detector.....	page 104
Instrument Setup – Virtual Detector	page 105




Fig. 12 UV Detector S 2600 tab

Autosampler Tabs

Different autosamplers can be controlled by the ChromGate® software. Any method window will include the tab for the configured one, to setup the needed injection parameters. You will find detailed setup information for the autosamplers starting on page 111.

Fraction Collector Tabs

Several fraction collectors including the Knauer multi valve FC and a virtual one can be controlled by the ChromGate® software. Any method window will include the tab for the configured one, to setup the needed control parameters. You will find detailed setup information for the fraction collectors on page 235.

Miscellaneous Instrument Tabs

Further instruments to be controlled by this software are:

Manager 5000/5050, IF2 (analog/digital outputs of the Knauer A/D converters)

Knauer Switching Valves

Valco Switching Valves

Column oven 4050, Jetstream oven

Flowmeter

For each of them a tab will be shown in the method window as far as the instrument has been configured. If you have included more than one of the Knauer or Valco Switching Valves from one configuration window, even if they are of different type (K-6, K-12, or K-16) they are represented together in one tab. The detailed description is given on page 131.

Baseline Check tab

The *Baseline Check* tab is presented if *Baseline Check* is selected on the *Configuration Options* dialog box (see section Setting up the Baseline Check on page 136). It is used to ensure that the baseline is satisfactory. The Baseline check must be enabled in the Options menu of the instrument configuration.



To use this feature, select the **Baseline Check** option in the **Single Run** dialog box or include **Baseline Check** in the current **Sequence line**.

When the baseline feature is selected, the initial conditions for the run will be used to acquire baseline data. If the data is satisfactory for all channels, the run will be performed; if the baseline data is not satisfactory, the run will be aborted and a message to that effect will be presented.

Trigger tab

The *Trigger* tab is used to indicate the action to start a chromatogram. The definition of the various options is presented on the tab. For details see page 135.

Quantification of the Compound(s) of Interest

Overview

The change of the signal from a detector is related to the concentration of the eluting compound. Under ideal circumstances, the peak for each compound of interest will be well-resolved from other peaks in the sample and will be a Gaussian peak. In this case, the area of the peak is directly related to the concentration of the compound of interest. The analyst could generate a calibration plot using standards of known concentration to determine the concentration/area relationship, and then use the calibration plot to determine the concentration of unknowns.

The calibration data is incorporated into the method, so that when you select a method to use, all of the information that is required for quantization is loaded.

Integration

Integration of the chromatogram requires two basic parameters, the **width** and the **threshold**, which define the starting point and ending point for a peak and to distinguish a peak from noise. When a method is established, the default values are entered, and these can be edited via the Integration Events window which is accessed by selecting **Method/Integration Events**.

#	Event	Start Time	Stop Time	Value
1	Width	0.000	0.000	0.2
2	Threshold	0.000	0.000	50
3				

Fig. 13 Integration Events window

The **Width** value is used to calculate a value for bunching or smoothing the data points before the integration algorithm is applied. Integration works best when 20 points are sampled across a peak. The width should be selected for the narrowest peak in the chromatogram.

The **Threshold** value is the first derivative of the chromatogram, and is used to distinguish the peak from noise and/or drift.

To set the noise and drift, select the desired event field, then select the start and stop time and value.

ChromGate® includes a large number of additional parameters that can be used to refine the chromatogram. These can be applied on a programmed basis or on a manual basis for data that is already collected. A detailed discussion on integration is presented in the ChromGate® Chromatography Data System, Reference Manual.

Generating and Using a Calibration Curve

A calibration curve is normally generated via a sequence, which is discussed in the chapter Quantitation and Using Sequence Files below.

Method File Commands

A method file can be saved, retrieved and edited as desired in the same way as with other Windows applications. As an example, if the Method dialog box is open and you want to save a new method, select File/Method/Save as and enter the desired name.




Fig. 14 Save Method File as... Window

Similarly, if you want to open a method, select File/Method/Open... and enter the file name.

A printed report describing all aspects of a method can be obtained by selecting File/Method/Print.

Quantification and Using Sequence Files

Overview

Calibration involves the generation of the relationship between the concentration of a compound and the area of the corresponding peak in the chromatogram. Once calibration is established, it can be used to determine the concentration of the compound of interest in unknowns.

A single level calibration involves a single sample to calculate the calibration constant; the curve is generated from that point and the origin. A multi-level calibration uses several standards (e.g. 10 µg/ml, 20 µg/ml, 50 µg/ml, 100 µg/ml).

Setting up a calibration involves the following steps:

1. Run the chromatograms containing the standards and save them
2. Identify the peaks from the standards
3. Generate a peak table
4. Generate the standard curve

The standards used to generate a calibration should be well defined and should be in a matrix that is similar to that containing the samples that you will analyze. The standards undergo exactly the same sample preparation process that is used for real samples.

There are two general ways that a standard can be employed:

- Internal Standard - A known amount of the compound used as a standard is added to the sample, and the area for the compound of interest is compared to the area of the standard.
- External Standard - The standard is a separate sample and the observed intensity is used to generate the calibration curve directly.

The application program allows you to generate calibration curves for a number of compounds in a single operation and an analytical report can be generated, printed and exported as desired.

Single Level Calibration

To set up a single level calibration:

1. Collect a chromatogram containing the standard and save it. The concentration of the compound in the standard should be similar to that in the samples that you want to analyze.
2. Open the stored data file
3. Press the **Analyze** button  to integrate the chromatogram and show the baseline.
4. Press the **Define Single Peak** button  to present the *Define Single Peak* dialog box. The retention time in the upper right corner will refer to the first peak in the chromatogram.




Fig. 15 Single peak definition

5. Press the **Next>>** button until the peak corresponding to the standard is selected. As an alternative, you can simply click on the peak in the chromatogram with the cursor.
6. Complete the dialog box for the first standard:
 - Enter the peak name
 - Indicate the concentration level. For a single level concentration, enter 1 for the level and indicate the concentration of the standard.
 - Indicate the desired unit (e.g. µg/ml)
 - Enter the internal standard ID# if you are using the internal standard method. This number is obtained from the peak table. Leave it blank until you have generated the peak table (see below), then enter the value directly into the **Peaks / Groups table**.


- Enter the Reference ID#, which is the Peak ID#s from the peak table. This number is obtained from the Peak Table. Leave it blank until you have generated the peak table, then enter the value directly into the **Peaks / Groups table**.
 - Indicate the desired basis for the calculation of the retention time window and the desired value.
7. Repeat the above step for each compound in the standard. When you have entered information for all standards, press **Done**.
 8. Press the **Peaks / Groups Table** button  or select the menu sequence *Method – Peaks / Groups ...* to present the calibration peak table. Each peak that you have defined will appear as a row in the table, along with the retention time and other parameters that you have entered.



#	Name	ID	Ret. Time	Window	Ref. ID #	ISTD. ID #
1	Compound 1		5.70958	0.285479	0	0
2	Compound 2		8.26361	0.413181	0	0
3						

Fig. 16 Peak table

9. Enter all values needed for the calibration as the concentration for all levels of your compounds.
10. If the table is complete, close Peak Table dialog box by clicking on the x in the upper right corner.
11. Save the method.
12. Select **Analysis/Analysis Single Level Calibration** from the menu bar to present the following dialog box.



Analysis / Single Level Calibration

Analysis information

Sample ID: PNA-STD5

Method: C:\EZChrom Elite\Enterprise\Projects\Default\M

Data path: C:\EZChrom Elite\Enterprise\Projects\Default\D

Data file: multi calibration level 2.dat

Print method report

Amount values

Sample amount: 1

Internal standard amount: 1

Multiplication factors: 1 1 1

Dilution factors: 1 1 1

Calibrate

Calibration level: 1

Clear all calibration

Clear calibration for level

Print calibration report

Clean replicates

Average replicates

Start

Cancel

Help

Fig. 17 Analysis Single Level Calibration Window

13. Enter the sample identification used for the standard.
14. Enter the method that you want to calibrate and the data path name for the data file.
15. Enter the name of the **calibration** data file.
16. Click on the Calibration check box and enter a 1 for Calibration level. If the method is presently not calibrated, it is not necessary to check any of the boxes (but if you want to be certain, select the Clear all calibration check box).
17. Repeat the steps 12. to 16. for all levels you want to calibrate. Enter the appropriate calibration level into the Calibration level field

18. Click **Start**. The peaks will be integrated, and calibration curves will be generated.

The method is now calibrated and can be used to analyze samples.

Creating a Sample Sequence

If you are using an autosampler, the sequence is used to define the samples, indicate how they are to be injected, the method used to separate them and how the data is to be calculated. A sample sequence can be used to acquire both calibration data as well as data from unknowns (and both can be performed in a single sequence). In addition, a sequence is used to generate a calibration curve and analyze unknowns from stored data.

To setup a sequence:

1. Click **File/Sequence/Sequence Wizard** or select **Create a Sequence** on the *Instrument Wizard* to access the Method page of the Sequence wizard.




Fig. 18 Sequence Wizard – Method

2. Enter the method name
3. Select the mode of data acquisition. If the data is to be collected as part of the sequence, the Amount Values fields will be activated. When you are beginning to develop a sequence, ignore them and press **Next** to present the Sequence Wizard –Unknowns. If the From existing data files option is selected, a dialog box will be presented to select the desired files.




Fig. 19 Sequence Wizard –Unknowns

4. Press the blue arrow on the Sample ID line to access a list of identification options for the samples. Typically, Line Number and Method Name are used, so that each sample will be identified with the sequence line number and method name.
5. Indicate where the data should be stored using the Data Path field.
6. Click on the blue arrow on the Data File line to access how the data file should be named. Typically Sample ID is used, this will assure that the sample is uniquely identified and easily related to the sample.
7. Enter the appropriate number of unknowns the number of repetitions.
8. Press next to present the Autosampler screen if you have indicated that the data files are to be acquired (initial screen in the wizard). (this screen is not presented if existing data files are used).



Fig. 20 Sequence Wizard – Autosampler

9. Enter the appropriate information and press **Next** to present:




Fig. 21 Sequence Wizard – Calibration

10. Enter the desired information and press **Next** to access the Reports page.



Fig. 22 Sequence Wizard – Reports

11. Enter the desired information and press **Finish** to present a spreadsheet that describes the overall sequence.

Run #	Status	Run Type	Level	Conc Override	Reps	Sample ID	Method	Filename
1		Summary Begin	0	n/a	1	<D>	test.met	001.d
2		Summary Run	0	n/a	1	<D>	test.met	002
3		Summary Run	0	n/a	1	<D>	test.met	003
4		Summary Run	0	n/a	1	<D>	test.met	004
5		Summary Run	0	n/a	1	<D>	test.met	005
6		Summary End	0	n/a	1	<D>	test.met	006
7								

Fig. 23 Sequence Wizard –Description of a finished sequence

12. Save the sequence.

Running a Sequence

To acquire data using the sequence file you just created, click **the Sequence Run** button on the command ribbon, or do a right-hand mouse click in the sequence spreadsheet, and select Run Sequence. The following dialog box will appear:




Fig. 24 Run Sequence dialog box

Enter the name of your sequence file by typing the name, along with path, in the **Sequence Name** field. You can also select it from a list of sequence files on your disk by clicking the **File** button next to the field. Leave the other parameters as their defaults.

Prepare your autosampler to inject your standard sample, followed by 3 unknown samples. When you are ready to inject your first sample, click Start. When the sequence is completed, you will have acquired and saved the data files for one standard and three unknown runs, and generated a simple result report for each unknown sample and a summary report for the sequence.



Since you have not yet defined a sequence summary report, do not check the box to print sequence reports.

Reporting

ChromGate® Client/Server comes with a complete suite of report templates that can be used without modification to generate reports. To see an example of one of these reports, use the **Reports/View/Area%** command from the menu bar. (Make sure your current chromatogram has been analyzed first.) The standard report will appear in a window on your screen.




Fig. 25 Area% Report

If your method contains no defined custom report, the system will use the standard report formats to print reports when required.

Knauer has added three reports in a fresh design (**Knauer Area% Report, Knauer ESTD Report, Knauer ISTD Report**), which also can be found as standard reports in the **Report** menu.


HPLC · SMB · Osmometry



Date File: C:\EZChrom Elite\Enterprise\Projects\Data\multi calibration level 2.dat
Sample Name: PNA-STDS

Acq. Operator: demo Location: N/A
Acq. Instrument: Knauer HPLC (Offline) Inj Volume: 0
Injection Date: 26.11.1990 20:49:09

Analysis Method: C:\EZChrom Elite\Methods\multilevel calibration.met
Last changed: 14.07.1999 21:14:22
Description: {Data Description}



KNAUER Area% Report

Analysis Time: 28.02.2012 10:48:35
Multiplier Factor1: 1
Dilution Factor1: 1
Sample Amount: 1

Channel A Results									
Ret. Time	k'	Area	Area %	Height	Height %	Asymm.	Width	Plates	Resolution
5,71	0,00	170464	26,30	20886	36,67	0,00	0,651		
6,55	0,00	139959	21,60	15331	26,92	0,00	0,651		
6,84	0,00	306	0,05	59	0,10	0,00	0,168		
8,53	0,00	337346	52,05	20673	36,30	0,00	0,927		
Totals		648075	100,00	56949	100,00				

End of Report

Knauer HPLC (Offline) 22.06.2012 10:50:18 mathias

Page 1 of 1

Fig. 26 KNAUER Area% Report

If you wish, you can modify the standard report templates, or create entirely new reports using the *Custom Report* capability of ChromGate® Client/Server. Modified standard reports can be stored as templates, *Custom Reports* will be available as part of the method. You can create custom method reports and / or custom Sequence reports.

To view the custom report template in the multilevel calibration.met file, first open the file if it is not already open. (Use the **File Open** button, followed by **Method**; then select the template from the file list.) Click the **Edit Custom Report** button on the command toolbar to access the method custom report editor. The current method custom report template will appear.

Collecting Data

General Data Collection Instructions with ChromGate



The user should review the operating manual for each component in the system for warm up requirements, preparing the component for operation, etc.

This section includes requirements and suggestions relative to the computer and e-line network when ChromGate® is used to operate the HPLC system.

Before starting the chromatographic runs:

- Make sure that the energy saving feature and the screen saver on the computer are off.
- Check that all modules are connected to the e-line cable
- To power up the system, first turn on the PC and finally power up the various modules.



If a unit loses power or if the communication is broken, turn off all components and start over.

- If a module has been operating via the local mode (UI interface), make sure that it is not in use when connecting the cable.
- If you are using a manual injector, make sure that the system is in the waiting for injection mode waiting for trigger before injecting.
- If you are using an autosampler with a single rack, the input numbers should be integers. If you are using multiple racks, input vial numbers, using a letter and number (e.g. A1, B3, etc.) should be used. If an incorrect vial number or invalid volume is noted, the sequence will stop.

Collecting Data

When you want to collect data, click on the desired instrument to present the Instrument Wizard dialog box (Fig. 9 on page 16) and select the desired operation.

If you select **Run one sample**, the Single Run Acquisition dialog box will be presented. Enter the desired Sample ID, method and name (e.g. sample.dat) and then press **Start**.

The chromatogram will be collected as described in the method. During the data acquisition, the chromatogram will appear in the window as shown below.




Fig. 27 Single run data acquisition window

If you have selected **Run a sequence of Samples**, the following dialog box will be presented:




Fig. 28 Sequence Run Acquisition dialog box

Enter the *Sequence name* and press **Start**. The chromatogram for the first run will be collected when the conditions have been met.

Instrument Status of a running Control Method

During Data acquisition, you can access the status of the system by selecting Control/Instrument Status on the main menu. The window that will be presented is shown below, and the format and tabs is dependent on the configuration of the system.

Open the Instrument Status Window using the menu sequence **Control – Instrument Status**.

The Instrument Status window contains the status tabs for all configured instruments, and as shown in the figure below the system status tab, providing an overview for the whole system.




Fig. 29 Instrument Status window, system status tab of a running method

The running method is visualized by the dark dotted lines between the displayed instruments.

For each pump, identified by name and serial number, the actual flow, pressure is displayed. If a gradient system is configured the solvent composition also will be shown.

For each detector channel the wavelength and output values are shown. Similarly for all included instruments the relevant data will be displayed, e.g. the actual switching positions in case of the valves.

All these information you also find on the single instrument tabs of this window, for details see pages 139 - 169. On these tabs you get the possibility of direct control interaction with the individual instruments. This is even possible during a method is running; however this option must be enabled while configuring the system (refer to Knauer instrument control method options on page 170).

Shutting Down the System

The system can be automatically shut down at the end of a sequence by including a shutdown method at the end of a sequence. A discussion of sequence types is presented in Chapter 4 of the reference manual.

The shutdown method should include the following;

- Column Oven temperature – value should be below ambient
- Pump conditions, the flow for all pumps should be 0.0 ml/min and 100% A and no events. For new pumps as Smartline pump 1050 also a **Standby** option is available.
- Deselect any timed events and Mix methods.
- Set the Trigger to **None**. Now an autosampler will not be controlled and will not perform an injection.
- Detector – activate the **Lamp Off** option or, if supported, the **Standby** option.

Setup and Control of Knauer HPLC Systems

This manual contains detailed information on the use of the ChromGate[®] instrument control software option for the Knauer Instruments.

It contains information on how to configure, set up and run samples automatically from the ChromGate[®] Chromatography Data System.

This manual is intended to be a supplement to the general ChromGate[®] Data System documentation and contains information specific to the use of the optional instrument control software.

Overview of Instrument Control

The instrument control software for the Knauer HPLC System enables you to enter instrument control parameters that become part of the method. The parameters are then executed in real time when samples are acquired using the data system. The software allows you to set and control parameters as pressures, flows and gradients for pumps, wavelengths and scans for detectors, injection parameters for autosamplers on the Knauer HPLC System, and monitor these in real-time.

Configuration – device communication port

Some of the devices have two communication ports on the rear panel for controlling the device by the computer. One is an RS-232 port, also called serial or COM port; the other is an Ethernet port, also called LAN (Local Area Network). Beside the Knauer valve drives, for all the devices the desired communication port must be defined in the device. Please pay attention to the corresponding notes in this manual and refer to the device's manual for more information, how the desired communication port must be configured. Please note, that for some devices the functionality depends on the selected communication port.

Configuring the Interface

Beside all detectors directly supported in ChromGate[®], any detector or sensor producing an analog signal which can be digitized with an A/D converter or Interface Box can be used for data acquisition. You must configure your data acquisition interface before you can acquire analog data using the data system. The Knauer interfaces can also be used to control other devices by analog or digital output.

Interface configuration is accessed through the **Tools** command on the ChromGate Main Menu.

Click the **Tools - Interface Configuration** command from the ChromGate[®] Main Menu. A window will appear displaying several possible interface devices. To configure a device, click on the icon to select it, then click the <Properties> button or double click on the instrument icon.



Fig. 30 Selection window for interface configurations

Knauer Interface Configuration

If you are using either a Knauer Interface Box/IF2 or a Manager 5000/5050, click on the corresponding icon, and then click **Properties**. A dialog box will appear where you can configure the Interface Box (see Fig. 32).



In regards to Knauer software, the Manager 5050 interface module is identical to the Manager 5000. To configure the Manager 5050, please select „Manager 5000“.

If you configure the interface newly, the message like Fig. 31 will remind you to enter the serial number of the interface device.



Fig. 31 Error message for missing serial number

Click OK and perform the required settings. The same message will appear if you try to close the setup window without an entered serial number.




Fig. 32 Dialog box for configuring the Manager 5000/5050

The Knauer Interface Box and the Manager 5000/5050 require the same parameters to be set up. Therefore, the configuration dialogs for both devices are nearly identical.

Number of boxes

Select the number of boxes used in the system. ChromGate® supports up to 4 Knauer Interface Boxes.

Configure box

Select the box you wish to configure. Each box has to be configured separately.

Serial Port

Manager 5000/5050 / Interface Box IF2

Select the number for the communication port on your PC to which the Manager 5000/5050 or the Interface Box IF2 is connected.

Interface Box Model 96

The Knauer Interface Box supports the binary Knauer Net. The networks 1, 2,... correspond to the number of the utilized communication port (COM1, COM2,...). Select the number for the communication port on your PC to which the Interface Box is connected.



You can use the Interface Box or the Manager 5000/5050 to acquire signals from different instruments at the same time. It is necessary to connect the interfaces to separate communication ports on a PC. See the Typical Cabling / Wiring chapter, Fig. 389 on page 285.

Serial Number

Type in the serial number of the Interface Box / Manager 5000/5050 you are using. The serial number is used for the device identification and addressing during serial communications. This serial number should match the serial number of the Interface Box. For Knauer HPLC-Boxes with an S/N higher than 65000, only the last 4 digits of the S/N must be entered (e.g. 5567 instead of 65567).

Box type

The only choice is Model 96. No others are selectable.

Channel

Four data acquisition channels are normally available. However, only channel one will be available if you have checked the **100 Hz** sampling rate option in model 96 or **50 Hz** in case of Manager 5000/5050 / IF2 (see Sampling rate on page 91).




Fig. 33 Dialog box for configuring the Model 96 Interface Box with activated 100 Hz option

Range

Select the Signal Range for the channel chosen from the channel drop-down list. The choices are **+/-10 V**, **+/-1 V**, **+/-0.1 V**, **+/-0.01 V** for the Interface Box and **+/-2.56 V**, **+/-1.28 V**, **+/-0.64 V**, **+/-0.16 V** for the Manager 5000/5050 / IF2. The Signal Range specifies the maximum

analog input voltage to be digitized by an A/D converter. Use the Signal Range to optimize the signal-to-noise ratio of the measured signal.

Manufacturer and Model

This area is for information only. It contains manufacturing information about the box which may be required in service situations.

Kontron Interface Configuration

If you are using a Kontron interface card click on the corresponding icon, and then click **Properties**. A dialog box will appear where you can configure the interface (see Fig. 34).

If you configure the interface newly the message like Fig. 31 on page 32 will remind you to enter the serial number of the interface device.



Fig. 34 Dialog box for configuring the Kontron interface card PCI AD

Number of PC cards

Select the number of PC cards used in the system. ChromGate[®] supports up to 4 interfaces.

Configure card

Select the card you wish to configure. Each card has to be configured separately.

Card type:

Select either ISA AD or PCI AD. In case of selected PCI AD the setup window is changed to the following appearance.




Fig. 35 Dialog box for configuring the Kontron interface card ISA AD

K-2700 / 2800 Interface Card and Driver Installation

ChromGate® communicates with the Diode Array Detectors K-2700 and 2800 via either a PCI interface card or LAN. Installation of the necessary PCI card is described in the detector manual. Please note, that the driver is not released for Windows Vista or Windows 7. For other PDA detectors (S 2800, S 2850...) only the LAN interface will work.

Configuring the KNAUER HPLC System

In order to access commands for control of Knauer instruments and autosamplers, you must configure them appropriately. To access the configuration screens, locate the icon on the main menu for the instrument you want to configure.



You must have system administration privileges in order to configure an instrument.




Fig. 36 ChromGate® main screen

Select File – New to create new instruments or groups:




Fig. 37 Creating new Instrument or Location/Group

If you have selected *Instrument*, the Instrument Configuration dialog (see Fig. 38) appears. For already existing instruments right- mouse click on the instrument, and select the **Configure...** command.

The Instrument Configuration dialog will appear.



Fig. 38 Instrument configuration dialog box

Instrument name

Enter a descriptive name by which the instrument is identified. This name will appear under the instrument icon, and when the instrument window is displayed.

Instrument type

Select the **Instrument type** from the drop-down list. Knauer devices are only available in the **Knauer HPLC system** type

Server name (in Client/Server only)

Select the **Server name**. This is the name of the CG server to which the instrument is attached.

Finally, click on **Configure** to proceed to the instrument configuration where one can define the detectors, pumps, and further devices for the instrument.

Instrument Configuration

Configuration – KNAUER HPLC System

When you click the **Configure** button, the following dialog will appear.




Fig. 39 Instrument selection window

Several icons will be displayed in the **Available Modules** box on the left. Add modules to be configured by double-clicking on each, or by clicking once on the icon, followed by the green arrow.




Fig. 40 Instrument selection window with selected items

Auto Configuration

Using the **Auto Configuration** button, all devices, connected via LAN and switched on, will be automatically added and configured.



Fig. 41 Auto Configuration start window

All Knauer devices use a single IP port for LAN communication. Only devices with the same IP port as selected in the Auto Configuration window will be found. The default IP address is 10001. It can be changed manually in the device's setup.

We recommend reviewing the configuration of all automatically added devices. Devices connected via RS-232 (serial connection) will not be found and must be added manually.

Now configure or review **each** module (detector, pump, autosampler, or event configuration) separately. Not configured instruments may be indicated by a question mark on the corresponding icons. If you try to exit the window with **OK** and one or more instruments are still not configured, a message like Fig. 42 will appear.




Fig. 42 Example of an operation message to complete configuration

Options

The **Options** button is used to initiate any optional analysis software for this instrument.

A list of available analysis options will be presented.




Fig. 43 analysis options window

Click the box(es) to select the optional software you want to enable for this instrument. The additional licenses are required beside the Baseline Check. For more information regarding the option, please check the **System Options** part above on page 14. The Knauer Fraction Collectors Control option (FRC) must not be enabled here. Whenever a configuration with a fraction collector is started, an installed FRC license will be activated.

Configuration – Knauer Pumps

For each pump in the Configured Modules window, double-click the icon and complete the configuration dialog.



For configuration information regarding the PLATINblue pump P-1, please refer to the PLATINblue system manual.



If the pump S 1050 should be controlled by RS-232, in the pump's own setup the "Interface" must be set to "RS232 115200". Please refer to the pump's manual, chapter "Setup menu" for more information.



The pumps K-1000, K-1001, S 1000 and K-1800 (preparative pump 1800) only can be controlled with ChromGate, if the pump's internal communication in the menu option "CONTROL" is set to "NET baud rate 9600".



To control a pump K-501, in the pump's internal setup the menu option "CONTROL" must be set to "Ext:RS232". The firmware version must be v. 1.23 or higher.



If you are operating a K-1001 pump fitted with a firmware version < 5.x, you must select the K-1000 icon for configuration. This is due to the status information, which the older firmware does not support. Otherwise you will get error messages when the pump starts running. For those pumps no advanced status information and no pressure trace recording are available.

Pump S 1000 Configuration

Name: S 1000 Pump Type: S 1000

Gradient Mode: Isocratic

Interface: Network2 Serial Number: 67890

Head: 10 ml Pressure Units: MPa

Material: n/a
Max. Head flow: 10.0 ml/min
Max. Head pressure: 40.0 MPa

Add.Info...

OK Cancel Help

Fig. 44 Pump S 1000 configuration window

Pump 10P / 20P / P2.1S / P4.1S Configuration

Name / Type: Pump 10P Gradient Mode: Isocratic

Serial Number: Head: 10 ml

Interface: LAN IP Address: Obtained using S/N

IP Port: 10001

Use S/N to identify the instrument Pressure Sensor Add.Info...

Material: n/a
Max. Head Flow: 10.0 ml/min
Max. Head Pressure: 40.0 MPa

OK Cancel Help

Fig. 45 Pump 10P/20P/P2.1S/P4.1S configuration window

Pump S 1050 Configuration

Name: S 1050 Gradient Mode: Isocratic

Serial Number: Head: 10 ml

Interface: LAN IP Address: Obtained using S/N

IP Port: 10001

Use S/N to identify the instrument Pressure Sensor Add.Info...

Material: n/a
Max. Head Flow: 10.0 ml/min
Max. Head Pressure: 70.0 MPa

Config.Service... Add.Info...

OK Cancel Help

Fig. 46 Pump S 1050 configuration window

Pump P2.1L Configuration

Name: Pump P2.1L Gradient Mode: Isocratic

Serial Number: Head: 100 ml

Interface: LAN IP Address: Obtained using S/N

IP Port: 10001

Use S/N to identify the instrument Pressure Sensor Add.Info...

Max. Head Flow: 100.0 ml/min
Max. Head Pressure: 40.0 MPa
Leak Sensor Sensitivity: Medium

Config.Service... Add.Info...

OK Cancel Help

Fig. 47 Pump P2.1L configuration window

The configuration windows for the other Knauer pumps are nearly identical.

Name

The name of the selected pump appears automatically. You can enter any other descriptive name for the pump. This name will appear on the pump tab when the instrument setup or status windows are displayed.



Use a name which is unique within the instrument (configuration). E.g. for a high-pressure gradient system, all pumps must have different names.

Pump Type (accessible for Pump K-120/501 only)

Select the type of pump you want to configure from the drop-down list. The choices are **K-120** and **K-501**.



Fig. 48 selection option for K-120, K-501 pumps

Name / Type (accessible for Pumps 10P/20P/P2.1S/P4.1S only)

Select the type of pump you want to configure from the drop-down list. The choices are **10P**, **20P**, **P2.1S** and **P4.1S**. The 10P/P2.1S is the device without the pressure sensor. As already described in the Name section above, you must use different names for all of the devices in an instrument (configuration).




Fig. 49 selection option for pumps 10P, 20P, P2.1S and P4.1S

Gradient Mode




Fig. 50 selection options for gradient modes

Select the desired gradient mode from the drop-down list. Please note, that for pumps K-1000/1, S 1000, S 1050, preparative pump 1800/K-1800 and P2.1L the gradient mode must correspond with the gradient mode set in the pump. The LPG option is not available for S 100/10P/20P/P2.1S/P4.1S, K-120 and K-501 pumps. The **Isobar** mode is only available for the pump P2.1L.

Interface / Serial Port

Select from the drop-down list the serial port or network number for the communication port on your PC where the instrument is connected.



To ensure an error-free operation, each serial-controlled device should be connected with an own COM port, also for devices that allow to connect more than one device onto the same COM port.

For pumps K-1000, K-1001, S 1000 and K-1800 (preparative pump 1800) the interface "Network" means the serial port number.



The pumps K-1000, K-1001, S 1000 and K-1800 (preparative pump 1800) only are controlled with ChromGate, if the pump's internal

communication in the menu option “CONTROL” is set to “NET baud rate 9600”.

The pumps S 1050, S 100 allow for communicate via RS-232 or LAN. Select the interface you want to use for communication with the device.



Fig. 51 Available interfaces pumps S 1050 / S 100 / 10P/20P/P2.1S/P4.1S

For pump P2.1L only communication via LAN is supported.

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the pump is connected.



If the pump S 1050 should be controlled by RS-232, in the pump's own setup the “Interface” must be set to “RS232 115200”. Please refer to the pump's manual, chapter “Setup menu” for more information.

LAN




Fig. 52 LAN interface dialog

All LAN-controlled devices can be added and configured automatically using the **Auto Configuration** button on the configuration main window. Beside this, you can let the software search for your already connected and switched-on device from the devices configuration window. The option “Use S/N to identify the instrument” must be enabled. Click on the  button to start the search for all connected LAN devices of the currently configured type. Please refer the picture below for the information which will be displayed for the found devices. If more than one will be found, all hits will be shown in the window. Select the desired device from the shown list. The software will read-out the information as serial number, pump head and gradient mode from the pump and will fill it automatically into the configuration screen. For pumps P2.1L and S 1050 the gradient mode can be changed using ChromGate, if necessary. Please see the corresponding section below. For all Knauer LAN devices port 10001 is selected by default. Don't change this setting, otherwise no communication can be established. To search for an instrument, the device must be connected, switched on and must have a valid IP address. As an option, you can enter the IP address of your device, if the option “Use S/N to identify the instrument” is disabled. Please note, that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.



We highly recommend using the S/N and not the IP address to identify the instrument because devices may automatically receive a new IP address when reset.




Fig. 53 LAN interface Device Information window

Serial Number

If you have read-out the device via LAN, the serial number has been added automatically. Otherwise type in the serial number of the pump you are using. The serial number of network controlled pumps is used for the device identification and addressing during communication. This serial number must match the serial number of the defined pump. If not, you will get an appropriate error message if you later try to start in run.

Head

Select the type of pump head you are using from the drop-down list. The pump head must also be selected in the pump's own menu. For the pumps K-1001, K-1000, K-501 and K-120 pumps the choices are **10 ml** and **50 ml** of different materials (steel, PEEK, ceramics, titan).




Fig. 54 Available pump heads for K-1001, K-1000, K-501, and K-120 pumps

For the pumps S 1000, 1050, S 100 and 10P/20P/P2.1S/P4.1S only **10 ml** or **50 ml** is to select, not considering the material. The pump P2.1L and the preparative pump 1800 can be equipped with four different pump heads with 100, 250, 500 or 1000 ml respectively. In case of titanium inlays also select the only available steel pump head for the preparative pump 1800.




Fig. 55 Available pump heads S 1000 / 1050 / 100




Fig. 56 Available pump heads P2.1L (A); preparative pump 1800 (B)

Once you have chosen the pump head, the related information about material, maximal flow and maximal pressure is displayed automatically. If you have read-out the device via LAN, the pump head type has been selected automatically.

Pressure Units

Select the pressure units you want to use for the pump from the drop-down list. The choices are **MPa**, **bar**, and **psi**. All control pressure fields in the pump tab of the instrument method window as well as the instrument status window will be displayed in the selected pressure units. The pressure unit on the pump's display can be changed for the P2.1L in the pump's own menu; all others will always show MPa.

Leak Sensor Sensitivity (P2.1L only)

Select the leak sensor's sensitivity in accordance with your requirements. The choices are **Off**, **Low**, **Medium** and **High**.

Config. Service (P2.1L / S1050 only)

Click the button **<Config. Service...>** to enter the Pump Configuration Setup window of the pump. This will only work, if the pump is connected via **LAN**.



Fig. 57 Pump Configuration Setup window of pump S 1050

Using the **<Read>** button, you can read-out the pump head and gradient configuration. If this is done, also the **<Write>** button will be active. You can change the settings for pump head and gradient and send these settings to the pump using the **<Write>** button. Note, that, if the auto detection for the pump head (RFID-based pump head recognition) is activated in the pump, the pump head cannot be changed.

Leaving the configuration window with the **<OK>** or **<Cancel>** button will leave the window with the last sent setting.



**Pressing the <OK> button will not send any changed settings.
Pressing the <Cancel> button will not recover the original settings.**

Pressure Sensor (S 100 only)

If your pump S 100 is equipped with a pressure sensor, please enable the corresponding check box. For pumps 20P/P4.1S the box is automatically enabled, because these models always have a pressure sensor, while models 10P/P2.1S do not have this component; here the box is disabled.

Add. Info

Click the button **Add. Info...** to enter more detailed description of the pump.




Fig. 58 Additional information box

Enter additional information or comments to be printed out when you print the instrument configuration.



Without the correct serial number, the communication with the pumps cannot be established.

When complete, click OK to exit the dialog and return to the instrument configuration icon list.

Configuration – Kontron Pumps

For each pump (320/322/325, 420/422/425, or 520/522/525) in the Configured Modules window, double-click the icon and complete the configuration dialog. You have to make sure that all settings correspond to the settings at the pumps themselves. The configuration windows of all Kontron pumps (Fig. 59) look very similar that they can be described together.




Fig. 59 Kontron Pumps configuration window

Name

The name of the selected pump appears automatically. You can enter any other descriptive name for the pump. This name will appear on the pump tab when the instrument setup or status windows are displayed. Use a name which is **unique** within the instrument.

Pump Type

Select the type of pump you want to configure from the drop-down list.

Gradient Mode

In case of the 3xx pumps it is automatically set to isocratic and for the 5xx pumps to LPG. For the 4xx pumps the choices are HPG A, HPG B, or HPG C.

Head

Select the type of pump head you are using from the drop-down list. The choices are **Analytical**, **Micro**, and **Semi Prep** for all Kontron pumps.




Fig. 60 Available pump heads for Kontron pumps

Once you have chosen the pump head, the related maximum pressure and maximum flow rate are displayed automatically.

Serial Number

Type in the serial number of the pump you are using. This serial number should match the serial number of the defined pump.

Pressure Units

Select the pressure units you want to use for the pump from the drop-down list. The choices are **bar**, **MPa**, and **psi**.

Serial Port

Select the number for the communication port on your PC where the instrument is connected. The selection <none> is for simulations. No real control of the instrument will be performed.

Baud rate

Select from the drop-down list the baud rate you want to use for the pump. The choices are **4800** and **9600**. The selection must be the same as at the pump itself. The default value is 4800 for all Kontron pumps.

Add. Info

Click the button **Add. Info...** to enter more detailed description of the pump. The entered additional information or comments will also be printed out when you print the instrument configuration.

When complete, click OK to exit the dialog and return to the instrument configuration icon list.

Configuration – Knauer Detectors

The configuration dialogs of the supported detectors differ only slightly from each other. The appearance of the corresponding detector configuration windows is nearly identical (Fig. 62 through Fig. 71).



For configuration information regarding the PLATINblue detectors PDA-1 and MW-1, please refer to the PLATINblue system manual.




Fig. 61 Fixed wave length detector configuration window




Fig. 62 Older Knauer single channel detector configuration windows




Fig. 63 UVD2.1L configuration window




Fig. 64 10D/UVD2.1S configuration window




Fig. 65 S 2520 configuration window




Fig. 66 S 2550 configuration window




Fig. 67 K-2600 configuration window




Fig. 68 S 2600 configuration window




Fig. 69 DAD configuration windows




Fig. 70 RF-10AxI/ RF-20A/AXS detector configuration window




Fig. 71 Alltech 650 Conductivity detector configuration window

Name

Enter a descriptive name for the detector. This name will appear on the detector tab when the instrument setup or status windows are displayed. The default detector name will be suggested automatically by selecting the detector type.



Use a name which is unique within the instrument.

Detector Type

For some detectors, the correct model type can be selected in the configuration window. Select the type of detector you want to configure from the drop-down list.




Fig. 72 Detector type selection

For some detectors this field is not accessible because no choices are possible.

Interface / Serial Port

Select from the drop-down list the serial port or network number for the communication port on your PC where the instrument is connected.



To ensure an error-free operation, each serial-controlled device should be connected with an own COM port, also for devices that allow to connect more than one device onto the same COM port.

Interface (DAD 2700 / 2800 / 2850 only)

For the diode array detectors, the choice is **PCI** and **LAN**. If you control the DAD via the LAN, you must either press the ...-button to search for the detector or enter the IP Address of the detector. The access to the IP Address field is blocked if PCI was selected.

Interface (UVD 2.1L / S 2550 / S 2520 / 10D/UVD2.1S only)

Select the Interface you want to use for communication with the device.



Fig. 73 Available interfaces detector S 2550 / S 2520 / 10D/UVD2.1S

For detector UVD2.1L only communication via LAN is supported.

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the detector is connected.



If the detectors S 2520 or S2550 should be controlled by RS-232, in the detector's own setup the "Interface" must be set to "RS232 115200". Please refer to the detector's manual, chapter "Setup menu" for more information.



The detectors 2550, 2520 and 10D/UDV2.1S have a limited data rate (max. 10 Hz) with a serial (RS-232) connection. The maximum data rate for a detector is only available via LAN connection for data acquisition, for the detector S 2550 only on one channel.



The detector S 2600 must be set in its own setup in "CONTROL" to "RS232 binary baud rate 115200".

LAN




Fig. 74 LAN interface dialog

All LAN-controlled devices can be added and configured automatically using the **Auto Configuration** button on the configuration main window. Beside this, you can let the software search for your already connected and switched-on device from the devices configuration window. The option "Use S/N to identify the instrument" must be enabled. Click on the  button to start the search for all connected LAN devices of the currently configured type. If more than one will be found, all hits will be shown in the window. Select the desired device from the shown list. The software will read-out the information as serial number and will fill it automatically into the configuration screen. For all Knauer LAN devices port 10001 is selected by default. Don't change this setting, otherwise no communication can be established. To search for an instrument, the device must be connected, switched on and must have a valid IP address.

As an option, you can enter the IP address of your device, if the option "Use S/N to identify the instrument" is disabled. Please note, that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.

Channels (available not for all detectors)

For multi-wavelength detectors more than one channel can be activated. According to the selected number, the channel fields become accessible for editing. The channel name will be shown later in the method window and in the result reports. Only channels, which have been activated in the detector configuration, will be available in the method setup later for data acquisition.

Serial Number

If you have read-out the device via LAN, the serial number has been added automatically. Otherwise type in the serial number of the detector you are using. The serial number of network controlled detectors is used for the device identification and addressing during communication. This serial number must match the serial number of the defined pump. If not, you will get an appropriate error message if you later try to start in run.

Y-Axis Units

Select the type of units to be displayed on the y-axis of your chromatogram from the drop-down list. For the UV and diode array detectors the choices are **AU** (Absorption Units), **mAU** (0.001 AU) or **µAU** (0.000001 AU). For the RI detectors, the only choice is **µRIU** and for the spectrofluorometric detector the units are given as **a.u.** (arbitrary units).



For DAD's you have additionally also the choice ADC counts. Selecting this, the signal mode in the Instrument setup (page 102) will be indicated as 'Intensity' without access for changing it. This change you can perform in the configuration window only.

Flow Cell

Enter a description for the flow cell you are using. The description will be printed out when you print the instrument configuration.

Add. Info

Click the button **Add. Info...** to enter more detailed description of the detector.




Fig. 75 Additional information box

Enter additional information or comments to be printed out when you print the instrument configuration.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Halogen Lamp (available not for all detectors)

Check this box when the detector is equipped with a halogen lamp. Also if you have read-out a detector equipped with halogen lamp, the option will not be enabled automatically. The detectors S 2600 and S 2520 can be equipped either with a deuterium or a halogen lamp. Please be sure that the correct lamp type is selected.

Leak Sensor Sensitivity (UVD2.1L only)

Select the leak sensor's sensitivity in accordance with your requirements. The choices are **Off**, **Low**, **Medium** and **High**.

Shutter Control (DAD, S 2600 only)

For a shutter-equipped detector, check this box if you want dark current correction to be carried out. If this option is chosen, the detector will close a shutter before every run, measure a dark current, and then subtract it from all measured spectra. Additionally the shutter will prevent the fiber optics for going blind.

Port Settings (Spectrofluorometric detector RF-10AxI / RF-20A/Axs only)



Fig. 76 Port setting dialog box

Select the settings from the pull down menus as shown in Fig. 76:

Baud Rate	9600
Data Bits	8
Parity	Even parity or No parity
Stop Bits	1



The default value for parity in ChromGate® is "Even". For RF-20A/Axs the parity default setting in the device currently is "No". Please check in the device for the correct setting. Refer to the description below how to do this.



A proper connection normally is only possible via onboard COM ports. The port settings must coincide with the RS-232 parameters on the detector. To control the detector via RS-232, the detector must be set into serial mode as described below (also refer the detector's manual).

Before you open the ChromGate instrument you must switch the detector in to serial mode.

RF-10 Axl: Press the function key <func> on the detector's keypad as often as "RS232C" will be displayed. Confirm with "E" key. The display will show "CONNECT", also confirm with "E" key.

If you press the "function" while "CONNECT" will be displayed, you will enter to the COM connection table. Changed values must be confirmed by pressing the key "E", to move to the next parameter, press the function key.

Parameter	Value
Baud	6 (means 9600)
Data Bits	8
Parity	2 (means even parity)
Stop Bit	1

The COM settings here and in the software must be identical. If, e.g., the parity is set to "0", means no parity, it must be set in the detector's setup as well in the software.

If all parameters have been entered, the display shows the message "CLOSE KEY". Press key "E" to leave the setup. To switch the device to serial control mode, please follow the steps above

To leave the serial mode, press the shift and CE/del key at the same time.

RF-20A/Axs: Press the function key <func> until the display shows "SYSTEM _ Enter to Select". Press the <enter> key. Press function key until the display shows "RS-232 PARAM _ Enter to select". Press the <enter> key. Now the settings for the serial interface can be reviewed and modified if necessary. To move to the next parameter, press the function key

Parameter	Value
Baud	9600
Data Bits	8
Parity	0 (means no parity)
Stop Bit	1

The COM settings here and in the software must be identical. If, e.g., the parity is set to "2", means even parity, it must be set in the detector's setup as well in the software.

If all parameters have been reviewed, the display shows "RS-232C _ Enter to Enable". Press the key <enter>. The display shows "RS-232C _ Enable" for only about 2 seconds, then the default screen will be shown. Now the detector is in serial control mode and can be controlled by software.

If the detector is switched Off and On, it will start again in RS-232 control mode. During the boot procedure the RS-232 control mode will be shown very shortly on the display.

To switch off the RS-232 control mode, the keys <CE> und must be pressed at the same time.



The driver for the model RF-20A/Axs only supports the firmware version 0.9 . This firmware is only available through Knauer. The control is compatible with model RF-10AxI and will not support the advanced features of the model RF20A/Axs.



The maximum data rate for the RF-10AxI using the serial port is 0.5 Hz, for the RF-20A/Axs 5 Hz. For higher data rates the 1 V Integrator Output with a Knauer interface must be used.



Both detectors may not work on additionally installed serial (COM) ports. It is recommended to connect it with a computer's onboard COM port.

Configuration – User Defined Detectors

The user defined detector can be assigned to any detector or sensor producing an analog signal which can be digitized with A/D converters or Interface Boxes. Such a detector could be any UV, RI, or fluorescence detector, or some type of pressure, flow, or temperature sensor, etc.

For each user defined detector in the Configured Modules window, double-click the icon and complete the configuration dialog.




Fig. 77 User defined detector configuration window

Detector Name

Enter a descriptive name for the detector. This name will appear on the detector tab when the instrument setup or status windows are displayed. Use a name which is **unique** within the instrument.

Detector Model

Select the model of the detector from the drop-down list (planned future option).

Interface

Select the A/D interface for the detector you have already configured (please see chapter "Knauer Interface Configuration" on page 32) and click the button to set up the acquisition channel.



Fig. 78 Acquisition channel configuration window

Select the **Device number** to which the detector is connected. The available choices depend on the global interface configuration as described on page 32.

Select the acquisition **Input** of the interface device.



It is not possible to use a channel twice in a system. A data acquisition channel cannot be used for control commands such as wavelength changes.

Click the **Add. Info...** button to enter additional information.

Y-Axis Units

Fill in the units of your signal to be displayed on the y-axis, e.g. microvolt, AU, MPa, etc., depending on the detector type. An interface device acquires your analog signal and stores it in microvolt. If you want to store and to display the signal in different units, you must use the correct multiplier.

Y-Axis Multiplier

ChromGate by default stores a signal in microvolt. If you wish to display another unit, enter the conversion factor to be applied. The following table shows commonly used y-axis labels and their corresponding multipliers.

Y-Axis Label	Y-Axis Multiplier
Volts	0.000001
Millivolt	0.001
Microvolt	1
Other	Enter a number such that when multiplied by microvolt, gives you the desired units.

For example, if you want to label the units, and each unit corresponds to 5 microvolt of signal, enter $1/5 = 0.2$. Please also note that the Multiplier also will be used for the Y-Axis Offset.

Y-Axis Offset

Type in the correction offset which will be used for the signal recalculation. The Y-Axis offset **must** be given in the same units as the Y-axis is labeled. Consider the formula:

$$DU = MS * Y-M + Y-O$$

DU = displayed units
 MS = measured signal in micro Volt
 Y-M = Y-Axis Multiplier
 Y-O = Y-Axis Offset



Please note that the Offset also will be calculated using the Y-Axis Multiplier. As an example: To change the Offset by 2 mV if the Multiplier is set on 0.001, a value of 2000 must be entered for the Offset. Negative Multipliers are allowed.

Add. Info

Click the button **Add. Info...** to enter more detailed description of the detector.

Enter additional information or comments to be printed out when you print the instrument configuration.

When complete, click OK to exit the dialog and return to the instrument configuration icon list.

Configuration – Kontron UV-Detectors




Fig. 79 Kontron detector configuration windows

The Kontron detectors are controlled similar to the Knauer detectors via the configured RS232 port. For the data acquisition, the analog output with an A/D card or box must be used.

Name

Enter a descriptive name for the detector. This name will appear on the detector tab when the instrument setup or status windows are displayed. The default detector name will be suggested automatically by selecting the detector type. The name must be unique within the instrument.

Type

Select the type of your detector.

Serial Port

Select the number for the communication port on your PC where the instrument is connected. The selection <none> is for simulations. No real control of the instrument will be performed.

Serial Number

Type in the serial number of the detector you are using.

Flow Cell

Enter a description for the flow cell you are using. The description will be printed out when you print the instrument configuration.

Data acquisition interface

The data acquisition is performed as it was described for the user defined detectors via the interface box or the Manager 5000/5050. Select the connected interface from the pull down menu. Most of the Kontron detectors do not support data acquisition via serial port; therefore the A/D converter is required.

Channel 1 Input:

Select the connected data channel of your interface.

Channel 2 Input:

Select the second connected data channel of your interface. It can be deactivated by the selection <none>. The Kontron x32 detectors are automatically set to <none>.

Y-Axis Units

Select the type of units to be displayed on the y-axis of your chromatogram from the drop-down list. The choices are **AU** (Absorption Units), **mAU** (0.001 AU) or **µAU** (0.000001 AU) and **ADC counts**.


Configuration – Kontron Diode Array Detectors

Fig. 80 Kontron DAD configuration window

The Kontron PDA detectors are controlled similar to the Knauer detectors via the configured RS232 port.

Name

Enter a descriptive name for the detector. This name will appear on the detector tab when the instrument setup or status windows are displayed. The name must be unique within the instrument.

Serial Number

Type in the serial number of the detector you are using.

Interface

For the diode array detectors, the choice is **PCI** and **LAN**. If you control the DAD via the LAN, you must enter the IP Address. The access to the IP Address field is blocked if PCI was selected.

Y-Axis Units

Select the type of units to be displayed on the y-axis of your chromatogram from the drop-down list. The choices are **AU** (Absorption Units), **mAU** (0.001 AU) or **µAU** (0.000001 AU) and **ADC counts**.



Selecting the ADC counts, the signal mode in the Instrument setup (page 102) will be indicated as 'Intensity' without access for changing it. This change you can perform in the configuration window only.

Flow Cell

Enter a description for the flow cell you are using. The description will be printed out when you print the instrument configuration.

Channels

Up to four channels can be activated. According to the selected number the fields **Chn.Name #1-4** become accessible for editing.

Halogen Lamp

Check this box when the detector is equipped with a halogen lamp.

Configuration – Virtual Detector

To configure a virtual detector you have only to open and to close the configuration window without any additional entering. The Y-axis multiplier is only to change if very weak signals are to simulate.




Fig. 81 Virtual detector configuration window



To avoid possible signal distortion, the Y-Axis multiplier in the configuration of the virtual detector ***must*** be equal to that of the selected data file.

Detector Connections

Up to four detectors can be connected to the Interface Box or Manager 5000/5050/IF2. The detector connections are principally identical. The following description for the Interface Box can also be used for the Manager 5000/5050/IF2, only the connector positions may differ.

Always use a shielded twisted pair cable to connect each detector to the




Fig. 82 Cable connections to the Knauer Interface Box, Model 96



The integrator inputs 5 and 6 as well as the I/O channels (WAGO connectors) 5 and 6 of the Interface Box are not accessible by ChromGate®.




Fig. 83 Cable connections to the Manager 5000/5050 or IF2

If you want your instrument to trigger the start of data acquisition or to perform an autozero, you must take care for making the trigger signal and autozero connections. The digital inputs and outputs used for the connections are situated on WAGO connectors adjacent to the Integrator Inputs. Each of the WAGO connectors corresponds to one of the data acquisition channels.

Trigger Signal Connections

If you want your instrument to trigger the start of data acquisition on your ChromGate system using the Interface Box, you must connect a trigger signal cable to the Interface Box. The trigger signal should be connected to one of the digital input connectors on the rear panel of the Interface Box. The digital input must correspond to that one selected in your method.

Connect one pole of the trigger cable to one of the connectors labeled D.IN1, and the other pole to the connector labeled D.GND (digital ground).

Autozero Signal Connections

If your detector has an autozero option and you want to use the feature to zero the baseline signal at the start of data acquisition, you must connect an autozero signal cable to the Interface Box. The autozero signal is connected to one of the digital output connectors on the rear panel of the Interface Box.

Connect the ground connector of the instrument to the D.GND (digital ground) labeled connector of the Interface Box and the AUTOZERO connector of the instrument to the connector labeled D.OUT3. The connectors D.GND and D.OUT3 must be situated on the same connector strip.




Fig. 85




Fig. 85 Wiring example Manager 5000/5050/IF2 and detector

These examples show the connection lines between the Interface Box or Manager 5000/5050 / IF2 and a detector (as an example the Knauer K-2501 detector). If a Knauer detector is connected via serial or LAN port, the autozero and the measurement signal will be transferred via this connection. The pictures should only show the principle of connection.

Configuration – Assistant ASM2.1L

The Azura assistant ASM2.1L is a modular instrument that allows combining up to 3 devices. The following devices can be included:

- pumps P2.1S, P4.1S
- detector UVD2.1S
- Knauer valve drives with 2, 6, 12 and 16 positions
- Valco valve drives with 2, 6, 8, 10, 12 and 16 positions

Configuration rules:

- Two pumps are only supported as HPG, both pumps must have the same pump head, three pumps are not supported.
- Only one UV detector is allowed.
- One valve drive can be used as a fraction collector, if the Knauer Fraction Collector Control option is installed and the appropriate license option is used. Cascading fraction valves is not supported. This is also applicable, if the valves are installed in different ASM2.1L housings.
- All the devices will be controlled by only one LAN port; a serial control via RS-232 is not implemented.
- Pumps cannot run in a HPG, if they are installed in different ASM2.1L housings.




Fig. 86 Assistant ASM2.1L configuration window, tab General

If one of the modules should not be controlled, it must be removed from the ASM 2.1L configuration. There is no option to disable a configured module later on in the Method Setup.

The configuration window has four tabs, one for general settings and one for each of the three device positions, left, middle and right. While the **General** tab allows for settings for the whole assistant, the device tabs give the option to configure the included devices.

General

Assistant Name

The name of the assistant ASM2.1L appears automatically. You can enter any other descriptive name. This name will appear on the assistant tab when the instrument setup or instrument status windows are displayed.



Use a name which is unique within the instrument (configuration).

Serial Number

If you have read-out the device via LAN, the serial number has been added automatically. Otherwise type in the serial number of the instrument.

Interface

The Assistant ASM2.1L supports LAN connection; only one LAN port for all 3 devices is required.



Fig. 87 LAN interface dialog

As all other LAN-controlled devices, the assistant ASM2.1L can be added and configured automatically using the **Auto Configuration** button on the configuration main window. All included device will be recognized and configured. Beside this, you can let the software search for your already connected and switched-on ASM2.1L from the devices configuration window. The option "Use S/N to identify the instrument" must be enabled. Click on the  button to start the search for all ASM2.1L instruments.

Leak Sensor Sensitivity

Select the leak sensor's sensitivity in accordance with your requirements. The choices are **Off**, **Low**, **Medium** and **High**.

Add. Info

Click the button **Add. Info...** to enter more detailed description of the pump.




Fig. 88 Additional information box



If you have already read-out the assistant instrument using <Auto Configuration> or the  -button, you may not configure your assistant instrument manually. This is the recommended approach. The automatic configuration eliminates the risk of a faulty configuration. Nevertheless, the manual configuration will be described below.

Left Device

Due to that for all three device positions the same devices and settings can be selected, the options will only be described for the Left Device tab.




Fig. 89 Configuration windows left device tab, device type selection

Device Type

Select the device type that is installed on the left side in the assistant housing, from the drop-down menu.

Pump P2.1S

If selected, you can choose the device name (please refer to the device name section below), the gradient type (none, HPG A, HPG B) and the pump head (10 ml, 50 ml).

Pump P4.1S

If selected, you can choose the device name (please refer to the device name section below), the gradient type (none, HPG A, HPG B) and the pump head (10 ml, 50 ml). Due to that the pump P4.1S is equipped with a pressure sensor, you can also select the unit for the pressure monitoring (MPa, bar, psi).

Detector UVD2.1S

If selected, you can choose the device name (please refer to the device name section below) and the Y-Axis Units (μ AU, mAU, AU). The "Flow Cell" field allows for entering information regarding the used detector flow cell.

Knauer Valve

If selected, you can choose the device name (please refer to the device name section below) and the number of position. Please be sure to select the correct number of position, if you select it manually. 2, 6, 12 or 16 position can be selected. For one valve in the assistant you can check the "**Use as a Fraction Collector**" box. If selected, there are additional options: You can enter the **Vial Volume** (please be sure to enter the correct volume) and setup **Tubing** and **Solvent/Peak Recycling**. Please refer to the chapter **Multi Valve Fraction Collector Configuration** for detailed information. The **Number of Vials** is fixed and inaccessible, depending on the selected number of positions of the used valve and the Solvent/Peak Recycling setup. If for another valve in the assistant the option "**Use as a Fraction Collector**" is already selected, for the next valve an error message "The fraction collection option is already selected for another valve" will be shown.

Valco Valve

If selected, you can choose the device name (please refer to the device name section below) and the number of position. Please be sure to select the correct number of position, if you select it manually. 2, 6, 8, 10, 12 or 16 position can be selected. For one valve in the assistant you can check the **“Use as a Fraction Collector”** box. If selected, there are additional options: You can enter the **Vial Volume** (please be sure to enter the correct volume) and setup **Tubing** and **Solvent/Peak Recycling**. Please refer to the chapter **Multi Valve Fraction Collector Configuration** for detailed information. The **Number of Vials** is fixed and inaccessible, depending on the selected number of positions of the used valve and the Solvent/Peak Recycling setup. If for another valve in the assistant the option **“Use as a Fraction Collector”** is already selected, for the next valve an error message “The fraction collection option is already selected for another valve” will be shown.

Device Name

The default name corresponds to the selected device type. All device names must be unique in an instrument (configuration). Therefore, enter different names for the same device types, e.g., if a high pressure gradient is used. If the configuration has been read-out by **Auto Configuration** or using the  -button, the software will name devices from the same type automatically with different names. For valves a letter in regard of the valve’s position in the assistant will be added, e.g. a valve in left position will be named as “ValveL”. Pumps will get an A or B increment for HPG A or B pump. However, you can edit the automatically given name as well as the HPG position for the pumps.

Middle Device

Due to that for all three device positions the same devices and settings can be selected, the options will only be described for the Left Device tab.

Right Device

Due to that for all three device positions the same devices and settings can be selected, the options will only be described for the Left Device tab.



If one of the modules should not be controlled, it must be removed from the ASM 2.1L configuration. There is no option to disable a configured module later on in the Method Setup.

Configuration – Autosampler

For each system only one autosampler will be accepted. Select the autosampler in the Configured Modules window, double-click the icon and complete the configuration dialog. For configuration of the **Injection Module**, please refer the FRC section, chapter “Stacked Injections”.



For configuration information regarding the PLATINblue autosampler AS-1, please refer to the PLATINblue system manual.

Configuration – Autosampler 3950

Smartline Autosampler 3950 (Alias)




Fig. 90 Configuration window for the autosampler 3950

Name

Enter a descriptive name for the sampler. This name will appear on the autosampler tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

You must type in the serial number from the manufacturer's label on the rear panel of the autosampler (**not** the Knauer S/N).

Interface

Select the Interface you want to use for communication with the device. Please note, that the autosampler 3950 can only have RS-232 or LAN interface. The communication port on the autosamplers rear panel is labeled as "Communication".



Fig. 91 Available interfaces Smartline autosampler 3950

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the autosampler is connected.

LAN




Fig. 92 LAN interface dialog

All LAN-controlled devices can be added and configured automatically using the **Auto Configuration** button on the configuration main window. Beside this, you can let the software search for your already connected and switched-on device from the devices configuration window. The option "Use S/N to identify the instrument" must be enabled. Click on the  button to start the search for all connected LAN devices of the currently configured type. If more than one will be found, all hits will be shown in the window. Select the desired device from the shown list. Configuration

information as syringe, loop and tubing volume as well as the tray will be added with default settings. Please check, if these settings match with the installed option and change it, if necessary. All configuration data will be downloaded into the autosampler, if the instrument will be open. As an option, you can enter the IP address of your device, if the option "Use S/N to identify the instrument" is disabled. Please note, that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.

Loop Volume

Enter the volume of the installed loop in μl . If the Prep mode option is checked, the loop volume is fixed on **10000 μl** .

Tubing Volume

Type in volume of the installed needle and tubing in μl . The valid range is **0 - 999 μl** . If the original tube is installed, you can find the tubing volume on a label on the tube.

Syringe Volume

The volume of installed syringe can be set on **250, 500 or 1000 μl** . If the Prep mode option is checked, the syringe volume is fixed on **2500 μl** .

Options

Tray Cooling

Check this box if you want to control the tray cooling. The appropriate corresponding temperature is to be entered later as part of the instrument setup.

SSV (Solvent Selection Valve)

Check this box if you have a solvent selection valve installed.

Prep mode

Check this box if you have the preparative instrument installed. The loop and syringe volumes then are set to a fixed volume (10,000 μl loop, 2,500 μl syringe) and the tray type to the 12 vials type.

Tray...

Click the Tray... button to set up the autosampler tray:




Fig. 93 Tray configuration for the autosampler 3950




Fig. 94 Tray configuration 84+3 vials for the autosampler 3950



The “84+3 vial” tray for the autosampler 3950 is supported in firmware version 1.14 or higher. A firmware upgrade from former firmware versions is not available.

Left (Right) tray type

Select the type of the installed tray for both sides. The choices are **96 low, 96 high, 384 low, 48 vials, and 12 vials**. The selection can be different for each tray side. In case of activated prep mode it is fixed on 12 vials. If the tray type **84+3 vials** is enabled, this choice is for both side; this tray uses the left and right tray side. Only for the tray type **84+3 vials** it is possible to select a source vial for the transport liquid in the method setup later. The 84+3 vials tray is only available for autosamplers with **firmware version 1.14** or higher. Autosamplers with a previous firmware version cannot be updated.

Plate process order

The vials of the two trays will be processed either in columns or rows order.

First Destination Vial

Enable this option if you want to use destination vials. Enter the first destination vial position. For each injection a destination vial will be used for mixing a diluting.

Reagent A(B) Vial

Enable these options if you want to use reagent vials. Enter the corresponding vial position. The selected vial will be used for **Reagent A** and/or **Reagent B** for all injections.

Use Mix & Dilute

Enable this option to select for each injection one destination, reagent A and reagent B vial (position). This allows a higher number of injections if a higher volume of the reagent(s) is required.

The Mix & Dilute option is only available if you have selected the same tray type for left and right. It is not available for the tray type **84+3 vials**.




Fig. 95 Tray configuration for the autosampler 3950 with option Use Mix & Dilute

For μl pick-up no transport vial must be chosen. For the tray type **84+3 vials** the transport vial can be selected in the method setup, for all other tray types the transport liquid will be taken from the bottle with wash solution.



All configuration settings will be transmitted to the autosampler if the instrument will be opened. In difference to the other Spark-manufactured autosamplers, the configuration data will not be stored permanently into the device and compared with your settings. However, please make sure, that the configuration settings are correct due to that the device will take over your settings without verification.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Configuration – Autosampler Knauer Optimas




Fig. 96 Configuration window for the autosampler Knauer Optimas

Name

Enter a descriptive name for the sampler. This name will appear on the autosampler tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the detector you are using.

Device ID

Select the device identifier for the autosampler from the drop-down list. This must match the settings of the autosampler being used. To check the device identifier for the autosampler Knauer Optimas, press the SYSTEM, MENU and then SERIAL buttons on the autosampler keyboard. The default identifier is 60.

Serial Port

Select the number for the communication port on your PC where the autosampler is connected. On the autosamplers rear panel use the port S2 IN, to connect the autosampler with the computer. It is recommended to switch all DIP switches for the port to OFF position for a proper communication.

Loop Volume

Enter the volume of the installed loop in μl . If the Prep mode is enabled, the loop volume is fixed on 10000 μl .

Tubing Volume

Type in volume of the installed needle and tubing in μl . The valid range is 0 - 999 μl . The volume is noted on a label on the tube.

Syringe Volume

Select the volume of installed syringe from the drop-down list. The choices are 250 μl , 1000 μl , and 2500 μl . If the 2500 μl syringe is selected, automatically the **Prep mode** check box is enabled. In the Prep mode the values for the syringe (2500 μl) and loop (10000 μl) and the 24 vial tray type are fixed.

Options (Tray cooling / Prep mode)

Check the appropriate boxes if you want to control the **tray cooling** or enable the **Prep mode**. The desired tray temperature can be entered later as part of the instrument setup. In the Prep mode the values for syringe (2500 μl) and loop (10000 μl) and the 24 vial tray type are fixed. The tubing volume is set to 45 μl . Please check the correct volume of the tube between the needle and the injection valve on the tube's label.

Tray...

Click the Tray button to set up the autosampler tray:




Fig. 97 Tray configuration for the autosampler Knauer Optimas

Tray Type

Select the type of the installed tray. The choices are **84+3**, **96** and **24 vials**. In Prep mode the tray type is fixed on 24 vials tray.

Vial Type

Select the type of the vials used. The choices are **Standard** and **2.5 ml**. The last is used only when <84+3> or <96> tray type is selected.

First Destination Vial

Enable this option if you want to use destination vials.

First Destination Vial Position

Enter the first destination vial position.

When complete, click **OK** to exit the dialog and return to the instrument configuration window.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.



The autosampler Knauer Optimas must be set in to serial mode to be controlled by software. Press the <Serial> button to switch the autosampler into serial mode.

Configuration – Autosampler 3800

Smartline Autosampler 3800 and K-3800 (Basic Marathon)




Fig. 98 Configuration window for the autosampler 3800

Name

Enter a descriptive name for the sampler. This name will appear on the autosampler tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Device ID

Select the device identifier for the autosampler from the drop-down list. This must match the settings of the autosampler being used. To check the device identifier for the Smartline Autosampler 3800 or the K-3800 (Basic Marathon) press the PROG/END button on the autosampler keyboard once in serial mode.

Serial Port

Select the number for the communication port on your PC where the autosampler is connected.

Serial Number

Type in the serial number of the detector you are using.

Loop Volume

Enter the volume of the installed loop in μl .

Sample Tray

Select the type of the installed sample tray. The choices are **96 vials of 1.5 ml** and **48 vials of 5.0 ml**.

Tray cooling

Select this box if you have a built-in Peltier cooling option installed and want to cool down your samples. The temperature setting will be done within the instrument setup.

Partial loop fill

Check this box to select the partial loop fill injection mode (for the Basic Marathon Plus only).

Add. Info

Click the button **Add. Info...** to enter a more detailed description of the autosampler. Enter additional information or comments to be printed out when you print the instrument configuration.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.



The autosampler 3800 must be set in to serial mode to be controlled by software. Press the buttons <F> and <4> on the autosamplers keypad, the display will show then “Serial Mode Press F0 to exit”.

Configuration – Autosampler 3900

Smartline Autosampler 3900 and K-3900 (Midas)




Fig. 99 Configuration window for the autosampler 3900

Name

Enter a descriptive name for the sampler. This name will appear on the autosampler tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the detector you are using.

Device ID

Select the device identifier for the autosampler from the drop-down list. This must match the settings of the autosampler being used. To check the device identifier for the autosampler 3900, press the SYSTEM, MENU and then SERIAL <SETTINGS> buttons on the autosampler keyboard.

Serial Port

Select the number for the communication port on your PC where the autosampler is connected.

Loop Volume

Enter the volume of the installed loop in μl .

Tubing Volume

Type in volume of the installed needle and tubing in μl . The valid range is 0 - 999 μl . If the original tube is installed, you can find the tubing volume on a label on the tube.

Syringe Volume

Select the volume of installed syringe from the drop-down list. The choices are 250 μl , 1000 μl , and 2500 μl .

Options (Column oven / Tray cooling)

Check this box if you want to control the column thermostat and/or the tray cooling. The appropriate corresponding temperatures are to be entered later as part of the instrument setup.

Tray...

Click the Tray button to set up the autosampler tray:




Fig. 100 Tray configuration for the autosampler 3900

Tray Type

Select the type of the installed tray. The choices are **84+3**, **96** and **24 vials**.

Vial Type

Select the type of the vials used. The choices are **Standard** and **2.5 ml**. The last is used only when <84+3> or <96> tray type is selected.

First Destination Vial

Enable this option if you want to use destination vials.

First Destination Vial Position

Enter the first destination vial position.

When complete, click **OK** to exit the dialog and return to the instrument configuration window.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Configuration – Autosampler Triathlon / Endurance

The configuration windows of the autosamplers Triathlon and Endurance are identical. The only difference is in the tray setup window, accessed by pressing the **Tray...** button.




Fig. 101 Configuration windows for the autosamplers Triathlon and Endurance

Name

Enter a descriptive name for the sampler. This name will appear on the autosampler tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the autosampler you are using.

Device ID

Select the device identifier for the autosampler from the drop-down list. This must match the settings of the autosampler being used. To check the device identifier for the Triathlon or Endurance autosampler, press the SYSTEM, MENU and then <COMMUNICATION SETTINGS> buttons on the autosampler keyboard.

Serial Port

Select the number for the communication port on your PC where the autosampler is connected.

Loop Volume

Enter the volume of the installed loop in μl .

Tubing Volume

Enter the volume of the installed needle and tubing in μl . If the original tube is installed, you can find the tubing volume on a label on the tube.

Syringe Volume

Select the volume of the installed syringe from the drop-down list. The choices are **100 μl** , **250 μl** , **500 μl** , **1000 μl** , and **10000 μl** . The 10000 μl syringe can only be selected when a user defined injection method is used.

Solvent Selection Valve (SSV) Option

Check this box if you have a solvent selection valve installed.

Tray Cooling

Check this box if you want to control the tray cooling. The setting of the appropriate temperature will be done later as part of the instrument setup.

Tray...

Click the Tray button to set up the autosampler tray. The tray configuration dialog will depend on the autosampler type (see below).

The Triathlon Trays

Since the Triathlon can be equipped with a mix of four available types of tray segments, it is necessary to define the tray configuration.




Fig. 102 Tray configuration for the autosampler Triathlon

Segment

The Triathlon tray consists of eight segments of different types.

Segment Type

Select the type of the corresponding installed segment. The choices are **Standard (A)** segment, **LSV (B)** segment, **Super LSV (C)** segment and **Micro vial (D)** segment. The selected type will be shown in the tray configuration window.




Fig. 103 Tray configuration, displaying different tray types

Segment Number

Select the number (from 1 to 8) of the corresponding installed segment. Duplicate segment numbers of one type are not allowed. However, if different segment types are involved, each of them shall be counted starting with 1.

Transport Vials

Enable this box if you want to use μ l pick-up injections with transport solvent from as many vials as necessary.

<input checked="" type="checkbox"/>	Transport vials		
		1	Standard (A)
		3	
<input checked="" type="checkbox"/>	Reagent A	4	Standard (A)
<input checked="" type="checkbox"/>	Reagent B	3	Standard (A)
<input type="checkbox"/>	Reagent C	1	Standard (A)
<input type="checkbox"/>	Reagent D	1	Standard (A)
<input checked="" type="checkbox"/>	Destination	1	LSV (B)

Fig. 104 Tray configuration; transport, reagent, and destination vials

First Transport Vial

Enter the first transport vial position.

Last Transport Vial

Enter the last transport vial position. All the vials from the first up to the last transport vial must be of the same type.

Transport Vial Segment

Select the type of the segment where the transport vials are installed. The choices are **Standard (A)** segment, **LSV (B)** segment, **Super LSV (C)** segment and **Micro vial (D)** segment. All the vials with transport solvent must be from one segment type, corresponding to the segment definition in the upper part of the dialog.

If all conditions for the transport vial definition are met, the selected vial positions will appear in Fig. 102 highlighted in yellow.

Reagent A/B/C/D / Destination Vial

Enable the corresponding boxes if you want to use the destination vial or vials with Reagent A/B/C/D.

Reagent A/B/C/D / Destination Vial Position

Select the Reagent A/B/C/D / destination vial position.

Reagent A/B/C/D / Destination Vial Segment

Select the type of the segment where the corresponding vials are installed. The choices are **Standard (A)** segment, **LSV (B)** segment, **Super LSV (C)** segment and **Micro vial (D)** segment.

The Endurance Trays

Tray configuration		1	2	3	4	5	6	7	8	
Plate type:	48 vials									
Plate process order:	Rows									
First transport vial:	No									
Last transport vial:	No									
First destination vial:										
Reagent A vial:	No									
Reagent B vial:	No									
Reagent C vial:	No									
Reagent D vial:	No									
		A	1	2	3	4	5	6	7	8
		B	9	10	11	12	13	14	15	16
		C	17	18	19	20	21	22	23	24
		D	25	26	27	28	29	30	31	32
		E	33	34	35	36	37	38	39	40
		F	41	42	43	44	45	46	47	48
		OK Cancel								

Fig. 105 Tray configuration for the autosampler Endurance

Plate Type

Select the type of the installed plate. The choices are **96-low wells**, **96-high wells**, **384-low wells** and **48 vials**.

Plate Process Order

Select the order of vial processing. The choices are in **Rows** and in **Columns**.

First Transport Vial

Select the first transport vial position. The vial can be placed in any of the positions 1 through 4. Select **No** if the transport vials are not used.

Last Transport Vial

Select the last transport vial position (≤ 4 and \geq first transport vial position). The transport vials must be placed in a continuous row. Select **No** if the transport vials are not used.

First Destination Vial

Select the first destination vial position. Leave the field blank if no destination vial is used.

Reagent A/B/C/D Vial

Select vial position 1 of the vial positions 1 through 4 for the Reagent A/B/C/D. Select **No** if the reagent vial is not used.

Configuration – Kontron Autosamplers

The configuration windows of the Kontron autosamplers 360, 460, 465, 560, and 565 are identical. The only difference is the vial number (see below).




Fig. 106 Configuration windows for the Kontron autosamplers

Name

Enter a descriptive name for the sampler. This name will appear on the autosampler tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the autosampler you are using.

Type

Select the AS type 460 or 465 respectively 560 or 565. In the 300 series only 360 is possible.

Serial Port

Select the number for the communication port on your PC where the autosampler is connected.

Port settings

In the dialog box for the port settings only the baud rate will be active for selection of 4800 Or 9600.




Fig. 107 Port settings

Vial number

For the AS 360, 560, and 565 a number in the range from 1 through 65 can be selected. In a method only the number of vials entered here are available. A wrong number will lead to an error message:




Fig. 108

Vial type

Select the volume of the used vials from the drop-down list. The choices are **1.8 ml**, **0.5 ml**, and **0.3 ml**.

Syringe Volume

Select the volume of the installed syringe from the drop-down list. The choices are **250 µl**, **500 µl**, **1000 µl**, and **2500 µl**.

Loop Volume

Select the volume of the installed loop from the drop-down list. The choices are **20 µl**, **100 µl**, **500 µl**, and **1000 µl**.

Configuration – Miscellaneous Instruments

Configuration – Switching Valves

The switching valves icon refers to a group of the Knauer switching valves; the Valco/Vici valves icon to the latest Vici valve drivers. Each of them must be configured before using in an instrument method.

Double-click the icon and complete the configuration dialog.




Fig. 109 Configuration window for the Switching Valves

Valves Number

Enter the number of valves used in a group. In maximum 32 valves can be configured in this configuration window.

Configure Valve

Select the number of the valve you would like to configure. You must configure all of the valves [Valves Number] before leaving the dialog.

Valve Name

Enter a descriptive name for the selected valve. This name will appear on the valves tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Valve Type

Select the type of the valve from the drop-down list. The choices are **2-**, **6-**, **12-** and **16-positions**, for the Valco/Vici valves additionally **8-** and **10-positions**.

Interface

Select the Interface you want to use for communication with the device.



Fig. 110 Available interfaces Knauer Switching Valves

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the valve is connected.

LAN (Knauer Switching Valves only)




Fig. 111 LAN interface dialog

To let the software search for your already connected and switched-on device, click on the  button to start the search for all connected LAN devices of the currently configured type. The option “Use S/N to identify the instrument” must be enabled. Select the desired device from the shown list. The software will read-out the information as serial number and will fill it automatically into the configuration screen. To search for an instrument, the device must be connected, switched on and must have a valid IP address.

As an option, you can enter the IP address of your device, if the option “Use S/N to identify the instrument” is disabled. Please note that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.

Add Valve Info...

Click the button **Add Valve Info** to add the valve specific description or comments for the selected valve.

Add. Info

Click the button **Add. Info...** to enter more detailed description of the valve group. Enter additional information or comments to be printed out when you print the instrument configuration.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Configuration – Manager 5000/5050, IF2

As with user defined detectors, the Manager 5000 / IF2 device must be configured to control any user defined instrument (necessary for programmed wave length control, pump control, switching valve control). Please note, that the Manager 5050 must be configured as a Manager 5000.




Fig. 112 Configuration window for the Manager 5000/5050 / IF2

I/O Name

Enter a descriptive name for I/O of the manager. This name will appear on the manager tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Interface

Select the used interface type Manager 5000 / IF2 (Manager 5000 if you own a Manager 5050).

Device number

Select the device number. The available choices depend on the global interface configuration as described on page 32.

Output

Select the channel of the interface in use.



It is not possible to use the cinch connector and the WAGO connector from the same channel for different instruments. They are not independent of each other.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Configuration – Column Oven 4050, Column Oven Jetstream

Select the Smartline column oven 4050 or the Jetstream in the Configured Modules window, double-click the icon and complete the configuration dialog.



Fig. 113 Configuration window for the column oven 4050

Name

Enter a descriptive name for the oven. This name will appear on the oven tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the device you are using.

Temperature units (column oven 4050 only)

Select from the drop down list either °C or °F.

Interface

Select the Interface you want to use for communication with the device.



Fig. 114 Available interfaces column oven 4050

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the column oven is connected.

LAN




Fig. 115 LAN interface dialog

To let the software search for your already connected and switched-on device, click on the  button to start the search for all connected LAN devices of the currently configured type. The option “Use S/N to identify the instrument” must be enabled. Select the desired device from the shown list. The software will read-out the information as serial number and will fill it automatically into the configuration

screen. To search for an instrument, the device must be connected, switched on and must have a valid IP address.

As an option, you can enter the IP address of your device, if the option "Use S/N to identify the instrument" is disabled. Please note that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Configuration – Flowmeter

Select the flowmeter in the Configured Modules window, double-click the icon and complete the configuration dialog.




Fig. 116 Configuration window for the Knauer and GJC flowmeter

Name

Enter a descriptive name for the flowmeter. This name will appear on the flowmeter tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the flowmeter you are using.

Serial Port

Select the number for the communication port on your PC to which the flowmeter is connected.

Y-Axis Units

The y-axis units are given as ml/min.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.




Fig. 117 Configuration window for the Bronkhorst flowmeter

Name

Enter a descriptive name for the flowmeter. This name will appear on the flowmeter tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the flowmeter you are using.

Serial Port

Select the number for the communication port on your PC to which the flowmeter is connected.

Capacity Mode

Select the desired capacity mode from the drop-down menu. The choices are "Volume Flow" and "Mass Flow". The Y-Axis unit will change dependently from the selected capacity mode.

Capacity Limit

Set the capacity limit which depends from the flowmeter model and the selected capacity mode / y-axis units.

Y-Axis Units

The available y-axis units depend from the selected capacity mode.

When complete, click **OK** to exit the dialog and return to the instrument configuration icon list.

Creating an Instrument Control Method

If you have the Instrument Control option installed and your instrument is configured as a KNAUER HPLC System, the instrument control functions appear as part of the Instrument Setup area of your instrument window.

There are two ways to access the Instrument window in which you can set up your devices and methods:

1. Double-click on the instrument icon you wish to start from the main screen. The window will be opened in the online mode.
2. Right-click on that icon and select the online or offline opening mode.

It may be required to log-in before you can access the instrument window functions.



Fig. 118 Login dialog box

Enter your name and your password and select the project from the pull down menu and then press the **Login** button.



Make sure that all devices of the selected instrument are switched on and connected to the described communication ports of the computer before opening the instrument window. The keypad of the column oven Jetstream must display “choose function”.

If a device is not switched on or connected to the computer if opening the instrument window please close the instrument window, switch or connect the device and reopen the instrument window. Otherwise the correct communication of the devices with the software can be affected.

All of the functions and procedures described in this chapter can be performed in the off line operation mode. The instrument wizard will appear the first time you open the instrument window. Its shape depends on whether you are working in offline or online mode.




Fig. 119 Instrument wizard, **A** offline, **B** online

For details on other aspects of method development (e.g. calibration, integration) and sequence operations, please refer to the ChromGate® Data System Reference Manual.

Instrument Setup

The instrument setup portion of your method displays how the devices (pumps, detectors, autosamplers, switching valves) are to be controlled and how data is to be acquired. This information is entered in the **Instrument Setup** dialog. Click the **Instrument Setup** button , or select *Method – Instrument Setup* from the menu. A window will appear which displays the parameters required for the equipment configured for this instrument.




Fig. 120 Instrument setup window (Smartline Pump 1050)

For each of the single configured devices you get a separate tab in the instrument setup window. Each of these contains several areas where you can enter parameters related to the individual device.

In any case you find at the bottom a unique graph part displaying the wavelength flow, gradient and possibly other related profiles.

With a right-hand mouse click you get access to multiple setup possibilities of the profile. For details see the ChromGate® Data System User's Manual.

The definition and setups of triggers (see page 135) and auxiliary traces (see page 135) must be performed in a separate tab for the whole instrument, not in the tabs of the single devices.

Instrument Setup – Pumps

The pump setup window shows several areas. Their appearance may depend as well on the selected pump as on the settings while configuring it.

Working Mode




Fig. 121 Working mode in pump setup window

Using the Kontron pump 325 additionally a CAM correction can be set. The correction value can be selected from 1 through 5.




Fig. 122 Working mode in Kontron 325 pump setup window

In case of isocratic or high pressure gradient no setting can be made. In LPG mode (except the Kontron x25 pumps) you can define the cycle time of the valves via the drop down list.

Control Pressure Limits




Fig. 123 Pressure limits in pump setup window

This section allows you to enter the control pressure limits in the units which have been chosen in the pump configuration dialog. The pump will automatically shut down its operation when the measured pressure is out of the limits. This option is not accessible in the case of the pump K-120, 10P, P2.1S or S 100 without pressure transducer.

Min Pressure

Enter the minimum pressure allowable for your pump. The zero setting disables the shut down due to minimum pressure.

Max Pressure

Enter the maximum pressure allowable for your pump. The default value always corresponds to the configured pump head. A setting above this value is not allowed.

If for the min and max limit "0" is set, the software will not check for the pressure and the limits will only be monitored by the pump internally.




Fig. 124 Error message when pressure limits are out of range




Fig. 125 Solvent description in pump setup window

The number of displayed rows depends on the configuration. Isocratic systems only show the component A. High pressure gradients show up to four components corresponding to the system configuration. In case of low pressure gradient pumps three (P2.1L, preparative pump 1800 and Kontron 42x) or four (S 1050, S 1000, K-1000, K-1001) components can be defined.



ChromGate® does not distinguish the two different LPG valve blocks for the preparative pumps P2.1L and 1800 (refer the pumps manual). The working limits of these blocks are: three solvents with a maximum summarized flow of 220 ml/min or two solvents with a maximum summarized flow of 1000 ml/min. Be sure to select the correct gradient mode for the connected valves, otherwise the pump or valves may be damaged.



In the K-501 pump's internal setup screen in the CONTROL menu, you should set the "Ext:RS232" mode to be able to control the pump from the ChromGate® software.

Using the Kontron pump 322 additionally the Compression stage can be selected. The choices are Low, Medium, and High.




Fig. 126 Solvent description in Kontron 322 pump setup window

Flow OFF at the end of run

Checking this option the flow will be set to zero after elapsing the runtime.

[Action] at the end of run (S 1050 / P2.1L only)

The pumps S 1050 and P2.1L allow not only setting the flow to zero, but also to switch the pump into Standby after the run.



Fig. 127 Possible action at the end of run for pump S 1050 / P2.1L

Pretreatment

This option enables several pretreatment operations of pumps and valves before injection and starting the data acquisition. For setting up see below on page 89.

Solvent Type... (S 1050 / P2.1L only)

Select in the Solvent Type menu the used solvent for a proper compression and flow correction. Note that for an LPG pump only one solvent type can be selected.




Fig. 128 Solvent type settings for pump S 1050 / P2.1L

Pump Program

Accordingly the spreadsheet for the pump program is displayed. Up to four columns for the composition will be displayed. Additionally, for pumps with digital outputs, the option for digital outs is given in a separate column. The spreadsheets of the S 100 / 10P/20P/P2.1S/P4.1S, K-120, K-501, and the Kontron pumps do not include this column.




Fig. 129 Pump program in pump setup window: **A** Pump 1800 LPG, **B** K-501 isocratic



For all time programs, the minimum time step between two lines is 0.02 min, except for the preparative pumps P2.1L / 1800; here the minimum time step is 0.06 min. For the pump P2.1L in isobar mode, the smallest time step is 0.25 minutes.

Use the spreadsheet in the pump setup view to enter your time (gradient) program. The number of columns shown depends on the defined system (isocratic, LPG, or HPG) and the number of included components.

In the first row (at time 0.00 min), enter the initial flow (in ml/min), solvent composition (in %), and digital events. Then for (time in minutes), you can define the pump parameters. For KNET pumps up to 50 steps are allowed with an accepted minimum time difference of 0.02 minutes between two lines.

The complete program is transferred or “uploaded” to the pump when you start a run.

All time programs are displayed in a spreadsheet similar to the one shown below. Although the information in the fields will vary, the spreadsheets always support certain basic features.

#	Time (min)	Flow (ml/min)	Comp. A (%)	Comp. B (%)	Comp. C (%)	Comp. D (%)	Events (0/1/P/_)	Comments
1	0.00	1.000	75	18	7	0	00000000	
2	2.00	1.000	75	18	7	0	00000000	
3	2.02	1.000	65	28	7	0	00000000	
4	10.00	1.000	100	0	0	0	00000000	

Fill Down
 Insert Line
 Delete Line
 Delete Lines

Fig. 130 Editing tools for the time (gradient) program

Each row is assigned a program line **number** and **time** in minutes, followed by columns for information for each time in the program. Rows and field information can be cut, copied, pasted, and cleared. To open the menu with these commands, right-click anywhere within the spreadsheet. To select a field, click on that field to highlight

To select a row, click on a number in the # column to highlight the entire row.

Certain spreadsheet fields will have choices available for you to select. Fields with selections available will display a combo-box button when the field or row is selected. Click the button to display the available choices.

Fill Down

This enables you to automatically copy spreadsheet information from one field or row down through the rest of the spreadsheet.

Insert Line

This command inserts a copy of the line in the spreadsheet below where the cursor is located.

Delete Line

This command deletes the line in the spreadsheet where the cursor is located.

Delete Lines

This command deletes the selected lines in the spreadsheet.

The flow and gradient composition profile is displayed in the profile window area below the spreadsheet.



For the HPG mode, the changes made for a gradient program will be automatically transferred to the gradient programs of all the pumps which form this HPG system. This does not affect the event configuration which can be setup for all pumps independently.

Configuring Events / Digital Outputs

This option is given only for the S 1050, P2.1L, S 1000, K-1001, and K-1800 pumps. The spreadsheets of the S 100 / 10P/20P/P2.1S/P4.1S, K-120, K-501, and the Kontron pumps do not include this column. Configure the events by pressing the  button which is displayed when you place the cursor on the **Events** cell. The following dialog will appear.



Fig. 131 Configuring Digital Outputs

The Events (Digital Outputs) correspond to the outputs available on the rear panel of the pump.

Select the state for each Event you want to set from the drop-down list. The choices are **Off** (0), **On** (1), and **Pulse** (P). Refer to the pump's user manual for details about digital outputs parameters.

Use the **All On** or **All Off** buttons to switch all of the outputs either on or off.

Constant pressure mode (Isobar) (P2.1L only)

If in the instrument configuration the gradient mode for the pump P2.1L is set to **Isobar**, the setup windows looks different:




Fig. 132 Instrument setup window P2.1L constant pressure mode (isobar)

Constant Pressure Mode: Flow Limits (ml/min): Set the allowed flow range for this operating mode. The allowed absolute minimum is 0 ml/min, the maximum depend on the selected pump head.

Pump program – Pressure (unit): Enter the desired system pressure. The pump will vary the flow, if the pressure in the system will change. The pressure unit depends on the unit selected in the pump's instrument configuration. The minimum time step between two time lines is 0.25 min.

Pretreatment (not for Kontron pumps)

Checking this option the **Setup...** button becomes active. Clicking on it the following dialog window will be opened:



Fig. 133 Pretreatment setting dialog window

The appearance of the pump program corresponds to that of the main pump setup window. It must be edited in the same way as described above. To add a new line, make a right mouse click into the table and select "Insert line" from the menu. The only difference is that **negative** times must be entered, because they are relative to the zero injection time.

During a pretreatment, no data acquisition is possible. However, the run status during a pretreatment is "Equilibrating".



If a pretreatment is to define for more than one instrument (pump and/or valve drive) in a method, it is required to set for all of them the same earliest time. All pretreatment procedures will start simultaneously.



If a pump should exactly switch off after 5 minutes of pretreatment, include a line with 0 ml/min flow into the pretreatment's pump program. Otherwise the pump may run for a longer time as it need some seconds to load the new flow rate from the method's "normal" pump program

An example for both :

A pump should run for exactly 5 minutes before the injection and a valve should switch 2 minutes before injection from position 2 to position 5, means, 3 minutes after the pump's pretreatment is started. Here are the setup tables create a time table for the valve, which starts with -5 minutes:

Pump

-5.04 min 1 ml/min

-0.04 min 1 ml/min

-0.02 min 0 ml/min (smallest time step is 0.06 minutes for pump 1800 / P2.1L or 0.02 min for all other pumps)

0.00 min 0 ml/min

Valve

-5.04 min pos. 2

-2.00 min pos. 5

0.00 min pos. 5

Instrument Setup – Detectors

All Detectors

In most cases you find at the bottom a unique graph part displaying the wavelength flow, gradient or other related profiles.

With a right-hand mouse click you get access to multiple setup possibilities of the profile. For details see the ChromGate® Data System User's Manual.



If the detector lamp is disabled, the event outputs cannot be controlled anymore; they will stay in the last set status.

The detector tab contains several areas where you can define the acquisition. Some of them are independent on type and model of the detector:

Acquisition channel on

Time constant: 0,1 Sec

Sampling rate: 1 Hz

Suitable for minimum peak width at base: 0,3330 Min

Run time: 10 Min

Acquisition delay: 0 Min

Fig. 134 Detector type independent settings

Acquisition channel on

Check the acquisition channel on box to enable data acquisition for the detector. This function is unavailable when using diode array detectors, the S 2550 or the K-2600.

Time Constant

Select the time constant for the detector from the drop-down list. The choices are **0.1, 0.2, 0.5, 1.0, 2.0, 5.0, and 10.0 sec**. The detectors UVD 2.1L, UVD 2.1S, S 2550, S2520 allow additional values (**0.000, 0.005, 0.01, 0.02, 0.05**, for LAN connection only). Smaller time constants are required for high speed analyses. A good rule is to select a time constant, which is the reciprocal of the data rate in Hz (e.g. 2 Hz data rate -> 0.5 sec time constant).

Sampling rate

This is the rate at which data will be sampled by the system. For any selected sampling rate the narrowest peak width will be shown, for which this sampling rate will be adequate. This selection is in Hz (data points per second). This is the selection for most chromatography applications. We recommend having more than 20 data points per peak (30 – 40). This means that for a peak with a peak width of 20 seconds (from peak start to peak end) a data rate of 2 Hz is recommended. Click on the down-arrow to get a list of the frequencies available for the configuration of your system. The choices are: **1, 2, 5, and 10 Hz**. The detectors UVD 2.1L, UVD 2.1S, S 2550, S2520 also support higher data rates **20, 50 and 80 Hz** (for LAN connection only, for S 2550 for one channel only). In case of the Interface Box Model 96 with activated 100 Hz option or the Manager 5000/5050 or IF2 with activated 50 Hz option (see page 33) you can additionally select **20, 50, or (Model 96 only) 100 Hz** for the configured User Defined Detector. This option is restricted to channel No. 1 if the interface while the other channels are not accessible.



It is recommended that you use the peak width's (see **Graphical Events Programming** or if included peak report) to determine the optimal sampling rate for your chromatography.

Run Time

Run time determines the length of time data will be sampled

Acquisition Delay

Acquisition delay is the interval between the start of the run (Trigger) and the time when sampling starts for this detector.

Instrument Setup – RI Detectors

(S 23[4]00, K-23[4]00/1)




Fig. 135 RI Detector setup window

The only additional setting is the optional:

Autozero at start

Check the **Autozero** field if you wish the signal to be zeroed at the moment a measurement is started (to define the current signal as the baseline).

Instrument Setup – UV Detectors (S 200, K-200)

The detectors K-200 and S 200 support settings that for all of the detectors described on page 91. The only additional setting is the optional

Autozero at start

Check the **Autozero** field if you wish the signal to be zeroed at the moment a measurement is started (to define the current signal as the baseline).

Because these detectors have a fixed wave length of 254 nm this setting is shown without access for changing.

Instrument Setup – UV Detectors (K-2000, K-2500)



In the detector's setup screen / EXT.CONTROL menu, you should set the keyboard mode to be able to control the wavelength from the ChromGate® software.

The detectors K-2000 and K-2500 do not support additional commands. Therefore, the detector tab gives no access to further settings than that for all of the detectors described on page 91.

Instrument Setup – UV Detectors (UVD 2.1L, UVD 2.1S, S 2520, S 2500, K-2001, K-2501)



For the detectors K-2001/2501 you should set in the detector setup screen / EXT.CONTROL menu the keyboard mode to be able to control the wavelength from the ChromGate® software.

The detector basic settings, as time constant, sampling rate etc. are described on page 91.




Fig. 136 UV Detector setup window

Autozero at start

Check the **Autozero** field if you wish the signal to be zeroed at the moment a measurement is started (to define the current signal as the baseline).

Autozero at WL change

No access to this option. The grayed area is only a reminder, that the autozero will take place at any wavelength change during a run.

Lamp OFF at the end of run

Check the lamp OFF box to switch off the detector lamp after a measurement is finished. This option may be used for the shutdown methods.



Detectors UVD 2.1L, S 2520 and 2550: If the lamp was switched off by the method and you switch off the detector then, on the next detector start the lamp is still off; it must be switched on manually in the detector's own setup or with the software via Instrument Status – Direct Control.

[Action] at the end of run (UVD 2.1L, UVD 2.1S, S2520 only)

The detectors S UVD 2.1L, UVD 2.1S and S 2520 allow not only to switch off the lamp at the end of run, but also to switch the device into Standby.



You should switch the lamp on manually before you start the next working session. After power up, allow 30 minutes for the lamp to stabilize before running samples that may require sensitive measurements.

Wavelength table

The wavelength table is a spreadsheet to enter your time (detector) program. In the first row (at time 0.00 min), enter the initial wavelength(s) in nm. The editing principle is the same as for the time (gradient) program, described in the chapter *Instrument Setup – Pumps, Pump Program*.

Configuring Events / Digital Outputs (UVD 2.1L, S2520 only)

This option is given only for the detectors UVD 2.1L and S 2550. Configure the events by pressing the  button which is displayed when you place the cursor on the **Events** cell. The following dialog will appear.




Fig. 137 Configuring Digital Outputs

The Events (Digital Outputs) correspond to the outputs available on the rear panel of the pump.

Select the state for each Event you want to set from the drop-down list. The choices are **Off** (0), **On** (1), and **Pulse** (P). Refer to the pump's user manual for details about digital outputs parameters.

Use the **All On** or **All Off** buttons to switch all of the outputs either on or off.



The standard version of the K-2001 detector is delivered with a deuterium lamp as well as with 4 filters for the wavelengths 200, 220, 254, and 280 nm. For the K-2001 you can only use these wavelengths in the time program. An optional tungsten-halogen lamp is also available. Two additional wavelength filters of your choice can also be built in. In this case you can use those corresponding wavelengths in the time program.



The Start Input of the detectors K-2501 and S 2500 will also be used to start an internal wave length program. If the detector is used for synchronization input when working under software control, the internal wave length program also will be started after receiving the synchronization signal (trigger). As a result of this fact, the wave length for the measurement may be different from the wave length programmed in the detector's instrument setup menu. Therefore, it is highly recommend to delete all program lines of the detector's internal wave length program (please refer to the detector's manual) and for the first line (time 0), which cannot be deleted, to set the wave length to the same value as in the software.

Refer to the section setting up a trigger on page 135 for details.

Instrument Setup – Kontron Detectors (3xx, 4xx, and 5xx)




Fig. 138 Kontron Detector setup window

Except the not available option Acquisition Channel On/Off the general detector settings as described on page 91 can be made. Additionally the range can be set to define the range of the analog detector outputs. The choice is **0.100, 0.200, 0.500, 1.000, and 2.000 Absorption Units Full Scale (AUFS)**.

The Kontron x30 and x35 detectors can be used in a dual channel mode. As far as the option Used Channel #2 is checked, a second wavelength column will appear in the spreadsheet below. Also the otherwise not accessible cycle time becomes active for changing. This defines the alternating measuring time for the two channels. The data acquisition is not performed truly simultaneously. The cycle time you can set within the range of 0.40 through 1.00 seconds with 0.05 s increments.

For the x32 detectors the option Used Channel #2 is not available and the cycle time is not displayed at all.

The wavelength table is a spreadsheet to enter your time (detector) program. In the first row (at time 0.00 min), enter the initial wavelength(s) in nm. The editing principle is the same as for the time (gradient) program, described in the chapter *Instrument Setup – Pumps, Pump Program*.

Instrument Setup – K-2600 Detector

The detector K-2600 supports the general detector settings described on page 91.

The appearance of some of the details of the K-2600 detector tab depends on the availability of a PDA option.




Fig. 139 Detector K-2600 setup window with PDA option

Used Channels:

Check the box(es) to turn the acquisition for the selected channel(s) on. If the box is not selected, no data will be acquired and no time program will be executed for this channel. As many channels are activated (at least one) as many wavelength columns appear in the spreadsheet below. You can only enable the number of channels you have activated in the instrument configuration setup of the device.

The wavelength table is a spreadsheet to enter your time (detector) program. In the first row (at time 0.00 min), enter the initial wavelength(s) in nm. For up to 50 steps (time in minutes), you can then define the detector parameters. The editing principle including programming the digital outs is the same as for the time (gradient) program, described in the chapter *Instrument Setup – Pumps, Pump Program*.

Autozero at start

Check the **Autozero** field, if you wish the signal to be zeroed at measurement start (to define the current signal as the baseline). You cannot disable this option if you are acquiring 3D data in the advanced mode of the detector.

Autozero at WL change

Check this option if you wish the signal to be zeroed at the moment the wavelength is changed (to define the current signal as the baseline).

Lamp OFF at the end of run

Check the Lamp OFF box to switch off the detector lamp after a measurement is finished. This option may be used for the shutdown methods.



You should switch the lamp on either manually or from the direct control window before you start the next working session. After power up, allow 30 minutes for the lamp to stabilize before running samples that may require sensitive measurements.

Scan, nm [Range 190 – 740 nm]

This section allows you to specify the scan spectral interval, in nm.

Start

Enter the start wavelength for the scan in nm.

End

Enter the end wavelength for the scan in nm.



The scan section limits Start and End are important even if no PDA option is available. All measuring wavelengths in the used channels 1 to 4 must be within these limits.

Acquire 3D data

You have access to this selection field only with the PDA option. The **Advanced...** button becomes accessible if the Acquire 3D data check mark is set.




Fig. 140 Detector K-2600 setup window with PDA option

A click on the **Advanced...** button will open the following dialog box.




Fig. 141 Detector K-2600 Advanced setup window

Interpolate peak profile after scan

While scanning, the data acquisition of the detector is interrupted. This leads to a plateau shape for that time and thus to a false quantitative result. When this option is activated the software interpolates a reasonable course of the chromatogram. The option is limited for runs with a time constant $\leq 1s$.



Warning! This option leads to a better shaped chromatogram; however the quantitative results are still not fully reliable. For accurate and reliable qualitative and quantitative results the run should be performed twice, with and without scanning.

Scan if a peak is detected

At any detected peak the software will perform a scan automatically if this option is activated. All further settings in this area are valid for the whole method. You may decide whether the scan shall be performed **At the peak apex** or **At the beginning of the peak** with the given **Delay [min]**.

With the **Slope [$\mu AU/sec$]:** and **at Height [μAU]:** settings you can define the thresholds, necessary for evaluating the signal as a peak. The unit [μAU] depends on the configuration setting.

Use channel defines the data channel where the analysis is performed. You cannot change this selection during a run. However you can change the measuring wavelength of this channel as often as necessary.

“Scan” column in time table

If the option “Acquire 3D data” is enabled, a new column “Scan” will be available in the time table for the detector K-2600. If you want to perform a time-defined scan, enter a new line into the table with the desired time and enable the scan by clicking the check box. Please note, that during a scan the measurement signal(s) cannot be acquired. In the chromatogram you will find a plateau for the scanning time. It is not allowed to enable a scan in the first and last line of the time table. The time between two scan must be 0.33 min in minimum.

#	Time (min)	Wave-length 1	Wave-length 2	Scan	Dig. Outs (0/1/P/_)	Comments
1	0,00	254	254	<input type="checkbox"/> On	00	
2	3,00	254	254	<input checked="" type="checkbox"/> On	00	
3	10,00	254	254	<input type="checkbox"/> On	00	

Fig. 142 Detector K-2600 time based scan

Instrument Setup – S 2550 Detector

The detector S 2550 supports the general detector settings described on page 91.




Fig. 143 Detector S 2550 setup window

Used Channel

Check the box(es) to turn the acquisition for the selected channel(s) on. Only channels that have been selected in the instrument configuration of the detector are available. If the box is not selected, no data will be acquired and no time program will be executed for this channel. As many channels are activated (at least one) as many wavelength columns appear in the spreadsheet below.

The wavelength table is a spreadsheet to enter your time (detector) program. In the first row (at time 0.00 min), enter the initial wavelength(s) in nm. The editing principle including programming the digital outs is the same as for the time (gradient) program, described in the chapter *Instrument Setup – Pumps, Pump Program*.

Autozero at WL change

Check this option if you wish the signal to be zeroed at the moment the wavelength is changed (to define the current signal as the baseline).

[Action] at the end of run

The detector S 2550 allows not only to switch off the lamp at the end of run, but also to switch the device into Standby.



You should switch on the lamp either manually or from the direct control window before you start the next working session. After power up, allow 20 minutes for the lamp to stabilize before running samples that may require sensitive measurements.

Scan, nm [Range xxx – yyy nm]

This section allows you to specify the scan spectral interval, in nm.

Start

Enter the start wavelength for the scan in nm.

End

Enter the end wavelength for the scan in nm.

Acquire 3D data

You have access to this selection field only with the PDA option. The displayed range depending on the detector type and the installed lamps will only be updated, if you open that window online. The **Advanced...** button becomes accessible if the Acquire 3D data check mark is set.




Fig. 144 Detector S 2550 Scan setup window with PDA option

A click on the **Advanced...** button will open the following dialog box.




Fig. 145 Detector S 2550 Advanced setup window

Scan if a peak is detected

At any detected peak the software will perform a scan automatically if this option is activated. All further settings in this area are valid for the whole method. You may decide whether the scan shall be performed **At the peak apex** or **At the beginning of the peak** with the given **Delay [min]**.

With the **Slope [μ AU/sec]**: and **at Height [μ AU]**: settings you can define the thresholds, necessary for evaluating the signal as a peak. The unit [μ AU] depends on the configuration setting.

Use channel defines the data channel where the analysis is performed. You cannot change this selection during a run. However you can change the measuring wavelength of this channel as often as necessary.

The option **Continuous quick scan** cannot be enabled. This option is not available at this time.



Warning! While scanning, the data acquisition of the detector is interrupted. This leads to a plateau shape for that time and thus to a false quantitative result. For accurate and reliable qualitative and quantitative results the run should be performed twice, with and without scanning.

“Scan” column in time table

If the option “Acquire 3D data” is enabled, a new column “Scan” will be available in the time table for the detector S 2550. If you want to perform a time-defined scan, enter a new line into the table with the desired time and enable the scan by clicking the check box. Please note, that during a scan the measurement signal(s) cannot be acquired. In the chromatogram you will find a plateau for the scanning time. It is not allowed to enable a scan in the first and last line of the time table. The time between two scan must be 0.1 min in minimum.

#	Time (min)	Wave-length 1	Wave-length 2	Scan	Events (O/I/P/L)	Comments
1	0,00	254	254	<input type="checkbox"/> On	0000	
2	3,00	254	254	<input checked="" type="checkbox"/> On	0000	
3	10,00	254	254	<input type="checkbox"/> On	0000	

Fig. 146 Detector S 2550 time based scan

Instrument Setup – Diode Array Detectors

(S 2600, DAD 2850, DAD 2800, K-2700 and Kontron DAD 540/545)

The S 2600, S 2850, S 2800, K-2700, and K-2800 detectors are fully-functional DAD (diode array detector) or PDA (photo diode array) detectors with the capability to simultaneously collect data over the entire UV and visible spectra.

These detectors support the general detector settings described on page 91.




Fig. 147 Detector S 2600 setup window with PDA option

Options

An. output (no access with Kontron DAD 540)

You can activate the option to select the scaling of the analogue output signal. This option is not compatible with a set sampling rate of 10 Hz. The following error message will be displayed when the window will be closed:




Fig. 148 Error message

Scale

The selected value represents the full scale. Values are selectable from 0.001 AU up to 10 AU. For the Kontron DAD 540 it is fixed on 10 AU

Autozero at WL change

Check this option if you wish the signal to be zeroed at the moment the wavelength is changed (to define the current signal as the baseline). At the run start always an Autozero is performed.

Lamp OFF at the end of run

The option is disabled for the K-2700 detector. Check the Lamp OFF box to switch the detector lamp off after a measurement is finished. This option may be used for the shutdown methods.



You should switch the lamp on from the direct control window before you start the next working session. After power up, allow 30 minutes for the lamp to stabilize before running samples that may require sensitive measurements.

Advanced

Pressing the **Advanced** button opens the window for reference correction:



Fig. 149 reference correction window

Referring to a wavelength with practically zero absorption the drift due for instance to a too short heating time of the lamp can be eliminated.

Reference Correction

According to the configuration up to 4 wavelengths can be monitored from the 3D data stream. Select the checkbox to activate the associated channel.

Sample

The sample wavelength is defined in the **Acquisition** tab.

Bw

The bandwidth for the sample channels are also defined in the **Acquisition** tab.

Ref

Reference is the wavelength at which a reference absorbance is measured as a function of time. The reference wavelength compensates for fluctuations caused by changes in baseline absorbance, for example due to changes in solvent composition during gradient elution. The selected wavelength must be higher than the sample wavelength and be inside of the scan range for this detector.

Bw

Enter the bandwidth of the reference wavelength. Limits are 1 to 32 nm in 1 nm steps.

Scan, nm [Range xxx – yyy nm]

The displayed range depending on the detector type and the installed lamps will only be updated, if you open that window online.

Fig. 150 DAD setup, scan range

This section allows you to specify the scan spectral interval, in nm. The given range limits are read out from the used instrument. Therefore they may differ from the shown values.

Start

Enter the start wavelength for the scan, in nm.

End

Enter the end wavelength for the scan, in nm.

Bandwidth

The Bandwidth is the wave length range that will be used to measure a data point. For the spectral acquisition it is highly recommended to use a bandwidth of 1 nm.



The scan section limits Start and End are important even if no PDA option is available. All of the measuring wavelengths in the Wavelength Program for the used channels 1 to 4 must be within these limits.

Signal Mode

The active signal mode is already selected in the instrument configuration (see page 50). The choices are **Absorption** and **Intensity**. For the Kontron DAD 540 it is fixed on **Absorption**.

Acquire 3D data

You have access to this selection field only with the PDA option.

Wavelength Program

#	Time (min)	Wave-length 1	Wave-length 2	Dig. Out (O/I/P)	Comments
1	0.00	254	254	00	
2	10.00	254	254	00	

Fig. 151 Detector S 2600 wavelength program table

Used Channel

Check the box(es) to turn the acquisition for the selected channel(s) on. If the box is not selected, no data will be acquired and no time program will be executed for this channel. As many channels are activated (at least one) as many wavelength columns appear in the spreadsheet below.

Bandwidth [nm]

The bandwidth is a measure of the range over which the calculations and measurements are used to calculate the intensity for a particular wavelength. For example, with a 10 nm bandwidth, the intensity for 250 nm would be calculated from 245 nm to 255 nm.

A wider bandwidth results in an increased signal to noise ratio because the multiple measurements cause the random noise to be cancelled out. That is, if the noise is truly random, more measurements will decrease

the noise because for each positive noise excursion there is a good likelihood that it will be cancelled by a negative excursion. But since the signal is not random, it will remain the same over multiple measurements.

A high bandwidth value will decrease the noise level but also decrease the possible resolution which can be measured. The default value is 8 nm but if your peaks have substantially wider UV spectra you can increase the bandwidth value in ChromGate to obtain a lower noise level.

Each discrete channel can use its own wavelength program and can be measured with its own bandwidth, so you can always obtain the lowest noise conditions while still having a low enough bandwidth to not degrade the resolution measured.

This means that the diode array detectors can always be optimized for maximum performance. The scan parameters can be set to the values which are optimal for doing the scans, and yet the discrete channels can be measured independently at their optimal values.

The wavelength table is a spreadsheet to enter your time (detector) program. In the first row (at time 0.00 min), enter the initial wavelength(s) in nm. The editing principle including programming the digital outs (S 2600 only) is the same as for the time (gradient) program, described in the chapter *Instrument Setup – Pumps, Pump Program*.

Instrument Setup – Fluorescence Detector RF-10Axl / RF-20A

The Fluorescence Detectors RF-10Axl and RF-20A/Axs support the general detector settings described on page 91.



The detector RF-10Axl only supports a (fixed) sampling rate of 0.5 Hz. For higher sampling rates the data acquisition via the 1V integrator output and an A/D converter (as Knauer interface box IF2) must be used.



The detector RF-20A/Axs supports sampling rates up to 0.5 Hz. For higher sampling rates the data acquisition via the 1V integrator output and an A/D converter (as Knauer interface box IF2) must be used. The serial control only allows the same features as for the detector RF-10Axl. All advanced features cannot be supported. Please be sure that the detector has a firmware version that allows for serial control.




Fig. 152 Fluorescence Detector RF-10Axl setup window

Lamp OFF at the end of run

Check the lamp OFF box to switch off the detector lamp after a measurement is finished. This option may be used for the shutdown methods. This option is supported only for detector firmware versions 3.32 and higher.

All additional settings are integrated as additional columns into the spreadsheet for the detector time program. The creation of the detector

time program is similar in principal to that described above in the chapter *Instrument Setup – Pumps, Pump Program*.

The **Time**, **Excitation** wavelength, and **Emission** wavelength values must be entered manually, whereas the options **Autozero**, **Event Out**, and **Event Mark** can be checked separately for each time line. Event Out short circuits the event out terminal for approximately 1 second. Event Mark adds an event mark to the analog output.

The values for the **Sensitivity**, **Gain**, **Recorder Range**, and the **Output Mode** can be selected from the pull-down menus. Fields with selections available will display a combo-box button when the field or row is selected. Click the button to display the available choices.

The available choices for the sensitivity are HIGH, MED, and LOW. For the Output Mode the available choices are Excitation Energy, Light Source Correction, and Emission NOT corrected.




Fig. 153 pull down menus in the Fluorescence Detector RF-10Ax1 setup

Instrument Setup – Alltech 650 Conductivity Detector

The Alltech 650 conductivity detector supports the general detector settings described on page 91.




Fig. 154 Alltech 650 conductivity detector setup window

Temperature Control

Enter the desired temperature of the measuring cell, the allowed temperature tolerance and the waiting time for conditioning.

Program table

The program table can be edited in the usual way. For any time window you can select the range and the polarity of the signal via drop-down lists. The choice for the range is 0.01, 0.05, 0.10, 0.50, 1.00, 5.00, 10, 50, 100, 500, 1000, and 5000 μS .

Instrument Setup – User defined Detector

For any user defined detector only the general detector settings, described on page 91, can be set.

Instrument Setup – Virtual Detector

A virtual detector is not a real instrument but a software tool having special functions. For instance, it is possible to carry out simulations or corrections of chromatograms with this “detector”. The result will be stored as normally created chromatogram.

The virtual detector requires the same general settings as any other detector as described on page 91 (left-hand side of the setup window). They have to be set also for the virtual detector.




Fig. 155 Virtual detector setup window

Ext. Trigger On in 0.00 Min after ready

Check this option to use a logical “trigger signal” for starting the system. This allows a synchronous system start also without an external trigger source (e.g. autosampler or manual injection valve). With “0.00 Min” the system will be started if all connected devices have sent a “ready” to the software. For a delayed start, enter the desired time.

Argument X / Y

X and Y are optional arguments; use of these can be switched off. A result formula can be calculated without either X nor Y, with only X or with both X and Y (single Y is forbidden). If using X or Y is switched on, “X” and “Y” signs must be used in the result formula expression.

After checking this field, the pull-down menu of the **Source** becomes accessible. The options are Time, Data File and External Source.

Time

With this selection no further settings for Argument X are possible. The argument value is equal to the current time value.

Data File

Load an already existing chromatogram (data file) that should be shown starting a run with the virtual detector. The value of this argument is calculated using the corresponding point in the selected data file. The “speed” of the data “acquisition” can be set with the sampling rate. If the original chromatogram was acquired for 10 min with 1 Hz, the chromatogram will be finished within 5 min if the sampling rate for the virtual detector is set to 2 Hz. If no data points left of the original chromatogram, the virtual detector will show a value of 0 (“zero”).

The Settings button will become accessible and clicking on it will bring up the following dialog box:



Fig. 156 Virtual detector setup, data file selection

The very last file selection will be displayed. You can either enter a new file or you can use the **Browse...** button.



To avoid possible signal distortion, the Y-Axis multiplier in the configuration of the virtual detector must be equal to that of the selected data file.

The dialog box which appears includes in the upper part the typical Windows searching facilities. Additionally you can search using the Sample ID and/or the Analyst.




Fig. 157 Virtual detector setup, data file selection

Use the "Trace" column fields to select the required trace if more than one is available in the data file.

After checking the **Shift Time (min)** option you can enter the desired time difference e.g. for correcting time delays of serial arranged detectors.

External Source

The **Settings** button also becomes accessible. Clicking on it the following dialog box appears.



Fig. 158 Virtual detector setup, channel selection

You can select from the pull-down menu one of the detectors or detector channels accessible in the configured system.



To avoid possible signal distortion, the Y-Axis multiplier in the configuration of the virtual detector ***must*** be equal to that of the selected User Detector.

Argument Y

This option becomes accessible, if the Argument X field is activated. The setting dialogs for Y are than the same as for X.

Formula (X,Y)

You can either fill-in the desired formula or you can select it from the pull-down menu, providing all of the formulas used before.

Formula possibilities

Following functions, constants, and operations will be supported:

Functions: "SIN", "COS", "TG", "ASIN", "ATG", "EXP", "LG" (log10), "LN" (ln), "SH" (hyperbolic sine), "CH" (hyperbolic cosine), "TH" (hyperbolic tangent), "SQRT", "ABS", "RAND" (randomizer).

Constants: "PI" ($\pi = 3.141592653$), "E" ($e = 2.71828183$).

Operations: Addition "+", Subtraction "-", Multiplication "*", Division "/", Power "^", unary plus and minus.

Advanced...

Pressing this button a dialog box appears to create a time table for changing the formula during a run.




Fig. 159 Virtual detector setup, time table

Instrument Setup – Assistant ASM2.1L

Due to the modular system, there are a lot of possible configurations for the ASM2.1L. The setup window will look different, depending on the configured modules. For all possible modules, the setup is described below.




Fig. 160 ASM 2.1L setup window for pump, valve, detector

Run time

Run time determines the run time of all devices of the Assistant.

Acquisition delay (only for UVD2.1S)

Acquisition delay is the interval between the start of the run (Trigger) and the time when sampling starts for this detector. This option will only be shown, if a detector is configured in the ASM 2.1L.

UV (only for UVD2.1S)

This section only appears, if the detector UVD2.1S is configured.

Time Constant

Select the time constant for the detector from the drop-down list. The values are 0.1, 0.2, 0.5, 1.0, 2.0, 5.0, and 10.0 sec. The additional values 0.000, 0.005, 0.01, 0.02, 0.05 will only be available, if the detector is connected by LAN. Smaller time constants are required for high speed analyses. A good rule is to select a time constant, which is the reciprocal of the data rate in Hz (e.g. 2 Hz data rate -> 0.5 sec time constant).

Sampling rate

The sampling rate is the rate at which data will be sampled by the system. Click on the down-arrow to get a list of the frequencies available for the configuration of your system. The choices are: 1, 2, 5, and 10 Hz and, if the detector is connected by LAN, also 20, 50 and 80 Hz. For any selected sampling rate the narrowest peak width will be shown, for which this sampling rate will be adequate. This selection is in Hz (data points per second). This is the selection for most chromatography applications. We recommend having more than 20 data points per peak (30 – 40). This means that for a peak with a peak width of 20 seconds (from peak start to peak end) a data rate of 2 Hz is recommended.

Autozero at start

Check the Autozero field if you wish the signal to be zeroed at the moment a measurement is started (to define the current signal as the baseline).

Autozero at WL change

No access to this option. The grayed area is only a reminder, that the autozero will take place at any wavelength change during a run.

Pump/HPG: Pressure Limits, [unit] (only for pump P4.1S)

This section is only accessible, if a pump P4.1S (**Pump**) or two pumps as a HPG with including a pump P4.1S (**HPG**) have been configured. The shown pressure unit depends on the unit selected in the pump configuration (MPa, bar, psi)

Min:

Enter the minimum pressure allowable for your pump. The zero setting disables the shut down due to minimum pressure.

Max:

Enter the maximum pressure allowable for your pump. The default value always corresponds to the configured pump head. A setting above this value is not allowed.

If for the min and max limit "0" is set, the software will not check for the pressure and the limits will only be monitored by the pump internally.

Options

The settings accessible in this section depend on the configured modules.

Standby at the End of Run

If the standby option is checked, the Assistant will switch into Standby if a run is finished. The pump motor and, the detector lamp will be switched off as well as the Assistant's display will show the Standby status of the Assistant. If this option is selected, all other options of the Options section become inaccessible.

Flow OFF at the end of run

Checking this option the flow will be set to zero after elapsing the runtime.

Lamp(s) OFF at the End of Run

Check the lamp OFF box to switch off the detector lamp after a measurement is finished. This option may be used for the shutdown methods.

Time program

The time program show the program lines for all configured modules. Please refer to the setup information of the separate instruments (pump, detector, valve, MultiValve FC) for the specific setup options.

If a HPG is configured, there will only be shown a line for Pump A. The header of the 3rd column will show the unit for the module selected in the table. For pumps the header shows "Flow (ml/min)", for the detector "WL (nm)" and for valves "Pos". The 4th column only shows "%A", if a HPG is configured. The smallest allowed time step in the time program is 0.02 minutes.

Show All Lines

The drop down menu allows for showing either the lines for all configured modules or for just one of the modules. The available modules depend on the configured modules for the Assistant. The option "Show ASM2.1L lines" allows for programming the 3 events of the Assistant.

Time (min)

Add a new time line with the time you want to change a module's setting. If you have added a line, select the module you want to program in the device column. If you want to change settings for several devices at the same time, you must add separate lines for each module.

Device

The device column allows for selecting one of the configured modules or the Assistant. If a HPG is configured, only the HPG A pump can be selected.

[3rd column]

The 3rd column allows for entering a value for the selected device. If no line is selected the header shows

If a pump is selected, the column header show "Flow (ml/min)" and the flow for the pump or HPG system can be selected. The allowed flow depends on the pump head. For a 10 ml pump head the flow rate range is 0.00 – 9.99 ml/min, for a 50 ml pump head 0.0 – 50.0 ml/min. Please note, that for a HPG both pumps must have the same pump head.

If the detector is selected, the column header shows "WL (nm)". Following the specifications of the UVD2.1S the allowed wavelength range is 190 – 500 nm.

If a valve is selected, the column header shows "Pos.". A drop-down menu will be accessible to select the desired valve position.

If the ASM2.1L is selected as a device, the column header shows "Events". The field shows now the selected events. "-" will be shown for events no option is selected. In practice this means, that the current status will not be changed. The setup will open by clicking the small arrow key. The choices for the events are Off, On and Pulse. If selected, "0" represents event Off, "1" event On and "P" Pulse.

[4th column]

The 4th column is only to show "%A", if as a device a pump of a HPG is selected. This allows for selecting the %A part of a gradient, while %B is calculated by the software (100% - %A). If no pump is selected the header shows

Pretreatment

This option enables several pretreatment operations of pumps and valves before injection and starting the data acquisition. For setting up refer the **Pretreatment** section on page 89.

Assistant Setup / FRC Setup

If a valve is configured as a **fractionation valve**, the Assistant setup includes two tabs, one for the Assistant and one for the fractionation valve. Please refer to the section **Fraction Collector Setup** for detailed information for the fraction valve setup (MultiValve FC).




Fig. 161 ASM 2.1L setup window for fractionation valve

Instrument Setup – Autosamplers

In any system you can include only one autosampler. Therefore all supported samplers will be explained separately without cross references.

Instrument Setup – Autosampler 3800

Smartline Autosampler 3800 and K-3800 (Basic Marathon)

The method parameters required for automatic operation of your autosampler are contained on the *Method – Instrument Setup – AS 3800* tab. The control parameters for the autosampler become part of your method and sequence files.




Fig. 162 Autosampler 3800 setup window

Injection Type

The AS 3800 can perform injections based on the flushed loop or partial loop filling (for the Plus only) principles. Loop filling can be done on a volume basis, controlled by the digital dispenser system, or on a time

basis based on the headspace pressure of the sample. You can choose between the **dispenser controlled** or **pressure controlled** injection type.

Dispenser Controlled Injection

For the dispenser controlled injection, specify the **Preflush Volume** (30-999 μ l).

Pressure Controlled Injection

For the pressure controlled injection, specify the **Flush Time** (1-999 seconds).

Check the **Needle Wash** button to select needle wash.

Temperature Control

Check the **Enable tray cooling** field if your sampler is equipped with the tray cooling facility and you want to use it.

Aux(iliaries)

This section allows you to program auxiliary output(s) for the AS 3800.

Enable Aux1

Select this box if time programmable auxiliary output 1 is to be used.

(Aux1) Time ON

Enter the time at which auxiliary 1 must switch ON (0 - 99 min 59 sec).

(Aux1) Time OFF

Enter the time at which auxiliary 1 must switch OFF (0 - 99 min 59 sec).

Enable Aux2

Select this box if time programmable auxiliary output 2 is to be used.

(Aux2) Time ON

Enter the time at which auxiliary 2 must switch ON (0 - 99 min 59 sec).

(Aux2) Time OFF

Enter the time at which auxiliary 2 must switch OFF (0 - 99 min 59 sec).

Instrument Setup – Autosampler Knauer Optimas / 3900 (Midas)

The method parameters required for automatic operation of your autosampler are contained in the *Method – Instrument Setup – Optimas/AS3900* tab. The control parameters for the autosampler become part of your method and sequence files.




Fig. 163 Autosampler 3900 setup window

Injection Method

Select the injection method by clicking on one of the available options: **Full loop**, **Partial loop fill**, **µl Pick-up**, or **None**.



The maximum programmable injection volumes are:

Flushed loop: not programmable, is equal to the loop volume but needs more sample to fill the loop completely (3 x loop volume for loop volumes < 100 µl and 2 x loop volume for loop volumes from 100 µl to 500 µl 1.5 loop volume for loop volumes > 500 µl)

Partial loop fill: 50 % of the programmed loop volume

µl Pick-up: Injection volume = (loop volume – 3 x needle volume) / 2

Syringe Speed

The aspirating speed of the syringe used in injection methods can be adapted depending on the viscosity of the samples (the choices are **Low**, **Normal**, and **High**). The speed of the syringe during washing or rinsing procedures is not affected by this setting.

Flush Volume

Enter the amount of sample in µl taken from a vial to preflush the needle before the loop is filled with the sample. The default value is 30 µl for a tube volume of 15µl. In the flushed loop mode this option is deactivated.



Flush volumes lower than twice the volume of the needle and tubing will result in a decrease of performance of the autosampler.

Wash times

Enter the wash volume in number of syringe volumes.

Air Segment

Check the Air Segment box to enable the use of an air segment between the sample and the wash solvent to reduce the amount of flush volume.

Headspace Pressure

Check the Headspace Pressure box to enable the use of headspace pressure in combination with the syringe to transport sample to the loop. The compressor will always be used during a wash procedure.



The accuracy and reproducibility of the autosampler may decrease if headspace pressure is switched off.

Needle Wash

Check the Needle Wash box to enable a wash procedure.

Alarm Buzzer

Check the Alarm Buzzer box to enable the alarm beep which is emitted when an error occurs during a run.

Run Screens

Check the Run Screens box to display the actual run information on the autosamplers display during a run.

Temperature control (

If your autosampler is equipped with a tray cooling and/or a facility thermo stated column oven (AS 3900 only) this options can be checked and the desired temperatures can be entered.



The tray cooling allows to enter a temperature between + 4 °C and + 15 °C.



The AS 3900 oven can only control the temperature at ambient temperature + 5 °C and up. The maximum temperature is + 60 °C.

Timed Events

These types of events or methods enable you to control other devices via the two auxiliaries (contact closures) or to activate a new oven set point. Each event can be programmed on a time base, which starts at the moment of injection.

Click the **Timed events...** button to set up the time base methods. The time base methods set up dialog will appear.




Fig. 164 Autosampler timed events

Click the **New** button to add one more lines to the programmed methods.

Method type:

Select one of the available **AUX-1** and **AUX-2**. Every auxiliary **AUX-1** through **AUX-2** can be set to **ON** (contact closure) four times and four times to **OFF** at a programmed **Time**. It is not possible to program the same auxiliary to ON and OFF at the same time.

The **maximal** programmable time for the time base methods is **9 h 59 min 59 sec**.

Click the **Commit** button to actualize your selection.

Click the **Delete** button to remove the selected line from the programmed methods.

The **Stop Events at Time** allows setting the end time of the timed events. If the **Stop Events at Time** box is not checked, the end time is equal to the analysis time (the run time). The **maximal** Stop events time is **9 h 59 min 59 sec**.



If the END time exceeds the programmed run time, this END time overrules the analysis time for the autosampler. The next injection will be started as soon as the timed events program has finished.

Mix Methods...

Click the Mix Methods button to create a program of sample preparation. The Mix Methods set up dialog will appear.



Fig. 165 Autosampler 3900/3950 Mix Methods

Clicking the New button first the method type drop-down list becomes active. The mix method can contain a maximum of 15 programmable lines including the **END OF MIX** command.

Click the **New** button to add one more line to the programmed methods.

Select the **Method Type** from the drop-down list. The choices are **Add**, **Mix** and **Wait**.




Fig. 166 Autosampler 3900/3950 Mix Step Add

Specify the **Amount** of a liquid in μl (max = the syringe volume) to be aspirated from a source (choices are **Sample**, **Reagent A**, **Reagent B** and **Wash**) and dispensed to a target (choices are **Sample** and **Destination**).



The aspirate and dispense speed depends on the selected syringe and syringe speed.

The reagents and transport solvent are located at the following positions on the tray:

Transport solvent:	vial 85
Reagent A:	vial 86
Reagent B:	vial 87



Fig. 167 Autosampler 3900/3950 Mix Step Mix

According to the number of **Mix Times** (1-9 times) the programmed **Amount** (μl , max = the syringe volume) will be aspirated and dispensed. If no destination vial is used the mixing step is performed in the sample vial.




Fig. 168 Autosampler 3900/3950 Mix Step Wait

Wait a programmed period of time **x:xx:xx** (0:00:00 – 9:59:59) before continuing with the next step (reaction time).

Click the **Commit** button to complete the operation.

Click the **Delete** button to delete the selected step.

Instrument Setup – Autosampler 3950 (Alias)




Fig. 169 Autosampler 3950 setup window

Injection Method

Select the injection method by clicking on one of the available options: **Full loop**, **Partial loop fill**, or **µl pick-up**. If µl pick-up is selected, a drop-down menu for the Transport vial position (**Transport pos.**) is available.



The maximum programmable injection volumes are:

Flushed loop: not programmable, is equal to the loop volume but needs more sample to fill the loop completely
(3 x loop volume for loop volumes < 100 µl and
2 x loop volume for loop volumes from 100 µl to 500 µl
1.5 loop volume for loop volumes > 500 µl)

Partial loop fill: 50 % of the programmed loop volume

µl pick-up: Injection volume = (loop volume – 3 x needle volume) / 2

Transport pos.: (only available for µl pick-up with tray type **84+3 vials**)
Select the position of the vial with the transport liquid. Available are vials 85, 86 and 87. The transport plug is 2.5 x needle volume (needle volume = tubing volume, please refer to the autosampler configuration chapter). The autosampler itself calculates the liquid level in the transport vials. It will be reset to 8000 µl every time if the autosampler is switched off and on. Additionally the level can be reset from the autosamplers direct control tab. The level cannot be read-out from the software. If the autosampler calculates, that the vial is empty, the method/sequence will stopped. Dependent from the transport liquid level, the autosampler will adapt the needle immersion depth. Therefore it is important to fill the transport vials whenever the autosampler was switched off or the transport liquid level was reset manually. If a vial is selected as the transport vial, the configuration of the vial as a Reagent A or B vial will be ignored.

Syringe Speed

The aspirating speed of the syringe used in injection methods can be adapted depending on the viscosity of the samples (the choices are **Low**, **Normal**, and **High**). The speed of the syringe during washing or rinsing procedures is not affected by this setting.

Flush Volume

Enter the amount of sample in µl taken from a vial to preflush the needle before the loop is filled with the sample. The default value is 30 µl for a tube volume of 15µl. The preflush will only be done for the injection modes **Full loop** and **Partial loop fill**. In the flushed loop mode this option is deactivated.



Flush volumes lower than twice the volume of the needle and tubing will result in a decrease of performance of the autosampler.

Wash times

Select the Wash times (in number of syringe volumes). The valid interval is between 1 and 9.

Needle height

It defines the penetration depth of the needle or the distance between the vial bottom and the needle tip. The valid settings are between 2 and 6 mm with the step 1 mm.

Air Segment

Check the Air Segment box to enable the use of an air segment between the sample and the wash solvent to reduce the amount of flush volume.

Headspace Pressure

Check the Headspace Pressure box to enable the use of headspace pressure in combination with the syringe to transport sample to the loop. The compressor will always be used during a wash procedure.



The accuracy and reproducibility of the autosampler may decrease if headspace pressure is switched off.

Needle Wash

Check the Needle Wash box to enable a wash procedure.

Temperature control

If your autosampler 3950 is equipped with a tray cooling facility, this option can be checked and the desired temperature (4 – 22 °C) can be entered.

Wait with tolerance (\pm °C)

The system becomes ready only if the actual tray temperature is inside the set range.

Wait time (min)

The waiting time is the equilibrating time in minutes after the tray temperature reached the temperature set range before the injection takes place.

Mix Methods...

Click the Mix Methods button to create a program of sample preparation. The Mix Methods set up dialog will appear. It is the same dialog as already described for the AS 3900, starting with Fig. 165 on page 115.

Please note, that for the autosampler 3950 for all tray types except the 83+3 vials tray no transport vial for the μ l pick-up can be defined. The transport liquid will be escaped from the wash solution. If the 84+3 vials tray is selected, the vials 85, 86 or 87 can be selected as a transport vial. The transport plug is 2.5 x needle volume (needle volume = tubing volume, please refer to the autosampler configuration chapter). The autosampler itself calculates the liquid level in the transport vials. It will be reset to 8000 μ l every time if the autosampler is switched off and on. Additionally the level can be reset from the autosamplers direct control tab. The level cannot be read-out from the software. If the autosampler calculates, that the vial is empty, the method/sequence will stopped. Dependent from the transport liquid level, the autosampler will adapt the needle immersion depth. Therefore it is important to fill the transport vials whenever the autosampler was switch off or the transport liquid level was reset manually. If a vial is selected as the transport vial, the configuration of the vial as a Reagent A or B vial will be ignored.

Stacked injections...

The option Stacked injections allow to define additional injections during a run. Such injections can only be programmed either for the **autosampler 3950** or the **injection module** (please refer the FRC section, chapter "Stacked Injections" of this manual). Click on the Stacked injections – button to open the Stacked Injections table dialog.

Enable the check box "Use stacked injection".

Stacked Injection Table

Use stacked injection

Injection method: Full loop (max inj. volume 100 µl)

Start sample preparation: 0,25 min before stacked injection

Do not wash the needle between injections from same vial

#	Time (min)	Vial	Inj. Volume (µl)	Comments
1	0,00	<INJ_POS>	<INJ_VOL>	Initial injection

The actual time(s) of the stacked injection(s) may deviate from that programmed.

Fig. 170 Autosampler 3950 Stacked injection table

Injection method

The injection method is the same defined in the instrument setup of the autosampler. It cannot be selected separately for the stacked injections.

Start sample preparation X.XX min before stacked injection

Define the time before the injection time; the autosampler should start with the sample/injection preparation. Please make sure, that the given time is sufficient to prepare the injection before the injection time arrives, otherwise the injection may be delayed. If the autosampler is ready to inject, it will wait until the programmed time for the stacked injection. If the autosampler needs a longer time than programmed, it will inject immediately after the injection preparation is finished. In this case the next stacked injection will be delayed by this “extra time”.

Do not wash the needle between injections from same vial

Checking this option prevents the needle wash, if the stacked injection is made from the same vial as the previous injection. This option may be useful for saving time or increasing reproducibility of injections.

Stacked injections time table

Make a right mouse click into the time table and select *Fill table* from the corresponding menu.

Stacked Injections

Number of Injections: 2

Sample Vial: Next

Injection Volume (µl): 100

Sample Preparation Time (min): 0,25

Time between Injections (min): 1,00

Do not wash needle between injections from same vial

Fig. 171 Autosampler 3950 Stacked injections table setup

Number of Injections

Enter the number of stacked injections. If you, e.g., enter “2”, two stacked injections will be made; in addition of the initial injection 3 injections at all

will be made per run. Please make sure that you have enough samples for the entered number of injection.

Sample Vial

Select the vials the stacked injection should be made from. Select **Same** to injection from the same vial as the previous injection, **Next** to inject from the vial number <n+1> from the previous injection or **Custom** to select free the vial number.

Injection Volume (µl)

Select the desired injection volume. Note that the volume depend from the selected injection mode of the autosampler instrument setup.

Sample Preparation Time (min)

Define the time before the injection time; the autosampler should start with the sample/injection preparation. Please make sure, that the given time is sufficient to prepare the injection before the injection time arrives, otherwise the injection may be delayed. If the autosampler is ready to inject, it will wait until the programmed time for the stacked injection. If the autosampler needs a longer time than programmed, it will inject immediately after the injection preparation is finished. In this case the next stacked injection will be delayed by this "extra time".

Time between Injections (min)

Enter the time between the injections.

Do not wash needle between injections from same vial

Enable this option if you wish to prohibit the needle wash after each injection, programmed in the autosamplers instrument setup, for the stacked injections from the same vial.

If you close the Stacked Injections window with **Apply**, the selected settings will be copied into the Stacked injections table. Any existing entries from previous programming will be deleted.

If you make a right mouse click on the table, you can select

- **Fill Down**, to fill the table with the same settings as in the selected line
- **Insert Line**, to add a new line
- **Delete Line**, to delete the selected line
- **Delete Lines**, to delete all selected lines
- **Fill Table...**, will open the Stacked injection setup again. The new table will overwrite the existing one. If you select the command **Fill table** to open the setup dialog for the stacked injection table, the window will show **default values** for all parameters, beside the sample preparation time.

Click on **OK** to close the Stacked Injections table window; the autosampler S 3950 instrument setup window will be shown.

In the **Aux Traces** tab a trace "Stacked Inj. [AS 3950]" can be enabled. The trace will show the start of the sample preparation and the stacked injection.




Fig. 172 Autosampler 3950 Stacked injections auxiliary trace




Fig. 173 Autosampler 3950 Stacked injections trace

For every stacked injection an entry for the sample preparation and the injection will be made.

Instrument Activity Log				
User	Logged	Source	Activity	
System	09.09.2010 16:20:06	PC-PM-LAB-TEST AS 3950:	Stacked injection: trigger marker is detected in 6.02 min after start.	
System	09.09.2010 16:19:50	PC-PM-LAB-TEST AS 3950:	Stacked injection 5 µl from vial 1C1 will be done in 5.75 min after start.	
System	09.09.2010 16:17:06	PC-PM-LAB-TEST AS 3950:	Stacked injection: trigger marker is detected in 3.02 min after start.	
System	09.09.2010 16:16:50	PC-PM-LAB-TEST AS 3950:	Stacked injection 5 µl from vial 1B1 will be done in 2.75 min after start.	
System	09.09.2010 16:13:41	PC-PM-LAB-TEST	Run Queue - Start Single Run - C:\EZChrom Elite\Enterprise\Projects\Default\Data\	
System	09.09.2010 16:13:40	PC-PM-LAB-TEST	Run Queue - Add Single Run - C:\EZChrom Elite\Enterprise\Projects\Default\Data\	

Fig. 174 Autosampler 3950 Stacked injections instrument activity log entry

Instrument Setup – Triathlon/Endurance Autosampler

The method parameters required for automatic operation of your autosampler are contained in the **Method/Instrument Setup/Triathlon (Endurance)** tab. The control parameters for the autosampler become part of your method and sequence files.




Fig. 175 Autosampler Triathlon setup window

Standard Injection

In the **Standard Injection** section you can enter parameters for the standard actions that can be performed with the autosampler.

Injection Method

Click one of the **Flushed**, **Partial** or **µl Pick-up** fields to select the injection method.



The maximum programmable injection volumes are:

Flushed loop: not programmable, is equal to the loop volume but needs more sample to fill the loop completely (3 x loop volume for loop volumes < 100 µl and 2 x loop volume for loop volumes > 100 µl)

Partial loop fill: 50 % of the programmed loop volume

µl Pick-up: Injection volume = (loop volume – 3 x needle volume) / 2

Syringe Speed and Scale Factor

The aspirating speed of the syringe used in injection methods can be adapted depending on the viscosity of the samples (the choices are **Low**, **Normal**, and **High**). Alternatively, syringe speeds can be reduced by entering a scale factor (**1** to **10**). The syringe speed will be the scale factor multiplied by the syringe speed. The syringe speed during the washing or the rinsing procedures is not affected by this setting.

Air Segment

Check the Air Segment box to enable the use of an air segment between the sample and the wash solvent to reduce the amount of flush volume.

Headspace Pressure

Check the Headspace Pressure box to enable the use of headspace pressure in combination with the syringe to transport sample to the loop. The compressor will always be used during a wash procedure.



The accuracy and reproducibility of the autosampler may decrease if headspace pressure is switched off.

Needle Wash

Check the Needle Wash box to enable a wash procedure.

Flush Volume

Enter the amount of sample in µl taken from a vial to preflush the needle before the loop is filled with the sample. The default value is 30 µl for an installed tube volume of 15 µl.



Flush volumes smaller than twice the volume of the needle and tubing will result in a decrease of performance of the autosampler.

Wash Volume

Enter the volume of wash solvent in μl to be used. The minimum programmable volume is 300 μl .

Needle Height

Enter the distance between the needlepoint and the plate holder in mm (default value is 5 mm).

Temperature control

If your Triathlon is equipped with a tray cooling facility, this option can be checked and the desired temperature can be entered.

Alarm Buzzer

Check the alarm buzzer box to turn ON the alarm beep, given when an error occurs during a run.

Run Screens

Check the run screens box to display the actual run information during a run.

Timed events...

These type of events or methods enable you to control optional ISS valves and other devices via auxiliary (contact closures) or a binary output (BCD or HEX).

Click the **Timed events...** button to set up the time base methods. The time base methods setup dialog will appear.




Fig. 176 Autosampler Triathlon Time base Methods

Click the **New** button to add one more line to the programmed methods.

Select the **Method Type** from the drop-down list. The choices are AUX-1, AUX-2, AUX-3, AUX-4, ISS-A, ISS-B, CODE-Out, and SSV Port.

Every auxiliary **AUX-1** through **AUX-4** can be set to **ON** (contact closure) four times and four times to **OFF** at a programmed **Time**. It is not possible to program the same auxiliary to ON and OFF at the same time.

The ISS valve can only be programmed if the optional ISS valves are installed. The valves **ISS-A** and **ISS-B** can be set to positions **6-1** four times and four times to **1-2** at a programmed **Time**.

The SSV port can only be programmed if the SSV option is installed. Set the time and the **SSV port** number, a value between **1** and **6**.

Eight SSV port switch commands are possible in the time base program.

The **Code** event needs a **time** and a value between **0** and **15 (F)**, hexadecimal output. Eight code output commands are possible in the time base program.

The **maximal** programmable time for the time base methods is **9 h 59 min 59 sec**.

Click the **Commit** button to actualize your selection.

Click the **Delete** button to remove the selected line from the programmed methods.

The **Stop Events at Time** allows setting the end time of the timed events. If the **Stop Events at Time** box is not checked, the end time is equal to the analysis time (the run time). The **maximal** Stop events time is **9 h 59 min 59 sec**.



If the END time exceeds the programmed run time, this END time overrules the analysis time for the autosampler. The next injection will be started as soon as the timed events program has finished.

User Defined Program (UDP)

The **User Defined Method** is a user defined combination of all possible programming actions of the Triathlon or the Endurance like: aspirating from or dispensing to a certain position, moving the syringe, switching the syringe valve, switching the compressor ON or OFF, switching the injection valve and the optional ISS valves, wait and programming all inputs and outputs.

You can build your own sample handling sequence, within the possibilities of the Triathlon or the Endurance.

Enable UDP

Check the **Enable UDP** box to use the user defined method. This selection sets the standard injection parameters inactive.

UDP Setup

Click the **Setup UDP...** button to set up the user defined method. The User Defined Program dialog will appear.




Fig. 177 Autosampler Triathlon edited user defined example program

Click the **Add** button to add one more line to the programmed steps.



The total number of steps for the user program cannot exceed 240.

The **Add New Step** dialog will appear.




Fig. 178 Autosampler Triathlon editing user defined programs, new step

Select the **Action** from the drop-down list. The appearance of the dialog will be changed, depending on the selected action.




Fig. 179 Autosampler Triathlon UDP action choices

Aspirate

Aspirates a programmed volume from a **sample** vial, a **destination** vial, a **wash** vial, or one of the **reagent vials** into the buffer tubing. Speed and height of the syringe can be entered. The maximum volume that can be aspirated is the total volume of the syringe. If more than one consecutive aspirating action (from different sources) is programmed, the total volume of the consecutive aspirating actions cannot exceed the volume of the syringe.

Dispense

Dispenses a programmed volume from the buffer tubing into the **sample** vial, a **destination** vial, the **waste**, a **wash** vial, or one of the **reagent vials**. Speed and height of the syringe can be entered.

Syringe Valve

This action controls the inter-connections of the syringe to one of its three ports:

- NEEDLE:** connection to the sample needle.
- WASH:** connection to the wash solvent bottle.
- WASTE:** connection to the syringe waste tubing.

Syringe

This action controls the movement of the syringe

- LOAD:** loads the syringe with the programmed volume.
- UNLOAD:** unloads the syringe with the programmed volume.
- HOME:** initializes the syringe to its HOME position. The volume previously aspirated will be dispensed to the last programmed position.

Needle Wash

With this action a needle wash can be performed. The content of the buffer tubing will not be rinsed to waste before the washing procedure

starts. The programmed volume of wash solvent is used to wash the needle at the wash position.



The wash position may be contaminated with the contents of the buffer tubing, which may generate cross-contamination. To prevent contamination of the wash position, program a “dispense to waste” action before programming a wash action.

Valve

The positions of the installed high pressure valves are programmed in this screen. The injector valve has two positions: **INJECT** and **LOAD**. If the ISS option is installed, the two valves have the two positions: **1-6** and **1-2**.

Wait

The WAIT action can be used to let the autosampler wait for a programmed time with maximum of **9 h 59 min 59 sec**.



During WAIT the needle turns back to home position if the next action is not at the same position. Should the needle be required to remain in the current position, it is necessary to program a dummy step after the WAIT in which an aspirate or dispense action of 0 µl is programmed at the actual position.

Compressor

With the compressor action the compressor is activated to transfer air pressure to a vial. The compressor will stay active until it is switched OFF. The compressor will be automatically switched off at the end of the needle wash routine if a needle wash is used.

Auxiliaries

With the Auxiliaries action the four standard auxiliaries (contact closures) can be controlled (see time base events).

Wait for Input

The autosampler will wait for one of the four inputs to become **HIGH** or **LOW** before continuing with the next step.

Prog. Output

The Prog. Output action is used to define two programmable outputs (contact closures). These are similar to the auxiliaries.

Code Output

The Code Output action is used to program the output to the connector P3 TIMED OUTPUTS on the rear panel of the autosampler. This is a HEX output in the range **0** to **15** (F).

Markers

The markers normally generated by the autosampler are not active in the user defined method, but can be programmed in this screen. The choices are: **Inject**, **Vial**, and **Labeled**.



Markers must be programmed, because they are not automatically generated as in the regular injection method. Refer to the Triathlon/Endurance User Manual for the cable connections.

SSV Port

Program the solvent selection valve (SSV) port position, range **1** to **6**.

Repeat

The last **x** steps will be repeated **y** times.

Button functions:

Click the **Add** button to complete the operation.

Click the **Edit** button to modify the selected step.

Click the **Remove** button to delete the selected step.

Click the **Up** button to move the selected step up.

Click the **Down** button to move the selected step down.

Click the **Save** button to save the user defined program or its part in to a file.




Fig. 180 Autosampler Triathlon, saving user defined programs

You can save **Entire Program**, **Selected** steps, and Steps in the specified **Range**. Enter the name of the file and, if wanted, your comments to the program. Click the **Save** button in the dialog to complete the operation.

Click the **Open** button to read the saved user defined program.




Fig. 181 Autosampler Triathlon, opening a user defined program

Select the file with the user defined program and the **Open Mode**. The choices are: **Overwrite existing program**, **Insert to begin**, **Append to end**, and **Insert to focused position**. Click the **Open** button in the dialog to complete the operation and to return to the main autosampler method set up dialog.

Click the **Store UDP** button to store the user defined method in the battery backup memory of the autosampler with the entered number of the local program.

Large Injection

Check the **Enable large injections** box to use this option for injections larger than 1 ml. This selection sets the standard injection parameters inactive. It cannot be combined with any user defined injection program.

Instrument Setup – Kontron Autosamplers




Fig. 182 Kontron Autosamplers setup window

Injection type

Needle wash

The standard procedure involves two needle wash cycles after each injection. Checking this option additionally a needle wash is performed before each injection.

Injection volume

The injection volume has to be entered for the AS 360, 560, and 565 in fact as a set of three parameters. The possible values depend on the configured syringe. They are summarized in the following table:

Syringe	Air Gap	Sample / Into Loop	Sum ^{*)}
250	5 to 120	0 to 230	240
500	5 to 240	0 to 470	480
1000	5 to 480	0 to 950	960
2500	5 to 800	0 to 1590	1600

*) Includes the air gap twice

Air Gap: A defined Air segment is used to separate the sample plug from the wash liquid in the sampler lines.

Sample: Defines the actual sample volume of sample taken from the sample vial.

Into Loop: Defines the volume required to transport the sample into the sample loop.

Aux

Two auxiliary outputs (four for AS 460 and 465) can be programmed. At the given times a pulse will be sent to the selected output number.

Mix Methods...

This button appears for the AS 460 and 465 only. Click the Mix Methods button to create a program of sample preparation. The Mix Methods set up dialog will appear.




Fig. 183 Kontron AS 460/465 Mix Methods

Clicking the **New** button first the method type drop-down list becomes active.

Click the **New** button to add one more line to the programmed methods.

Select the **Method Type** from the drop-down list. The choices are **Add**, **Mix** and **Wait**.




Fig. 184 Kontron AS 460/465 Mix Step Add

Specify the **Amount** of a liquid in μl (max = the syringe volume) to be aspirated from a source (vial 1 to 110) and dispensed to a target (vial 1 to 110). Additionally the air gap can be defined.

To continue click the **Commit** button and then the **New** button to program the next step.




Fig. 185 Kontron AS 460/465 Mix Step Mix

The programmed **Amount** (μl , max = the syringe volume) will be aspirated and dispensed three times. Again the air gap can be defined. The mixing step is performed in that vial which is specified in the last adding step by **add to**.



Fig. 186 Kontron AS 460/465 Mix Step Wait

Wait a programmed period of time **x:xx:xx** (0:00:00 – 9:59:59) before continuing with the next step (reaction time).

Click the **Commit** button to complete the operation.

Click the **Delete** button to delete the selected step.

Instrument Setup – Miscellaneous Instruments

Instrument Setup – Manager 5000/5050/IF2 I/O




Fig. 187 Setup window, Manager 5000/5050 / IF2

Edit the spreadsheet as already described for other time programs, e.g. as described in the chapter *Instrument Setup – Pumps, Pump Program*.

Enter the maximum output voltage in the **voltage** column, e.g. for maximum delivery of a pump and the percentage values into the **AOF** (**A**nalog **O**utput **F**actor) column. The output voltage, produced by the Manager 5000/5050/IF2 I/O, will be the result of

$$\text{Voltage} * \text{AOF} / 100$$

The influence of the **Step Mode** option is illustrated in the next figure:




Fig. 188 Analog output profile without and with activated Step Mode

Instrument Setup – Switching Valves




Fig. 189 Setup window, Knauer Switching Valves

Use this spreadsheet to enter your time program. The number of columns corresponds to the number of configured valves in the system. In the first row (at time 0.00), enter the initial position for each valve. Then you can define the states of the valves at the given time.

The choices for a Knauer 2-position valve are **BLANK** and **L** (Load) and **I** (Inject), for a Vici 2-position valve **BLANK** and **A** (Load) and **B** (Inject), **1** through **6** for a 6-position valve, **1** through **8** for a 8-position valve, **1** through **10** for a 10-position valve, **1** through **12** for a 12-position valve and **1** through **16** for a 16-position valve. The BLANK space means that the valve will not change its state at the given time.

You can enter the settings either by typing the desired values into the corresponding field or by selecting your choice from the pull down menu which becomes accessible via the  button, appearing in the selected field.

For details on the spreadsheet operations, refer to the section *Creating a Time (Gradient) Program* as described in the chapter *Instrument Setup – Pumps, Pump Program*.

The valves program profile will be adapted automatically to the spreadsheet inputs when you save the method.

If you select any field in the comment column the  button will appear, to open a comment window.




Fig. 190 Setup window, Knauer Switching Valves, Comments

Pretreatment

Checking this option the **Setup...** button becomes active. Clicking on it the following dialog window will be opened:




Fig. 191 Pretreatment setting dialog window

The appearance of the valve program corresponds to that of the main valve setup window. It must be edited in the same way as described above. To add a new line, make a right mouse click into the table and select "Insert line" from the menu. The only difference is that **negative** times must be entered, because they are relative to the zero injection time.



If a pretreatment is to define for more than one instrument (pump and/or valve drive) in a method, it is required to set for all of them the same earliest time. All pretreatment procedures will start simultaneously.



If a pump should exactly switch off after 5 minutes of pretreatment, include a line with 0 ml/min flow into the pretreatment's pump program. Otherwise the pump may run for a longer time as it need some seconds to load the new flow rate from the method's "normal" pump program

An example for both:

A pump should run for exactly 5 minutes before the injection and a valve should switch 2 minutes before injection from position 2 to position 5, means, 3 minutes after the pump's pretreatment is started. Here are the setup tables create a time table for the valve, which starts with -5 minutes:

Pump

-5.04 min 1 ml/min
 -0.04 min 1 ml/min
 -0.02 min 0 ml/min
 0.00 min 0 ml/min

Valve

-5.04 min pos. 2

-2.00 min pos. 5

0.00 min pos. 5


Instrument Setup – Column Oven 4050

Fig. 192 Setup window, column oven 4050

For details on the spreadsheet operations, refer to the section *Creating a Time (Gradient) Program* as described in the chapter *Instrument Setup – Pumps, Pump Program*.

The oven program profile will be adapted automatically to the spreadsheet inputs.

Wait before injection

If you are working with a temperature program check this option. The setting **Wait with tolerance (+/- °C)** defines the accuracy of the oven temperature allowing a new injection.

A temperature independent delay between succeeding injections can be set in the field **Wait time (min)**. The system will wait for that time after having reached the initial temperature and before the next injection is performed.

Cooler Temperature (°C)

Enter the desired temperature.

Leakage Sensor Alarm (%)

Enter the desired sensitivity level.

Instrument Setup – Column Oven Jetstream

Make sure that the column oven Jetstream is switched on before opening the instrument window. Its keypad must display “choose function”. Otherwise the oven will not be controlled.




Fig. 193 Setup window, column oven Jetstream

For details on the spreadsheet operations, refer to the section *Creating a Time (Gradient) Program* as described in the chapter *Instrument Setup – Pumps, Pump Program*.

The oven program profile will be adapted automatically to the spreadsheet inputs (compare Fig. 192).

Wait before injection

If you are working with a temperature program check this option. The setting **Wait with tolerance (+/- °C)** defines the accuracy of the oven temperature allowing a new injection.

A temperature independent delay between succeeding injections can be set in the field **Wait time (min)**. The system will wait for that time after having reached the initial temperature and before the next injection is performed.

Instrument Setup – Flowmeter

The setup of the flowmeters is the same for all supported models is very simple.




Fig. 194 Setup window, Flowmeter

Run time

Enter the desired run time.

Acquisition delay

Here you can define a delay time between starting the run and the data acquisition.

Setting up Auxiliary Traces




Fig. 195 Setup window, Auxiliary Traces Any configured HPLC system with at least one pump with pressure transducer, a column oven or a flowmeter will be completed by an additional tab for auxiliary traces.

Only the selected traces will be stored with the run data.

Setting up a Trigger

Any configured HPLC system with at least one detector acquisition channel will be completed by an additional tab for the trigger.

Click on the **Trigger** tab to designate the trigger type and to setup the synchronization. The trigger type determines how the data sampling and the gradient program(s) are started.




Fig. 196 Setup window, trigger tab

Trigger Type

Select the type of synchronization from the drop-down list. The choices are **None**, **Manual**, and **External**.

- None** Time program(s) and Sampling start immediately after clicking on Start. Sequence acquisitions do not pause between runs.
- Manual** Operator must press **Enter** to start the run. Sequence acquisitions pause for confirmation between runs.
- External** Program execution and data sampling are started from an external trigger, if this option is selected. Further settings are not necessary. Sequence acquisitions pause for external trigger signal between runs.

Setting up the Baseline Check

This check is only available if you have activated it according to the next figure. For configuration of the system, refer to page 38.




Fig. 197 Activation of the baseline check option

Any configured HPLC system with activated baseline check option will be completed by an additional tab for this check. Click on the **baseline check** tab to define the check conditions.




Fig. 198 Defining the baseline check conditions

General

These selections will be valid for all active detectors or detector channels. Set the time for which the baseline must meet the set criteria **to pass** the check.

The check and the method will be stopped, if the system has not passed the check within the set time. This time must be at least equal to the set pass time.

Noise Test Method

The noise test will be performed either according to the root of the mean of squares (**rms**) method or to **ASTM**. Select your choice.

For any available detector channel of the configured system a separate line will appear in the spreadsheet. Here you can set the thresholds for the noise and drift for each channel separately. Any test (noise **and** drift) must be enabled separately for any channel.



To use this feature, select the **Perform Baseline Check** option in the **Single Run** dialog box or include **Baseline Check** in the current **Sequence** line.




Fig. 199 Activating the baseline check for a single run

Instrument Status of a (running) Control Method

If you have the Instrument Control option installed, and your instrument configured as a KNAUER HPLC System, the instrument control functions appear as part of the Instrument Status area of your instrument window.

To access the Instrument window where you can set up your devices and methods, double-click on the instrument icon you wish to start in the Main Menu. It may be required that you log-in before you can access the instrument window functions.

For details on other aspects of method development (calibration, integration) and sequence operations, see the ChromGate® Data System Reference Manual.

System Status

Open the Instrument Status Window using the menu sequence **Control – Instrument Status**.

The Instrument Status window contains the status tabs for all configured instruments, and as shown in the figure below the system status tab, providing an overview for the whole system.




Fig. 200 Instrument Status window, system status tab, method not running




Fig. 201 Instrument Status window, system status tab, method running

For some devices the software reads out the data only during a running method. Therefore, if no method is running, some data are not available (n/a) or might not be the actual data but rather data from a previous run. The running method is visualized by the dark dotted lines between the running instruments.

For each pump, identified by its name, the actual flow and pressure is displayed. If a gradient system is configured, the solvent composition will also be shown.

For each detector channel, the output values and as far as available the wavelength is shown. Similarly, the relevant data for all included instruments will be displayed, e.g. the actual switching positions in the case of the valves.

The single instrument tabs of this window provide more detailed information. They also enable direct system control. On these tabs you have the possibility of directly controlling the individual instruments. This is even possible while a method is running; however this option must be enabled when configuring the system (refer to the Knauer instrument control method options on page 170).




Fig. 202 Instrument status tab, example: Smartline Pump 1000

All of the instrument tabs provide communication status information at the bottom. The KNAUER Net status and the communication port assigned for communication is displayed.



Fig. 203 KNAUER Net instrument status

For instruments not controlled via Knauer net the serial port number is shown:



Fig. 204 KNAUER serial port instrument status

For instruments that will be controlled via LAN, the IP address is shown:



Fig. 205 KNAUER LAN instrument status

SN xxxxx

This field displays the serial number of the instrument which was entered in the configuration window. It should be the same as in the list in Fig. 206 where the serial numbers are read out from the instrument. The serial numbers serve as instrument identifiers; an error message will be generated in the case that the method is started with incorrectly entered serial numbers. In this case, the serial number must be corrected in the instrument configuration.

The active state of the communication with an instrument is represented by the green  symbol. Only in this case you will have access to the **Device Information**. Click the  button to display instrument(s) connected to that COM port. Information about the serial number, device name, and instrument software versions will appear:




Fig. 206 Information about the connected instrument(s)

Instrument Status – Pumps

The tabs for all types of pumps appear identical, independent of the features of the pump. The embedding of the pump into an HPG- or LPG-system will also not change the appearance of this window tab. Unavailable features and functions will only be disabled and inaccessible. In case of the Kontron pumps the GLP data and the Events are not displayed.




Fig. 207




Fig. 208

Instrument status tab, example: Kontron Pump 320

The tab consists of two main parts, the Monitor and the Direct Control areas.

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run. During a Pretreatment, the pump status is equilibrating.

Flow (ml/min)

The actual flow is indicated in this field.

Pressure Control (MPa)

The highlighted actual pressure is shown framed by the settings for the allowed minimum and maximum pressures.

Components (%)

In case of gradient systems (HPG or LPG) the percentage of components will be shown during the run. The components are also indicated by the (multi) colored bar. In the case of isocratic systems the bar is completely red colored, showing 100%.

Events (Digital Outputs)

The current state for each Event (Digital Output) will be displayed in this control. The **On** state is represented by a red active Light Emission Diode LED symbol (●), the **Off** state, including waiting for pulse, is represented by an inactive grey LED symbol (◐).

GLP

This section, only available with Knauer networking pumps, displays GLP (Good Laboratory Practice) related information transferred to or “downloaded” from the pump. Pressing the <Refresh>-button will actualize this information. The extent of the information depends on the given pump.

Direct Control

The direct control area will appear disabled as long as the direct control option is not activated. Except during a run it can be activated via the menu sequence METHOD – OPTIONS – RUNTIME SETTINGS..., (see Fig. 249 on page 170).

The direct control enables someone to control the instrument directly, even during a running method. Changes made to the settings and applied during a run will be automatically stored in the method. The corresponding new lines in the spreadsheets are marked with the comment DC Op (Direct Control Operation). Applied changes in a not-running method will not be saved. However, it will be recorded in the Activity Log.

Flow (ml/min)

Specify the flow in ml/min you would like to have at the moment and click the <Flow Apply> button to send the flow to the pump. If an LPG pump is used, also the gradient setting will be send.

<Flow Apply>

The Flow Apply button will send the flow, entered into the Flow field, to the pump.

<Stop>

The <Stop> button will interrupt delivery while the run time continues. The pump remains in slave mode. This button is also active during a run to stop the pump e.g. in an emergency situation. This will not stop the whole method, only the pump's flow! Due to, that this is a method change; the current method will be changed. The Stop will be added to the currently running method, if the option “Save changes in time table” of the “Direct control during a run” section of the **Runtime Settings** is enabled. Please refer the section “Runtime Settings” in the “Knauer Instrument Control Options” for more information.

<Flow & Gradient Apply>

Pressing this button will send the current flow and gradient settings for execution to the pump.

Gradient Components (%)

Specify the composition (% A, B ...), you would like to have at the moment and click the <Gradient Apply> button. The components field is applicable for the HPG (High Pressure Gradient) or LPG (Low Pressure Gradient) mode. Depending on how many components are configured, the corresponding number of fields (% A, B ...) will be accessible.

The <Gradient Modify> button is active only during a run. Clicking on it opens the following window.




Fig. 209 Gradient modify window

In the gradient spread sheet the elapsed lines are grayed. The red line shows the actual gradient status. The yellow ones are that which are open for modifying.

The Gradient Modify allows for changing the gradient table by adding new lines of freezing the gradient for a defined time.

To modify the gradient program, enter the time for the new line and the desired values in the fields on the top of the window. To change an already existing time line, enter this time in the "Time" field and the new values. Clicks on <Commit Line> insert the new or modified line into the gradient program. The changes by one or more committed lines becomes only active (send for execution to the pump) by a mouse-click on the button <Send Gradient>. You will see then the new line(s) in the gradient spread sheet.

Due to that the complete gradient in the pump will be replaced with a new one, the run time on the pump's display will start with time 0.00 min. This will not touch the method run time.




Fig. 210 Gradient table modify window

To delete any line of the gradient program, enter the time of the line to be deleted and press the <Delete Line> and then the <Send Gradient> button.

You can also **freeze** the current gradient status for a time you have entered by pressing the <Freeze Gradient> button.



Fig. 211 Gradient modify window freeze gradient

This will extend the pump table and therefore also the pump's run time.

In the following example the gradient (100% A / 0% B – 0% A / 100% B in 10 minutes) will be freeze for 2 minutes (from 3.34 – 5.34 minutes: at 3.34 minutes the freeze-command was sent). The gradient freeze will be sent and executed, if you click on the <Freeze Gradient>-button.




Fig. 212 Gradient modify window – freeze gradient for 2 minutes entered




Fig. 213 Gradient modify window – freeze gradient for 2 minutes executed

The pump table now has two new lines (3.34 and 5.34 minutes) with the frozen gradient settings. The gradient graphics on the right-handed side also shows the new gradient. The pump table now ends after 12 minutes instead of 10.

In this case you should prolong the set run time via the menu sequence **Control – extend run...**. In our example it is required to extend the run time by 2 minutes. You must close the *Set New Gradient* window to modify the method's run time.




Fig. 214 Extend run time window

All direct control settings will also be stored in the Instrument Activity Log. To see the log, click on “File – Instrument Activity Log – Display Log”.

User	Logged	Source	Activity
tester1	04/09/2012 10:56:28	PC-PM-SWTEST-1	Extending run by 2.00 minutes.
tester1	04/09/2012 10:54:30	PC-PM-SWTEST-1	Direct Control: S 1050: Gradient has been modified in 6.58 min aft
tester1	04/09/2012 10:51:15	PC-PM-SWTEST-1	Direct Control: S 1050: Gradient has been modified in 3.34 min aft
tester1	04/09/2012 10:47:25	PC-PM-SWTEST-1	Run Queue - Start Single Run - \\pc-pm-laboric\EZChrom Elite\Ent
tester1	04/09/2012 10:47:24	PC-PM-SWTEST-1	Run Queue - Complete Preview Run - \\pc-pm-laboric\EZChrom El
tester1	04/09/2012 10:47:21	PC-PM-SWTEST-1	Run Queue - Add Single Run - \\pc-pm-laboric\EZChrom Elite\Ente
tester1	04/09/2012 10:42:38	PC-PM-SWTEST-1	Run Queue - Start Preview Run - \\pc-pm-laboric\EZChrom Elite\E
tester1	04/09/2012 10:42:37	PC-PM-SWTEST-1	Run Queue - Add Preview Run - \\pc-pm-laboric\EZChrom Elite\E
tester1	04/09/2012 10:41:36	PC-PM-SWTEST-1	User tester1' logged in.

Fig. 215 Extend run time window

Control Pressure Limits

This section allows you to enter the control pressure limits in the units chosen in the pump configuration dialog. The pump will automatically shut down its operation when the measured pressure is outside of the limits.

Min Pressure

Enter the minimum pressure allowable for your pump.

Max Pressure

Enter the maximum pressure allowable for your pump.

Apply

Click the button to actualize the selected Min/Max pressure limits.

Events

Click the  button to set up the event outputs. Refer to the chapter *Instrument Setup – Pumps*, section *Configuring Events* for details. Click the **Events Apply** button to accept the settings. This status area is only shown for pumps offering this feature.

Purge




Fig. 216 Purge settings window

Click the **Purge** button to open the purge set up window. In **Purge flow (ml/min)** enter the desired flow rate. Click on the **Start** button to start the pump. In the Gradient components section you can select each solvent channel directly or a “Mix” (ABCD). You can set a time limit for the purge. Enable the **Stop after** option and enter the desired purging time.

Instrument Status – Detectors

Instrument Status – RI Detectors (S 23[4]00, K-23[4]00/1)

The tabs for all RI detectors appear identical, independent of the features of the detector. Unavailable features or functions are only grayed and inaccessible.




Fig. 217 RI-Detector status tab, example: S 2300 / S 2400

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Signal (µRIU)

The actual signal is indicated in this field. The unit corresponds to the setting in the configuration.

Direct Control

The direct control option must be activated (see Fig. 249 on page 170).

The buttons **Flush** and **Autozero** can be pressed to perform the corresponding functions. The flush function (flushing for 1 min) is deactivated during a run. The Autozero defines the current signal as a zero level.

Instrument Status – UV Detectors

(S 2520, S 2500, S 200, K-200, K-2000/1, K-2500/1)

The tabs for most of the UV detectors appear identical, independent of the features of the detector. Unavailable features or functions are only grayed and inaccessible.



Fig. 218 UV-Detector status tab, example: S 2500

The only exception is the Direct Control Area of the fixed wave detectors.



Fig. 219 UV-Detector status tab, Direct Control: S 200

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Signal (μ AU)

The actual signal is indicated in this field.

Wavelength (nm)

The actual wavelength is indicated in this field.

Direct Control

The direct control option must be activated (see Fig. 249 on page 170).

D2 Lamp

The On state of the lamp is represented by an active green LED symbol (●), the Off state is represented by an inactive grey LED symbol (○). Clicking on the **Lamp Off** button will (after inquiry) switch the lamp off and change the button label to **Lamp ON**.

Autozero

Click this button to define the current signal as the baseline zero level.

Wavelength (nm)

Specify the wavelength you would like to have at the moment and click the **Apply** button.




Fig. 220

Channel 1 / 2

The actual signal in the units which have been chosen in the detector configuration dialog at the actual wavelength in nm is displayed for each active channel.

Communication – instrument

This field displays the serial port of the detector which is used to control the instrument. The active state of the communication with an instrument is represented by the green  symbol. If the serial port was selected as <None>, then no real control of the detector takes place.

Communication - A/D interface

The active state of the communication with an A/D interface is represented by the green  symbol. Only in this case you will have access to the **Device Information**. Click the  button to display the instrument(s) connected to that COM port. Information about the serial number, device name, and instrument software versions will appear:




Fig. 221 Information about the connected instrument(s)

Direct Control

Acquisition Channels

The number of channels that will be activated depends on how many possibilities are given to set new values for the measuring wavelengths. Pressing the Apply button forces the wavelength change on the detector channels.

Acquisition Parameters

Select the **Resp. time** from the drop-down list. Smaller Resp. time (or time constants) are required for high speed analyses. A good rule is to

select a time constant, which is at least one tenth of the peak width at the base of the narrowest peak of interest.

The **range** can be set to define the range of the analog detector outputs. The choice is 0.100, 0.200, 0.500, 1.000, and 2.000 **Absorption Units Full Scale (AUFS)**.

The Kontron x30 and x35 detectors can be used in a dual channel mode. As far as the option Used Channel #2 is checked, the otherwise not accessible **cycle time** becomes active for changing. The cycle time you can set within the range of 0.40 through 1.00 seconds with 0.05 s increments.

For the x32 detectors the option Used Channel #2 is not available and the cycle time is not displayed at all.




Fig. 222

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Digital Outs (Events)

The current state for each digital output will be displayed in this control. The On state is represented by an active red Light Emission Diode LED symbol (●), the Off state, including waiting for pulse, is represented by an inactive grey LED symbol (◐).

Channel 1 / 2 / 3 / 4

The actual signal in the units which have been chosen in the detector configuration dialog at the actual wavelength in nm is displayed for each active channel.

GLP

The GLP data of the detector are displayed in this area. The data will be actualized by pressing the **Refresh** button.

Direct Control

The direct control option must be activated (see Fig. 249 on page 170).

Lamp

The **On** state of the lamp is represented by an active green LED symbol (●), the **Off** state is represented by an inactive gray LED symbol (◐). Clicking on the **Lamp Off** button will (after inquiry) switch the lamp off and change the button label to **Lamp ON**.

Autozero

Click this button to define the current signal as a zero level.

Scan, nm

This section is disabled for changes. It only provides information on how the scan spectral interval is specified. Changes must be made from the instrument setup window (Fig. 150 on page 102). It is not possible to modify this function during a run.

Digital Channels (nm)

The number of channels that will be activated depends on how many possibilities are given to set new values for the measuring wavelengths. Pressing the WL Apply button forces the wavelength change on the detector channels.

Digital Outs (Events)

Click the  button to set up the digital outputs. Refer to the chapter *Instrument Setup – Pumps, Pump* section *Configuring Digital Outputs* for details. Click the **DO Apply** button to accept the settings. For detectors without the feature of digital outputs this area is grayed and inaccessible.



The D₂ lamp should be switched on manually from the direct control window or using the Wake-Up option from the K-2600 before starting the next working session. After power up, allow 30 minutes for the lamp to stabilize before running samples that may require sensitive measurements.

Instrument Status – Diode Array Detectors

(S 2600, DAD 2850, DAD 2800, K-2700 and Kontron DAD 540/545)




Fig. 223 DAD-Detector status tab, example: DAD 2600

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Channel 1 / 2 / 3 / 4

The actual signal in the units which have been chosen in the detector configuration dialog at the actual wavelength in nm is displayed for each active channel.

Direct Control

The direct control option must be activated (see Fig. 249 on page 170).

D2 Lamp (HAL Lamp)

The **On** state of the lamp(s) is represented by a green active LED symbol (●), the **Off** state is represented by an inactive gray LED symbol (●). Clicking on the **Lamp Off** button will (after inquiry) switch the lamp off and change the button label to **Lamp ON**. The appearance of the HAL lamp control depends on the individual detector and its configuration.

Autozero

Click this button to define the current signal (over the specified spectrum) as the baseline zero level.

Digital Channels (nm)

The number of channels that will be activated depends on how many possibilities are given to set new values for the measuring wavelengths. Pressing the WL Apply button forces the wavelength change on the detector channels.

Events (Digital Outs)

Click the  button to set up the digital outputs. Refer to the chapter *Instrument Setup – Pumps, Pump* section *Configuring Digital Outputs* for details. Click the **Events Apply** button to accept the settings. For detectors without the feature of digital outputs this area is grayed and inaccessible.



The D₂ lamp should be switched on manually from the Direct Control window before starting the next working session. After power up, allow 30 minutes for the lamp to stabilize before running samples that may require sensitive measurements.

Info

Click this button to display the device specific data in the detector information tab.




Fig. 224 Detector information window

Diagnostics

Click this button to use the diagnostic features of the ChromGate[®] software. The Diagnostics window appears which allows you to access information and to control important parameters and modules of the device.




Fig. 225 Diagnostics window of Diode Array Detectors

Integration Time [Range 12 –1000] msec

The integration time is the time interval when the detector “accumulates” light to receive one spectrum. This value is a measure of the intensity of light originating from a lamp. The more light that passes the optical scheme and the cell, the shorter the integration time will be. The integration time depends also on the specific spectral interval (Scan) selected for the measurement, since the maximum of intensity may be varied. The integration time range is defined by the device parameters (see K-2700 / 2800 device specific data). Click the button to determine the optimal integration time. You can use this measurement to check quality of the lamp and used optical elements. If the software cannot find an optimum integration time, the lowest time of the range (i.e. 12) will be shown.

The integration time will be checked before every run. If there not enough light can be detected, a single run will be aborted with an appropriate error message. In a sequence run, the sequence will not be aborted and the data will be acquired, but ChromGate will enter a message into the instrument activity log, that there was not enough light.

Lamp

Select the lamp you want to check. This option is not enabled for the K-2700 detector. The choices are Deuterium and Halogen lamps.

Scan, nm [Range 189 – 1024 nm]

This area is used to specify spectral parameters used for your test measurements. The range displays the minimal and maximal wavelengths available for the detector. The minimal and maximal wavelengths are the limits which are defined by the configuration of the detector.

The **Start of Scan** is the minimal wavelength with a default value of 200 nm. The **End of Scan** is the upper limit of the scans to be measured with a default value of 500 nm.

The **Bandwidth** is a measure of the range over which the calculations and measurements are used to calculate the intensity for a particular wavelength. For example, with an 11 nm bandwidth, the intensity for 250 nm would be calculated from 245 nm to 255 nm. The wider bandwidth

results in an increased signal to noise ratio because the multiple measurements cause the random noise to be cancelled out. That is, if the noise is truly random, more measurements will decrease the noise because for each positive excursion in the noise will have a good likelihood of being cancelled by a negative excursion. But the signal, since it is not random, will remain the same over multiple measurements. Bandwidth is a tradeoff, with a high value decreasing the noise but also decreasing the possible resolution which can be measured. 5 nm is the default value, but if your peaks have substantially wider UV spectra you can increase the bandwidth value in ChromGate® to obtain a lower noise level.

Mode

Select the signal mode for your measurements from the drop-down list. The choices are **Absorption** and **Intensity**.

Open/Close Shutter

A click on this button will open/close the shutter. The label will be changed accordingly.

Filter On/Off (S 2600 only)

You can click on this button depending on whether or not you are using an edge filter.

Dark Signal

Click the button to measure the dark current. This operation is enabled only for devices equipped with a shutter.

Set Reference

Click the button to measure a spectrum to be used as a reference spectrum.

Spectrum

Click the button to measure an absorption or intensity spectrum, depending on the signal mode selected for the measurement. The measured spectrum will be shown in the spectrum view.

Instrument Status – Fluorescence Detector RF-10AxI / RF-20A




Fig. 226 Detector Status tab, Fluorescence Detector RF-10AxI

Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

The actual excitation wavelength, emission wavelength, signal value, sensitivity, gain, and the selected output mode will be displayed.

Direct Control

All settings you have entered in this area will become effective after clicking on the **Set** button. To achieve changes during a run, the direct control option must be activated (see Fig. 249 on page 170).

Excitation wavelength, Emission wavelength

You can enter the values of the excitation and emission wavelengths.

Sensitivity, Gain, Recorder Range, Output Mode

The time constant, sensitivity level (HIGH, MED, or LOW), the gain (1, 4, 16), recorder range, and the output mode (Emission NOT Corrected, Excitation Energy, Light Source Corrected) can be selected from the corresponding pull down menus. Press the **Set** button to send the parameters to the detector.

Autozero

Activate this option to perform an autozero at the start of the diagnostic measurement. Once the Autozero button has been pressed the autozero is performed immediately.

Event Out, Event Mark

Event Out short circuits the event out terminal for approximately 1 second. Event Mark adds an event mark to the analog output.

Diagnostics

Click this button to use the diagnostic features of the ChromGate[®] software. The Diagnostics window appears which allows you to access information and control of important parameters and modules of the device.




Fig. 227 Diagnostics window of Fluorescence Detector RF-10AxI with an anthracene spectrum

Scan, nm [max range 200 – 900 nm]

This area is used to specify spectral parameters used for your test measurements. Set the **Start** and **End** values within the given limits according to your needs.

Spectrum type

Select the excitation or emission spectrum from the pull-down menu. Depending on your choice (the field description will change with the selection) you can enter either the **Emission wavelength** or the **Excitation wavelength**.

Scanning speed

Select a value between 24 and 3000 nm/min from the pull-down menu.

Sensitivity, Gain, Output Mode

Select the sensitivity level (HIGH, MED, or LOW), the gain (1, 4, 16), and the output mode (Emission NOT Corrected, Excitation Energy, Light Source Corrected) from the corresponding pull-down menus.

Analog output

An analog output can be activated. The **Recorder range** is combined with this option. You can define this range by selecting one from the pull-down menu.

Scan

A click on this button will start the diagnostic measurement. It will take some time depending on your scan range and speed settings. During the scan, the info "please wait" will blink in the head bar. The spectrum will be displayed after finishing the scan.



The obtained diagnostic spectrum cannot be saved directly or exported! For later comparisons you must take a screenshot and save this in an appropriate way.

Instrument Status – Conductivity Detector Alltech 650




Fig. 228 Detector Status tab, Conductivity Detector Alltech 650

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

The actual cell temperature, the conductance and total conductance as well as the range will be displayed.

Direct Control

All settings you have entered in this area will become effective after clicking on the **Apply** buttons. To achieve changes during a run, the direct control option must be activated (see Fig. 249 on page 170). Also an autozero can be initialized.

Instrument Status – User Defined Detector

The detector status tab of a user defined detector displays only the monitor area. No direct control is possible.

Monitor




Fig. 229 Detector status tab, user defined detector

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Communication - A/D interface

The active state of the communication with an instrument again is represented by the green  symbol. Only in this case you will have access to the **KNAUER Net Device Information**. Click the  button to display a list of Knauer networking instruments which are connected to that COM port. A list containing serial numbers, device names, and instrument software versions for Knauer networking instruments will appear:




Fig. 230 Device information for connected Manager 5000 interface module

Instrument Status – Virtual Detector




Fig. 231 Detector status tab Virtual detector

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Signal (mA)

The actual signal is indicated in this field. The units correspond to the setting in the configuration of the user defined or virtual detector.

Direct Control

Signal / Apply

You can set the detector to a desired signal value for the time being.

External Trigger ON

Pressing this button will give the trigger signal to your system if it is in the position waiting for trigger.

Instrument Status – Assistant ASM2.1L

The Assistant status window allows for checking the current status of all configured modules of the Assistant and the direct control.

The status and direct control options on this window depends on the configured modules. Below the possible status options are described.




Fig. 232 ASM2.1L status tab – valve, valve, pump

Monitor

Run Status

The status **On** will be displayed when the ASM 2.1L is ready to run (pump is found, valves have been initialized, detector lamp is on). The status **On** is also shown, if a Pretreatment is running (please refer to Pretreatment section on page 89). During a run the status **Run** and the current run time in minutes will be shown. The **Off** status will be shown, if no connection with the ASM 2.1L can be established or if the detector lamp is switched off.

Events (Digital Outputs)

The current state for each Event (Digital Output) will be displayed in this control. The **On** state is represented by a red active Light Emission Diode LED symbol (●), the **Off** state, including waiting for pulse, is represented by an inactive grey LED symbol (◐).

Pump

The left-handed part of the **Pump** section shows the pump status as flow in ml/min and pressure (depending on the selected unit in MPa, bar or psi).

Flow [ml/min] (Status)

The flow shown here is the current pump flow set either by direct control or a method.

Pressure [unit] (only for P4.1S) (Status)

The shown pressure is the current system pressure, measured by the pump's pressure sensor. The unit depends on the pressure unit selected in the pump's configuration (MPa, bar or psi). Please note, that also a pressure may be measured if the pump is stopped, e.g. if another pump is running or if the flow path is blocked. If only a pump without a pressure sensor is configured, the pressure section is shown grayed out.

The right-handed part of the **Pump** section allows for direct control of the pump, independently from the method setup. During a run this part is not accessible.

Flow [ml/min] (Direct Control)

Enter the desired flow for the pump and click on the <Apply> button to start the pump or change the current flow. If the flow rate is out of the acceptable flow range, defined by the mounted pump head, an appropriate error message will be shown.

Click on the <Stop> button to set the flow rate to 0 (zero).

Pmin / Pmax (only for P4.1S) (Direct Control)

Enter the desired values for the allowed minimum and maximum system pressure. The pump will immediately stop, if the maximum pressure value is exceeded or if the pressure is below the minimum value for 30 seconds. The unit depends on the pressure unit selected in the pump's configuration (MPa, bar or psi). Please check for the pressure unit shown in the Status section of the pump.

Click on the <Apply> button to send the pressure limits to the pump. If the pressure limits are out of acceptable pressure range, defined by the mounted pump head, an appropriate error message will be shown.

HPG (only if an HPG is configured)

HPG			
Flow (ml/min):	1,000	Flow (ml/min): <input type="checkbox"/>	A (%) <input checked="" type="checkbox"/> <input type="button" value="Stop"/>
Components (%):	100 : 0	<input type="text" value="0,000"/>	<input type="text" value="60"/> <input type="button" value="Apply"/>
Pressure (MPa):	11.5	PMin: <input type="text" value="0,0"/>	PMax: <input type="text" value="0,0"/> <input type="button" value="Apply"/>

Fig. 233 ASM2.1L status tab – HPG

The left-handed part of the **HPG** section shows the pump status as flow in ml/min, the gradient components in % and pressure (depending on the selected unit in MPa, bar or psi).

Flow (Status)

The flow shown here is the current flow of both pumps, set either by direct control or a method.

Components (Status)

The Components show the current pump gradient, set either by direct control or a method.

Pressure (only one pump has a pressure sensor) (Status)

The shown pressure is the current system pressure, measured by the pump's pressure unit. Please note, that also a pressure may be measured if the pump is stopped, e.g. if another pump is running or if the flow path is blocked.

The right-handed part of the **HPG** section allows for direct control of the pumps, independently from the method setup. Please note, that only %A can be entered, the %B component of the gradient will be completed automatically by the software. Only a parameter, the check box is enabled for, can be edited and the changed value will be send by clicking the <Apply> -button.

Flow (*Direct Control*)

To enter a flow, enable the Flow check box. If enabled, you can enter a flow, which will be the resulting flow of both pumps. Click on the upper <Apply> button (same line as for the Flow and A(%) values) to start the pumps or change the current flow. If the flow rate is out of the acceptable flow range, defined by the mounted pump head, an appropriate error message will be shown. If the check boxes for **Flow** and **A(%)** have been enabled, the values for both will be send to the pumps while clicking the <Apply> button. If you wants to update only the flow or gradient, please disable the A(%) check box.

Click on the <Stop> button to set the flow rate to 0 (zero).

A(%) (*Direct Control*)

If the corresponding check box is enabled, you can change the gradient composition. Enter the desired value for HPG pump A and click on the upper <Apply> button (same line as for the Flow and A(%) values) to send the new gradient composition to the pumps. The program will complete automatically the value for HPG pump B. If the check boxes for **Flow** and **A(%)** have been enabled, the values for both will be send to the pumps while clicking the <Apply> button. If you wants to update only the A(%), please disable the flow check box.

Pmin / Pmax (only for P4.1S) (*Direct Control*)

Enter the desired values for the allowed minimum and maximum system pressure.

Click on the <Apply> button to send the pressure limits to the pump. If the pressure limits are out of acceptable pressure range, defined by the mounted pump head, an appropriate error message will be shown.

UV

The left-handed part of the **UV** section shows the detector status as Deuterium lamp status, wave length in nm and signal. The signals unit depends on the unit selected in the detector's configuration (AU, mAU, μ AU).




Fig. 234 ASM2.1L status tab – UV

D2 Lamp (*Status*)

The button shows the current lamp status. The **On** state is represented by a green active Light Emission Diode LED symbol (●), the **Off** state, including initializing (heating), is represented by an inactive grey LED symbol (◐).

WL: (*Status*)

WL shows two different values: the current wave length and the signal. The wave length's unit is nm; the signal unit depends on the unit that has been chosen in the ASM2.1L configuration for the detector.

The right-handed part of the **UV** section allows for direct control of the detector, independently from the method setup. During a run this part is not accessible beside the autozero (AZ) button.

WL: (*Direct Control*)

Enter the desired wave length for the detector click on the <Apply> button to change the wave length in the detector. If the wave length is out of the acceptable wave length range, an appropriate error message will be shown.




Fig. 235 Error message Wavelength out of range

Lamp OFF / Lamp On button (*Direct Control*)

This button allows for switching off or on the detector lamp. The button's labeling indicates, what will happen if the button is pressed; it does not show the lamp status! If the lamp is on, the button shows "Lamp OFF", because the lamp will be switched off, if you click on it. It shows "Lamp ON" if the detector lamp is off – clicking the button will switch on the lamp. While the deuterium lamp is heating, the (inaccessible) button shows "Heating".

AZ button (*Direct Control*)

If you click on this button, the detector will perform an autozero, which means, the current signal will be set to 0 (zero). This button will also be active during a run, please be careful then: Clicking this button if a peak elutes will make the peak information (area, maybe also retention time) unusable.

Valve(s)

The left-handed part of the **Valve(s)** section shows the valves current position.




Fig. 236 ASM2.1L status tab – Valve(s)

Position: (*Status*)

The current position of all configured valves will be shown. For each valve the name given in the configuration is shown. If the AMS2.1L is configured using the "Auto Configuration" option, the default name for a valve is "Valve", followed by an extension for the valve's position in the ASM2.1L case (L = left, M = middle, R = right). For Knauer two-position valves the position can be either I (injection) or L (load), for Valco two-position valves A or B. For multi-position valves the position can be 1 – 16, depending on the valves configuration.

The right-handed part of the **Valve(s)** section allows for direct control of the valve(s), independently from the method setup. During a run this part is not accessible.

Click on the small arrow key  of the drop-down menu. The menu shows a list of possible valve positions. For Knauer two-position valves the positions are I (injection) and L (load), for Valco two-position valves A and B and for multi-position valves 1 – 16, depending on the valves configuration. That means, for a 6-position valves the available positions are 1 – 6, for a 16-position valve 1 – 16.

GLP

The GLP section of the ASM 2.1L status window allows for reading-out the GLP data of all modules of the ASM 2.1L.




Fig. 237 ASM2.1L status tab – GLP (UV)

Click on the <Refresh> button to read-out the GLP information of the built-in modules. You can use the scroll bar to scroll through the GLP data. The following GLP information will be shown:

Instrument data

Serial number
Firmware version
Instrument installation date
Last service date
Last service code
Device working time (in hours)

Pumps

Pump head
Pump operation time

UV

Lamp operation time (in hours)
Number of lamp ignitions
Optical bandwidth (in nm)
Light intensity I-Sig at UV-max
Light intensity I-Ref at UV-max

Valve(s) (for Knauer and Valco valves)

Valve type (number of ports and positions)
Switching cycles

System

The System section of the ASM 2.1L status window allows for changing the Events of the ASM 2.1L and the access for the UV Diagnostics, if an UV detector is configured.




Fig. 238 ASM2.1L status tab – System with UV Diagnostics

Events: (Status)

The three digits show how the events have been set in the “*Events / Digital Output Control*” menu (This may not be the current status of the Events, if this settings have not been send by clicking the <Apply> button. The current status of the events will be represented by the LED symbols on the top of the ASM 2.1L status window.)

Click on the  button to setup the events (digital outputs). Refer to the chapter *Instrument Setup – Pumps, Pump section Configuring Digital Outputs* for details. Click on the <Apply> button to send the settings to the device.

Pump Purge ->

Click on the  button to open the purge setup window. Please refer to the chapter **Instrument Status – Pumps** for more information.

UV Diagnostics ->

Click on the  button to open the UV diagnostics window. Please refer to the chapter **Instrument Status – Detectors** for more information.

Instrument Status – Autosampler

Only one autosampler can be included in any system, therefore all supported samplers will be explained separately without cross references

Instrument Status – Autosampler 3800

Smartline Autosampler 3800 and K-3800 (Basic Marathon)




Fig. 239 Autosampler 3800 status tab

Monitor**AS Status**

The status **Ready** will be displayed when the autosampler is waiting to perform an injection. If an injection program is running, the actual program step will be indicated. After a completed injection, the status is described as Analysis timer running during the chromatographic run.

Current Vial

The actual position of the injection needle is shown.

Tray Temperature (°C)

If the cooling option is installed and configured the actual tray temperature will be shown.

Aux.

The auxiliary output(s) status of the AS 3800 is indicated. The **On** state is represented by an active red LED symbol (●) and the **Off** state is represented by an inactive LED symbol (○).

Direct Control**Tray Cooling**

The tray cooling can be switched **On** or **Off** as long as the cooling option has been included and configured. The AS 3800 does not support to set a temperature.

Advance Tray

Clicking briefly on this button will move the tray to the next position. By holding down the **Advance Tray** button, the tray will rotate until the button is released.


Instrument Status – Autosampler Optimas/3900

Fig. 240 Autosampler 3900 status tab

Monitor**AS Status**

The status **Ready** will be displayed when the autosampler is waiting to perform an injection. If an injection program is running, the actual program step will be indicated. After a completed injection, the status is described as Analysis timer running during the chromatographic run.

Current Vial

The actual position of the injection needle is shown.

Oven Temperature (°C)

If the column oven option has been configured the actual column oven temperature will be shown.

Tray Temperature (°C)

If the tray cooling option has been installed and configured the actual tray temperature will be shown.

Direct Control

The direct control features are not accessible during an injection.

Oven Temperature (°C)

If the column oven option is installed and configured, the column oven temperature can be set. The entered value will be transferred to the oven

by clicking the **Apply** button. The **On / Off** button can be used for enabling / disabling the heating control. The labeling of the button changes with the actual status.

Tray Temperature (°C)

If the cooling option is installed and configured the tray temperature can be set. The entered value will be transferred to the sampler by clicking the **Apply** button. The **On / Off** button can be used for enabling / disabling the cooling control. The labeling of the button changes with the actual status.

Start Needle Wash

Press this button to force the autosampler to perform the needle wash.

Advance Tray

Clicking briefly on this button will move the tray to the next position.

Instrument Status – Autosampler 3950




Fig. 241 Autosampler 3950 status tab

Monitor

AS Status

The status **Ready** will be displayed when the autosampler is waiting to perform an injection. If an injection program is running, the actual program step will be indicated. After a completed injection, the status is described as Analysis timer running during the chromatographic run.

Current Vial

The actual position of the injection needle is shown.

Tray Temperature (°C)

If the tray cooling option has been installed and configured the actual tray temperature will be shown.

Direct Control

The direct control features are not accessible during an injection.

Tray Temperature (°C)

If the cooling option is installed and configured the tray temperature can be set. The entered value will be transferred to the sampler by clicking the **Apply** button. The **Off** button can be used for disabling the cooling control. To switch it on again set the required temperature and press the **Apply** button.

Reset Transport Vials (only with 84+3 vials tray)

If the tray type **84+3 vials** is configured, the button Reset Transport Vials allows for reset the transport liquid level of all transport vials (vials 85, 86, 87) to the default value of 8000 µl. If you click on that button, the software will ask if you wants to reset the filling level of all transport vials to 8000 µl. Click <Yes> to do so or <No> to leave the window without a level reset. Please be sure, that you also fill the transport vial(s), because the autosampler adjusts the needle immersion depth dependent on the vials filling level. If the vials will not fill, the autosampler can aspirate air instead transport liquid from the transport vial. The software cannot read-out or set the filling level. The transport vial filling level will automatically be reset, if the autosampler will be switched on.

Start Needle Wash

Press this button to force the autosampler to perform the needle wash.

Move Tray

Clicking briefly on this button will change the tray position:
FRONT ←→ HOME.

[>>...]

The Stacked Injection table can be changed during a run. To have access, the option **Direct Control during a run** in *Method – Runtime settings* must be enabled. Pressing [>>...]-button during a run will open the known Stacked Injection setup with the default values. For the re-programming of the Stacked injection table, please refer the instrument setup chapter of the autosampler 3950. If you click on the **Apply**-button, the new stacked injection table will be executed immediately; the program will not check when the last injection was made.

If the **Direct Control during a run**-option **Save changes in time table** is enabled, all program lines for stacked injections, which have not been executed yet during the current run, will be replaced by the stacked injections from the direct control, whereas lines for stacked injections, which have already been executing during the current run, remain untouched.

If the **Direct control during a run**-option **Save changes in time table** is not active, the new table will be executed, but not stored in the instrument setup. However, in this case the method does not include the correct injection program.

[>]

Pressing the [>]-button opens a window, that allows you for perform several service operations.

Tray position

You can move the tray to several positions. The choices are Tray Home, Needle exchange position and Tray Front.

Syringe position

The syringe can be moved to Home, End or Exchange position.

Valve position

The valve can be switch to Inject and Load position.

Instrument Status – Triathlon/Endurance Autosampler




Fig. 242 Autosampler Triathlon status tab

Monitor

AS Status

The status **Ready** will be displayed when the autosampler is waiting to perform an injection. If an injection program is running, the actual program step will be indicated. After a completed injection, the status is described as Analysis timer running during the chromatographic run.

Current Vial

The actual position of the injection needle is shown.

Tray Temperature (°C)

If the tray cooling option has been installed and configured the actual tray temperature will be shown.

Direct Control

The direct control features are not accessible during an injection.

Tray Temperature (°C)

If the tray cooling option has been installed and configured the tray temperature can be set. The entered value will be transferred to the sampler by clicking the **Apply** button. The **On / Off** button can be used for enabling / disabling the cooling control. The labeling of the button changes with the actual status.

Start Initial Wash

Press this button to force the autosampler to perform the needle wash.

Move Plate (Endurance only)

Clicking briefly on this button will move the tray to the next position.

Advance Tray (Triathlon only)

Clicking briefly on this button will move the tray to the next position. By holding down the **Rotate Tray** button, the tray will rotate until the button is released.


Instrument Status – Kontron Autosamplers

Fig. 243 Kontron Autosamplers status tab

Monitor**AS Status**

The status **Ready** will be displayed when the autosampler is waiting to perform an injection. If an injection program is running, the actual program step will be indicated. After a completed injection, the status is described as Analysis timer running during the chromatographic run.

Current Vial

The actual position of the injection needle is shown.

Direct Control**Auxiliary**

The auxiliaries 1 or 2 (1-4 for 460/465 autosamplers) can be selected and activated.

Position

The tray can be moved directly to any entered position.

Initial Wash

Press this button to force the autosampler to perform the needle wash

Advance Tray

Clicking briefly on this button will move the tray to the next position.


Instrument Status – Miscellaneous Instruments**Instrument Status – Manager 5000/5050/IF2 I/O**

Fig. 244 Manager 5000/5050 / IF2 status tab

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed since the start of the run.

An. Output (mV)

The actual output signal is indicated in this field.

Digital Output (Event)

The current state for the digital output will be displayed in this control. The On state is represented by an active red Light Emission Diode LED symbol (●). The Off state is represented by an inactive LED symbol (○).

Direct Control

The output voltage **AO(mV)** and the percentage factor values of the **AOF(%)** can be set. The settings become effective by clicking the **Apply** button, resulting in an output voltage = $AO * AOF / 100$ at the given channel.

The digital output defines the event possibilities **0**, **1**, and **Pulse**, selectable from the pull-down menu. The setting becomes effective by clicking the **Apply** button.

Instrument Status – Knauer Switching Valves

Run Status

The Run Status displays a run time in minutes elapsed since the start of the run. While the method is not running the status HOLD will be displayed.

Valves Status, Direct Control

The actual position for each valve is displayed in the **Position** row.

To switch the valves specify the required positions and click the **Set Position** button. The choices are **BLANK** and **I** (Inject) and **L** (Load) for a 2-position valve, **1** through **6** for a 6-position valve, **1** through **12** for a 12-position valve, and **1** through **16** for a 16-position valve. The BLANK space means that the valve will not change its state after pressing the Set Position button.




Fig. 245 Status tab, switching valves

If there is for any reason the communication to a valve interrupted, the corresponding field for the valve position is highlighted and the position is indicated as zero.

Instrument Status – Column Oven 4050 and Jetstream




Fig. 246 Status tab, column oven 4050



Fig. 247 Status tab, column oven Jetstream



To achieve the communication with the column oven Jetstream its control unit must display *choose function* as it does after powering up. Otherwise communication can be established, but the oven will not execute commands sent by software.

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed since the start of the run.

State

The state **On**(Idle) or **Off** will be indicated if the oven will be controlled or not.

Temperature (°C)

The actual oven temperature will be shown.

Cooler Temperature (°C) (S 4050 only)

The actual cooler temperature will be shown.

Digital Outs (Events) (S 4050 only)

The current state for the digital outs will be displayed in this control. The On state is represented by an active red Light Emission Diode LED symbol (●). The Off state is represented by an inactive LED symbol (○).

Direct Control

Temperature (°C)

The oven temperature can be set. The entered value will be transferred to the oven by clicking the **Set** button.

Cooler Temperature (°C) (S 4050 only)

The oven temperature can be set. The entered value will be transferred to the oven by clicking the **Set** button.

Leakage Alarm Level (%) (S 4050 only)

The leakage alarm level can be set. The entered value will be transferred to the oven by clicking the **Set** button.

Digital Outs (Events) (S 4050 only)

Click on the  button to set up the digital outputs. Refer to the chapter described in the chapter *Instrument Setup – Pumps*, section *Configuring Digital Outputs* for details. Click the **DO Apply** button to accept the settings.

Set Temperature Control (S 4050 only)

Press the **Off** button to switch the temperature control off. It will be switched on again by setting the desired temperature and pushing the **Set** button.

GLP (S 4050 only)

The minimum and maximum reached temperatures, the total working time and the working time above 150°C will be monitored. The read out will be actualized by clicking on the **Refresh** button.

Instrument Status – Flowmeter



Fig. 248 Status tab, Knauer Flowmeter

Monitor

Run Status

The status HOLD will be displayed when no program is running; otherwise the Run Status displays a run time in minutes elapsed from the start of the run.

Flow (ml/min)

The actual flow will be shown in units as defined in the instrument setup.

Knauer Instrument Control Method Options

General Settings

Some general method settings belong exclusively to the Knauer instrument control. They are accessible by the METHOD menu.

Runtime Settings

The runtime setting options are general settings. The selection is valid for all ChromGate methods until it is changed.




Fig. 249 Knauer Instrument Control Method Option - Runtime Settings

Emergency during a run options

In the case that an emergency situation is detected, the software can be forced to stop all pumps and/or to switch off all lamps and/or switch off all ovens.

If enabled, this option is valid for the currently used computer.

Direct Control during a run

If this option is enabled, all devices can be controlled manually using the Direct control section of the device's instrument status tab. The direct control during a run option is a valuable tool for method development but can also lead to the possibility for unintended changes to the method. For this reason, it is recommended that this option not be **enabled**, especially for routine methods. All changes during a run will be added to the instrument activity log, but cannot be stored in the original/acquisition method, attached to the data file. **This may not meet your GLP/GMP requirements.**

Save Changes in time table

This option can only be enabled/disabled, if the "Enabled" checkbox for "Direct Control during a run" is activated. If the "Save changes..."-Option is checked, any method changes by direct control will be saved in the method. Otherwise the changes only can be monitored in the instrument activity log and in the audit trail, refer to the EZChrom reference manual. Due to that the changes cannot be stored in the method, which is attached as the original/acquisition method to the data file; we recommend reanalyzing the chromatogram with the changed method. All methods used to analyze a chromatogram, also will be added to the chromatogram data file and can be extracted as "From Results" in the *Open Data* dialog.

<Advanced...>

If you click the **<Advanced>** button, three additional options are available.




Fig. 250 Knauer Instrument Control Method Option - Advanced Runtime Settings

Check run time when saving method

ChromGate allows different run times for all devices. If you want to have the same run time for all devices and let the software check if, e.g. after a change, all devices have the same run time, enable this option. If there are different run times, the software displays an error message with the run time of all devices. The run times must be equalized before you can save the message.




Fig. 251 Error message Run time check failed



If enabled, this option is valid for all instruments of the Enterprise; it is not limited on the currently used instrument or computer. That means, that the run time will be checked for all instruments and therefore in all methods that will be stored on the mentioned Enterprise machine.

System wake-up and all lamps on

If the option “System wake-up and all lamps on” is enabled and a date and time is selected, all device, switched to Standby and the lamps of all connected detectors will be switched on at the selected time. The method must be opened and the detectors switched on or in Standby at this time. The wake-up date time you can set by direct editing the time in the date-and-time line. If you click on the arrow key , you can select a date from a monthly overview (refer to the picture below).

Please note, that Standby only works for all Azura and the newer Smartline devices pump S 1050 and detectors S 2550 und S 2520. Older Smartline devices as pump S 1000, detector S 2600 cannot wake-up from Standby by a software command; also the lamp of the detector cannot be switched on, if the detector is in Standby.




Fig. 252 Advanced Runtime Settings – System wake-up... - date setting

If enabled, this option is valid for the currently used computer.

Trace name

If the option “Trace name” is enabled, the wave length of the initial line of the wave length table becomes part of the trace name in the chromatogram window as well as in the report (chromatogram, channel name in run report) instead of the channel name. This belongs to all detectors with a wave length table, as UV, DAD and fluorescence detectors. The channel name in the channel selector will not be changed.




Fig. 253 Advanced Runtime Settings – Trace name changes

The option will only change the trace name for new chromatograms.



If enabled, this option is valid for all instruments of the Enterprise; it is not limited on the currently used instrument or computer. That means, that the trace name will include the wave length for all instruments and therefore for all methods that will be stored on the mentioned Enterprise machine.

Download Tab / Method

The actually set parameters (may be even not saved) of the current instrument setup tab as well as those of the whole method can be downloaded to the corresponding instruments. Select the menu

CONTROL – DOWNLOAD METHOD or

CONTROL – DOWNLOAD TAB




Fig. 254 Knauer Instrument Control – Control – Download Method

If the device(s) is(are) not properly connected, the Instrument Activity Log will open and show the appropriate error message(s).

Error: Switching Valves : Set position of the valve failed
 Error: JetStream Oven: Error: Error setting temperature
 S 1000: Cannot download device settings: The instrument with the given serial number is not found on the communication port chosen.

Fig. 255

The action starts operation of devices but not the data acquisition.



In case of Download Tab, only a single pump will be controlled even for an HPG system.

Solvent Control

The Solvent Control option is accessible only with the pumps S 1000, K-1001, and K-1800. It is accessible via the menu sequence

CONTROL – SOLVENT CONTROL.




Fig. 256 Solvent Control window (quaternary LPG)

To reach the dialog box for editing the settings, either click on the **Change** button or select the menu sequence

METHOD – CONFIGURE SOLVENT...




Fig. 257 Solvent Control dialog box (quaternary LPG)

If you are working with a low pressure gradient system, you can enter the values for all components in the same window. You can activate the option for each of the components separately. As many solvents are configured as many components will be accessible. Consequently, an isocratic system will give access only to component A.

Check the **On** option to get a continuously monitored solvent amount for each pump separately. With this option activated, you have access to entering the values for the bottle size, actual volume, the acoustic warning level, and the shutdown level. **Accept** the entered values and close the window. If you close the window without having accepted any changes, those will not become valid.

The actual volume of solvent in each bottle is displayed. The eluent volume of the symbolized bottles displayed depends on the actual volume and bottle size.




Fig. 258 Solvent Control window – binary HPG system

If you are working with a high pressure gradient system, again you can enter the values for all components in the same window. However, only that component will be accessible, to which the pump was selected.




Fig. 259 Solvent Control dialog box – HPG system



The amount of solvents used will be calculated correctly only under software control. Should the pumps have been controlled manually, or the bottles have been refilled, the actual solvent amount must be corrected.

Qualification Procedures

ChromGate® provides a performance qualification tool for checking your HPLC systems. This is not a qualification procedure for each individual device (pump, detector etc.) but for the whole configured chromatographic instrument. The procedure is at last a normal method used to determine the measurement reproducibility.

The PQ procedure is started from the main instrument screen by selecting the menu sequence

CONTROL – QUALIFICATION PROCEDURES




Fig. 260 PQ start screen

You have the choice to create a new PQ session method or to select earlier created ones. These may be completed or not. If you select a not completed one you can continue with it. Selecting a completed one the final Report will be displayed.

The delete button becomes active if one of the earlier sessions is selected.





Fig. 261 PQ configure screen

Simultaneously the instrument setup screen and the PQ configure screen will be opened. If more than one detector (channel) and/or pump are configured you have to select the active ones via the drop down menus. Next you can select either any customized method or the Knauer Test Kit.



If you have selected the Knauer Kit, no further entries are to do in this screen. You have no access for editing.

For custom tests you have to enter the parameters according to your method.

Clicking on Continue or selecting Settings on the left hand side will open the next dialog window.




Fig. 262 PQ setting screen

Here the injection configuration, the peak detection configuration and the result processing are to define. The peak detection configuration is predefined and not access is possible as far you have selected the Knauer Kit. The result processing parameters define the thresholds for passing the test.

Clicking on the continue button first a consistence test is performed. If this test failed you will get the hints what changes are necessary.




Fig. 263 PQ consistence check, failed




Fig. 264 PQ consistence check, passed

Clicking on the **Start>** button starts the automatic performance of the whole PQ procedure. It will run like a usual sequence. The direct control should not be used while running the PQ sequence.

While the data acquisition is running you can add traces of other runs as usual but with one exception. You cannot add traces of earlier runs of the same PQ session because they are at this time only temporarily stored and not accessible. The final storage takes place at the end of the session. This is indicated by monitoring the performance qualification report.




Fig. 265 Performance qualification report

The chromatograms of the single session runs can be inspected by selecting them on the left hand side. They are now permanently stored in
Installation folder \
Enterprise_OEM\Performance Qualification\Generic PQ Procedures\ ...

Validation of Integration

This tool is a hardware independent test of the installed software package that indicates how reliably the integration of the chromatograms will be performed. To carry out this test, select any method setup window, to gain access to the menu sequence.

METHOD – OPTIONS - VALIDATION OF INTEGRATION...




Fig. 266 Validation of Integration dialog window

The test can be carried out using either a Gaussian curve (as shown in Fig. 266) or a pre-selected real chromatogram. If you select the latter option, the usual chromatogram window is opened, enabling the software to integrate this chromatogram.




Fig. 267 Validation of Integration dialog window, real chromatogram

To perform the validation procedure, just press the **Analyze** button. All integration parameters are preset and cannot be altered. If you had previously changed the integration events, the validation will not be performed and the following information box will be displayed.




Fig. 268 Validation of Integration attempt with manual integration fixes

The only integration parameter which can be modified is the baseline drift when using the Gaussian sequence. Check this option to consider a baseline drift of 10%. Press the **Analyze** button to repeat the validation procedure.




Fig. 269 Validation of Integration with baseline drift

If you press the **View report** button first, you will get the following message:




Fig. 270 Error message

After a few moments calculation time, the Gaussian sequence will also appear in the background window. In both cases, the performed integration will now be visualized by the baseline.

Press the **View report** button, to open the automatically created validation report text file.

ChromGate Chromatography Data System 3.3.2.980

Validation report of the Integration Algorithms

Reference file: alkylbenzoate_reference.dat

Calculation Mode: USP

	Pk#	Ret. Time (min)	Height (µAU)	Area (µAU'sec)	USP Width (min)	Theor. Plates	Cap. Factor	Resolution	Asymmetry (10%)
Calculated	10	4,52	447078,00	2823826,00	0,1665	11787,0400	3,0357	2,4088	1,1216
Expected	10	4,52	447078,00	2823826,00	0,1665	11787,0400	3,0357	2,4088	1,1216
Difference %	10	0,00	0,00	0,00	0,02	0,00	0,00	0,00	0,00
Calculated	11	5,44	388441,00	2868281,00	0,1957	12381,2637	3,8601	5,0983	1,1037
Expected	11	5,44	388441,00	2868281,00	0,1957	12381,2637	3,8601	5,0983	1,1037
Difference %	11	0,06	0,00	0,00	0,01	0,00	0,00	0,00	0,00
Calculated	12	6,93	296108,00	2600684,00	0,2318	14299,5430	5,1875	6,9554	1,0469
Expected	12	6,93	296108,00	2600684,00	0,2318	14299,5430	5,1875	6,9554	1,0469
Difference %	12	0,00	0,00	0,00	0,00	0,00	0,00	0,00	0,00
Calculated	13	9,25	218154,00	2468432,00	0,3020	15020,6299	7,2619	8,7046	1,0118
Expected	13	9,25	218154,00	2468432,00	0,302	15020,6299	7,2619	8,7046	1,0118
Difference %	13	0,04	0,00	0,00	0,00	0,00	0,00	0,00	0,00
Calculated	17	12,74	149308,00	2268936,00	0,4052	15821,8555	10,3780	1,7915	0,9896
Expected	17	12,74	149308,00	2268936,00	0,4052	15821,8555	10,3780	1,7915	0,9896
Difference %	17	0,03	0,00	0,00	0,01	0,00	0,00	0,00	0,00

Date : Wed Dec 07 18:17:29 2011

Fig. 271 Validation report, reference chromatogram

ChromGate Chromatography Data System 3.3.2.980

Validation report of the Integration Algorithms

Reference file: Gaussian trace, baseline drift = 0.

Calculation Mode: USP

	Pk#	Ret. Time (min)	Height (µAU)	Area (µAU·sec)	USP Width (min)	Theor. Plates	Cap. Factor	Resolution	Asymmetry (10%)
Calculated	1	1,00	999996,00	15039481,00	0,4058	97,1724	-0,1071	0,0000	1,0000
Expected	1	1,00	1000000,00	15039769,65	0,4000	100,0000	-0,1071	0,0000	1,0000
Difference %	1	0,00	0,00	0,00	1,44	2,83	0,00	0,00	0,00
Calculated	2	2,00	499997,00	7519692,00	0,4058	388,6902	0,7857	2,4644	1,0000
Expected	2	2,00	500000,00	7519884,82	0,4000	400,0000	0,7857	2,5000	1,0000
Difference %	2	0,00	0,00	0,00	1,44	2,83	0,00	1,42	0,00
Calculated	3	3,00	249998,00	3759825,00	0,4058	874,5573	1,6786	2,4644	1,0000
Expected	3	3,00	250000,00	3759942,41	0,4000	900,0000	1,6786	2,5000	1,0000
Difference %	3	0,00	0,00	0,00	1,44	2,83	0,00	1,42	0,00
Calculated	4	4,00	125000,00	1879914,00	0,4058	1554,7704	2,5714	2,4644	1,0000
Expected	4	4,00	125000,00	1879971,21	0,4000	1600,0000	2,5714	2,5000	1,0000
Difference %	4	0,00	0,00	0,00	1,44	2,83	0,00	1,42	0,00
Calculated	5	5,00	62500,00	939956,00	0,4058	2429,3408	3,4643	2,4644	1,0000
Expected	5	5,00	62500,00	939985,60	0,4000	2500,0000	3,4643	2,5000	1,0000
Difference %	5	0,00	0,00	0,00	1,44	2,83	0,00	1,42	0,00
Calculated	6	6,00	31250,00	469964,00	0,4058	3498,2683	4,3571	2,4644	1,0000
Expected	6	6,00	31250,00	469992,80	0,4000	3600,0000	4,3571	2,5000	1,0000
Difference %	6	0,00	0,00	0,01	1,44	2,83	0,00	1,42	0,00
Calculated	7	7,00	15625,00	234969,00	0,4058	4761,6870	5,2500	2,4644	1,0000
Expected	7	7,00	15625,00	234996,40	0,4000	4900,0000	5,2500	2,5000	1,0000
Difference %	7	0,00	0,00	0,01	1,44	2,82	0,00	1,42	0,00
Calculated	8	8,00	7808,00	117288,00	0,4056	6224,9165	6,1429	2,4650	1,0009
Expected	8	8,00	7812,00	117490,68	0,4000	6400,0000	6,1429	2,5000	1,0000
Difference %	8	0,00	0,05	0,17	1,40	2,74	0,00	1,40	0,09
Calculated	9	9,00	3902,00	58550,00	0,4054	7886,8540	7,0357	2,4662	1,0019
Expected	9	9,00	3906,00	58745,34	0,4000	8100,0000	7,0357	2,5000	1,0000
Difference %	9	0,00	0,10	0,33	1,34	2,63	0,00	1,35	0,19
Calculated	10	10,00	1948,00	29140,00	0,4042	9791,8213	7,9286	2,4704	1,0053
Expected	10	10,00	1953,00	29372,67	0,4000	10000,0000	7,9286	2,5000	1,0000
Difference %	10	0,00	0,26	0,79	1,06	2,08	0,00	1,19	0,53

Date : Thu Dec 08 09:22:11 2011

Inspector: matthias Signature:

Fig. 272 Validation report, Gaussian trace

Generic Drivers

Generally spoken, any device that supports RS-232 ASCII communication can be controlled with ChromGate[®] software. No changes in source code of ChromGate[®] are required but for each device an additional driver is required.

The configuration of a generic device is to be done by Knauer. For each of those devices the new generated driver has to be installed.

If the driver is installed, its icon will be present in the list of instruments for configuration.




Fig. 273 Instrument selection window, example Jasco FP-920

For each item on the left hand side the corresponding parameters are shown and described right.




Fig. 274 Instrument configuration window, example Jasco FP-920

These configuration windows are strongly depended on the given type of instrument. Some settings are fixed. They will appear grayed with no access.

The setup programming is done using a script language, as demonstrated below.

The upper part of the window is for the programming inputs, whereas the lower one gives the instrument specific programming information e.g. the allowed ranges for different parameters. Before any program is saved it is checked for its correctness. If this check fails you will get a corresponding error message.




Fig. 275 Instrument setup programming window, example Jasco FP-920




Fig. 276 Instrument setup programming window, example Jasco FP-920

After completing the generic driver program an instrument tab as for any other instrument will be displayed in the instrument status tab.




Fig. 277 Instrument status tab, example Jasco FP-920

Depending on the device, a few most important parameters will be monitored and displayed in the status area.

Direct control is not supported for Generic Drivers.



Generic drivers only support basic functionality. Please ask Knauer for available new Generic drivers or for the requirements for the development of a new Generic driver at software@knauer.net.

ChromGate® System Suitability Setup

Before ChromGate® can make System Suitability calculations, you must enter the required acceptance ranges for the peaks of interest. This is done in using the **Method – System Suitability** command.

1. Click on the **Method – System Suitability** command to view the System Suitability Setup dialog.




Fig. 278

2. Select the first peak to be used for calculations by highlighting it with the mouse in the **Compounds** list.
3. In the adjacent spreadsheet, click on the first field in the **Parameter** column. A drop-down list of available parameters is displayed. Select a parameter from the list. If you do not wish to perform a suitability test on any given peak, simply leave the **Parameter** fields blank. Similarly, if you do not wish to have one of the test criteria used (for example, %RSD), leave it blank. Some parameters have a choice of calculation methods. These parameters have the calculation method displayed after the parameters in parenthesis (e.g. Plates/Meter (JP) indicates Plates/Meter calculated using the Japanese Pharmacopeia calculation). For details on these calculations, see the equations section.
4. For each parameter selected, enter a minimum value (**Min**), maximum value (**Max**), and maximum allowed percent relative standard deviation (**%RSD**).
5. If you want a noise test to be performed, click on the lower spreadsheet in the column labeled **Test**. Select one of the noise calculations **Noise (rms)/ASTM Noise (rms)**, or **6-sigma Noise (rms)** from the drop-down list. If you select this box, you must enter a **Start** and **Stop** time for the test, and the **Threshold** value to determine acceptable limit. The RMS noise value for the portion of chromatogram between the Start and Stop times will be calculated and compared to the Threshold value to determine whether the test Passed or Failed. Note that the times you enter for the noise test should be representative of a baseline area of your chromatogram where not peaks elute. Keep in mind, that the noise may depend on the detector type, wave length, flow rate, solvent and other parameters used in the test.

6. If you want a drift test to be performed, click on the lower spreadsheet in the column labeled **Test**. Select **Drift (uV/min)** from the drop-down list. If you select this box, you must enter a **Start** and **Stop** time for the test, and a **Threshold** value (in $\mu\text{V}/\text{Min}$) for the acceptable limit.
7. When you have completed the System Suitability Setup, close the box.

Copy & Paste

After highlighting sections of the spreadsheet that you wish to copy, push:

Ctrl C- to copy the highlighted section of the spreadsheet to the clipboard.

Shift + Insert- to paste the clipboard item to where the cell is highlighted.

Shift + Arrows- highlights the spreadsheet cells, and can move them.

Suitability Calculation Selection

System Suitability calculations require information about your chromatography column. These parameters are entered in the **Method/Advanced/Performance tab** dialog box.

1. Click on Method/Advanced/Performance tab.

Enter values for the column parameters shown, required for calculation of performance options.

The screenshot shows a dialog box titled "Advanced Method Options -- Detector 1". It has several tabs: "Export", "Custom Parameters", "Column / Performance" (which is selected), "Files", and "Advanced Reports".

Under the "Column / Performance" tab, there is a section for "Column Information" with the following fields:

- Unretained peak time: 0.85 Minutes
- Column length: 0.25 meters (radio button selected), cm (radio button unselected)
- Particle diameter: 0 microns
- Column serial number: [empty text box]
- Column installation date: [empty text box]
- Column description: [empty text box]

Below this section is a checked checkbox: "Calculate performance parameters for this channel".

Underneath this checkbox is a list of calculation methods:

- USP
- EMG
- DAB, BP, EP, ASTM
- AOH
- JP

Fig. 279

2. Select a method for calculation of the performance values by clicking on the **Calculation Method** down-arrow.
3. Choices include **USP**, **EMG** (Exponential Modified Gaussian), **DAB**, **JP**, and **AOH (Area/Height)** method. The calculations used are given in the equations section. For BP, EP, and ASTM methods, choose DAB.
4. Click on the **Calculate Performance Parameters** box to enable the suitability calculations for the method

Running a Suitability Test

Once you have completed the **System Suitability Setup** and **Performance** dialogs, create a Sequence or modify an existing Sequence to include one or more system suitability runs according to the needs and SOP's of your laboratory. If you want the System Suitability Report to print automatically at the end of the sequence, you must also select the **Print Sequence Reports** option from the **Sequence Properties** dialog box.

For example, following USP standards, five replicate standards are run at the beginning of the sequence. These are designated as system suitability standards in the sequence. At the end of the sequence, suitability calculations are made, and EZChrom Elite Client/Server generates a system suitability report.

Run #	Status	Run Type	Level	Sample ID	Method	File
1		Sys Suit Begin	0	Suitability1	Multilevel Calibration.met	calibratio
2		Sys Suit Std.	0	Suitability2	Multilevel Calibration.met	multi calibra
3		Sys Suit Std.	0	Suitability3	Multilevel Calibration.met	multi calibra
4		Sys Suit Std.	0	Suitability4	Multilevel Calibration.met	multi calibra
5		Sys Suit Std.	0	Suitability5	Multilevel Calibration.met	multi calibra
6		Sys Suit End	0	Suitability5	Multilevel Calibration.met	multi calibra

Fig. 280

Run Type is selected by clicking on the **Run Type** field, then selecting the appropriate sample type from the choices provided.




Fig. 281

The first sample in your suitability set should be designated as **Begin System Suitability**. Additional suitability standards should be designated as **System Suitability Standards**, and the final sample in your suitability set should be designated as **End System Suitability**. Multiple sample run types can be selected for a given sample.

A System Suitability Report will be generated at the end of the set of suitability standards when the sequence is run.

To view the Suitability Report on screen, click **Reports/View/Sequence Custom Report** choice. The System Suitability Report will appear in a list of available Sequence Custom Reports.

To print the Suitability Report, click on **Reports/Print/Sequence Custom Report**, and select the **System Suitability Report** choice. Note that in order to have your report printed automatically at the end of the sequence, you must have selected the **Print Sequence Reports** option in the **Sequence Properties** dialog.

Suitability Reports

EZChrom Elite Client/Server provides a default suitability report template. This report template can be used as-is, or it can be customized in the Sequence Custom Report editor.

System Suitability Report

Page 1 of 1

Sequence : C:\EZC\from Elite\SEQUENCE\Suitability.seq
 User : System
 Printed : 3/6/97 9:43:24 AM

SYSTEM IS NOT SUITABLE

Channel A	Compound	Parameter	Min	Max	% RSD		
	Acetaminophen	estd	8.000	11.000			
		tperm	46000.00	50000.00	5.000		
Sample ID	Compound	Parameter	Average	Low	High	% RSD	Status
	Acetaminophen	estd	18.009	10.000	50.046	99.444	
Suitability		10.000					Passed
Suitability		10.000					Passed
Suitability		10.000					Passed
Suitability		10.000					Passed
Suitability		50.046					FAILED
		tperm	49016.13	46025.98	49763.66	3.410	
Suitability		49763.66					Passed
Suitability		49763.66					Passed
Suitability		49763.66					Passed
Suitability		49763.66					Passed
Suitability		46025.98					Passed
Channel A	Sample ID	Test	Start(Min)	Stop(Min)	Threshold	Result	Status
	Suitability	Drift	0.000	0.012	1000.000	600.000	Passed
	Suitability	Drift	0.000	0.012	1000.000	600.000	Passed
	Suitability	Drift	0.000	0.012	1000.000	600.000	Passed
	Suitability	Drift	0.000	0.012	1000.000	600.000	Passed
	Suitability	Drift	0.000	0.012	1000.000	428.571	Passed
	Suitability	Noise	0.000	0.012	5.000	0.000	Passed
	Suitability	Noise	0.000	0.012	5.000	0.000	Passed
	Suitability	Noise	0.000	0.012	5.000	0.000	Passed
	Suitability	Noise	0.000	0.012	5.000	0.000	Passed
	Suitability	Noise	0.000	0.012	5.000	0.000	Passed

Fig. 282

ChromGate® PDA Option

The PDA Option enables you to view and analyze spectra generated with a diode array detector or a K-2600 on your LC. In order to use the PDA Option, you must have control software installed for your LC that contains PDA detector capabilities.

PDA Method Setup

In addition to instrument control parameters that control the operation of your PDA detector, there are setup options for the PDA available in the method. These include setting up libraries, spectral filtering, definition of multiple wavelengths for multi-wavelength display, and setup of two wavelengths to be used for the Ratio view. These options are available from the main menu by selecting Method/PDA Options.

PDA Options Library

Using this tab, you can enter library parameters for spectral library searches that will be saved as part of your method. When you do a spectral library search "using method parameters", these parameters will be used.




Fig. 283

Library/Enabled

Enter the spectral library to be searched or select from available libraries by clicking the Library field followed by clicking the file button. You can select more than one library. To enable the library for searching, click the **Enabled** box. If this list is empty, or no library is enabled, no hits will be returned when a search is performed.

Search parameters

Enter the parameters to use for the search.

Wavelength range

Enter the wavelength range to search

Wavelength step

Enter a step number for the search. Larger numbers will make the search faster, but if you use too large of a step, spectral details may not be picked up.

Max hits

Specify the number of hits that will be reported in the results of a library search. Note that this works in conjunction with the *Similarity Threshold* parameter to limit the number of hits reported.

Similarity threshold

Enter a number for threshold of similarity. The library search results will only display matches whose similarity to the unknown exceeds this value.

Pre-filters

The options in this group allow you to specify search pre-filters that will be performed on library spectra prior to the test for similarity. All pre-filters are optional.

Retention time range

When values are entered in these fields, PDA limits its search to library spectra obtained from peaks whose apex is within the specified retention time range. This pre-filter is optional and may be left blank.

Lambda max

When one or more of these values is specified, library search will be restricted to those library entries containing lambda max value(s) within 5 nm of all of the specified values. Entries without matching lambda max

value(s) are automatically excluded from the search (no similarity calculation is made). Entering values for lambda max is optional.

Compound name filter

If you are searching only for spectra that contain a certain name, enter the name here. Only spectra containing that name will be searched.

PDA Options Purity

The Purity tab is used to set the parameters necessary to perform on-demand peak purity calculations and peak purity calculations that occur as part of analysis.




Fig. 284

Purity Calculations

Wavelength range

Specify the wavelength range over which the purity calculations will be performed, for example, from 200 nm to 400 nm.

Wavelength step

Specify the wavelength spacing (in nm) to be used when purity calculations are performed.

Purity threshold

The Threshold setting is used to modify the impact of the noise spectrum on total peak purity calculations. It has no effect on the calculation of three-point purity. Higher values for the Threshold setting increase the weight of spectral noise on the purity calculation. This has the effect of raising the computed similarity values for spectra and thus increasing the total peak purity value.

Absorbance threshold

This value represents the percentage of peak height that spectra will be included in purity calculation. Spectra in sections of the peak that do not exceed this threshold will not be included in the purity calculation. This provides a method of eliminating spectra where the concentration of the compound is so low that the solvent spectrum interferes.

Background compensation

Checking this box will cause spectra to be corrected for background using the peak baseline prior to being used in the calculation of purity.

Per-peak spectrum calculations

Checking any of these boxes indicates that the indicated value will be calculated on a per-peak basis during analysis. The result of this

calculation will then be available in reports and as chromatogram annotations.

Disabling values that are not of interest will speed up analysis. If a box is unchecked and the field appears in a run report, it will be reported as zero

Total Purity

Checking this box will cause total peak purity to be calculated on a per-peak basis during analysis. The result of this calculation will then be available in reports and as chromatogram annotations. If this box is unchecked and a report contains the field, a zero will be reported.

3 Point Purity

Checking this box will cause 3-point peak purity to be calculated on a per-peak basis during analysis. The result of this calculation will then be available in reports and as chromatogram annotations. If this box is unchecked and a report contains the field, a zero will be reported.

PDA Options Spectrum

The Spectrum tab is used to specify the types of filtering and processing to be performed on spectra that are extracted from the 3D data during analysis as well as to spectra displayed in the PDA Spectrum View.



The processing specified on this page is performed prior to any use of the spectra in the software, including display, searching and reporting.




Fig. 285

Spectral filtering

In this area you designate how spectral filtering (if any) will be performed.

Filtering type

Choose the type of filtering for the spectral plot. The choices are None, Smooth, 1st Derivative, and 2nd Derivative. Selecting one of these smoothing algorithms can remove noise from the spectrum.

Background correction

If this box is checked, then, prior to display, a correction for spectral background is made, as follows:

1. The spectra from the baseline start and baseline stop times for the peak are extracted from the data.

2. For each spectrum in the peak, a corresponding background spectrum is generated by linear interpolation between the baseline start and baseline stop spectra.
3. This background spectrum is subtracted from the original spectrum.



The baseline start and stop times that are used in calculating background compensation are based on the detected peaks for the channel currently selected.

Interpolate spectra

Select this check box to automatically perform 10:1 interpolation of spectra displayed in the spectrum window using a cubic spline curve fit. This interpolation is performed after the applying any spectral filtering option (1st derivative, 2nd derivative or smooth) to the display spectrum.



Interpolated spectra may not be stored or added to a library.

The peak spectrum is defined as

This specifies if and how spectra extracted from the 3D data are averaged before being exported or used in the Spectrum Report.

Spectrum Averaging is the calculation of a spectrum whose absorbance values are the average of the absorbencies.

The peak apex spectrum

When this is selected, the spectrum at the apex of the peak will be used, without averaging.

The avg. of upslope, apex and down slope spectra

When this is selected, three spectra (the spectra at the upslope inflection point, the apex and the down slope inflection point) are averaged.

The avg. of every n spectra across the peak

When this is selected, enter a value for **n**. Every nth spectrum from the peak is extracted, beginning at the start of the peak and continuing to the end of the peak. These extracted spectra are averaged and become the spectrum that is reported or exported.

Per-peak spectrum calculations

Checking any of these boxes indicates that the indicated value will be calculated on a per-peak basis during analysis. The result of this calculation will then be available in reports and as chromatogram annotations.

Unchecking values that are not of interest will speed up analysis. If a box is unchecked and the field appears in a run report, it will be reported as zero

Similarity

Checking this box causes the Peak Apex Similarity to the reference spectrum to be calculated on a per-peak basis during analysis. The result of this calculation will then be available in reports and as chromatogram annotations. If this box is unchecked and a report contains the field, a zero will be reported.

Peak Apex Similarity applies only to named peaks and its use requires that a reference spectrum be included in the peak table.

Upslope Similarity

Checking this box causes upslope similarity to be calculated on a per-peak basis during analysis. Upslope similarity compares the peak apex spectrum to the spectrum at the peak inflection point on the upslope side

(to the left of the apex). The result of this calculation will then be available in reports and as chromatogram annotations. If this box is unchecked and a report contains the field, a zero will be reported.

Down slope Similarity

Checking this box causes down slope similarity to be calculated on a per-peak basis during analysis. Down slope similarity compares the peak apex spectrum to the spectrum at the peak inflection point on the down slope side (to the right of the apex). The result of this calculation will then be available in reports and as chromatogram annotations. If this box is unchecked and a report contains the field, a zero will be reported.

Lambda max

Checking this box indicates causes the lambda max (wavelength at which the highest absorbance occurs) to be computed for each peak of each multi-chromatogram extracted from the 3D data during analysis. The result of this calculation will then be available in reports and as chromatogram annotations. If this box is unchecked and a report contains the field, a zero will be reported.

Lambda max Calculation range

This allows the user to restrict the wavelength range over which the Spectrum Max Plot will be calculated. Default is the detector acquisition range.

PDA Options Multi-Chromatogram

The Multi-Chromatogram tab is used to specify data channels within the 3D data for integration and quantization. It allows you to select wavelengths to display on the Multi-Chromatogram plot. When you select an Enabled box, enter a wavelength and bandwidth to be used for one plot. After you have set up several wavelengths, you can disable them temporarily by clicking the Enabled box such that the check mark is not displayed. The Multi-Chromatogram display is selected from the View Gallery. The Multi-Chromatogram view is not functional unless one or more valid wavelengths are enabled in the Multi-Chromatogram Definition tab.




Fig. 286

Enabled

To enable a wavelength, click the checkbox. To disable a wavelength, click the checkbox again to remove the checkmark. Only enabled wavelengths will appear in the multi-chromatogram view.

Wavelength

Enter the wavelength you want to view.

Bandwidth

Specify the nm range to be averaged in generating the analogue signal for each channel.

The multi-chromatogram data is an average of the absorbencies monitored at each wavelength in the wavelength range. The wavelength range is equal to the selected wavelength +/- one half of the bandwidth.

For example, given a wavelength of 600 nm and a bandwidth of 4 nm, for each sample time, the data points from the chromatograms at 598, 599, 600, 601 and 602 nm are added, and the sum divided by 5.

If, in the example above, the detector's range is 190 - 600 nm (inclusive), the software only averages the points across chromatograms at 598, 599 and 600 nm (dividing by 3). The range would not be exceeded.

PDA Options Ratio

The X:Y Ratio tab is used to set the displays of the Ratio Plot View and the Max Plot. The Ratio View displays two PDA wavelength channels and the ratio of those two channels. These may be viewed during real-time acquisition as well as during post-run analysis. The flat tops on the ratio peaks are a preliminary indication of peak purity.




Fig. 287

Ratio plot

These controls specify the wavelength of the ratio multi-chromatogram channels. The extracted chromatograms will be centered about the specified wavelength.

Channel 1, Channel 2

Wavelength

Enter the wavelength of the chromatogram channel.

Bandwidth

Enter the nm range to be averaged in generating the chromatogram channel.

Configuration of user defined detectors

Enter the ratio threshold. The threshold is the minimum absorbance value required in the chromatograms of both channels for calculation of the ratio. It is a method of eliminating ratio calculation at minor absorbance

values, such as those that occur with noise. If the threshold value is not met, the ratio plot will be zero.

PDA Views

The PDA views provide a variety of ways to view PDA data. You can choose to display one view at a time, or a combination of views, such as 3D and Contour or Mixed with 3D. The view is selected using the **View/PDA View** menu, or from the PDA toolbar by selecting **Views**. A right mouse click within any of the view windows provides a unique menu of options for that view.

3D View

The 3D View provides a three-dimensional view of absorbance versus time and wavelength. Wavelengths of appreciable absorbance and interference, which may be invisible in a single wavelength plot, are easy to locate with the 3D View. The plot can be elevated and rotated around its axis for display from any angle.




Fig. 288



During a run, the user must manually refresh this display. It does not update automatically.

Right-click inside the 3D view to display the pop-up menu that gives you access to plot rotation and axis setup features.



Fig. 289

3D Plot movement selections

Click on the type of plot movement you wish to use (elevate, roll, rotate, rotate XYZ, or spin). A checkmark will appear next to the selected movement option. Once you have selected one of the options, you will be returned to the plot and the cursor will indicate the type of movement selected.

Hold down the left mouse button and move the cursor in the direction you wish to move the plot. The plot will move as you move the cursor. The movement option will remain in effect until you turn it off. When finished, click the right mouse button, and de-select the movement option. If you wish to return to the original view, select the **Reset** command.



These values can be viewed or changed using the 3DPlot/Properties dialog.

3D Properties

This dialog allows you to set the properties of the 3D plot.




Fig. 290

Style

Select **Color** or **Grayscale** for the plot. Checking the **Wire Mesh** box will cause the data to be rendered as a wire frame plot, rather than as a solid fill plot.

Colors

Range

Select how you want the coloration on the plot to appear. When **Light & Dark Range** is selected, alternating light and dark bands are used. When **Full Spectrum** is selected, then a continuous color spectrum is used.

Background

Select the color for the background of the plot.

Axes

Select the color for the plot axes.

Display

This specifies the relative quality of the displayed contour plot. Lower Display Detail will result in faster drawing of the plot.

Print

This specifies the relative quality of the printed contour plot. Lower Print Quality will result in faster printing.

Viewing

Use these fields to view current rotation settings, or to set them manually.

Light source

Selecting this option causes the plot to be shaded as if an external light source were shining upon it.

Rotation

This reports the current level of rotation (front to back) in the aspect position of the plot. When a new value is entered, the plot is redrawn to reflect the new value upon exiting the dialog.

Elevation

This reports the current level of tilt (forward or backward) in the aspect position of the plot. When a new value is entered, the plot is redrawn to reflect the new value upon exiting the dialog.

Roll

This reports the current level of roll (side to side) in the aspect position of the plot. When a new value is entered, the plot is redrawn to reflect the new value upon exiting the dialog.

Zoom

This reports the current level of magnification in the plot. When a new value is entered, the plot is redrawn to reflect the new value upon exiting the dialog.

Performance - Use Zoom/Rotate Bounding Box

When this box is checked, the plot will be temporarily replaced by a box during zoom and rotation operation. When the operation is completed, the plot will be redrawn. Checking this box will increase performance on computers with slower graphic subsystems.

3D Properties Axis

This dialog lets you set up axis limits for your 3D plot.




Fig. 291

Limits

Enter the limits for the 3D plot.

Time

Click **Auto scale** if you want to have the software automatically scale the time axis to the maximum values. To enter a manual range, de-select the **auto scale** box, and add your own limits, or click the **Get Limits** button to enter the limits displayed on the current 3D graph.

Wavelength

Click **Auto scale** if you want to have the software automatically scale the wavelength axis to the maximum values. To enter a manual range, de-select the **auto scale** box, and add your own limits, or click the **Get Limits** button to enter the limits displayed on the current 3D graph.

Absorbance

Click **Auto scale** if you want to have the software automatically scale the absorbance axis to the maximum values. To enter a manual range, de-select the **auto scale** box, and add your own limits, or click the **Get Limits** button to enter the limits displayed on the current 3D graph.

Labels

You can customize the labeling of the 3D plot using the parameters in this area. You can select font, size, color and style using the selections provided.

3D Plot Rotation Options

A right mouse click on the 3d plot provides a menu containing various rotation options. When you select one of these rotation options, a checkmark will appear next to the option and that rotation option will become active on the plot.



Fig. 292

The cursor will change to indicate the option is in effect.

While the option is active, you can move the plot in the designated way by holding down the left mouse button and moving the mouse in the desired direction on the plot. When you have the plot in the desired position, click the right mouse button and turn off the option by selecting it again (checkmark removed). If you want to return to the original plot view, click the right mouse button and select **Reset**.

Contour View

The Contour Plot (also referred to as an Isoabsorbance Plot) provides an aerial view of the absorbance of the sample at each wavelength versus time. The contour view supplies quick and easy-to-assimilate information about those wavelengths at which the sample exhibits appreciable absorbance. With contour view, it is also possible to generate a Chromatogram View for an individual wavelength and a Spectrum View for a given point in time.




Fig. 293

Right-click inside the window to display the pop-up menu. Select Properties to display the Contour Properties dialog box.

To generate a chromatogram view from contour view of Mixed View or Mixed View w/ 3D:

1. Select View /PDA Views /Mixed View to display the Contour Map, Chromatogram, and Spectrum.
2. Move the cursor to the triangle-shaped handle located on the left-hand wavelength axis of the Contour Map and press the left mouse button.
3. Drag the cursor up or down to the desired wavelength and release the mouse button.
4. The chromatogram associated with the specified wavelength is displayed in the Chromatogram View.

To generate a spectrum view from contour view of Mixed View or Mixed View w/ 3D:

1. Select View /PDA Views /Mixed View to display the Contour Map, Chromatogram, and Spectrum.
2. Move the cursor to the triangle-shaped handle located on the time axis of the Contour Map and press the left mouse button.
3. Drag the cursor to the desired peak and release the mouse button.
4. The Spectrum associated with the specified retention time value is displayed in the Spectrum View.

Contour Properties

A right mouse click anywhere on the contour plot, followed by selecting the Properties button will display a dialog where you can select the way the contour plot is displayed.




Fig. 294

Style

Select how you want the plot to appear, **Grayscale** or **Color**.

Colors

Range: Use this to select how the coloration of the plot is to be displayed. When **Light & Dark Range** is selected, alternating light and dark bands are used. When **Full Spectrum** is selected, then a continuous color spectrum is used.

Background: Select the color to be used for the background of the plot.

Display

This specifies the relative quality of the displayed contour plot. Less Display Detail will result in faster rendering of the plot.

Print Quality

This specifies the relative quality of the printed contour plot. Coarse Print Quality will result in faster printing.

Cursors

The values in these boxes reflect the current contour cursor positions. When you change the values, the cursors on the plot will be changed when you exit the dialog or when you click the Apply button.

Time

This specifies the time position (X Value) of the cursor. When a value is entered, the X cursor is updated to the new position upon exiting the dialog.

Wavelength

This specifies the wavelength position (Y Value) of the cursor. When a value is entered, the Y cursor is updated to the new position upon exiting the dialog.

Bandwidth

This control specifies the wavelength band that will be averaged when a chromatogram is extracted from the contour plot. The extracted chromatogram is an average of the absorbencies at each wavelength in the wavelength band. The wavelength band is equal to the selected wavelength (see above) +/- one half of the bandwidth.

Contour Axis Setup

You can change the setup of the axes for the contour plot using this dialog. If the Auto scale box is selected, the software will automatically set the axis parameter. If you want to enter the parameter manually, de-select the Auto scale box, then enter the range values, or click the Get Limits button to bring in the limits from the current view.




Fig. 295

Time

Click **Auto scale** if you want to have the software automatically scale the time axis to the maximum values. To enter a manual range, de-select the **auto scale** box, and add your own limits, or click the **Get Limits** button to enter the limits displayed on the current contour graph.

Wavelength

Click **Auto scale** if you want to have the software automatically scale the wavelength axis to the maximum values. To Contour Axis Setup

Absorbance

Click **Auto scale** if you want to have the software automatically scale the absorbance axis to the maximum values. To enter a manual range, de-select the **auto scale** box, and add your own limits, or click the **Get Limits** button to enter the limits displayed on the current contour graph.

Labels

You can customize the labeling of the contour plot using the parameters in this area. You can select font, size, color and style using the selections provided.

Mixed View

This view displays the contour, chromatogram and spectrum views along with the fourth pane displaying similarity, purity or peak profile.




Fig. 296

To zoom in on a portion of either the contour plot, chromatogram or spectrum, hold the left mouse button down, move the mouse over the plot until the area of interest is highlighted and release the mouse button.

To quickly move to the previous level of zoom, double-click on the plot.

To zoom out to the full plot after multiple zooming operations, use Ctrl+Z, or use Shift + double-click in the window, or click the right-mouse button anywhere in the window and select Full Unzoom from the pop-up menu.

Synchronize



You can synchronize the axis and zoom limits with the other two plots using the Synchronize button on the PDA toolbar. When the contour plot is zoomed, the time range of the chromatogram and the wavelength range of the spectrum plot are automatically adjusted to the same range as the contour. When synchronized zooming has been performed, the x axis plot limits for the chromatogram and spectrum plot are changed. To unzoom these plots to the limits of the data, you must select Full Unzoom from the contour plot context menu.

Mixed View Arrange Button

If you have changed the size of the 4 views, clicking on the Arrange button will arrange all the views in the same size again.

Mixed View Actions Button

This button presents the following actions:




Fig. 297

Search Library...

Opens the Spectral Library Search window.

Overlay Spectra...

Allows you to select a spectrum to overlay in the Spectrum view pane.

Select Peak for Similarity/Purity Display...

Enables you to select a peak for display in the Purity view pane.

Add to Multi-Chromatogram Table...

Enables you to select a chromatogram to add to the Multi-Chromatogram table in the PDA options.

Add Spectrum to Spectral Library

Selecting this menu item allows you to create a new library entry from the spectrum currently displayed in the spectrum pane. The spectrum will be added to the currently open library (opened by selecting File > Spectrum Library > Open). A name of the form "Spectrum at time xx.xx" will be used as the Compound Name for the new library entry.

Mixed View Options Button

This button presents options that enable you to toggle the view for the lower right pane of the mixed view:




Fig. 298

Show 3D Plot
 Show Peak Similarity and Threshold Plot
 Show Peak Purity Plot
 • Show Spectrum Similarity Table

Chromatogram View

The Chromatogram View may be used alone or as part of the Mixed View displays. Generate a chromatogram from the Contour plot as described in the section on contour plot.

Right-click within the chromatogram window to display the pop-up menu. This menu contains the same options as the basic chromatogram window for all detector types, and in addition enables you to overlay chromatograms from different wavelengths and change the Gallery view.

Overlay Chromatograms

When this option is available when you have a chromatogram view and a contour plot displayed simultaneously (Mixed View). Select **Overlay Chromatograms** to add chromatograms from different wavelengths to the view whenever you slide the wavelength selection cursor on the contour plot to a new wavelength.

Max Plot

A Max Plot is a chromatogram with each point plotted at its maximum absorbance. This plot gives an indication of the appearance of the chromatogram when the wavelengths are optimized for each peak.




Fig. 299

Peak Purity Plot

The Total Peak Purity view displays the purity profile for a chromatogram extracted from the 3D data. The view displays purity information for the Max Plot chromatogram that is displayed in the Chromatogram pane of the Mixed View. The pane will be blank until the data has been analyzed, and a peak has been selected for the purity calculation.




Fig. 300

Select Peak for Purity Similarity Display

To select the reference spectrum for the purity calculation, make sure the data has been analyzed. (If you are not sure, click the analyze button on the toolbar.) Then, click **Actions/Select Peak** from the Mixed View toolbar. A dialog will appear that instructs you to select a peak by holding down the Ctrl button and then clicking on a peak from the Max Plot pane. The retention time of the reference spectrum is displayed along with a value indicating peak purity.




Fig. 301

You can continue to select or change the peak. When you are finished, click the **Close** button.

Purity View Properties

To set up the properties for the Purity View, do a right mouse click, and select Properties.




Fig. 302

By default, the Similarity View will use the current method settings. To override the method settings, click the Override method settings box, then enter the values you wish to use.

Wavelength range

Specify the wavelength range over which the purity calculations will be performed, for example, from 200 nm to 400 nm.

Wavelength Step

Specify the wavelength spacing (in nm) to be used when purity calculations are performed.

Purity threshold

The comparison of two spectra gives a Similarity Index (SI), the closer to 1.00 number is the more similar the spectra are. The purity threshold is used to eliminate spectra that do not match. If the SI is greater than the purity threshold, the spectra will be considered pure. In general, a spectra with an SI greater than 0.9900 would have a high probability of being the same as the apex spectra. A SI greater than 0.9000 but less than 0.9900 shows some similarity but would need to be evaluated with care and caution. SI less than 0.9000 should be considered non-similar.

Absorbance threshold

This value represents the percentage of peak height that spectra will include in purity calculation. Spectra in sections of the peak that do not exceed this threshold will not be included in the purity calculation. This provides a method of eliminating spectra where the concentration of the compound is so low that the solvent spectrum interferes.

Background compensation

Checking this box will cause spectra to be corrected for background using the peak baseline prior to being used in the calculation of purity.

Similarity View Properties

To set up the properties for the Similarity View, do a right mouse click, and select **Properties**.




Fig. 303

By default, the Similarity View will use the current method settings. To override the method settings, click the Override method settings box, then enter the values you wish to use.

Wavelength range

Specify the wavelength range over which the purity calculations will be performed, for example, from 200 nm to 400 nm.

Wavelength Step

Specify the wavelength spacing (in nm) to be used when purity calculations are performed.

Purity threshold

The comparison of two spectra gives a Similarity Index (SI), the closer to 1.00 number is the more similar the spectra are. The purity threshold is used to eliminate spectra that do not match. If the SI is greater than the purity threshold, the spectra will be considered pure. In general, a spectra with an SI greater than 0.9900 would have a high probability of being the same as the apex spectra. A SI greater than 0.9000 but less than 0.9900 shows some similarity but would need to be evaluated with care and caution. SI less than 0.9000 should be considered non-similar.

Absorbance threshold

This value represents the percentage of peak height that spectra will include in purity calculation. Spectra in sections of the peak that do not exceed this threshold will not be included in the purity calculation. This provides a method of eliminating spectra where the concentration of the compound is so low that the solvent spectrum interferes.

Background compensation

Checking this box will cause spectra to be corrected for background using the peak baseline prior to being used in the calculation of purity.

Multi-Chromatogram View

This displays multiple chromatographic plots of absorbance versus time, each at a different wavelength. Maximize the PDA display and select **View / PDA Views / Multi-Chromatogram** (or use the PDA toolbar/Views) to display all of the chromatograms in their respective channels which were specified in the **Method / PDA Options / Multi-Chromatogram** tab.

Right-click within the window to display the pop-up context menu.




Fig. 304

Spectrum View

This view displays the spectrum associated with a time on the chromatogram. The time of the spectrum to be displayed can be changed from Mixed View by sliding the spectrum selection box (on the x-axis).




Fig. 305

Spectrum Properties

Select this menu item to access the trace properties dialog. This enables you to add another trace to the view, or change scaling. It also lets you selectively remove overlaid traces from the view.

Up to 3 wavelength maxima and -minima can be displayed. Use the Annotations... command from the right mouse click menu in the Spectrum view and select Lambda Max and/or Lambda Min.




Fig. 306

Spectrum Background Correction

A background correction can be performed if the current spectrum has been extracted from an integrated peak. Select the **right mouse click/Operations/Background Correction** command from the Spectrum window. Then click on the spectrum on which you wish to perform the calculation. A background correction for the displayed spectrum will be performed based on the working chromatogram. The background-corrected spectrum is then displayed along with the original spectrum. Please refer to the **PDA Analysis and Calculations** section for details on how this calculation is performed.

Background Correction may only be performed once on a spectrum. If a second background correction is attempted on an already corrected spectrum, a message box is displayed and the operation ignored.

If the current spectrum has not been extracted from an integrated peak of the working chromatogram, selecting this menu item will have no effect on the spectrum.

Background subtraction must always be the first operation performed on a spectrum. If a background correction is attempted on a spectrum after another operation has already been performed (including operation specified in the **Spectrum** tab of **PDA Options**), a message box is displayed and the background correction request ignored.

Interpolate Spectrum

Spectrum Interpolation can be performed if the current spectrum has been extracted from an integrated peak. Select the **right mouse click/Operations/Interpolate** command from the Spectrum window. Then click on the spectrum on which you wish to perform the calculation. A 10:1 interpolation will be performed on the spectrum. The new interpolated spectrum is then displayed along with the original spectrum. Please refer to the **PDA Analysis and Calculations** section for details on how this calculation is performed.

Refer to the **Data Analysis and Calculations** section for details on how the values are computed.

Export

Selecting this menu item will display a **File Save As** dialog. When a valid filename is entered and Ok is pressed, the currently displayed spectrum will be exported as an ASCII data file. Refer to the **Spectrum Export** section for details on the data file format.

Overlay Spectra

You can overlay spectra using the **Actions/Overlay Spectra...** command from the Mixed View toolbar. When you select this option, a dialog appears where you can enter the retention time from which you want to extract a spectrum. You can also overlay spectra quickly by moving the wavelength cursor in the Contour or Chromatogram view to the desired wavelength. To remove the spectra you have added to the view, use the **Clear overlays** command from the right mouse click menu in the Spectrum view. You can also clear selected spectra using the Trace Setup dialog accessed by doing a right mouse click in the view and then selecting Properties.




Fig. 307

Ratio View

The Ratio view displays two PDA wavelength channels and the ratio of those two channels. These may be viewed during real-time acquisition as well as during post-run analysis. The flat tops on the ratio peaks are a preliminary indication of peak purity. The ratio wavelengths and parameters are set in the **Method /PDA Options Ratio** tab.

The Y-axis of the ratio chromatogram is auto-scaled to 1.




Fig. 308

Right-click within the window to display the pop-up context menu. The context menu is the same as for a standard chromatogram graph window.

PDA Utilities

The Utilities menu allows displayed spectra to be printed, copied, saved or exported. The Utilities menu is available from the right mouse click menus of various PDA views.

Print

Select Print to automatically print the currently displayed spectra.

Copy to Clipboard


Select Copy to Clipboard to copy the displayed spectra to the Clipboard. The contents of the Clipboard may then be pasted into other software.

Save trace

Select Save trace to save the spectrum to a file with an .spc extension for later inclusion in a library or report.

Spectrum Similarity Table

The Spectrum Similarity table is displayed in the PDA Mixed View when the Options Mode is set to Show Spectrum Similarity Table. Spectra are automatically added to and deleted from this table whenever a spectrum is added to and deleted from the Spectrum pane of the Mixed View.



Spectrum Name	Similarity	Spectrum Source
✓ 5.80 Min	1.0000	(Current Data)
pna spe...	0.0510	C:\enterprise\Pr...
6.90 Min	0.4373	(Current Data)

Fig. 309

Spectrum Name

This column displays identifications for each of the spectra added to the table. For spectra extracted from the 3D data, the identification includes the time. For spectra extracted from a chromatogram based on a peak name, the identification includes the peak name. For spectra loaded from a file, the identification includes the retention time.

A check mark next to the colored line for a spectrum indicates that the spectrum will serve as the Reference Spectrum for the similarity calculation.

Spectra may be added to the table in the Mixed View by selecting the **Add to Similarity Table** from the **Actions** button menu.

Similarity

This column displays the similarity of the spectrum on that row relative to the Reference Spectrum. The Reference Spectrum is determined by double-clicking on a row of the table or by highlighting a row and pressing the 'Set Reference' button.

Until a reference spectrum has been selected the similarity for all spectra in the table will be zero.

Refer to the PDA Analysis and Calculations section for details on how this calculation is performed.

Source

This column displays the source of the spectrum. (Current Data) refers to working spectra extracted from the current 3D data. When spectra come from a stored data file, this column will display the source filename.

Set Reference

Pressing this button will designate the spectrum in the currently highlighted row of the table as the new reference spectrum. The similarities for all of the rows will be recalculated based on the new reference spectrum.

Print

Pressing this button will output a simple text report representation of the table to the default printer.

Properties

Pressing this button will display a dialog with parameters related to the calculation of similarity.

Similarity Table Properties

Enter the wavelength range and data point spacing (in nm) for similarity calculations.




Fig. 310



During analysis, it is possible that a portion of this wavelength range will be outside the range of the acquired data. In that case, the wavelength range will be truncated to the limits of the acquired data.

Spectral Library Definition

This window appears when you select the File/Spectral Library/New command and is used to define a spectral library to use with the PDA option. In order to define a library, you must have the spectra to be added saved on disk as spectra files (.spc) or you can create a library from the spectra obtained during the current run.




Fig. 311

Spectrum Data Source

Click on the arrow in this field to select either **Current Data** or **Spectrum File** to be used to select the spectrum to be added to the library.

To add a spectrum from a spectrum file, the spectrum must have been previously saved as a file with the .spc extension, typically, using the **Utilities /Save Trace** function of the pop-up menu in the Spectrum pane.

If you chose to enter a spectrum from the current data file, a dialog will appear where you must select the retention time at which you wish to select the spectrum.

Component Name

Enter the chemical name or description of the spectrum. This field is used for display and search pre-filter purposes.

Lambda Max From ... To

The wavelength range of the spectrum over which the software will calculate the lambda (absorbance) max value is displayed in these columns.

Retention time

When the spectrum was extracted from a chromatographic run, the software automatically identifies the retention time of the spectrum and enters the value in this field. When desired, the value in this field may be edited. This field can also be used as a search pre-filter.

Comment

Enter any descriptive information desired related to the spectrum. This field is used only for display and documentation purposes.

Library Notes

This field is used to document relevant information about the library as a whole, e.g., documentation about run conditions or general sample information.



Fields in existing library entries may be edited by selecting the fields with the cursor. Right-click on a row in the library table to display a pop-up menu used to cut, paste, copy, and insert and delete lines from the library as necessary.

Spectrum Information

Select the retention time at which you wish to use the spectrum, or the name of the peak whose spectrum is to be used.




Fig. 312

How to Collect Spectra for a Library

1. Perform an acquisition using a sample with known components or standards to be placed in the library.
2. Perform an analysis using the appropriate integration parameters. (Ideally, Spectral Background Correction would be checked in the Method/PDA Options/Spectrum tab.)
3. Using the contour plot in the Mixed view window, drag the vertical cursor to the apex of a peak to display the corresponding spectrum in the spectrum view. See also Extracting Spectra and Chromatograms from the Contour.
4. Right-click within the spectrum view to display the Spectrum pop-up menu, select Utilities/Save Trace and enter a file name in the dialog box. Repeat this procedure for each spectrum to be added to the library. Files are saved automatically with the .spc extension.



To add the 1st or 2nd Derivative of a spectrum to the library, select the appropriate filter in the Method/PDA Options/Spectrum tab and repeat steps 3 and 4.

How to Add Spectra to a Library

1. When all of the spectra are saved as .spc files, select File/Library/New to display the Library Definition dialog box and create a new library. Alternatively, select File/Library/Open to add spectra to an existing library.
2. Click in the Spectrum File cell of row 1 to display the Open dialog box. Double-click on the appropriate .spc file in the list box. The .spc

file name is entered into that cell and the associated spectrum is simultaneously displayed. A component name and comment may be entered in the appropriate columns. Repeat this procedure with subsequent rows in the table until all spectrum files are entered into the library.

3. Select File /Library /Save As and enter a name for the library. The .lib extension is automatically appended.

A new easy way to add spectra to an opened library is to select the option Add Spectrum to Spectrum Library... from the Actions... menu. Move the cursor in the chromatogram view to the appropriate position. The corresponding time is displayed in the dialog box. Simply click on Add to save the spectrum in the library.




Fig. 313

Library Search

How to Perform a Manual Library Search

1. In the Contour map of the Mixed View window, drag the vertical cursor to the apex of the peak of interest to display the corresponding spectrum in the spectrum pane.
2. Click on the Actions/Search library... Click either Method to use the library parameters from the method, or Quick to enable you to modify the search parameters. Before you do the search, make sure you have either selected a library in your method, or have opened a library to do a quick search.

Click Search Now to display the Library Search Results window, showing the three closest matches in the specified library. When appropriate, click on the >> or << button to display additional matches.

Spectral Library Search

This window appears when you click on the Actions/Search Library command from the Mixed View window. It allows you to perform a library search on a spectrum you have selected in the contour map, and also lets you select a spectrum from the current data or stored spectrum file to search.




Fig. 314

Hits Displayed

Select the number of hits you wish to display. Hits will be displayed in order of similarity.

Search

Select **Method** if you want to use the parameters from the current method for your search. Select **Quick** if you want to enter or change the parameters yourself before the search.

Display Parameters

Click this button if you want to display the parameters. If you always use the method parameters, you may not wish to display them when you search.

Search Now

When you click this button, a library search will be performed on the selected spectra using the parameters, and the number of hits requested will be displayed.

Search Spectrum

Click the arrow to select a new spectrum to search on. You can choose from the current data, or from a stored spectrum file. The current spectra source will be displayed.

Search Parameters

This tab contains search parameters. If you have elected to use the Method search parameters, the parameters displayed will be those in your current method in the PDA Options/Library tab. If you have elected to do a Quick search, the parameters will be default parameters.

Wavelength Range

These specify the wavelength range over which the library search will be performed.

Wavelength Step

This specifies the data point spacing to be used when a library search is performed.

Max Hits

This specifies the number of hits that will be reported in the results of a library search. Note that this works in conjunction with the *Similarity Threshold* parameter to limit the number of hits reported.

Similarity Threshold

Specifying a value for this control will cause the library search results to only display matches whose similarity to the unknown exceeds this value.

Note that this control works in conjunction with the *Max Hits* parameter to limit the number of hits reported.

Library

Displays the library from the method, or for Quick enables you to select a library to be used.

Pre-Filters tab

When Quick is selected as the Search Mode, the items on this tab allow the user to specify search pre-filters that will be performed on library spectra prior to the test for similarity. All pre-filters are optional. If Method is selected as the Search Mode, the items on this tab are read-only and reflect the parameters values on the Library tab of PDA Options.

Retention Time Range

When a Retention Time Range is specified, library search will be restricted to those library entries whose retention time is within the specified range. Entries outside this range are automatically excluded from the search (no similarity calculation is made). Entering a value for this is optional.

Lambda Max

When one or more of these values is specified, library search will be restricted to those library entries containing a lambda max within 5 nm of one of the specified values. Entries without a matching lambda max are automatically excluded from the search (no similarity calculation is made). Entering values for this is optional.

Compound Name Filter

When a Compound Name Filter is specified, library search will be restricted to those library entries whose name contains the specified string as a case-insensitive sub-string. Entries without a matching sub-string are automatically excluded from the search (no similarity calculation is made). Entering a value for this is optional.

Custom Report

A variety of PDA information can be placed in a custom report. These items are inserted in the report by placing the cursor at the location where you want to insert the item, then do a right mouse click to access the custom report menus.

PDA Insert Graph Items

The Insert Graph menu has the following additional items available when the PDA Option is enabled: 3D Data Graph, and Contour Graph.

3D Data Graph

The 3D View provides a three-dimensional view of absorbance versus time and wavelength. Wavelengths of appreciable absorbance and interferences, which may be invisible in a single wavelength plot, are easy to locate with the 3D View. Select this function to automatically enter the 3D map into the report. Click the right-mouse button within the 3D map and select Properties to display the 3D Properties dialog box and enter appropriate changes

The 3D map can be elevated and rotated around its axis for display from any angle. These functions work the same as in the 3D view window.

Insert Contour Graph

A data file containing PDA data may include a Contour Plot in a custom report. The contour view (also referred to as Isoabsorbance plot) provides an aerial view of the absorbance of the sample at each wavelength vs. time. The contour view supplies quick and easy-to-assimilate information about those wavelengths at which the sample exhibits appreciable absorbance.

Changing parameters for the contour graph in a custom report works exactly as in the Contour view.

PDA Insert Report Items

When the PDA option is enabled, the following report items will be available to insert into a custom report: Library Search Report, Purity Report, and Spectrum Report.

Library Search Report

This command will insert a table of library search results in your custom report.

To modify the parameters of the report, do a right mouse click in the report table, then select Properties. The following dialog will appear.




Fig. 315

Search Libraries for

This specifies what peaks should be searched for the report. When *Apex Spectra of Peaks* is selected, the *Specific Spectrum* button is disabled. When *A Specific Spectrum* is selected, the other controls of this group are disabled.

Detected on

This allows the user to select the channel or channels from which peaks will be detected.

Limit to named peaks only

This allows the user to restrict the peak selection (made above) to a time range. When checked, only peaks within the specified time range will be searched. When this box is unchecked, the Time Range edit fields are disabled.

Limit Peaks to time range

This allows the user to restrict the peak selection (made above) to a time range. When checked, only peaks within the specified time range will be searched. When this box is unchecked, the Time Range edit fields are disabled.

Time Range

This allows the user to restrict the peak selection (made above) to a time range. When the **Limit to time range** box is unchecked, these fields are disabled.

Specific Spectrum button

Pressing this button will display a menu allowing the user to specify the spectrum to be searched. The menu will allow the user to specify a spectrum based on one of the following choices: Current Data, Named File, or Spectrum Peak.

Spectral Display

The in this group specify the layout of the search results graphs.

Graph Height

This specifies the relative height of a search results graph. A value of 100% corresponds to a standard sized graph. A larger value may be selected to provide more a detailed plot.

Graph Properties

Pressing this button displays the graph properties dialog as found on the standard Spectrum graph.

Library Definition Report

Selecting this option will insert a Library Definition Report into your custom report. To specify parameters for the report, do a right mouse click in the report table, then select **Properties**. The following dialog will appear.



Fig. 316

Enter or select the library you wish to include in the report, then click OK.

Purity Report

Selecting this option will insert a Purity Report into your custom report. To specify parameters for the report, do a right mouse click in the report table, then select Properties. The following dialog will appear.




Fig. 317

Detected on

This allows the user to select the channel or channels from which peaks will be detected.

Limit to named peaks only

This allows the user to restrict the peak selection (made above) to a time range. When checked, only peaks within the specified time range will be searched. When this box is unchecked, the Time Range edit fields are disabled.

Limit Peaks to time range

This allows the user to restrict the peak selection (made above) to a time range. When checked, only peaks within the specified time range will be searched. When this box is unchecked, the Time Range edit fields are disabled.

Time Range

This allows the user to restrict the peak selection (made above) to a time range. When the *Limit to time range* box is unchecked, these fields are disabled.

Graph Height

This specifies the relative height of a search results graph. A value of 100% corresponds to a standard sized graph. A larger value may be selected to provide more a detailed plot.

Graph Properties

Pressing this button displays the graph properties dialog as found on the standard graph.

Spectrum Report

This will display a table showing spectra extracted from peaks of the chromatogram. Based on the selection on the Spectrum tab of PDA Options, the report will contain either the apex spectrum, an average of the upslope, apex and down slope spectra, or several averaged spectra. To set parameters for this report, do a right mouse click on the table, then select Properties. The following dialog will appear.




Fig. 318

Report peak spectra

Detected on

This allows the user to select the channel or channels from which peaks will be detected.

Limit to named peaks only

This allows the user to restrict the peak selection (made above) to a time range. When checked, only peaks within the specified time range will be searched. When this box is unchecked, the Time Range edit fields are disabled.

Limit Peaks to time range

This allows the user to restrict the peak selection (made above) to a time range. When checked, only peaks within the specified time range will be searched. When this box is unchecked, the Time Range edit fields are disabled.

Time Range

This allows the user to restrict the peak selection (made above) to a time range. When the *Limit to time range* box is unchecked, these fields are disabled.

Report information

Checking any of these controls indicates that the indicated value will be calculated and printed to the right of the graph. Unchecking values that are not of interest will speed up analysis.

Peak Area

Checking this indicates that the peak area should be printed to the right of the spectrum graph.

Lambda Max

Checking this indicates that the three largest lambda max values should be printed to the right of the spectrum graph.

Lambda Min

Checking this indicates that the three largest lambda min values should be printed to the right of the spectrum graph.

Total Purity

Checking this indicates that total peak purity should be printed to the right of the spectrum graph.

3 Point Purity

Checking this indicates that 3 point peak purity should be printed to the right of the spectrum graph.

Similarity

Checking this indicates that peak apex similarity to the reference spectrum should be printed to the right of the spectrum graph. The use of this value requires that a reference spectrum be included in the peak table. See the section on *Analysis and Calculations* for details.

Spectral Display

The controls in this group specify the content and labeling of the spectrum graphs.

Show Lambda Maxima

Checking this box will annotate each spectrum graph with the 3 largest lambda max absorbance.

Show Lambda Minima

Checking this box will annotate each spectrum graph with the 3 largest lambda minima absorbance.

Background Correction

Checking this will cause each apex spectrum to be corrected for background using the chromatographic baseline prior to being used elsewhere in analysis. Refer to the *PDA Analysis and Calculations* section for details on the formula used.

Filtering Type

This allows the user to specify a mathematical filtering function to be performed on all spectra extracted from the 3D data during analysis. Refer to the *PDA Analysis and Calculations* section for details on the formula used for each filter.

Graph Height

This specifies the relative height of the spectrum graph. A value of 100% corresponds to a standard sized graph. A larger value may be selected to provide more a detailed plot.

Graph Properties

Pressing this button displays the graph properties dialog as found on the standard Spectrum graph.

Peak Table

When using the PDA option, certain columns will be added to the peak table that enable you to analyze peaks from the PDA detector.

Detection

Select the basis for the identification of the peak. If you choose Ret Time, only the retention time will be used for identification of the peak. If you choose Ret Time with Spectral Confirm, the Similarity of the peak's spectrum to that of a designated reference spectrum will be used in addition to the retention time as the basis of peak identification.

Spectrum

If you want Similarity to be used as a basis for peak identification, then click on the arrow to the right of this field to specify the stored reference spectrum to be used for comparison. During identification, this reference spectrum is compared to the peak apex spectrum and a similarity index is

computed. A peak is considered identified if this calculated similarity index is at least the value specified in the Similarity column of the peak table.

If Similarity is not specified as a basis for peak identification, then this field is ignored.

Similarity

If Similarity is specified as a basis for peak identification, then this field specifies required minimum similarity index for a peak to be considered identified. During identification, the reference spectrum (see previous section) is compared to the peak apex spectrum and a similarity index is computed. A peak is considered identified if the calculated similarity index is at least the value in this field.

If Similarity is not specified as a basis for peak identification, then this field is ignored.

Analysis Channel

Specify which Diode Array wavelength channel is to be used for analysis of the peak. The choices will be those specified in **Instrument Setup/DAD**.

PDA Analysis and Calculations

This section describes the various types of calculations that are related to PDA data and how the data is analyzed.



In order to maintain accuracy during the application of multiple operations, all calculations are performed using double-precision floating-point numbers.

Chromatograms Extracted from the 3D Data

Two types of chromatograms may be extracted from the 3D data:

- **Multi-Chromatogram Channels:** One or more chromatograms defined on the *Multi-Chromatogram* tab of *PDA Options*
- **Working Chromatogram:** A single chromatogram extracted from the Contour plot on the *Mixed View* display. (This chromatogram is displayed in the *Chromatogram View*)

Specific Analysis-related capabilities for these chromatograms are detailed below.

Multi-Chromatogram Channels

The following apply to multi-chromatogram channels:

- In the Channel Selection drop-down list, each multi-chromatogram channel has a unique identifier that includes the wavelength and bandwidth of the channel.
- Each Multi-Chromatogram Channel has a separate Integration Events Table, Manual Integration Table, Export Table and Performance Table associated with it.
- All Multi-Chromatogram Channels share a common peak and group table. Within the peak and group tables, any channels may be selected as the analysis channel for quantitative information.
- When an analysis is performed, all multi-chromatogram channels are automatically analyzed.

- Each multi-chromatogram channel is an average of the absorbencies monitored at each wavelength in the wavelength range. The wavelength range is equal to the selected wavelength +/- one half of the bandwidth.

Working Chromatogram

The following apply to the working chromatogram:

- In the Channel Selection drop-down list, a single entry is added that identifies the channel as applying to the working chromatogram.
- This channel has a separate Integration Events Table, Manual Integration Table, Export Table and Performance Table associated with it. These are separate from the tables of the Multi-Chromatogram Channels (above).
- When an analysis is performed, the working chromatogram will automatically be analyzed.
- The working chromatogram is an average of the absorbencies monitored at each wavelength in the wavelength range. The wavelength range is equal to the selected wavelength +/- one half of the bandwidth.
- The working chromatogram is not available as a trace in Custom Reports.

Spectra Extracted from the 3D Data

Two types of spectra may be extracted from the 3D data:

- Analysis Spectra: One or more spectra automatically extracted from the 3D data during analysis and used for peak identification, similarity and/or library search.
- Working Spectrum: A spectrum extracted from the Contour plot on the Mixed View display. (This chromatogram is displayed in the Spectrum View)
- Specific Analysis-related capabilities for these spectra are detailed below.

Analysis Spectra

The following applies to the analysis spectra:

- Prior to being used elsewhere, all analysis spectra extracted from the 3D data are filtered according to the settings on the Spectrum tab of PDA Options.

Working Spectrum

The following apply to the working spectrum:

- Working spectra extracted from the 3D data are never filtered according to the settings on the Spectrum tab of PDA Options.
- The working Spectrum is not available as a trace in Custom Reports.

Background Correction

Background Correction may only be performed on a spectrum as a result of the settings on the *Spectrum* tab of *PDA Options*, the settings on the *Purity* tab of *PDA Options*, or by setting the properties of the Peak Purity Pane.

Background subtraction is always the first operation performed on a spectrum. If a background correction is attempted on a spectrum after another operation has already been performed (including an operation

specified in the *Spectrum* tab of *PDA Options*), a message box is displayed and the background correction request ignored.

Background subtraction may only be performed once on a spectrum. If a second background correction is attempted on a spectrum after the operation has already been performed (including the operation specified in the *Spectrum* tab of *PDA Options*), a message box is displayed and the background correction request ignored.

A background correction is performed as follows:

- The spectra from the baseline start and baseline stop times for the peak are extracted from the 3D data. The Max Plot is used to determine the peak that is used.
- For each spectrum in the peak, a corresponding background spectrum is generated by linear interpolation between the baseline start and baseline stop spectra.
- These background spectra are subtracted from the original spectra.

Spectrum Interpolation

Spectrum interpolation may be performed on a spectrum as a result of the settings on the *Spectrum* tab of *PDA Options* or by selecting *Operations > Interpolate* from the Spectrum View context menu.

Interpolation may only be performed once on a spectrum. If a second interpolation is attempted on a spectrum after the operation has already been performed (including the operation specified in the *Spectrum* tab of *PDA Options*), a message box is displayed and the interpolation request ignored.

Interpolated spectra may not be stored in spectrum libraries.

The calculation is performed by doing a 10:1 interpolation of the spectrum data points using a cubic spline curve fit. This interpolation is performed after the applying any spectral filtering option (1st derivative, 2nd derivative or smooth) to the spectrum.

Spectrum Smoothing

Spectrum smoothing may be performed on a spectrum as a result of the settings on the *Spectrum* tab of *PDA Options* or by selecting *Operations > Smooth* from the Spectrum View context menu.

Interpolation may be performed repeatedly on a spectrum.

The calculation is performed by doing a 9 point Savitsky-Golay smooth on the spectrum data points.

Spectrum Derivatives

Calculation of the 1st and 2nd derivatives of a spectrum may be performed on a spectrum as a result of the settings on the *Spectrum* tab of *PDA Options* or by selecting *Operations > Smooth* from the Spectrum View context menu.

Derivatives may be computed repeatedly on a spectrum.

The absorbance values of the 1st derivative of a spectrum are computed by calculating the differences between adjacent absorbance values to create a new spectrum. The 2nd derivative of a spectrum is defined as the 1st derivative of the 1st derivative of the spectrum.

Upslope and Down slope Spectra

The upslope and down slope spectra of a peak are identified by calculating the second derivative of the portion of the chromatogram containing the peak. The two times at which the second derivative plot crosses zero are known as the inflection points. For normal peaks (i.e. non-Negative peaks), the Upslope Spectrum is the spectrum at the time represented by the first inflection point, while the Down slope Spectrum is the spectrum at the time represented by the second inflection point.

Library Search Calculations

General

During analysis, if one or more search libraries are defined on the *Library* tab of *PDA Options*, then an automated library search is performed on every integrated peak of every PDA analysis channel.



Unless a *Library Search Results* object is part of the method custom report, no automated library searching will be done when analysis is performed.

In this section, a *Query Spectrum* is defined as the unknown spectrum that is being searched. A *Reference Spectrum* is defined as a spectrum from a spectrum library file.

During a search, the apex spectrum of the peak (the query spectrum) is compared to each spectrum contained in the libraries (reference spectra) to determine the similarity of the query spectrum to the reference spectrum. The similarity is quantified through the calculation of a Similarity Index for each query/reference pair. The Similarity Indices are used to generate a hit list of the best matching entries. A perfect match will have a Similarity Index of 1.0000. Similarity indices less than 1 indicate differences in the spectral patterns.

If the query spectrum and the reference spectrum have different wavelength ranges, then the intersection of the two ranges is used in the similarity calculation.

If the query spectrum and the reference spectrum have different wavelength steps, then the higher resolution spectrum is de-resolved to match the other spectrum before being used in the similarity calculation.

Pre-Filters

A pre-filter is a criterion on a reference spectrum that must be met before that spectrum is used in similarity calculations. One or more pre-filters may be specified in the search parameters. When multiple pre-filters are specified, a reference spectrum must meet all of the individual pre-filters in order to be considered for similarity calculations.

Reference spectra that do not meet all of the pre-filter criteria are automatically excluded from calculations and from being a candidate for a hit list. No similarity calculation is performed on these spectra.

The following pre-filters are supported:

Retention Time Range

When a Retention Time Range is specified, a reference spectrum is excluded unless the retention time for the spectrum in the library is within the specified range.

Lambda Max

Up to three Lambda Max values may be specified. When one or more Lambda Max values are specified, a reference spectrum is excluded

unless all of the specified maxima are within 5 nm of a maxima for the spectrum in the library.

Compound Name

When a Compound Name is specified, a reference spectrum is excluded unless the specified string is a case-insensitive sub-string of the compound name for the spectrum in the library.

Ratio Chromatogram Calculation

The ratio chromatogram plot consists of data points calculated as follows:

$$\text{Ratio Pt.} = \text{abs}_1 / \sqrt{(\text{abs}_1^2 + \text{abs}_2^2)}$$

where: abs_1 = the absorbance in chromatogram 1 at this wavelength
 abs_2 = the absorbance in chromatogram 2 at this wavelength

If the absorbance of any point in chromatogram 1 or chromatogram 2 is less than the threshold, the Ratio Pt. for that wavelength is set to 0.

Similarity Calculations

The Similarity Index (SI) compares two spectra across the wavelength range defined in *Method > PDA Setup > Library*. A Library Search is performed using the Similarity Index, determined as follows:

A spectrum is considered as an array of absorbencies at each wavelength.

$$(a(\lambda_1), a(\lambda_2), \dots, a(\lambda_n))$$

Where $a(\lambda_i)$ is the absorbance at λ_i

The spectrum can also be represented as an n-dimensional vector, with each absorbance corresponding to one dimension of the vector, as follows:

$$S = (a(\lambda_1), a(\lambda_2), \dots, a(\lambda_n))$$

When two spectra are obtained from the same compound, the ratio between corresponding elements in S1 and S2 is constant; consequently, the vectors have the same direction. In such a case, the angle between the two vectors would be 0. As a general rule, the greater the angle between the two vectors, the less similarity they have. This relationship is calculated as follows:

$$SI = \cos(\Theta)$$

where Θ = The angle between the two vectors

The closer the SI value is to 1, the more similar two spectra are considered to be. Two spectra are considered to be identical when the SI value is close to 1.

Lambda Max/Min Calculations

A lambda max (min) is defined as a local maximum (minimum) of the absorbance values of a spectrum. The user must define the wavelength range over which the calculation should be performed.

To compute n lambda max (min) values for a spectrum, the software will find the n local maxima (minima) with the largest (smallest) absorbance values.

The calculation range of Lambda Min/Max is based on the X-axis settings for the graph in the *Axis Setup...* context menu item of the graph.

Noise Spectrum Calculations

A noise spectrum is used in Purity calculations to compensate for solvent and detector effects. The spectrum is computed using a 3D data file and a user-supplied retention time range.

The noise spectrum is an array of data points over the spectrum range. Each data point in the noise spectrum is 3 standard deviations of the mean of the absorbance values for that wavelength over the user-supplied time range of the 3D data.

Peak Purity Calculations

Peak purity calculations use the spectra from 3D PDA data to assess the purity of a chromatographic peak. Spectra from the peak are compared to the apex spectrum to compute a measure of uniformity of the spectra.

Background Correction

Background Correction may only be performed on a spectrum as a result of the settings on the **Spectrum** tab of **PDA Options**, the settings on the **Purity** tab of **PDA Options**, or by setting the properties of the Peak Purity Pane.

Background subtraction is always the first operation performed on a spectrum. If a background correction is attempted on a spectrum after another operation has already been performed (including an operation specified in the **Spectrum** tab of **PDA Options**), a message box is displayed and the background correction request ignored.

Background subtraction may only be performed once on a spectrum. If a second background correction is attempted on a spectrum after the operation has already been performed (including the operation specified in the **Spectrum** tab of **PDA Options**), a message box is displayed and the background correction request ignored.

A background correction is performed as follows:

1. The spectra from the baseline start and baseline stop times for the peak are extracted from the 3D data. The Max Plot is used to determine the peak that is used.
2. For each spectrum in the peak, a corresponding background spectrum is generated by linear interpolation between the baseline start and baseline stop spectra.
3. These background spectra are subtracted from the original spectra.



While compensation for background will provide more precise purity calculations, it can slow down re-analysis of large data files.

Calculating Total Purity

Total purity for an integrated peak is calculated as follows:

1. The purity calculation range of the peak is determined by starting from the peak apex and working toward peak start and peak end. Each chromatographic point in the peak is compared to the absorbance threshold. The start and stop values of the range are defined as the point at which the absorbance drops below the absorbance threshold.
2. For each point in the calculation range, a similarity index is calculated by comparing the spectrum for that point to the peak apex spectrum.

3. All spectra whose similarity index meets or exceeds the Purity Threshold will be considered similar to the apex spectra.
4. A sum is calculated of the chromatographic area represented spectrally similar points.
5. The Purity Index is calculated by dividing this sum by the total area over the calculation range of the peak.

The purity index will have a range of between 0.000000 and 1.000000.

Three Point Purity

Three Point Purity is calculated by comparing the spectrum at the apex of a peak (**Point 1**) with the spectra at the up-slope (**Point 2**) and down-slope (**Point 3**) of the peak. The up-slope spectrum is defined as the spectrum located 20% of the peak width from the peak start. The down-slope spectrum is defined as the spectrum located 20% of the peak width from the peak end.

$$\text{Time}_{\text{up}} = \text{Time}_{\text{start}} + (\text{Time}_{\text{apex}} - \text{Time}_{\text{start}}) * 1/5$$

$$\text{Time}_{\text{down}} = \text{Time}_{\text{apex}} + (\text{Time}_{\text{end}} - \text{Time}_{\text{apex}}) * 4/5$$

Similarity indices are generated for the up-slope and down-slope spectra, relative to the apex. Since spectra taken at different points along a pure peak will look identical, they will have high similarity indices. The closer the Similarity Index is to 1.0000, the more similar or more pure the peak is deemed to be.

While it is easy to decide that a peak with a very low Purity value (0.0000-0.8900) is indeed not pure, it is more difficult to decide whether peaks with purity values between 0.9000 and 0.9500 are actually contaminated. The user's judgment is required to make this determination and may necessitate using other methods to determine purity, such as overlaying the spectra or comparing first and second derivatives.

Spectrum Export

The format for spectrum export is as follows:

Export Type:	Spectrum
Version:	2
Sample ID:	pna_6 (1:2)
Data File:	D:\Data\PNA_Level_2_vp.dat
Method:	D:\Methods\test.met
User Name:	Tontala
Spectrum Type:	Peak Apex
Operations Applied:	Smooth
X Axis Title:	nm
Y Axis Title:	mAU
X Axis Multiplier:	1.0
Y Axis Multiplier:	1.0
Number of Spectra:	23
Peak Ret. Time:	13.440
Number of Points:	201
	200 193
	201 30
	202 -7

PDA Data Export

To export the PDA data, at first a multichromatogram must be enabled in the PDA options. Then elect the new multichromatogram in the channel selector and open *Method – Advanced...* In the *Data Export* tab enable the export and select *PDA Data* from the drop-down menu. Add *ASCII 3D Data* to the export parameters. Select the Field separator and the path for the export files. Store the method. If you analyze PDA data with the modified method, the export file will be created. If you open the file in Microsoft Excel, note that for each wavelength a separate column is required. Excel supports up to Excel version 2003 256 columns. If the scan range includes more than 256 wavelength steps, not all will shown.

ChromGate® Preparative Option

To use the ChromGate® Preparative Option it has to be installed separately as described in KNAUER FRC control option on page 8. Then, the implemented fraction collectors will appear in the device list for creating instruments.




Fig. 319 Fraction collector selection window

Fraction Collector Configuration

The configuration windows are the same for different fraction collectors. Therefore, the procedure is the same.




Fig. 320 Fraction collector configuration window

FC Name

Enter a descriptive name for the collector. This name will appear on the fraction collector tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Serial Number

Type in the serial number of the detector you are using.

Rack Type

Select from the drop down list the required rack type. The possible choices depend on the Type of fraction collector. The **Vials number** will be displayed automatically.

Vial Volume [ml]

The rack depended value will be displayed automatically. However, you are allowed to overwrite it by self-defined values. If any fraction is larger than the set volume the next free collection position will be used for continued collection.

Serial Port

Select from the drop-down list the serial port number for the communication port on your PC where the collector is connected.

Interface (S 3050 / MultiValve fraction collector / Foxy R1/R2 only)

Select the Interface you want to use for communication with the device.

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the autosampler is connected.

LAN (S 3050, MultiValve fraction collector)

To let the software search for your already connected and switched-on device, click on the  button to start the search for all connected LAN devices of the currently configured type. The option "Use S/N to identify the instrument" must be enabled. Select the desired device from the shown list. The software will read-out the information as serial number and will fill it automatically into the configuration screen. To search for an instrument, the device must be connected, switched on and must have a valid IP address.

As an option, you can enter the IP address of your device, if the option "Use S/N to identify the instrument" is disabled. Please note that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.

LAN (Foxy FC R1/R2)

There is no option for the Foxy fraction collectors to search via LAN or read-out the configuration. If selected, please enter the IP address of the FC. The IP address must match with the IP address in the FC; the Foxy FC R1/R2 do not support LAN with DHCP. The LAN communication is supported in **firmware version 1.01.03** and higher.



Please note that the touch pad of the Foxy R1/R2 FC will NOT be blocked if the device will be controlled by software. The software cannot recognize inputs via the touch pad. Therefore it is highly recommended, not to use the touch pad under software control. Otherwise the vial position of fractions and other information may not be stored correctly in the „Device Monitor - Show Rack“.



The communication port of the fraction collector S 3050 must be select using a switch on the collector's rear panel.

Tubing Parameters




Fig. 321 Tubing parameters setup

These settings enable someone to consider the time delay for transporting the sample from the flow cell to the collector valve. If the option **Use recycling valve** is checked, the capillary from the flow cell via the recycling valve to the fractionation valve must be taken into account.

If you activate the option **Fixed time** you can enter this time delay directly. This value will be valid, even if the flow rate is changed!

Entering the **Capillary Volume** directly or calculating it by entering the **Capillary Length** and **Inner Diameter**, the delay time will be adapted to the programmed flow rate of the method.

Solvent and/or Peak Recycling

For an optional solvent and/or peak recycling, for most fraction collectors an additional recycling valve is necessary. As a recycling valve only a Knauer 6-port valve can be used. The Knauer fraction collector 3050 includes a valve for recycling. For a multi valve fraction collector some ports of the fractionation valve can be defined for the recycling.

To define a separate recycling valve, check the option **Use recycling valve** and click on the  button. For the separate recycling valve a COM port or the IP address is to define. The settings for the volumetric delay will be displayed without access for changing. The settings are valid as well for the recycling valve as for the fractionation valve. Therefore, the capillary between the recycling valve and the fraction collector should be as short as possible.

Only for the Knauer fraction collector 3050 and the multi valve fraction collectors the option **Use own ports for recycling** will be displayed. If you want to use this option, click into the check box for this option. For the multi valve fraction collector you will find, that the number of available ports for fractionation will be reduced by 2. Additionally the ports numbers for solvent and peak recycling will be shown. It is not possible only to enable one of the recycling options.




Fig. 322 Recycling valve setup

The dialog box for the port settings will be active only for fraction collectors with different options.




Fig. 323 Port settings

Virtual Fraction Collector Configuration




Fig. 324 Virtual fraction collector configuration window

The virtual fraction collector is a tool to optimize the fractionation without need for sample or eluents. Therefore, only opening and closing the configuration window would be sufficient. However you have two areas for individual settings. .

Rack Type

Select from the drop down list that rack type best in comparison with your real one. The choice is the same as for the Büchi fraction collector: Micro with 240 vials, standard with 120, and preparative with 48 positions. The **Vials number** will be displayed automatically.

Vial Volume [ml]

The rack dependent value will be displayed automatically. However, you are allowed to overwrite it by self-defined values. If any fraction is larger than the set volume the next free collection position will be used for continued collection.

Tubing Parameters

These settings enable someone to consider the time delay for transporting the sample from the flow cell to the collector valve. Here you can adapt it to the real conditions of your system.

The solvent and peak recycling options you cannot simulate.

Multi Valve Fraction Collector Configuration




Fig. 325 Multi valve fraction collector configuration window

Name

Enter a descriptive name for the collector. This name will appear on the fraction collector tab when the instrument setup or status windows are displayed. Use a name which is unique within the instrument.

Tubing Parameters

The same setup procedure as with any other fraction collector must be done. However you cannot define the tubing parameters for each included valve. One setting will be valid for all cascaded valves!

Solvent and/or Peak Recycling

With the Knauer MultiValve fraction collector the optional solvent and/or peak recycling an additional recycling valve can be used but it is not necessary.



Fig. 326 Additional recycling valve setup

For an additional valve its communication port (serial or LAN) must be configured.




Fig. 327 Own port recycling valve setup for 12 position valve

For use of own ports no further settings are required.

Valves number

Enter the number of valves which build up the multi valve fraction collector. The number is in practical patterns unlimited.

Configure Valve

The drop down list offers as many positions as the valve number was entered above. The following settings have to be made for each valve separately.

Interface

Select the Interface you want to use for communication with the device.

RS-232

Select the number for the RS-232 communication port (**COM Port**) on your PC to which the autosampler is connected.

LAN

To let the software search for your already connected and switched-on device, click on the  button to start the search for all connected LAN devices of the currently configured type. The option "Use S/N to identify the instrument" must be enabled. Select the desired device from the shown list. The software will read-out the information as serial number and will fill it automatically into the configuration screen. To search for an instrument, the device must be connected, switched on and must have a valid IP address.

As an option, you can enter the IP address of your device, if the option "Use S/N to identify the instrument" is disabled. Please note that in case the device receives the IP address from a DHCP server, e.g. a router, the IP address may change if the device is switched off and switched on.

Valve type

The drop down list offers 2-, 6-, 12-, and 16-position valves. Select for each valve the desired type. If you read-out the valve via LAN, the appropriate valve type has been selected automatically.

Vial Volume [ml]

Enter any desired vial volume. If any fraction is larger than the set volume the next free collection position will be used for continued collection.

Vial number

The Vial number will be displayed automatically.

Fraction Collector Setup

In any system you can include only one fraction collector. However, the setup window is independent on the configured fraction collector. Therefore they will be explained together.

This window differs in one respect from other setup windows. It is not only important for setting up the fraction collector. You need it also during the runs. However, those aspects will be described later in the chapter instrument status of a (running) control method (page 137).




Fig. 328 Fraction collector setup window (upper part)

Detection Channel

Select from the drop down list of active detector channels that one to be responsible for fractionations.

Collector Reset

This option must be checked for cumulative collection from repeating runs.

Advanced collector settings

With a click on the  button the **Advanced collector settings** window will open.




Fig. 329 Advanced collector settings

Rack Settings

If the option **Exchange rack if full** is enabled, automatically a new, empty rack will be created if the current one is full. A copy of the full rack will be stored; the data from the fractions will not be lost. This allows changing the rack manually during the sequence without starting the sequence again. If this option is disabled, a sequence will stop if no empty vial is left for the next fraction.

The **Empty vial check** tests before the injection, if enough empty vials for the fractions are left. Enter the number of vials needed for one run. The software will check then before the injection, if this number of vials is left. If not, the system will not start the next run.

Stacked Injections

Restart from first vial of current run

If stacked injections are used, check this box if you want the collector to start for every stacked injection with the same vial as for the previous runs for the fraction collection. This is an option, if you inject the same sample several times and want to collect the fractions of the same substances from every run in the same fraction vials.

Adjust FRC table after stacked injection

This option allows you to copy the already executed part of the fraction table for the next injection, if a stacked injection is performed. In this case you will have the same collection conditions for every injection without any manual intervention. This makes sense, if you use Single Event parameters, which otherwise should have to be programmed separately for every injection. If you use the peak recognition mode, just program the collection condition for the whole run time, the **Adjust FRC table** option leave unchecked.

Pump

If more than one pump is configured in the system, select the used one from the drop down list for volume and delay calculations.

Update Flow

The set flow of the selected pump after the last actualization will be displayed. If the flow is not constant for the whole method, the value of 0.000 will be shown. Please enter then the flow rate for the time the fractions will be collected. The flow is important to calculate the time delay between detecting a component in the detector and the approach to the fraction valve as well as the filling level of the vials.

Fraction Collection Program

The editing principles are again the same as it was described for the pump gradient programming. Some fields are free for entering data whereas others are either only informative without access or you can select the possible options via a drop down list.

Exp. Time (min)

Enter the desired time values.

Mode

#	Exp. time (min)	Mode	Event
1		Single Event	Unconditional
2	10.00	Peak Recognition	Unconditional

Fig. 330 Fractionation mode selection

The single event option is to select for programming a detailed exact fractionation program, whereas the peak recognition is preferably for use in method development. Every new detected peak according to the defined thresholds (see Parameters on page 238) will be collected in the next empty vial. The selected mode will influence the further setting possibilities.

Event

The possible options depend on the selected mode:

#	Exp. time (min)	Mode	Event	Exp. init (min)
1	0.00	Single Event	Unconditional	
2	10.00	Single Event	Unconditional	

#	Exp. time (min)	Mode	Event	Exp. init (min)
1	0.00	Peak Recognition	Signal Level	
2	10.00	Peak Recognition	Signal Level	

Fig. 331 Event definition

Unconditional

Signal independent action as selected in the parameter window Fig. 332.

Signal Level

The collection takes place according to set level thresholds.

Signal Slope

The collection takes place according to set slope thresholds.

Signal Level/Slope

The collection takes place according to set level and slope thresholds which can be combined in logical **and** and **or** combinations.

Local Maximum

The collection will be started by a local maximum.

Local Minimum

The collection will be started by a local minimum.

Spectral Similarity

The collection takes place according to set similarity thresholds.

Exp. interv.

These fields are only informative according to the settings in the parameter window for the expected possible time shift for the retention.

Rel. time shift

These fields are only informative according to the settings in the parameter window if Use as reference for correction next events is activated, refer to page 242.

Parameters

Activating any field in this column, an arrow appears. Clicking on that a dialog box for the FC parameters at the given time will be opened. Its appearance depends as well on the mode as on the selected event.


Unconditional

Fig. 332 FC parameters unconditional (only single event)

Action

This selection, as far as accessible, is independent on the set mode (unconditional, signal level ...). The choices are Waste, Collect To, Peak Recycling, Solvent Recycling, None, Collect Slices.

Waste

The eluate will be wasted.

Collect To

The eluate will be led to the set position.

Peak Recycling

The eluate will be led back to the column (peak shaving).

Solvent Recycling

The eluate will be collected in a recycling container because its purity is better than the set limits.

None

Neither the recycling valve nor the fractionation valve will be switched. The previous action is continued, only a marker is set.

Collect Slices

The eluate will be collected in set volume fractions into the **next** free positions.

You have only to select the type of action and in case of collect to the destination vial. After closing the window the parameter will be designed as **Timed**.

The options *Relative to previous event* and *Use as reference for correction next events* will be explained generally at the end of the section parameters on page 242.

Signal Level

FC Parameters at 4.00 min

Signal

Ascending Descending

Signal Level, Units: 25

Expected interval (min): +/-: 0.5

Please notice that the parameters will be used in the Units selected for the Detection Channel

Relative to previous event (min): 0

Use as reference for correction next events

Action

Action: Collect To

next position: 11

OK Cancel

Fig. 333 FC parameters Signal Level (single event)

Enter the level threshold value and define it to be on the ascending edge of the peak. Then enter the expected possible shift (interval) for the retention time, select the action and the destination vial. After closing the window the parameter will be designed as **+L**.

FC Parameters at 4.50 min

Signal

Ascending Descending

Signal Level, Units: 25

Expected interval (min): +/-: 0.5

Please notice that the parameters will be used in the Units selected for the Detection Channel

Relative to previous event (min): 0

Use as reference for correction next events

Action

Action: Waste

next position: 11

Warning: this event may overlap with one or more events

OK Cancel

Fig. 334 FC parameters Signal Level (single event)

Enter the level threshold value and define it to be on the descending flange of the peak. Then enter the expected possible

shift (interval) for the retention time, select the action waste. After closing the window the parameter will be designed as **-L**.

The red displayed warning message for possible event overlapping must not be a reason to change the event programming but at least to check it. It may happen that one event cannot become active. In the shown case (Fig. 333 and Fig. 334) the time difference between both events is 0.5 min; however both allow a time interval of also 0.5 min leading to this possible overlapping.

In case of **peak recognition mode** and event **signal level** the parameter window only enables you to enter the threshold values for the peak at beginning and end as well as to check the solvent recycling option. No further settings are possible. The action is set to collect to next. After closing the window the parameter will be designed as **L** without leading sign.




Fig. 335 FC Parameters Signal Level (Peak Recognition)

Signal Slope

When using the slope instead of the level, the settings are made similarly. The only difference is that the slope values, defined as signal difference per 1 sec, are to be used as thresholds. After closing the window the parameter will be designed as **+S**, **-S** or **S**.

Signal Level/Slope

The settings using the slope and the level are again quite similar to perform. The only difference is that both the level and slope values are to be used as thresholds.




Fig. 336 FC parameters Signal Level/Slope (single event)

Additionally you have to select either the **and** option, meaning that both conditions must be exceeded, or the **or** option, meaning that only one of these conditions must be exceeded. This can be defined differently for the ascending and descending flange of the peak.




Fig. 337 FC parameters Signal Level/Slope (Peak Recognition)

After closing the window the parameter will be designed for example as **+L and S**, **-L or S** for single events and **Lb or Sb – Le or Se** for peak recognition.

Spectral Similarity

The settings using the spectral similarity are quite similar to perform. The difference is that the similarity values are to be used as thresholds. The maximum value is 1.000 for absolute identical spectra.




Fig. 338 FC parameters Spectral Similarity (single event)

Additionally you must select the desired reference spectrum. Clicking on ... will open the selection window.




Fig. 339 Spectrum selection window

After closing the window Fig. 338 the parameter will be designed as **+Sim**, **-Sim** or **Sim** for peak recognition.




Fig. 340 FC parameters Spectral Similarity (Peak Recognition)

Local Maximum / Minimum

If the mode single event is selected you can also use local maxima or minima to define action conditions.




Fig. 341 FC parameters Local Maximum (single event)

Special settings in parameter window

These settings are not available in the peak recognition mode.

Relative to previous event

If you check this option you have to enter the relative time shift from the previous event. In the Exp. time column this time difference is automatically considered. The time in this row is displayed blue colored indicating that there is no further access for free editing.




Fig. 342 FC parameters Relative to previous event (single event)

#	Exp. time (min)	Mode	Event	Exp. interv. (min)	Rel. time shift (min)	Parameters	Fraction Position	Fraction Name	Fraction Check	Comments
2	1.00	Single Event	Unconditional	0.00	0.00	Timed	next		<input type="checkbox"/>	
3	1.30	Single Event	Local Maximum	0.20	0.30	Max	next		<input type="checkbox"/>	
4	1.30	Single Event	Local Maximum	0.00	0.00	Max	next		<input type="checkbox"/>	

Fig. 343 FC parameters Relative to previous event (single event)

Use as reference for correction next events

Normally any shift in retention times is increasing with the magnitude of the retention. Using the reference option ChromGate adapts the shift for the calculation of the expected (allowed) interval. It is possible to set more than one event as reference. According to the elapsed run time always the last one will be used for calculation.




Fig. 344 FC parameters Use as reference (single event)

After closing the window the parameter indication will be supplemented by **(Ref)**.

Action selection: None

Selecting this action option neither the recycling valve nor the fractionation valve will be switched. The previous action is continued, only a marker is set. This is useful especially in combination with the option *relative to the previous event*. Thus for instance a collection can be started at a defined time after the set conditions are met.

Fraction Position

This field is only informative according to the settings in the parameter window. Either the defined vial position is displayed or next for the first free vial position. Further possibilities are Waste, Peak recycling, Solvent Recycling, and in case of action None 0.00 is displayed. In case of **Collect Slices** always next is shown.

Fraction Check

If this option is checked any run will be interrupted if the defined fractionation conditions will not met. It is not possible, directly to continue such an interrupted run. In case of peak recognition mode this option cannot be used.

Rack configuration path

In any method the fraction collector program must be connected to a rack configuration. The rack configuration file contains information about collected fractions in one or several runs. Opening the fraction collector tab of the setup window first time, no rack configuration path will be displayed.




Fig. 345 Fraction collector setup window

You must either create a new configuration or you can load an already existing one. Trying to leave this window without doing so you will get the error message:



Fig. 346

Press the **Create/Open...** button to open the next window. Here you can select an existing file or type in a new file name which should be descriptive for the configured rack.




Fig. 347 Selecting / creating the rack configuration file

In both cases press the open button. The window will be closed and the path is displayed in setup window (Fig. 348).

Chromatogram window




Fig. 348 Fraction collector setup window

This area is important as well for setting up the method as well for the instrument status during the run. Here only the setup aspects shall be described.

If you have performed at least one run with the method, the last taken chromatogram will be displayed. Independent on this you can load any already stored chromatogram in this area to use it for the fractionation. Select Open Data... and then the desired chromatogram.

This will open the usual chromatogram view:




Fig. 349 Opened chromatogram

Going back to the setup window the chromatogram will there also appear. Moving the cursor inside the graph area it will become a cross hair and additionally a vertical line appears. Depending on the position the actual values Time, Amplitude (level), and Slope will be displayed directly above the graph.




Fig. 350 Opened chromatogram in FRC setup window

Move the cursor to the desired position and double click with pressed Control (<Ctrl>) key. In the fraction collection program a new line is

created and the parameter window for this line is opened. The time, level, and slope values are automatically inserted. Complete the parameter settings as necessary close the parameter window. Repeat this as often as different events has to be programmed. For all lines, added in this way, as a comment "DC operation" will be displayed.




Fig. 351 Programming the fraction table

In Fig. 352 an example for the fraction program of the above chromatogram is shown.

#	Exp. time (min)	Mode	Event	Exp. interv. (min)	Rel. time shift (min)	Parameters	Fraction Position	Fraction Check	Comments
1	0,00	Single Event	Unconditional	0,00	0,00	Timed	Waste	<input type="checkbox"/>	
2	5,51	Single Event	Signal Level/Slope	0,10	0,00	+L and S	1A1	<input type="checkbox"/>	
3	5,85	Single Event	Signal Level/Slope	0,10	0,00	-L and S	Waste	<input type="checkbox"/>	
4	6,30	Single Event	Signal Level/Slope	0,10	0,00	+L and S	1A2	<input type="checkbox"/>	
5	6,70	Single Event	Signal Level/Slope	0,10	0,00	-L and S	Waste	<input type="checkbox"/>	
6	7,98	Single Event	Signal Level/Slope	0,10	0,00	+L and S	1A3	<input type="checkbox"/>	
7	8,38	Single Event	Local Minimum	0,10	0,00	Min	1A4	<input type="checkbox"/>	
8	8,75	Single Event	Signal Level/Slope	0,10	0,00	-L and S	Waste	<input type="checkbox"/>	
9	10,00	Single Event	Unconditional	0,00	0,00	Timed	Waste	<input type="checkbox"/>	

Fig. 352 Example for a fraction table

The colored squares below the graph are to indicate the status in running methods and for fraction tracing. However, those aspects will be described later in the chapter instrument status of a (running) control method (page 137).

Pretreatment

The pretreatment option is only shown, if a recycling valve has been selected in the fraction collector configuration dialog. If you enable the pretreatment check box, the button **Setup...** is accessible. If you click on the button, the pretreatment configuration window will open. Here you can select the position of the recycling valve, if any other device runs a pretreatment. The choices are **Waste** and **Slv. Recycling**.

Fraction Collector Instrument Status

The fraction collector status is shown in two separated tabs, one for monitoring the instrument status and direct control options (Fig. 353) and the other for the rack view (Fig. 354).




Fig. 353 Fraction Collector status tab - Monitor / Direct control




Fig. 354 Fraction Collector status tab - Rack view

In case of the Knauer MultiValve fraction collector both tabs are combined to one (Fig. 355).

Monitor

Status

The actual operation is displayed. It is additionally indicated by the LED color (waste → grey, fraction → green, peak recycling → red, or solvent recycling → blue).

Signal

The actual value of the controlling detector channel is displayed.

Vial

The actual fraction vial number is displayed.

Fr. Time [min]

The collection time for the actual fraction is displayed.

Fr. Volume [ml]

The collected volume of the actual fraction is displayed.

Direct Control

Operation

The operation which shall be performed can be selected from the pull down menu. The choices are waste, peak recycling, or solvent recycling. It will be executed by clicking on the **Set** button.

Collect to

The operation fraction collection can be initiated here. The destination vial is either defined as next, or in case of collect to by its number. The collection starts by clicking on the **Go** button.

Run Time Settings

All operations started in the direct control mode are performed unconditioned as timed events. For consecutive conditioned runs the values for level, slope and the expected interval will be stored according the settings. However, in any case it is recommended to check the automatically stored fractionation program before using it.

Rack View

The displayed rack corresponds to the configured one. All already used vials will appear green highlighted. Moving the cursor across the rack will change the upper descriptive line. There the vial number, the begin and end of collection and the collected volume is displayed (see the arrows in Fig. 354).

Multi valve fraction collector

For this fraction collector the rack view is integrated into the instrument status tab. The shape of the rack view again depends on the configured valves. If different valves (6, 12, or 16 positions) are combined the graphical separation may not be shown. The number of displayed positions also may succeed the actual existing number of collection positions, due to the rectangular arrangement.




Fig. 355 Multi Valve fraction collector status tab with addtl. recycling valve

Differing from other instrument setup windows, that of the fraction collector is also valuable during the run. In the lower part the measured chromatogram will be displayed as it is in progress. Additionally the areas of fractions, peak recycling, and solvent recycling are shown by highlighting in the corresponding colors. Refer to the chapter *Instrument Setup – Fraction Collector* beginning page 235 and Fig. 348 on page 245.




Fig. 356 Fraction collector instrument setup window (Running method)

elp, press F1

Fraction Annotations

The fractions as well as the fraction and vial numbers can be shown in the chromatogram. Open the chromatogram's annotations (right mouse click on the chromatogram – *Annotation*). Select "Fractions" from the drop-down menu and add **Fraction** and/or **Vial Number** using the green arrow key. In the menu "Other" on the bottom of the menu window, please enable the option **FRC Actual** to show the selected annotations. Be sure to have selected the correct Trace for your settings.




Fig. 357 Chromatogram - FRC Annotation

To change the settings for colors, font type etc., open the chromatogram's appearance (right mouse click on the chromatogram – *Appearance*). Select the correct detector channel in *Trace* and select the annotation you want to change the settings for.




Fig. 358 Chromatogram - FRC Appearance

Please note, that for report previews as well as for pdf prints, the fraction color will not be shown correctly (always pink/yellow). If the report is been printed with a “physical” printer, the colors are as selected in the *Appearances* menu.

Stacked Injection

The option Stacked Injection allows to perform injection during a run. As an injection device the autosampler 3950 or the Injection Module can be used. For the configuration and instrument setup of the autosampler 3950, please refer the appropriate chapters.

Injection module configuration

The injection module will be installed with the Knauer FRC AddOn. In the instrument configuration dialog, select the injection module and add it to the Configured modules. Note that it is not allowed to add an autosampler and the injection module.




Fig. 359 Injection module configuration window

The injection module consists of a pump and a Knauer 2-position switching valve. Optional an injection loop can be used. As a pump a pump P 2.1L, P 2.1S, P 4.1S, S 100, S 1050 or 1800 can be selected. If

selected, click on the -button. The pump configuration window will open. Please configure the pump accordingly the pump configuration described in the appropriate chapter. Note that you cannot configure the same pump as for the solvents. However, it is not possible to use e.g. one channel of an LPG solvent pump for the injection module. Also configure the valve as described in the appropriate chapter. If you want to use a sample loop, enable the check box *Use Sample Loop*.

The injection module also allows for using the AZURA ASM 2.1L as an injection module, if a pump and a valve have been installed. The installed modules in the ASM 2.1L have own IP ports, which allows for separate control. The configuration does to allow to select an ASM 2.L or search for the modules by using the -button – the pump and valve must be configured in the Injection Module configuration manually. To use an ASM-pump or -valve in the injection module, select the appropriate pump or valve model in the Injection Module configuration, disable the option “Use S/N to identify the instrument” and enter manually the serial number and IP address of the ASM 2.1L device. Now change the IP port for the modules (pump and valve): If the ASM 2.1L has IP port 10001, the module installed on the left side has IP port 10002, the middle module 10003 and the module installed on the right-handed side the IP port 10004. If the IP port of the ASM 2.1L has been changed in the ASM's own setup, the IP ports of the build-in modules will be changed automatically; they will always have the next possible IP ports. However, before you enter the IP ports of the pump and valve, please check for the IP port of the used ASM 2.1L device.

Injection module instrument setup

If you open the configured instrument and the instrument setup, you will find an own tab for the injection module.




Fig. 360 Injection module instrument setup window

The injection module will be used for the initial and the stacked injections.

Injection program

The settings in this section will only be used for the initial injection, that means, for the first injection at the start of the run.

Injection Volume

Enter the desired injection volume in ml.

Injection Flow

Enter the flow rate of the injection pump during the injection. If no injection is prepared, the pump is stopped; during a run the pump and valve cannot be controlled beside the programmed stacked injection.

Injection Duration

The duration of the injection is calculated from the injection volume and the injection flow. The reduced flow during start and stop of the pump is taken into account with 6 seconds in sum.

Pump Settings

Control Pressure Limits (MPa)

Enter the values for the desired pressure limit. 0 will disable the software control of the pressure. Keep in mind that the pressure for the sample loading may be different from the system pressure during a run if a different flow rate is used or if you use a sample loop. Note that for a pump S 100 without pressure transducer this option is not accessible.

Solvent

Enter the name of the solvent.

Stacked Injections [>>...]

The option Stacked injections allow to define additional injections during a run. Such injections can be programmed either for the autosampler 3950 or the injection module (please refer the Instrument Setup for the autosampler 3950 for details of the stacked Injection setup of the autosampler. Click on the [>>...] – button to open the Stacked Injections table dialog.




Fig. 361 Injection module Stacked injection table

Stacked injections time table

Make a right mouse click into the time table and select *Fill table* from the corresponding menu.




Fig. 362 Injection module Stacked injections table setup

Number of Injections

Enter the number of stacked injections. If you, e.g., enter “5”, two stacked injections will be made; in addition of the initial injection 6 injections at all will be made per run. Please make sure that you have enough samples for the entered number of injection.

Injection Volume

Enter the desired injection volume in ml.

Injection Flow

Enter the flow rate of the injection pump during the injection. If no injection is prepared, the pump is stopped. During a run the injection module's pump can only be controlled using the option “Direct control during a run” (please refer the manual's appropriate chapter Runtime Settings, page 141).

Injection Duration

The duration of the injection is calculated from the injection volume and the injection flow. The reduced flow during start and stop of the pump is taken into account with 6 seconds in sum.

Time between Injections (min)

Enter the time between the injections. The value must be higher than the time calculated for the injection duration.

Click the **Apply**-button to apply the settings into the Stacked injections table.

#	Time (min)	Inj. Volume (ml)	Inj. Flow (ml/min)	Comments
1	0,00	1,00	1,000	Initial injection
2	15,00	5,00	9,500	
3	30,00	5,00	9,500	
4	45,00	5,00	9,500	
5	60,00	5,00	9,500	
6	75,00	5,00	9,500	

Fig. 363 Injection module Stacked injection table

If you make a right mouse click on the table, you can select

- **Fill Down**, to fill the table with the same settings as in the selected line
- **Insert Line**, to add a new line
- **Delete Line**, to delete the selected line
- **Delete Lines**, to delete all selected lines
- **Fill Table...**, will open the Stacked injection setup again. The new table will overwrite the existing one. If you select the command **Fill table** to open the setup dialog for the stacked injection table, the window will show **default values** for all parameters, beside the sample preparation time.

Click on **OK** to close the Stacked Injections table window; the injection module instrument setup window will be shown.

In the **Aux Traces** tab a trace “Stacked Inj. [Inj. Module]” can be enabled. The trace will show the start of the sample preparation and the stacked injection.




Fig. 364 Injection Module Stacked injections auxiliary trace




Fig. 365 Injection Module Stacked injections trace

For every stacked injection an entry for the sample preparation and the injection will be made.

User	Logged	Source	Activity
System	08.10.2010 18:26:32	PC-PM-LAB-TEST	Run Queue - Complete Single Run - C:\EZChrom Elite\Enterprise\Projects\Default\Data\Stacked Injection with InjModule-3_1_001.dat
System	08.10.2010 18:23:23	PC-PM-LAB-TEST	Inj. Module: Stacked injection: trigger marker is detected in 5.04 min after run start.
System	08.10.2010 18:22:38	PC-PM-LAB-TEST	Inj. Module: Stacked injection with inj. volume 2.00 ml, flow 3.000 ml/min (duration 0.72 min) will begin in 5.29 min after run start.
System	08.10.2010 18:20:24	PC-PM-LAB-TEST	Inj. Module: Stacked injection: trigger marker is detected in 3.06 min after run start.
System	08.10.2010 18:19:38	PC-PM-LAB-TEST	Inj. Module: Stacked injection with inj. volume 2.00 ml, flow 3.000 ml/min (duration 0.72 min) will begin in 2.29 min after run start.
System	08.10.2010 18:16:02	PC-PM-LAB-TEST	Run Queue - Start Single Run - C:\EZChrom Elite\Enterprise\Projects\Default\Data\Stacked Injection with InjModule-3_1_001.dat
System	08.10.2010 18:16:01	PC-PM-LAB-TEST	Run Queue - Add Single Run - C:\EZChrom Elite\Enterprise\Projects\Default\Data\Stacked Injection with InjModule-3_1_001.dat

Fig. 366 Injection Module Stacked injections instrument activity log entry

Injection module instrument status

If you open the instrument status, you will find an own tab for the injection module.




Fig. 367 Injection Module status window (initial injection)

Monitor

Run Status

Shows the current run status of the module. If no run is performed, the status is *Hold*. During the initial injection the status *Injecting* with the time of the current injection duration is shown. During the run the run status shows the run time.

Pump Flow

Shows the current flow rate of the injection module's pump.

Pump Pressure

Shows the current pressure of the injection module's pump. This option is only shown if the pump is equipped with a pressure transducer.

Vial Position

The vial position can be either L (Load) or I (Inject).

If no sample loop is used, that means, that the sample will be injected directly on the column, the vial position during the injection (while the injection module's pump is running) is I (Inject). If the injection is made, the valve will switch into L (Load) position.

If a sample loop is used, the valve will switch to L (Load) position if the pump is running and fill the sample loop. For the injection, the valve will switch to I (Inject) position.

Direct Control

Pump Flow

During a run, the pump cannot be controlled using the pump flow option. The buttons Stop and Flow Apply will be inaccessible.

[>>...]

If the option **Direct control during a run** is enabled, after the initial injection the [>>...]-button is accessible. If you click on the button, the known stacked injection setup with the default values will open. This is allow for create a modified stacked injection table. The setup works as described for the instrument setup of the injection module.




Fig. 368 Injection module Stacked injections table setup

If you click the **Apply**-button, the new table will be executed immediately; the program will not check when the last injection was made.

If the **Direct Control during a run**-option **Save changes in time table** is enabled, all program lines for stacked injections, which have not been executed yet during the current run, will be replaced by the stacked injections from the direct control, whereas lines for stacked injections, which have already been executing during the current run, remain untouched.

If the **Direct control during a run**-option **Save changes in time table** is not active, the new table will be executed, but not stored in the instrument setup. However, in this case the method does not include the correct injection program.

SEC Option

Overview

Size Exclusion Chromatography (SEC) is a chromatography technique used to determine the Molecular Mass distribution of large molecules and polymers. SEC (also referred to as GPC) is performed by injecting a sample onto a column comprised usually of rigid polymer gels with known pore-sizes. As the sample migrates through the column, the smaller molecules penetrate the pores of the column and are retained longer than the larger molecules. There is no column adsorption involved in SEC as there is in classical partition chromatography. Since the size of the molecule is related to its Molecular Mass, the elution time (or volume) can therefore be used to approximate the Molecular Mass of the molecule; larger molecules eluting first and smaller molecules eluting last.

Due to the nature of the types of compounds analyzed using SEC (generally high Molecular Mass polymers), it is relatively rare to encounter a chromatogram with sharp, baseline-resolved peaks as is found in partition chromatography. The calculations involved are therefore generally designed to produce Molecular Mass distribution (MWD) or average Molecular Mass numbers rather than discrete Molecular Mass values.

Because the distribution of Molecular Mass is of interest here, results often include not only the average Molecular Mass values, but also a list of area slices for the sample, along with their associated average Molecular Masses. This type of report is called a "Slice Report".

SEC Calibration

ChromGate® SEC supports three methods for SEC Calibration: Narrow Standards calibration, Broad Range calibration (two versions) and Universal calibration.



The time is converted to seconds for calibration and Molecular Mass calculations.

Narrow Standards Calibration

A calibration curve is created by running a standard sample consisting of compounds of known Molecular Masses. The retention times (or volumes) of these "narrow" standards are plotted against the logarithm of the corresponding Molecular Masses (log M).

Broad Range Calibration

Using the Hamielec method: there are two methods available for calculation of Broad Range calibrations.

Broad Range 1 is a linear calibration and requires a SEC column that can provide a linear calibration curve. It uses only average MW values of the polymer standard but assumes a linear approximation of GPC calibration curve. (Mw and Mn values must be known for the polymer standard that can be determined by light scattering and osmotic pressure techniques). This method needs one broad-MWD standard of the same structure as the unknown sample. It should also be possible to use two different Molecular Mass standards with two known MW values in any combination of Mn and Mw.

The Broad Range 2 method is a non-linear calibration utilizing a broad-range standard.

Universal Calibration

Universal calibration allows a narrow range polystyrene standard calibration curve to be used to calibrate for a broader range of polymers. Using the relationship between the hydrodynamic volume, (the product of the Molecular Mass M and the intrinsic viscosity $[\eta]$) of the polymer and its elution time), the narrow range calibration curve can be adjusted to more closely fit the sample being analyzed.

The intrinsic viscosity $[\eta]$ of a polymer is experimentally determined from directly measuring the viscosity of the polymer solution. It is related to the polymer Molecular Mass through the Mark-Houwink equation

$$[\eta]=Km^a$$

where K and a are constants that vary with polymer type, solvent, and temperature. Tables have been published that contain K and a constants for a wide variety of polymers. Thus, using published values for K and a , or the measured viscosity of the polymer, one can determine Molecular Mass information for the sample.

Applications of SEC

SEC is used commonly for quality control of high Molecular Mass polymers. Molecular Mass distribution is indicative of the physical properties of the polymer. Such properties include strength, flexibility, and tackiness, among others. Because the shape of the broad SEC peaks can vary, in many cases simply examining the average Molecular Mass number is not adequate to characterize the polymer, as two polymers can have the same average Molecular Mass but very different weight distribution. It is therefore important to evaluate a range of different Molecular Mass calculations, as well as examine the actual peak shape (overlay comparisons are helpful here) to get a good indication of polymer quality.

The various different Molecular Mass numbers are used to characterize different properties of polymers, as follows.

- M_N** Number-average Molecular Mass, is used in determination of flexibility and tackiness of samples, which is a function of the proportion of low-Molecular Mass material.
- M_W** Weight-average Molecular Mass, is indicative of the strength of the polymer, as it gives an indication of the proportion of high-Molecular Mass material in the sample.
- M_Z** Z-average Molecular Mass, is used to determine brittleness, and indicates the proportion of very high Molecular Mass material in the sample.
- M_V** Viscosity-average Molecular Mass, is used to relate the average viscosity of the sample to its Molecular Mass.
- M_W/M_N** Polydispersity, is an indication of how homogeneous the sample is. A low number indicates more homogeneity, or narrow Molecular Mass distributions. Higher numbers indicate more complex polymers.

Using ChromGate® SEC Software

The ChromGate® Size Exclusion Software is an optional package that enables ChromGate® Chromatography Data System to perform SEC calculations. Once the SEC option is enabled, it is possible for the user to configure any one (or more) of the instruments connected to the system as an SEC instrument. It is also possible to perform normal partition chromatography calculations on some instruments, while performing SEC on other instruments. (Note that SEC and partition

chromatography are distinctly different techniques, and cannot therefore be performed simultaneously on the same instrument.)

In general, performing an SEC analysis involves the following steps.

1. Run SEC calibration standard(s) and use it to set up the calibration information for the SEC samples (Method/SEC Setup).
2. Create a custom report to include the types of SEC reports and plots you want for your analyses.
3. Create a sequence (if appropriate) and run your SEC samples.

Running an SEC Calibration Standard

Before you run your calibration standard(s) you must first set the acquisition parameters so that all the data are correctly sampled and saved on your disk. Do this using **Instrument Setup**. From the instrument window for the chromatograph you intend to use for your SEC run, click on the **Instrument Setup** button or select **Method/Instrument Setup** from the menu.

Instrument Setup

The instrument setup portion of your method tells how you want to acquire the data coming from your chromatograph.



The default acquisition data rate of 1 Hz is adequate for SEC runs. You can adjust the "slice width" for slice reports independently in the SEC setup part of the method

Detector 1...x

For each detector channel (x) configured on the instrument, define the following acquisition information.

Acquisition Channel On

Select this box to turn the acquisition for this channel ON. If this box is not selected, no data will be acquired for this channel.

Sampling

This is the rate at which data will be sampled by the system. You can choose how you want to specify the sampling rate. When you select a sampling rate, a prompt will appear indicating the narrowest peak width for which this sampling rate will be adequate.



It is recommended that you use the Graphical Events Programming to determine the optimum sampling rate for your chromatography.

Frequency

This selection is in Hz (samples per second). This is the selection for most chromatography applications. Click on the down-arrow to get a list of the frequencies available for the configuration of your system.

Period

When you select this type of sampling, you must select the number of seconds (or milliseconds) between data points. Enter the value, then select whether the period is in milliseconds (msec) or Seconds.

Run Time

Run Time determines the length of time data will be sampled

Acquisition Delay

Acquisition Delay is the interval between the start of run (Trigger) and the time when sampling starts for this channel.

Trigger

Select the Trigger tab to select the type of trigger for the instrument.




Fig. 369

Trigger Type

Determines how the data sampling is started.

- None** Sampling starts immediately after clicking on Start. Sequence acquisitions do not pause between runs.
- Manual** Operator has to press **Enter** to start the run. Sequence acquisitions pause for confirmation between runs.
- External** If the data sampling is started from an external trigger, select this option. The type of trigger is designated when the instrument is configured.

When you have completed the acquisition parameters, click the **X** box in the upper right corner of the dialog box to exit the dialog.

Single Run Acquisition

There are two ways you can acquire data using ChromGate® Client/Server. One way is with a sequence (for multiple runs), and the other way is to make a single run. To make a single data acquisition run, you need to specify the method to be used for analysis, and a file name for data storage.

Note: In order to use a method for data acquisition, its instrument setup should have the acquisition channel turned On, and a sampling rate and run time designated.

To make a single run, click the **Single Run** button, or select the **Control/Single Run** command from the menu. The following dialog will appear.

Single Run Acquisition

Run information:

Sample ID: Test Run (010)

Method: EZChrom Elite\Methods\multilevel calibration.met

Data path: c:\enterprise\Projects\Default\Data

Data file: Test Run (010).dat

Print method report

Amount values:

Sample amount: 1

Internal standard amount: 1

Multiplication factor: 1

Calibrate

Calibration level: 1

Clear all calibration

Clear calibration for level

Print calibration report

Clear replicates

Average replicates

Start

Cancel

Help

Description...

Fig. 370

Run Information

This section allows you to specify files for the run.

Sample ID

Enter a Sample ID for the run. This can contain text and numbers, and is saved with the data file.

Method

Enter the name of the method to be used for data acquisition and processing. Include the entire path name if the method is not in your default method directory. You can select the method from a list of methods available on your disk by clicking the File button adjacent to the field.

Data Path

Enter a path name where the data acquired for this run will be stored. Click the File button to select a path from a list of those on your disk.

Data File

Enter a file name to be used to save the data on disk. You can select a path by clicking the File button adjacent to the field. It is not possible to use an existing file name, unless the file exists in a directory whose path contains the term "public". For example, if your data files are saved in a directory entitled "C:\Public\Data", the ChromGate[®] Client/Server files saved in this directory can be overwritten.

Print Method Report

When this box is checked, the method report (or reports) will be printed at the end of the run.

Amount Values

In this section, you can enter values that affect how the concentrations are calculated. If you are making a single data acquisition prior to calibrating your method, simply leave these values at the default level.

Sample Amount

The Sample Amount value is used as a divisor during calculation of concentrations. It is intended to compensate for differences between samples due to weighing and when percentages of the total sample are being calculated rather than the amount detected in an injection.

Internal Standard Amount

For calibration runs, the Internal Standard Amount is taken from the method Peak Table. For unknown runs, enter the amount of the Internal Standard in your unknown sample.

Multiplication Factor

Enter a multiplication factor to be used for this run. All analyzed peaks will be multiplied by this factor.

Calibrate

This box is only used for calibration of partition chromatography samples. It is **not** used in SEC runs.

When you have completed the Single Acquisition Run dialog box, click **Start** to begin the acquisition. The current data will appear in the chromatogram window as it is acquired and stored on disk. At the end of the run, the chromatogram will be analyzed according to the method parameters, and a report generated if specified. If the sample is not analyzed at the end of the acquisition, click the **Analyze** button if you wish to view the results.

Sample Description

When you click the **Description** button, you can enter a text description for the sample that is saved in the data file. The description can be viewed from the Open File dialog, or when the file is open as the current data file using the **Data/Properties** command.

Submit

This button appears when data is currently being acquired using a sequence or single run. The **Submit** button allows you to submit a single acquisition to be run at the completion of the current run. The run is entered at the end of the Run Queue if you are currently running a sequence of runs.

Submit Priority

This button will appear when data is currently being acquired using a sequence. When you click **Submit Priority...** you can submit the single run to be executed immediately after the **current** sequence run record in the Run Queue. After this sample is finished, the sequence will resume.



If the chromatogram is not integrated at the end of the run, or if you were expecting a report and none was printed, check the Method/Properties section of your method to make sure data analysis and reporting are turned on for this method.

Defining SEC Baseline

You can use your stored SEC calibration standard chromatogram to establish the SEC baseline for the analysis.

With the chromatogram displayed in the chromatogram window, click on the SEC baseline button from the toolbar. Then click the mouse on the chromatogram at the beginning and then the end of the desired SEC baseline.




Fig. 371

The range you selected will be displayed in the dialog. Click Add to add this baseline range to your SEC method setup.

Alternatively, this information can be entered manually in the SEC Setup screen. An SEC baseline is required before SEC analysis can be performed.

Defining SEC Result Ranges

You can use your stored SEC calibration standard chromatogram to establish the SEC result range(s) for the analysis.

With the chromatogram displayed in the chromatogram window, click on the SEC Range button on the toolbar. Then click the mouse at the beginning and end of each desired result range. Molecular Mass calculations will be made for each range you select.




Fig. 372

The specified range will be displayed in the dialog box. If the result range falls within a baseline range already entered, that baseline range will be displayed. Click the Add button to enter the designated result range into your SEC setup.

Alternatively, this information can be entered manually in the SEC Setup screen. One or more result ranges are required before an SEC analysis can be performed.

Define SEC Peaks

You can use your stored calibration standard chromatogram(s) to enter peaks and their retention times into your SEC calibration setup.

1. Make sure your calibration standard chromatogram is displayed as the current chromatogram, and that it has been analyzed. (If it has not been analyzed, click on the **Analyze** button.)
2. Select on the **Method/ Define SEC Peaks** command, or click the Define SEC Peaks button on the tool bar of the chromatogram window. You will be prompted to click your mouse before the first

SEC peak and after the last SEC peak you wish to enter in your SEC setup peak list.

Once you have clicked to define the peak region, a box will appear where you select options for the defined peaks.




Fig. 373

Add

Choose this option to add the selected peaks to the current SEC calibration peak table. Existing peaks in the table will be maintained.

Replace

Choose this option to selectively replace peaks in the current SEC calibration table. The tolerance is used to determine if an existing time should be replaced with a new time. If the tolerance is 0 the numbers must be an exact match. The larger the tolerance percentage is the larger the variance can be between the new value and the existing value and still be considered a match.

Average

Choose this option to selectively average new times in the current SEC calibration table with existing ones. The tolerance is used to determine if an existing time should be average with the new time. If the tolerance is 0 the numbers must be an exact match. The larger the tolerance percentage is the larger the variance can be between the new value and the existing value and still be considered a match. Only five trial values will be used in the average. If you try to average a sixth time value, the oldest trial time will be thrown out.

Tolerance

Tolerance % is used to determine “matching” values when using the replace or average options. If the tolerance is 0 the numbers must be an exact match. The larger the tolerance percentage is the larger the variance can be between the new value and the existing value and still be considered a match.

Click OK to enter the defined peaks into the SEC Setup table.

Annotating SEC Chromatograms

When you view an SEC chromatogram in the chromatogram window, there are a number of annotation options available to you that are specific to SEC samples. These are available for you to set when you do a right-hand mouse click in the chromatogram window, and select **Annotations**.




Fig. 374

A dialog box will appear where you can select both SEC-specific and/or normal chromatogram annotations. To view the SEC-specific parameters for annotation, select SEC from the drop down list (peaks/SEC). A list of SEC calculated values will appear that you can choose for annotation on your chromatogram. To select a parameter, double-click on it to move it to the right-hand box.




Fig. 375

Note that the standard **Baseline** is not used in SEC, so it should not be selected. To annotate the SEC Baseline, select the checkbox for this option. Once the annotations are selected, click OK and they will appear on your chromatogram, as shown below.

Note: If you have not entered an SEC Baseline or Result Range into your SEC method yet, it will not appear on the chromatogram.




Fig. 376

You can change the colors and fonts of these annotations from the **Appearance** dialog. (Right-hand mouse click and select **Appearance**.)

Slice vs. Mol. Wt. Trace (Molecular Mass Distribution)

To create a Slice vs. Molecular Mass (Molecular Mass Distribution) trace, open a data file and analyze it. Select the **View/Overlay** command from the menu bar. Click the right-hand mouse button on the chromatogram, and select **Add Trace**. Select the **Current Data** for the Data Source (or another source if desired), and then select a **Mol Wt Dist** for the Trace. Then **Scale** to the largest peak, and click **OK**. (Note: This can also be added using the Rt. Mouse click / **Properties** command.)




Fig. 377

Because the scales are totally different, you should view your current chromatogram trace in a different window than the Mol Wt. Dist trace. To do this, select the **View/Tiled** command from the menu bar. The Mol Wt.

Dist. Trace will appear in a separate window. Notice the axes are Molecular Wt. Vs Slice Area for this trace.




Fig. 378

SEC Setup

SEC Setup involves entering Molecular Mass information about your standard, which will subsequently be used to calculate Molecular Mass numbers for your unknown samples. SEC baselines, result ranges, and peaks can be entered automatically with graphic programming using a stored calibration standard. These values can also be entered manually. You must then complete the SEC Setup by entering the information required for the selected SEC calibration method and calculations.

The information displayed in the SEC Setup screen will vary depending on the Calibration Type selected (Narrow, Broad1, Broad2, or Universal). For details, see the section describing the calibration method of choice.

When you have completed the SEC Setup, you are ready to acquire and process your SEC runs. You can acquire data either manually (single injection), or by batch acquisition (auto injector). See the ChromGate® User's Guide for information on creating a batch sequence.

SEC Setup Narrow Standard Method

When you select the **Narrow** Standard method for Calibration Type, the SEC Setup screen will look similar to the one displayed below.




Fig. 379

Enable Analysis

Click this box to enable SEC calculations for the current channel. Each channel can have its own analysis conditions.

Molecular Mass

For each calibration peak, enter a Molecular Mass that is representative of the calibration standard peak.

Average Time (min)

Enter the time in minutes of the calibration standard peak. If you have used the graphical **Define SEC Peaks** command, the peak times will be entered already for you. Note: the time is converted to seconds for calibration and Molecular Mass calculations.

Time Slice

Select a time slice width (in seconds) for the calculations. The slice width will determine the frequency of Molecular Mass fractions reported. If a Slice Report is generated, it will be based on the slice width you select here. The software will force this number to be an integral of the data collection rate during the analysis.

Baseline Start and Stop Times

The baseline start and stop times (in minutes) can be entered manually or selected graphically using the SEC buttons on the chromatogram toolbar.

Results Range Start and Stop

One or more ranges of the chromatogram must be defined to be used in calculation and reporting of the Molecular Mass distribution. Each range is determined by the results range start and stop values. The times (in minutes) can be entered manually or graphically using a stored chromatogram.

Horizontal Baseline

Select this button if you wish the SEC baseline to be drawn as a horizontal line forward from the Baseline Start Time to the Stop Time. The horizontal baseline stops when it crosses the data trace. A threshold value in microvolt must be specified to ignore any detector baseline noise in the trace that would cause the horizontal projection to stop when it encounters the noise in the trace. If the Horizontal Baseline reaches the end of the detected peak before the programmed Stop Time, the range will stop where the Horizontal Baseline crosses the data trace.

Sample Data

If you have selected Narrow standards calibration, this field will appear empty.

Reference Peak Time (min)

Select this box if you are using a retention time reference peak to adjust for changes in flow rate. All elution times/volumes used in calculations will be offset by the time difference between the value entered here and the time of the reference peak in the sample. Enter the expected retention time for the reference peak, in minutes.

Reference Window (%)

If you are using a reference peak, enter a value for the % window for the reference peak. A peak falling within +/- this % of the expected retention time for the reference peak will be considered the reference peak. (If two peaks fall within the window, the peak closest to the center will be chosen as the retention time reference.)

Fit Type

Select a calibration curve fit type from the drop-down list. The list includes: Point-to-point, Linear, Quadratic, Cubic, 4th Degree Polynomial, and Spline.

Goodness of Fit

The calibration curve for the data in the SEC Setup will be displayed in the calibration curve window. The r^2 value, which represents goodness of fit of the calibration points to the curve selected, will be displayed above the curve. A value of 1.00 indicates perfect fit. Select the **Time** or **Volume** button at the bottom of the calibration curve to display the x-axis of the curve as either retention time in minutes or retention volume. If you select volume, enter the flow rate in ml/min.

Remove calibration curve outliers

To temporarily remove an outlier from the calibration curve, click on the point with your mouse. It will change color and the curve will be recalculated with the point removed.

Time/Volume

Select the mode of viewing for your calibration curve: Time or Volume. If you select Volume, you must enter the flow rate so that the curve is displayed correctly.

Print

Click this button to print your current SEC method information.

Export settings

Click this button to select export options and path.

Zooming on the Calibration curve

You can zoom on the SEC Calibration curve by clicking and dragging with the mouse. Double-click the mouse in the graph area to unzoom the curve view.

SEC Setup Universal Calibration

If you select **Universal** for the Calibration Type, the following screen will appear.




Fig. 380

Enable Analysis

Click this box to enable SEC calculations for the current channel. Each channel can have its own analysis conditions.

Molecular Mass

For each calibration peak, enter a Molecular Mass that is representative of the calibration standard peak.

Average Time (min)

Enter the time in minutes of the calibration standard peak. If you have used the graphical **Define SEC Peaks** command, the peak times will be entered already for you. Note: the time is converted to seconds for calibration and Molecular Mass calculations.

Use Constants

Select the constants to be used in the calibration - either Mark-Houwink or Viscosity from the drop-down selection. Based on your selection, the appropriate columns in your spreadsheet will become active.

Std Mark-Houwink k

If you have selected Mark-Houwink for the Use Constants field, enter the k value for the standard in this field.

Std Mark-Houwink a

If you have selected Mark-Houwink for the Use Constants field, enter the a value for the standard in this field.

Intrinsic Viscosity

If you have selected Intrinsic Viscosity for the Use Constants field, enter the value here. Otherwise this field is not active.

Time Slice

Select a time slice width (in seconds) for the calculations. The slice width will determine the frequency of Molecular Mass fractions reported. If a Slice Report is generated, it will be based on the slice width you select here. The software will force this number to be an integral of the data collection rate during the analysis.

Baseline Start and Stop Times

The baseline start and stop times (in minutes) can be entered manually or selected graphically using the SEC buttons on the chromatogram toolbar.

Results Range Start and Stop

One or more ranges of the chromatogram must be defined to be used in calculation and reporting of the Molecular Mass distribution. Each range is determined by the results range start and stop values. The times (in minutes) can be entered manually or graphically using a stored chromatogram.

Horizontal Baseline

Select this button if you wish the SEC baseline to be drawn as a horizontal line forward from the Baseline Start Time to the Stop Time. The horizontal baseline stops when it crosses the data trace. A threshold value in microvolt must be specified to ignore any detector baseline noise in the trace that would cause the horizontal projection to stop when it encounters the noise in the trace. If the Horizontal Baseline reaches the end of the detected peak before the programmed Stop Time, the range will stop where the Horizontal Baseline crosses the data trace.

Calibration

Select Universal for the calibration type.



If Universal Calibration is selected, the Mark-Houwink constants or intrinsic viscosity of the sample(s) must be entered to display an adjusted calibration curve.

Sample Data

Use this area to input **K** and **a** or measured **intrinsic viscosity** values for your samples. These values in conjunction with the Mark-Houwink constants/intrinsic viscosity specified in the calibration table are used to adjust the calibration curve to better represent the sample being analyzed.



All samples analyzed with a Universal Calibration must have the Mark-Houwink constants or intrinsic viscosity specified here. If you analyze samples with different values, you must create a new method with the values for those samples.

Reference Peak Time (min)

Select this box if you are using a retention time reference peak to adjust for changes in flow rate. All elution times/volumes used in calculations will be offset by the time difference between the value entered here and the time of the reference peak in the sample. Enter the expected retention time for the reference peak, in minutes.

Reference Window (%)

If you are using a reference peak, enter a value for the % window for the reference peak. A peak falling within +/- this % of the expected retention time for the reference peak will be considered the reference peak. (If two peaks fall within the window, the peak closest to the center will be chosen as the retention time reference.)

Fit Type

Select a calibration curve fit type from the drop-down list. The choices include: Point-to-point, Linear, Quadratic, Cubic, 4th Degree Polynomial, and Spline.

Goodness of Fit

The calibration curve for the data in the SEC Setup will be displayed in the calibration curve window. The r^2 value, which represents goodness of fit of the calibration points to the curve selected, will be displayed above the curve. A value of 1.00 indicates perfect fit. Select the **Time** or **Volume** button at the bottom of the calibration curve to display the x-axis of the curve as either retention time in minutes or retention volume. If you select volume, enter the flow rate in ml/min.

Remove calibration curve outliers

To temporarily remove an outlier from the calibration curve, click on the point with your mouse. It will change color and the curve will be recalculated with the point removed.

Time/Volume

Select the mode of viewing for your calibration curve: Time or Volume. If you select Volume, you must enter the flow rate so that the curve is displayed correctly.

Print

Click this button to print your current SEC method information.

Export settings

Click this button to select export options and path.

When you have completed the SEC Setup, you are ready to acquire and process your SEC runs. You can acquire data either manually (single injection), or by batch acquisition (auto injector). See the ChromGate[®] User's Guide for information on creating a batch sequence.

SEC Setup Broad Range1 Calibration

When you select Broad Range1 as the calibration method, the SEC Setup will appear as shown below.




Fig. 381

Enable Analysis

Click this box to enable SEC calculations for the current channel. Each channel can have its own analysis conditions.

Molecular Mass Table

The Broad Range 2 method is a linear calibration method that does not utilize the Molecular Mass Table. Enter values in the Calibration Data area.

Time Slice

Select a time slice width (in seconds) for the calculations. The slice width will determine the frequency of Molecular Mass fractions reported. If a Slice Report is generated, it will be based on the slice width you select here. The software will force this number to be an integral of the data collection rate during the analysis.

Baseline Start and Stop Times

The baseline start and stop times (in minutes) can be entered manually or selected graphically using the SEC buttons on the chromatogram toolbar.

Results Range Start and Stop

One or more ranges of the chromatogram must be defined to be used in calculation and reporting of the Molecular Mass distribution. Each range is determined by the results range start and stop values. The times (in minutes) can be entered manually or graphically using a stored chromatogram.

Horizontal Baseline

Select this button if you wish the SEC baseline to be drawn as a horizontal line forward from the Baseline Start Time to the Stop Time. The horizontal baseline stops when it crosses the data trace. A threshold value in microvolt must be specified to ignore any detector baseline noise in the trace that would cause the horizontal projection to stop when it encounters the noise in the trace. If the Horizontal Baseline reaches the end of the detected peak before the programmed Stop Time, the range will stop where the Horizontal Baseline crosses the data trace.

Calibration Data

Enter values for M_n and M_x in the fields provided or D1 and D2, then click the Calibrate button. The system will calculate D1 and D2 or M_n and M_x values to be used for the calibration.

Reference Peak Time (min)

Select this box if you are using a retention time reference peak to adjust for changes in flow rate. All elution times/volumes used in calculations will be offset by the time difference between the value entered here and the time of the reference peak in the sample. Enter the expected retention time for the reference peak, in minutes.

Reference Window (%)

If you are using a reference peak, enter a value for the % window for the reference peak. A peak falling within +/- this % of the expected retention time for the reference peak will be considered the reference peak. (If two peaks fall within the window, the peak closest to the center will be chosen as the retention time reference.)

Remove calibration curve outliers

To temporarily remove an outlier from the calibration curve, click on the point with your mouse. It will change color and the curve will be recalculated with the point removed.

Time/Volume

Select the mode of viewing for your calibration curve: Time or Volume. If you select Volume, you must enter the flow rate so that the curve is displayed correctly.

Print

Click this button to print your current SEC method information.

Export settings

Click this button to select export options and path.

SEC Setup Broad Range 2 Calibration

If you select **Broad Range 2** for the Calibration Type, the following screen will appear.




Fig. 382

Enable Analysis

Click this box to enable SEC calculations for the current channel. Each channel can have its own analysis conditions.

Molecular Mass

Enter a series of Molecular Masses that will be used in determining the broad range calibration curve.

Average Time (min)

Enter the time in minutes for the Molecular Mass specified. If you have used the graphical **Define SEC Peaks** command, the peak times will be entered already for you. Note: the time is converted to seconds for calibration and Molecular Mass calculations.

Use Constants

Select the constants to be used in the calibration - either Mark-Houwink or Viscosity from the drop-down selection. Based on your selection, the appropriate columns in your spreadsheet will become active.

Std Mark-Houwink k

If you have selected Mark-Houwink for the Use Constants field, enter the k value for the standard in this field.

Std Mark-Houwink a

If you have selected Mark-Houwink for the Use Constants field, enter the a value for the standard in this field.

Intrinsic Viscosity

If you have selected Intrinsic Viscosity for the Use Constants field, enter the value here. Otherwise this field is not active.

Time Slice

Select a time slice width (in seconds) for the calculations. The slice width will determine the frequency of Molecular Mass fractions reported. If a Slice Report is generated, it will be based on the slice width you select here. The software will force this number to be an integral of the data collection rate during the analysis.

Baseline Start and Stop Times

The baseline start and stop times (in minutes) can be entered manually or selected graphically using the SEC buttons on the chromatogram toolbar.

Results Range Start and Stop

One or more ranges of the chromatogram must be defined to be used in calculation and reporting of the Molecular Mass distribution. Each range is determined by the results range start and stop values. The times (in minutes) can be entered manually or graphically using a stored chromatogram.

Horizontal Baseline

Select this button if you wish the SEC baseline to be drawn as a horizontal line forward from the Baseline Start Time to the Stop Time. The horizontal baseline stops when it crosses the data trace. A threshold value in microvolt must be specified to ignore any detector baseline noise in the trace that would cause the horizontal projection to stop when it encounters the noise in the trace. If the Horizontal Baseline reaches the end of the detected peak before the programmed Stop Time, the range will stop where the Horizontal Baseline crosses the data trace.

Calibration Data

Enter values for M_n and M_x in the fields provided or k and a , then click the **Calibrate** button. The system will calculate k and a or M_n and M_x values to be used for the calibration. You must have completed the Molecular Mass table before entering the calibration data.

Reference Peak Time (min)

Select this box if you are using a retention time reference peak to adjust for changes in flow rate. All elution times/volumes used in calculations will be offset by the time difference between the value entered here and the time of the reference peak in the sample. Enter the expected retention time for the reference peak, in minutes.

Reference Window (%)

If you are using a reference peak, enter a value for the % window for the reference peak. A peak falling within +/- this % of the expected retention time for the reference peak will be considered the reference peak. (If two peaks fall within the window, the peak closest to the center will be chosen as the retention time reference.)

Fit Type

Select a calibration curve fit type from the drop-down list. The choices include: Point-to-point, Linear, Quadratic, Cubic, 4th Degree Polynomial, and Spline.

Goodness of Fit

The calibration curve for the data in the SEC Setup will be displayed in the calibration curve window. The r^2 value, which represents goodness of fit of the calibration points to the curve selected, will be displayed above the curve. A value of 1.00 indicates perfect fit. Select the **Time** or **Volume** button at the bottom of the calibration curve to display the x-axis of the curve as either retention time in minutes or retention volume. If you select volume, enter the flow rate in ml/min.

Remove calibration curve outliers

To temporarily remove an outlier from the calibration curve, click on the point with your mouse. It will change color and the curve will be recalculated with the point removed.

Time/Volume

Select the mode of viewing for your calibration curve: Time or Volume. If you select Volume, you must enter the flow rate so that the curve is displayed correctly.

Print

Click this button to print your current SEC method information.

Export settings

Click this button to select export options and path.

When you have completed the SEC Setup, you are ready to acquire and process your SEC runs. You can acquire data either manually (single injection), or by batch acquisition (auto injector). See the ChromGate® User's Guide for information on creating a batch sequence.

Zooming on the Calibration curve

You can zoom on the SEC Calibration curve by clicking and dragging with the mouse. Double-click the mouse in the graph area to unzoom the curve view.

Once you have zoomed in on a chromatogram, you can scroll the chromatogram to the right or left without losing the zoom. This is done by pressing the <CTRL><SHIFT> keys down and moving the mouse until the cursor changes to a "hand" and dragging it to the left or right.

You can also scroll the X- or Y- axis to view features which may be out of the range. To do this, press the <CTRL><SHIFT> keys down while the mouse cursor is outside the graph area, yet near the axis of interest. The cursor will change to an up/down arrow near the Y-axis, or a left/right arrow near the X-axis. Moving the mouse in this mode will scroll the graph up/down or left/right on the axis.

To restore the original view, do a right-hand mouse click in the chromatogram window, followed by the **Full Unzoom** command.

Removing Calibration Outliers

You can temporarily remove an outlier from your calibration curve by clicking on the outlier point on the calibration curve. When you select an outlier point this way, the point will turn red, as will its Molecular Mass

information in the spreadsheet, and the calibration curve will be re-drawn with the point excluded. To restore the point to the curve, click it again. You can also double-click on the row number in the spreadsheet to do the same thing.

SEC Export

When you select the **Export settings** button in the SEC Setup window, the following dialog box appears. Click the **Export Enabled** box to turn export functions **on**. Then select the boxes for the data you wish to have exported.

Enter a path, including file name, where the export file should be placed.




Fig. 383

Click OK when you have completed the dialog. The data will be exported in HTML format (.htm) to make it easy to view with a browser or for transfer to other programs.



This export applies only to the current channel. If you want to export data from a different channel, you must select the channel, then set up the export for that channel.

SEC Custom Reports

SEC reports are created using ChromGate® Custom Report formatting. When you click on the **Custom** reports button, the Custom Report formatter window will appear. For details on using the custom report editor, see the main ChromGate® User's Guide. When using the SEC software, the following options will be available in the custom report editor that allow you to place SEC-specific items in your report.

Adding an SEC Chromatogram

To add an SEC trace to your report, do a right-hand mouse click on the report form and select **Insert Graph/Data Graph**. In the resulting dialog, select **Current Data** if you want the report to display whatever chromatogram is in your current chromatogram window. Select the **Trace** (channel, for example) you wish to display.

The chromatogram will appear on your report. You can adjust the ranges and add annotations by selecting the chromatogram, do a right-hand mouse click, and then select the appropriate command. See the main ChromGate® User's Guide for details on how to change the appearance of your custom report chromatogram.

Adding a Molecular Mass Distribution Trace

To add a Slice vs. Molecular Mass (Molecular Mass Distribution) trace to your custom report, Click the right-hand mouse button on the report, and select **Insert Graph/Data Graph**. Select the **Current Data** for the Data Source (or another source if desired), and then select **Mol Wt Dist** for the **Trace**.




Fig. 384

Click OK to add the trace to your custom report. Note that this trace is slice area vs. Molecular Mass. You can add annotations to this chromatogram in a similar fashion as you do with other chromatograms. It is not possible to overlay this type of trace with a standard chromatogram because the scaling is different.




Fig. 385

Adding a Slice Report

To add a slice report to your custom report, do a right hand mouse click and select **Insert Report/SEC Slice Report**. A dialog will appear.




Fig. 386

Select the trace for which you want to generate the slice report from the drop-down list.

Select the **Specified Range** button if you want to report slices within a certain time range. Enter a **Start Time** and **Stop Time** for the range you wish slices to be reported. Otherwise, select the **All Slices** button.

A table place-holder will appear on your report. To view the current results, do a right-hand mouse click and select the **Show Data at Design Time** option.

Adding an SEC Summary Report

To add a summary report to your custom report, do a right hand mouse click and select **Insert Report/SEC Summary Report**. A dialog will appear where you select the trace for the summary.




Fig. 387

Select the trace for the report, and click **OK**.

A table placeholder will appear on your report. To view the current report data, click the right mouse button in the table, and select **Show Data at Design Time**.

SEC Summary information

Channel A	
Processing Start Time (min)	12.831
Processing Stop Time (min)	25.859
Number of Slices	39
Weight Average Molecular Mass	552812
Number Average Molecular Mass	76133
Z Average Molecular Mass	1250191
Z+1 Average Molecular Mass	1847926
Polydispersity index	7.261
Peak Molecular Mass	183506
Z Average / Weight Average	2.262
Z+1 Average / Weight Average	3.343

Saving an SEC Report Template

If you plan to use the same report for all your SEC methods, you can save the custom report as a template, which can then be opened from

any SEC method. To save the template, select the **File/Template/Save As...** command. Then enter the name you wish to use for your template.

To use a stored template for a new method, open the template using the **File/Template/Open** command.

ChromGate® SEC Equations

$$\begin{aligned} \text{Weight Average Molecular Mass} &= M_w &= \sum (A_i M_i) / \sum A_i \\ \text{Number Average Molecular Mass} &= M_n &= \sum A_i / \sum (A_i / M_i) \\ \text{Z Average Molecular Mass} &= M_z &= \sum (A_i M_i^2) / \sum (A_i M_i) \\ \text{Z+1 Average Molecular Mass} &= M_{z+1} &= \sum A_i M_i^3 / \sum A_i M_i^2 \\ \text{Polydispersity Index} &&= M_w / M_n \\ \text{Viscosity Average Molecular Mass} &= M_v &= (\sum A_i M_i^a / \sum A_i)^{1/a} \end{aligned}$$

Where:

- MW = Weight Average Molecular Mass
- Ai = Area of the slice of interest
- Mi = Average Molecular Mass of the slice of interest
- MN = Number Average Molecular Mass
- MZ = Z Average Molecular Mass
- Mz+1 = Z + 1 Average Molecular Mass
- Mv = Viscosity Average Molecular Mass
- a = Mark-Houwink Empirical Constant



M_v is only calculated when using universal calibrations with Mark-Houwink Constants.

Equations Universal Calibration

Using Universal Calibration, Molecular Masses are determined by first adjusting the calibration points by using the following equations:

$$\text{Intrinsic Viscosity } [\eta] = K m^a$$

Given Sample viscosity:

$$\text{Adjusted Molecular Mass } M = [\eta]_{\text{standard}} / [\eta]_{\text{sample}}$$

Given Sample Mark-Houwink Constants:

$$\text{Adjusted Molecular Mass } M = 10^{\log([\eta]_{\text{standard}} / K_{\text{sample}}) / (a_{\text{sample}} + 1)}$$

Next, they are plotted to create a general calibration curve that is used to lookup sample Molecular Masses.

Equations Broad Standard Method 1

The goal of method 1 is to solve the following Molecular Mass calibration curve equation. To do this the analysis routines solve for the two unknowns D1 and D2.

$$M(t) = D_1 \text{Exp}(D_2 t)$$

To solve this equation a user must supply a molecular number average Mn and Molecular Mass average Mw for a single broad Molecular Mass distribution MWD standard. Also, an F(t) detector response is needed from a mass concentration detector. With these three known values a

single variable search can be used to find D2 using the following equation.

$$M_w / M_n = \left(\int_a^b F(t) \exp(-D_2 t) dt \right) \left(\int_a^b F(t) \exp(D_2 t) dt \right)$$

The usual range for D2 is between 0 and 1 inclusive. Therefore, golden section search is used to find D2s value. For example, with a known value for Mw/Mn the search routine would substitute 0.5 for D2. If the equation is valid, both sides of the equation equal each other, the search ends and D2 = 0.5. If the equation is not valid, then values of 0.25 and 0.75 are substituted for D2 to see if either can satisfy the equation. If neither can solve the equation then the D2 value that is closest to satisfying the equation is the direction the search should continue.

Once D2 is solved D1 can be solved directly using one of the following equations.

$$M_n = \frac{D_1 \int_a^b F(t) dt}{\int_a^b F(t) \exp(D_2 t) dt} \quad M_w = \frac{D_1 \int_a^b F(t) \exp(-D_2 t) dt}{\int_a^b F(t) dt}$$

With D1 and D2 solved a calibration curve can be created over time.

Equations Broad Standard Method 2

The second method is to solve for a non-linear Molecular Mass calibration curve. The curves equation has two unknown Mark-Houwink constants K and A. It is these constants that must be solved before this curve can be created. To implement method 2 a user must supply a molecular number average Mn and Molecular Mass average Mw for a single broad Molecular Mass distribution MWD standard. Also, the user must supply a universal calibration curve $\phi(t)$.

$$M = \left(\frac{1}{K} \right)^{\frac{1}{1+A}} \phi(t)^{\frac{1}{1+A}}$$

Given below is an equation with unknown η and β .

$$M = \left(\frac{1}{K} \right)^{\frac{1}{1+A}} \phi(t)^{\frac{1}{1+A}} = \eta \phi(t)^\beta$$

$$\beta = \frac{1}{1+A}$$

$$\eta = \left(\frac{1}{K} \right)^{\frac{1}{1+A}}$$

Now there are four unknowns K, A, η , and β .

A single variable search using the following equation is used to find β , and that in turn will allow a direct calculation to find η .

$$M_w / M_n = \left(\int_a^b F(t) \phi^\beta(t) dt \right) \left(\int_a^b F(t) \phi^\beta(t) dt \right)$$

With β solved one of the equations below is used to directly calculate η .

$$M_w = \frac{\eta \int_a^b F(t) \phi(t)^\beta dt}{\int_a^b F(t) dt}$$

$$M_n = \frac{\eta \int_a^b F(t) dt}{\int_a^b F(t) \phi(t)^{-\beta} dt}$$

With β and η known, the equations are solved for A and K to create the curve.

$$\beta = \frac{1}{1+A}$$

$$\eta = \left(\frac{1}{K} \right)^{\frac{1}{1+A}}$$

Typical Wiring Schemes

Pump S 1000 / Manager 5000 / Detector S 2500 / Detector S 2300 /
Autosampler 3900




Fig. 388

- | | | |
|---|-----------------------------|--|
| 1 | serial NULL-modem cable | (A0895) |
| 2 | serial Y-cable | (A0755) |
| 3 | connection cable DIN -> DIN | (G0649) |
| 4 | Trigger: | AS 3800 (Control I/O Connector) M 5000 |
| | Pin 4 (N.O.) | -> DIG. |
| | Pin 5 (COMM) | -> IN |

Manager 5000 / User Defined Pump / User Defined Detector / manual Injection Valve




Fig. 389

- | | | |
|---|---------------------------------------|------------------------|
| 1 | serial NULL-modem cable | (A0895) |
| 2 | Detector signal cable (CINCH – CINCH) | (A0868) |
| | (CINCH – wire ends) | (G1021) |
| 3 | Trigger: | Injection Valve K-2600 |
| | Contact 1 | -> START IN |
| | Contact 2 | -> GROUND |
| 4 | Autozero cable | |
| 5 | Pump control cable | |

Pump K-1001 / manual Injection Valve / Detector K-2600




Fig. 390

- 1 serial Y-cable (A0755)
- 2 serial one-to-one cable (A0884)
- 3 Trigger:
 - Injection Valve K-2600
 - Contact 1 -> START IN
 - Contact 2 -> GROUND

**2 Pumps K-1001 / Autosampler Triathlon / Detector K-2501 /
Fluorescence Detector RF-10AxI**




Fig. 391

- | | | |
|---|-------------------------|------------|
| 1 | serial Y-cable | (A0755) |
| 2 | serial one-to-one cable | (A0884) |
| 3 | serial NULL-modem cable | (A0895) |
| 4 | Injection Marker Cable | (0830.811) |

Pump K-501 / Autosampler 3800 / Detector K-2301




Fig. 392

1	serial NULL-modem cable	(A0895)
2	serial NULL-modem cable	(A0895)
3	serial 25 pin male D-cable	(A1004)
4	Trigger:	
	AS 3800 (Control I/O Connector)	K-2501
	Pin 4 (N.O.)	-> START
	Pin 2 (COMM)	-> GROUND

**Pump K-1001 / Autosampler 3800 / Detector K-2600 /
Interface Box / User Defined Detector**




Fig. 393

- | | | |
|---|---|-----------|
| 1 | serial Y-cable | (A0755) |
| 2 | serial Y-cable | (A0755) |
| 3 | serial 25 pin male D-cable | (A1004) |
| 4 | Trigger: AS 3800 (Control I/O Connector) K-2501 | |
| | Pin 4 (N.O.) | -> START |
| | Pin 2 (COMM) | -> GROUND |
| 5 | Detector signal cable (CINCH – CINCH) | (A0868) |
| | (CINCH – wire ends) | (G1021) |

**Pump K-1001 / Autosampler Triathlon / Detector DAD 2800 /
Detector K-2301**




Fig. 394

- | | | |
|---|--------------------------|-------------|
| 1 | serial Y-cable | (A0755) |
| 2 | serial one-to-one cable | (A0884) |
| 3 | serial one-to-one cable | (A1004) |
| 4 | DAD communication cable | (A0887) |
| 5 | Trigger: | |
| | Triathlon (P4 Connector) | K-1001 |
| | Pin 1 (N.O.) | -> START IN |
| | Pin 2 (COMM) | -> GROUND |

Index

3

3D Data	96, 98, 216
3D Properties	197
3D View	196

A

Acquisition	
50/100 Hz option	33
Channel	33, 91, 95, 96, 98, 102, 148
Configuration	53
Delay	92, 108, 109
Source	53
Actualize flow	237
Add. Info	50
Air Segment	113, 117, 122
Alias	<i>see</i> AS 3950
Alltech 650	<i>see</i> Conductivity detector
Analog output	154
Analysis Options	37
AS 3800	
Configuration	70, 111
Instrument Status	161
Setup	111
Wiring Scheme	288, 289
AS 3900	
Configuration	71
Instrument Status	162, 166
Setup	112
Wiring Scheme	284
AS 3950	
Configuration	64
Instrument Status	163
Setup	116
Autosampler	
Configuration	60, 63
Instrument status	156, 161
Large Injection	128
Setup	107, 111
Wiring Scheme	287, 288, 289, 290
Autosampler 3900	<i>see</i> AS 3900, <i>see</i> AS 3900, <i>see</i> AS 3900
Autosampler 3950	<i>see</i> AS 3950
Autozero	58, 96
Auxiliaries	112, 126
Axis Offset	54
Axis Units	54, 81, 82

B

Background Correction	227
Bandwidth	102, 151
Baseline Check	18, 37, 136
Basic Marathon	<i>see</i> AS 3800

C

Calculating Total Purity	227
--------------------------	-----

Calibration Curve	20
Chromatogram View	204
Collector reset	236
Column Oven	68, 71, 135
Configuration	80
Instrument Status	168
Setup	133
Communication port	33
Conductivity detector	
Configuration	48
Instrument status	154
Setup	104
Configuration	
System	14
Configuration	
Acquisition channel	53
Autosampler	60, 63
Column oven	80
Detector	45
Flowmeter	81
Fraction collector	230
Interface	31
KNAUER HPLC System	35
Kontron detector	55, 56
Manager 5000	79
Pumps	38
Switching Valves	78
User defined Detector	53
Virtual Detector	57
Contour View	199
Control Method Options	170

D

DAD	
Analysis	222
Calculations	222
Configuration	48, 50, 56
Custom Report	216
Diagnostics	150
Instrument Status	149
Method Setup	189
Option	95, 102, 185, 189, 230
Setup	100
Utilities	210
Views	196
Wiring Scheme	290
Dark Current	152
Data	
Collection	28
Destination Vial	66, 75
Detection channel	236
Detector	
Configuration	45
Connections	57
Instrument Status	145
Program	102
Setup	91
Synchronization	94
Type	48

- Device ID
 Autosampler 68, 70, 71, 73
- Diagnostics
 DAD 150
 RF-10Axl 153
- Digital Output 89, 94, 140, 148, 149, 150,
 156, 167, 168, 169
- Diode Array *see* DAD
- Direct Control
 AS 3800 162
 AS 3900 162, 163, 166
 Column Oven 169
 DAD 149
 Enabling 170
 Fraction collector 248
 K-2600 148
 Manager 5000 167
 Pumps 141
 RF-10Axl 153, 154, 156
 RI Detector 145
Save Changes 170
 Switching Valves 167
 Triathlon/Endurance 165
 UV Detector 146
- Download
 Tab 172
- E**
- Emergency during a run 170
- Emission wavelength 153
- Endurance
 Configuration 72
 Setup 121
 Status 165
- Ethernet 35, 49, 56
- Excitation wavelength 153
- Export 209
- Extend run time 143, 144
- F**
- Fixed wave length detector 45
- Flow Cell 50, 55, 56
- Flowmeter 135
 Configuration 81
 Instrument Status 169
 Setup 134
- Fluorescence Detector
 Configuration 45
 Diagnostics 153
 Instrument Status 145, 152, 154
 Setup 103
Wiring Scheme 287
- Flush Volume 113, 117, 122
- Fraction collector
 Configuration 230
 Direct control 248
 Instrument status 247
 Program 237, 245
 Setup 235
- Fractionation conditions 238
- G**
- Generic Drivers 183
- GLP 141, 148
- Gradient modify 142, 143
- H**
- Halogen Lamp 50, 51, 56, 57
- Headspace Pressure 113, 118, 122
- HPG System 174, 175
- I**
- Injection Method 113, 117, 122
- Injection Program 124
- Injection Type 111
- Injection Valve*
Wiring Scheme 285, 286
- Installation
 ChromGate 8
 instrument selection 9
 preparative option 11
- Instrument
 Configuration 36
 Control 31
 Name 36
 Setup 84
 Status 29, 137
 Autosampler 156, 161
 Column Oven 168
 DAD 149
 Detector 145
 Endurance 165
 Flowmeter 169
 Fraction collectors 247
 Manager 5000 166
 Pumps 139
 RF-10Axl 152, 154
 Switching Valves 167
 Triathlon 165
 User defined Detector 155
 UV Detectors 145, 147, 148
 Virtual Detector 155
 Type 36
 Wizard 83
- Instrument Configuration 14
- Integration
 Events 19
 Time 151
 Validation of - 178
- Integration 19
- Interface 40, 45, 49, 55, 56, 61, 231
 Configuration 31
- Interface Box
 Configuration 32
 Wiring 57, 289
- Interface Box IF2 33
- Interpolate Spectrum 209
- IP Address 49, 56

J

Jasco FP-920 183

K

K-1001
 Configuration 38
 Setup 85
 Status 139
Wiring Scheme 286, 287, 289, 290

K-120
 Setup 85
 Status 139

K-1800
 Configuration 38
 Setup 85
 Status 139

K-2600
 Direct Control 148
 Instrument Status 148
 Setup 95
Wiring Scheme 286, 289

K-2700/2800
 Interface Card Installation 35
 Setup 100

K-2x00/1
 Configuration 45
 Instrument Status 145
 Setup 92
Wiring Scheme 287, 288, 290

K-501
 Setup 85
 Status 139
Wiring Scheme 288

KNAUER Net 138

Kontron
 detector
 Configuration 55
 Setup 95
 Interface 34
 Pumps *see also* Pumps

Kontron Autosamplers
 Configuration 76
 Setup 128

Kontron DAD 540 100

Kontron detector
 Configuration 56

Kontron Pumps
 Configuration 44

L

Lambda Max/Min Calculations 226
 Lamp OFF 93, 96, 98, 101, 103, 109
 Library Definition Report 218
 Library Search 214, 225
 license dongle 12
 Local Maximum 242
 Login 83
 Loop Volume 65, 68, 70, 71, 73, 77

M

Manager 5000
 Configuration 79
 Instrument Status 166
 Interface Configuration 32
 Setup 130
 Wiring 57, 284

Marker 126, 153

Method 13

Window 16

Wizard 15

Midas *see* AS 3900, *see* AS 3900

Mix Methods 115, 118, 119

Mixed View 203

Multi valve fraction collector 62, 63, 234, 248

Multi-Chromatogram View 207

N

Needle Height 123

Needle Wash 112, 114, 118, 122, 125

Noise Test 136

O

Output profile 131

P

P1x0
 Setup 85
 Status 139

Partial loopfill 70

PDA *see* DAD

Peak

Purity calculations 227

Recognition 237, 241

Recycling 232, 234

Table 22, 221

Performance qualification 175

Report 178

Plate Type 76

Preflush 112

Preparative option

Installation 11

Pressure

Limits 85, 144

Units 43, 45, 51, 61, 62, 63

Pump

Configuration 38

Head 42, 44

Setup 85

Status 139

Type 40, 44

Purity Calculations 191

Q

Quantification 19

R

Rack configuration	243
Ratio Chromatogram	226
Ratio View	210
Recycling valve	232
Reference Signal	152
Report	26
Performance qualification	178
RF-10Axl <i>see</i> Fluorescence detector	
RI Detector	
Configuration	45
Setup	92
status	145
RS232 communications	68, 70, 71, 73
Run Screens	114, 123
Run Status	
Column Oven	168
DAD	149
Flowmeter	169
K-2600	147, 148
Manager 5000	167
Pumps	140
RF-10Axl	152, 154, 155
RI Detector	145
Switching Valves	167
UV Detector	146
Run time	
Extension	143, 144
Run Time	92, 108

S

S 1000	
Setup	85
Status	139
<i>Wiring Scheme</i>	284
S 2x00	
Configuration	45
Instrument Status	145
Setup	92
<i>Wiring Scheme</i>	284
Sample Sequence	23
Sample Tray	65, 68, 70, 71, 73
Sampling Frequency	91, 92
Sampling Rate	91, 108
Scale Factor	122
Scan	96, 98, 101, 151, 153
Scanning speed	154
Search parameters	190
SEC	
Acquisition Delay	260
Acquisition Setup	259
Amount Values	262
Applications of	258
Average Time (min)	268, 270, 275
Calibration	257, 259, 271
Channel Status	259
Custom Reports	278
Define SEC Peaks	263
Enable Analysis	268, 270, 273, 275
Equations	281

Export	278
Fit Type	269, 272, 277
Goodness of Fit	269, 272, 277
Horizontal Line	269, 271, 274, 276
<i>Instrument Setup</i>	259
Intrinsic viscosity	258, 271
Mark-Houwink equation	258
Molecular Mass	268, 270, 275
<i>Narrow Standards</i>	257
Number-average molecular mass	258
Period	259
Polydispersity	258
Reference Peak Time (min)	269, 272
Reference Window (%)	269, 272
Results Range	268, 271, 273, 276
Run Time	259
Sample Data	269, 271
<i>Sample Description</i>	262
Sampling Frequency	259
Sampling Rate	259
Setup	267
Single Run Acquisition	260
Start Time (min)	268, 271, 273, 276
Submit	262
Time Slice	268, 271, 273, 276
Trigger	260
Trigger Type	260
Universal calibration	258
Use Constants	270, 275
Viscosity-average molecular mass	258
Weight-average molecular mass	258
Z-average molecular mass	258
Sequence	13
Files	20
Run	25
Wizard	23
Serial Number	42, 61
Server name	36
Setup	
AS 3900	112
AS 3950	116
Autosampler	107, 111
Baseline Check	136
Column Oven	133
Conductivity detector	104
DAD	100
Flowmeter	134
Fraction collector	235
Instrument	84
Kontron detectors	95
Manager 5000	130
RF-10Axl	103
Switching Valves	131
Trigger	135
Short Guide	13
Shutter	152
Shutter Control	51
Signal Mode	102, 152
Similarity	222
Calculations	226
Threshold	190
View properties	206

Index		295
Single event	237	
Single Level Calibration	21	
Solvent		
Control	173	
Recycling	232, 234	
Selection Valve (SSV)	65, 73	
Spectral Analysis	185, 189, 230	
Spectral Library	212	
Spectral similarity	238, 241	
Spectrum	152	
Spectrum Background Correction	209	
Spectrum Export	228, 229	
Spectrum View	208	
SSV	65	
Step Mode	130	
Switching Valves		
Configuration	78	
Instrument Status	167	
Setup	131	
Syringe Speed	113, 117, 122	
Syringe Volume	65, 68, 71, 73, 77	
System		
Options	14	
Shutting Down	30	
Status	137	
T		
Temperature control	104	
Threshold	19	
Time Constant	91, 108	
Time Program	102	
Timed Events	114, 123	
Traces	135	
Transport Vial	74	
Tray Configuration		
Endurance	75	
Triathlon	74	
Tray Cooling	65, 68, 70, 71, 73	
Triathlon		
Configuration	72	
Setup	121	
Status	165	
<i>Wiring Scheme</i>	287, 290	
Trigger	58, 135	
Tubing parameters	231	
Tubing Volume	65, 68, 71, 73	
U		
USB license dongle	12	
User defined Detector		
Configuration	53	
Setup	104	
Status	155	
<i>Wiring Scheme</i>	285, 289	
User defined Injection Program	124	
UV Detector		
Configuration	45	
Instrument Status	145, 147	
Setup	92	
V		
Validation of Integration	178	
Vial Position	66, 69, 72, 75	
Virtual Detector		
Configuration	57	
Formula	107	
Setup	105	
Status	155	
Virtual fraction collector	233	
W		
Wake-Up option	149	
Wash Times	113, 117	
Wash Volume	123	
Wavelength table	94	
Width	19	
Wiring		
Interface box:	57, 289	
Manager 5000	57	
Wiring Schemes	284	
Y		
Y-Axis Units	50, 56	

© Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
All rights reserved.
The information in this document is subject to
change without prior notice. Translation of the
original German edition of this manual.
2012-10-02
Printed in Germany on environmentally friendly
paper from sustainable forests.

® ChromGate
are registered trademark of
Dr. Ing. Herbert Knauer GmbH

▶ See up-to-date manuals online:
www.knauer.net/downloads

www.knauer.net

HPLC · SMB · Osmometry

Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
Hegauer Weg 38
14163 Berlin, Germany

Phone: +49 30 809727-0
Telefax: +49 30 8015010
E-Mail: info@knauer.net
Internet: www.knauer.net



© KNAUER 2012 V7057-8/0.05/10.12/Koe