

Application Note

► Determination of mannose and mannoooligosaccharides with an improved RI detector



AZURA


Category	Food
Matrix	H ₂ O _{dd}
Method	HPLC
Keywords	Oligosaccharides, functional food, Eurokat Na ⁺ , AZURA RID 2.1L, Konjac root
Analytes	Mannose, mannobiose, mannotriose, mannotetraose, mannopentaose, mannohexaose
ID	VFD0148N

Summary

This application note describes the determination of saccharides with the new AZURA RID (refractive index detector) 2.1L. A mixture of saccharides was separated on a KNAUER Eurokat Na⁺ column. The new AZURA RID is approximately three times more sensitive than its predecessor Smartline RI detector S2300, achieving now a limit of quantification (LOQ) between 10 and 20 µg/ml for the tested saccharides.

Introduction

Glucomannans are branched polysaccharides of mannose and glucose. They are found in many plants as energy reserve. The Konjac root (*Amorphophallus konjac*) contains high amounts of glucomannans. Flour from the Konjac root showed nutritional benefits such as cholesterol lowering. Due to its stimulation of gut-friendly bacteria it is called prebiotic¹. A prebiotic is a selectively fermented ingredient that allows specific changes, both in the composition and/or activity in the gastrointestinal microflora that confers benefits upon host well-being and health². Prebiotics are a specific class of functional food. Oligosaccharides have the advantage that they survive the entire gastrointestinal transit and thus reach gastrointestinal bacteria³. It was shown that glucomannooligosaccharides, obtained by partial hydrolysis of glucomannans, showed a better prebiotic function than pure glucomannan. The prebiotic effect for glucomannooligosaccharides with a degree of polymerization (DP) of 5 was better than for those with a DP of 10⁴. Partial enzymatic hydrolysis is therefore of great interest⁵. Liquid chromatography is one possibility to analyze samples after lysis and to estimate the composition and degree of polymerization. The new AZURA RID 2.1L was tested for its LOQ (limit of quantification). Further the KNAUER Eurokat Na⁺ column was tested for its ability to separate saccharides with a degree of polymerization higher than 2. A mixture of mannose and mannoooligosaccharides was used as a test substrate, since these oligosaccharides can be used as reference for the products of the hydrolysis.

**Fig. 1**

Chemical structures of the tested saccharides (structures kindly provided by company Megazyme)

Experimental: Sample preparation

Standard solutions of 12 mg/ml of the six saccharides were prepared in $\text{H}_2\text{O}_{\text{dd}}$. A mixture of all six saccharides with a concentration of 2 mg/ml each was prepared and passed through a 45 μm filter. From this standard mixture eight 1:2 dilution steps were prepared for the calibration points. Final concentrations were: 2 mg/ml; 1 mg/ml; 0.5 mg/ml; 0.25 mg/ml; 0.125 mg/ml; 0.0625 mg/ml; 0.03125 mg/ml; 0.015625 mg/ml; 0.0078125 mg/ml.

Saccharides were purchased from company Megazyme

Method parameters

Column	Eurokat Na^+ , 300x8 mm, 10 μm
Eluent A	$\text{H}_2\text{O}_{\text{dd}}$
Flow rate	0.5 ml/min
Injection volume	10 μl
Column temperature	75°C
System pressure	50 bar
Detection	RI
Run time	16 min

The LOQ of the new KNAUER AZURA RID 2.1L was determined and compared with the KNAUER Smartline S2300 RI detector. A mixture of saccharides was used as test substrate and separated using an analytical KNAUER Eurokat Na⁺ column. A series of dilutions starting from 2 mg/ml to 0.0078125 mg/ml were prepared. Calibration curves were calculated for each saccharide of the mixture.

Results

The mixture of mannose, mannobiose, mannotriose, mannotetraose, mannopentaose and mannohexaose was separated using the Eurokat Na⁺ column (Fig.2). Mannose, mannobiose and mannotriose were baseline separated. Mannotetraose was nearly baseline separated while the peaks for mannopentaose and mannohexaose were not baseline separated (Fig. 2). The chromatogram of this mixture using the Smartline RI detector S2300 showed a similar separation (data not shown).





Fig. 2

Separation of the saccharide mixture, each 2 mg/ml, Eurokat Na⁺ 300 x 8mm

**Fig. 3**

Calibration graphs for all six tested saccharides with indication of linearity r^2

The LOQ for the six tested saccharides was determined and compared using the AZURA RID 2.1L and Smartline RI detector S2300. The noise was determined and the LOQ calculated according to the peak height using the calibration graphs (Fig.3).

The LOQs for all tested saccharides determined with the AZURA RID 2.1L were significantly lower than those detected with the Smartline RI detector S2300 (Fig 4). The LOQ range for the Smartline RI detector S2300 was between 65 $\mu\text{g/ml}$ and 35 $\mu\text{g/ml}$, while the LOQ range for the new AZURA RID 2.1L was between 20 $\mu\text{g/ml}$ and 8 $\mu\text{g/ml}$. The LOQ of the AZURA RID 2.1L was thus determined to be 20 % to 38 % lower than the LOQ of the Smartline RI detector S2300, which is a significantly higher sensitivity (Fig. 5).

The total LOQ values of all measured saccharides are shown in Tab. 1. The relative standard deviation (RSD) values were calculated for the peak height of four replicates at 15.6 $\mu\text{g/ml}$ of each saccharide. The % RSD values calculated for the tested saccharides measured with the RID 2.1L were between 0.52 % and 3.33 % and for the Smartline S2300 between 2.52 % and 7.61 %. These results underline the increased sensitivity of the AZURA RID 2.1L and more accurate measurement at low sample concentration.




Fig. 4

Comparison of AZURA RID 2.1L (black bars) and Smartline RI detector S2300 (white bars) for their LOQ values for the six tested saccharides




Fig. 5

LOQ ratio of Smartline RI detector S2300 and AZURA RID 2.1L values for all six saccharides

Tab. 1

Comparison of total calculated LOQs and %RSD (peak height) of 1:128 dilution step (15.6µg/ml) of all saccharides measured with AZURA RID 2.1L and Smartline RI S2300

Saccharide	RID 2.1L		S2300	
	LOQ in µg/ml	%RSD (height) for 15.6 µg/ml	LOQ in µg/ml	%RSD (height) for 15.6 µg/ml
Mannose	20.35	1.67	67.72	6.74
Mannobiose	15.82	3.33	50.39	7.61
Mannotriose	14.09	2.11	43.49	6.38
Mannotetraose	8.77	0.58	40.02	5.30
Mannopentaose	12.79	1.07	36.93	5.66
Mannohexaose	14.15	0.52	36.29	2.52

Conclusion

A mixture of six saccharides was separated using a KNAUER Eurokat Na⁺ column. Mannose, mannobiose, mannotriose and mannotetraose were baseline separated. The new KNAUER AZURA RID 2.1L was used for detection and had an LOQ depending on the saccharide between 10 and 20 µg/ml. In comparison with the Smartline RI detector S2300 the AZURA RID 2.1L is up to 4.5 times more sensitive.

The AZURA RID 2.1L in combination with the Eurokat Na⁺ column is an excellent tool for the separation and identification of oligosaccharides with a degree of polymerization of 2 to 5. Concentrations can be analyzed for certain saccharides down to 10 µg/ml.

Eurokat Na⁺ is applicable for the separation of sugar oligomers

Physical properties of recommended column



Stationary phase	Sulfonated cross-linked styrene-divinylbenzene co-polymer in sodium form
Particle size	10 µm
Form	spherical
Cross linkage %	6
Dimensions	300x8 mm
Order number	30GX210EKN

Recommended instrumentation



The analysis of saccharides requires an isocratic HPLC system equipped with degasser, autosampler, column oven and a refractive index detector. Other configurations are also available. Please contact KNAUER to configure a system that fits your needs.

Description	Order No.
AZURA P 6.1L Quaternary HPLC pump, stainless steel 10ml	APH34EA
AZURA RID 2.1L	ADD31
AZURA CT2.1 Column Thermostat	A05852
Autosampler 3950	A50070
AZURA Eluent Tray E 2.1L	AZC00
OpenLab CDS EZChrom Edition	A2600-1

References

1. Al-Ghazzewi F. H., Khanna S., Tester R.F., Piggot j. The potential use of hydrolysed konjac glucomannan as a prebiotic. *J Sci Food Agric* 87:1758-1766 (2007)
2. Roberfroid MB. Prebiotics: The Concept Revisited. *J Nutr.* 137 (3 Suppl 2): 830S–75S. PMID 17311983 (March 2007)
3. Albrecht S., et al Enzymatic production and characterization of Konjac glucomannan oligosaccharides *J. Agric. Food Chem.* 59: 12658-12666 (2011)
4. Chen J., Liu D., Shi B., Wang H., Cheng Y., Zhang W. Optimization of hydrolysis conditions for the production of glucomanno-oligosaccharides from konjac using β-mannase by response surface methodology *Carbohydrate Polymers* 93: 81-88 (2013)
5. van Zyl W. H. Rose S. H., Trollope K., Görgens J. F. Fungal β-mannanases: Mannan hydrolysis, heterologous production and biotechnological applications *Process Biochemistry* 45: 1203-1213 (2010)

Authors

Dr. Yannick Krauke, LC LAB Application specialist/Project Manager

Hagen Schlicke, LC LAB Application specialist/Project Manager

Contact information

KNAUER
Wissenschaftliche Geräte GmbH
Hegauer Weg 38
14163 Berlin, Germany

Tel: +49 30 809727-0
Fax: +49 30 8015010
E-Mail: info@knauer.net
Internet: www.knauer.net

This project has received funding from the European Union's Seventh Framework Program for research, technological development and demonstration under grant agreement no FP7-KBBE-2013-7-613802; <http://www.valorplus.eu/>

