

ALTERNATIVE XYLITOL EXTRACTION VIA HPLC PURIFICATION FROM FERMENTED BIOMASS

Yannick Krauke, Hagen Schlicke, Matthias Luebbert; Krauke@knauer.net
KNAUER Wissenschaftliche Geräte GmbH, Hegauer Weg 38, 14163 Berlin; www.knauer.net

SUMMARY

A latest approach in bioethanol generation is the usage of yeast and bacteria that uses C5 sugars for fermentation and the valorization of bio refinery by products. Here it is shown that a hemicellulose-like fermentation mash has a high content of the artificial sweetener xylitol and that its purification by HPLC can be accomplished using polymer based Eurokat columns. The product is soluble in water and can easily be used for further applications.

INTRODUCTION

The second generation of bio refinery uses biomass with lower contents of C6 glucose and higher contents of C5 sugars. Besides ethanol generation its goal is the full usage of biomass by valorizing by products. Fermentation of C5 sugars with microorganisms result in mash that could be used for further applications. Polymer based Eurokat columns were tested for their ability to separate fermentation mash and among them the Eurokat Ca column had the best separation profile. Analysis of the mash revealed high contents of xylitol. Purification of highly pure xylitol was established.

RESULTS

The fermentation mash was analyzed on different columns (Eurokat Na, H, and Ca) to determine the optimal stationary phase. The Eurokat Ca column showed the best separation profile for xylitol **Fig.1** even though it has the longest run with about 28 min compared to Eurokat Na with 18 min and Eurokat H with 12 min (not shown). A more detailed analysis of the fermentation mash identified five components: xylose, arabinose, glycerol, mannitol and xylitol **Fig.1**. Xylitol had the highest concentration with 80 mg/mL in the sample, followed

by glycerol with 20 mg/mL. The other three components had concentrations of 7-8 mg/mL **Fig.1**. The baseline separation of xylitol indicated promising batch purification. Overload studies with a semi-preparative Eurokat Ca column were performed. This column has a three times higher volume (50 mL) than the analytical column (15 mL) and larger particle size (25-56 μm) enabling higher sample loading and faster flow rates with lower back pressure. The collected fraction of xylitol **Fig.2** had a purity of 99 %, measured with RI **Fig.3**.

Sample analysis




Fig. 1 Chromatogram of 1:10 dilution of fermentation mash 10 μL injection on Eurokat Ca - 1 xylose (8.2 mg/mL), 2 arabinose (8.3 mg/mL), 3 glycerol (21.0 mg/mL), 4 mannitol (7.0 mg/mL), 5 xylitol (80.6 mg/mL)

Batch purification




Fig. 2 Fractionation of xylitol from 1000 μL injection; yellow fraction area (9.5 mL)

Fraction analysis




Fig. 3 Comparison of sample and fraction chromatograms; blue = sample, red = fraction from batch purification

MATERIALS AND METHOD

Sample preparation

Vogelbusch Biocommodities GmbH provided the fermentation mash that resulted from fermentation with yeasts of hemicellulose-like hydrolysate with high xylose content. The sample was filtered through 0.45 μm filter after centrifugation. A 1:10 dilution was prepared and analyzed. For calibration a mixture of xylose, arabinose, glycerol, mannitol and xylitol was prepared and six dilution steps from 15 mg/mL to 0.3 mg/mL prepared.

Method parameters

Analytical runs were performed with KNAUER analytical Eurokat columns (300 \times 8 mm) with integrated pre-columns (30 \times 8 mm) with 10 μm particles at 75 $^{\circ}\text{C}$ running at flow rates of 0.5 mL/min using $\text{H}_2\text{O}_{\text{dd}}$ as eluent. The KNAUER AZURA analytical HPLC system comprising of the AZURA P 6.1L HPG 10 mL pump, 3950 autosampler, AZURA DAD 2.1L diode array detector with high sensitivity KNAUER LightGuide cartridge flow cell, AZURA RID 2.1L refractive index detector, AZURA CT 2.1 column thermostat controlled by the OpenLAB[®] EZChrom Edition software was used. The purification of xylitol was performed with KNAUER Eurokat Ca columns (250 \times 16 mm) with 25-56 μm particles at 75 $^{\circ}\text{C}$ running at flow rates of 2.5 mL/min using $\text{H}_2\text{O}_{\text{dd}}$ as eluent. The KNAUER AZURA Preparative HPLC system comprising of the AZURA P 6.1L HPG 50 mL pump, 3950 autosampler (preparative version), AZURA RID 2.1L refractive index detector, AZURA CT 2.1 column thermostat controlled by the OpenLAB[®] EZChrom Edition software was used. The refractive index detector's Extended Dynamic Range (EDR) feature was used for preparative experiments.

CONCLUSION

The Eurokat Ca column was found to be the best column for analysis of fermentation mash among tested Eurokat columns. The used fermentation mash has a high content of xylitol (80 mg/ml). A semi-preparative batch purification of the xylitol resulted in high recovery (95 %) of xylitol with a purity of 99 %. Upscaling of the batch process or application of SMB (simulated moving bed) chromatography would be promising for xylitol production from fermentation mash.

Acknowledgement: This project has received funding from the European Union's Seventh Framework Program for research, technological development and demonstration under grant agreement no FP7-KBBE-2013-7-613802



ALTERNATIVE XYLITOL EXTRACTION VIA HPLC PURIFICATION FROM FERMENTED BIOMASS

ADDITIONAL RESULTS

The fermentation mash was separated on Eurokat Na and Eurokat H columns **Fig. A1**. On the Eurokat Na column only three peaks were detected. The the Eurokat H column xylitol was also not baseline separated from the other substances. For the overload studies 50 μL to 1500 μL of the 1:10 dilution of the fermentation mash were separated on the Eurokat Ca column. Overlays of all the chromatograms show a shift in the early eluting phase (10-14 min) due to volume overload but less for xylitol **Fig. A2**.

Column screening




Fig. A1 Chromatograms of 1:10 dilution of fermentation mash; *left* Eurokat Na; *right* Eurokat H; 1 xylitol; 10 μL injection

Overload experiments




Fig. A2 Overload studies with semi-preparative Eurokat Ca and fermentation mash; 50 μL , 100 μL , 200 μL , 500 μL , 1000 μL , 1500 μL injection

ADDITIONAL MATERIALS AND METHODS

Tab. A1 Comparison of properties and method parameters of applied Eurokat columns

Column	Column Dimensions	Particle (μm)	Eluent	Flow rate (mL/min)	Injection Volume (μL)	Temperature ($^{\circ}\text{C}$)	Column Volume (mL)
Eurokat H	300 \times 8 mm + 30 \times 8 mm	10	H ₂ O/5 mM H ₂ SO ₄	0.5	20	75	15
Eurokat Ca, Na	300 \times 8 mm + 30 \times 8 mm	10	H ₂ O	0.6	20	60	15
Eurokat Ca	250 \times 16 mm	25-56	H ₂ O	2.5	10000	75	50

AZURA Analytical system

Instrument	Description	Article No.
Pump	AZURA P 6.1L, HPG, 10mL, SSt	APH35EA
Autosampler	3950 analytical version	A50070
Detector 1	AZURA DAD 2.1L	ADC01
Flow Cell	High Sensitivity LightGuide 50 mm, 6 μL	AMD59
Detector 2	AZURA RID 2.1L	ADD31
Thermostat	AZURA CT 2.1	A05852
Software	OpenLAB [®] CDS EZChrom Edition	A2600-1

AZURA Preparative system

Instrument	Description	Article No.
Pump	AZURA P 6.1L, HPG; 50 ml, SSt	APH38FA
Autosampler	3950 preparative version	A50054-1
Detector	AZURA RID 2.1L	ADD31
Thermostat	AZURA CT 2.1	A05852
Fraction collector	Foxy R1	A59100
Software	OpenLAB [®] CDS EZChrom Edition	A2600-1

