

# Easy and fast isolation of rosmarinic acid from lemon balm with mass-directed purification

Jan Wendrich, Yannick Krauke, Lilit Avagyan, Kate Monks; applications@knauer.net  
KNAUER Wissenschaftliche Geräte GmbH, Hegauer Weg 38, 14163 Berlin; www.knauer.net

## SUMMARY


Rosmarinic acid is a natural product widely spread over different plant families. Preparative reversed-phase chromatography was used for the mass-directed purification of rosmarinic acid from a lemon balm extract. The AZURA® Prep System together with the 4000 MiD mass spectrometer was showed to be well suited for this application.

## INTRODUCTION


The ubiquitous natural product rosmarinic acid shows antiviral, antimicrobial and anti-inflammatory characteristics. It is used in different kinds of medicinal products for example in ointments for sports injuries. Leaves of lemon balm contain a high concentration of rosmarinic acid and are therefore an interesting source for the isolation of this compound. Here, we present an effective and time-saving method for the isolation of rosmarinic acid from a lemon balm extract based on the technique of mass-directed purification.

## RESULTS


A method for the isolation of rosmarinic acid from a lemon balm extract was developed on an analytical scale using an AZURA Analytical System and a Eurospher II C18 column (Fig. 1). The UV spectra from the analysis showed the presence of many compounds with the structural motif of a phenyl acrylic acids. For a time-saving isolation of the target compound, the developed method was then transferred directly to the AZURA Preparative System with the ability to fractionate via molecular mass. One fraction with a compound of the desired mass ( $m/z$  359.2; [M-H]) was collected (Fig. 2 & Fig. 3). The following analysis of this fraction with the AZURA Analytical System showed that it was possible to isolate the target compound rosmarinic acid with the technique of mass-directed purification in a purity of >95 % (Fig. 4).




**Fig. 1** Analytical chromatogram of the crude lemon balm extract at 280 nm; gradient separation 20 %-100 % acetonitrile, \*rosmarinic acid peak



**Fig. 3** UV-Chromatogram of a purification run for the crude lemon balm extract at 280 nm, \*rosmarinic acid peak



**Fig. 2** SIM (single ion monitoring) chromatogram of a purification run for the target mass of  $m/z$  359.2



**Fig. 4** Analytical chromatogram of the fraction containing rosmarinic acid at 280 nm, isocratic separation 50/50 water/acetonitrile


## MATERIALS AND METHODS

AZURA Analytical HPLC System was used for the method development consisting of a low-pressure gradient AZURA P6.1L pump, an AZURA AS 6.1L autosampler, an AZURA DAD 2.1L diode array detector equipped with a 10 mm PressureProof flow cell and an Eurosphere II 100-5 C18 150 x 4.6 mm column. The gradient method was run for 20 min at a flow rate of 1 mL/min starting with 80/20 % water/acetonitrile increasing to 100 % acetonitrile over 20 min. Both eluents contained 0.1 % formic acid as an additive. The wavelength of the detector was set to 280 nm at a data rate of 20 Hz. 10  $\mu$ L of the sample was injected.

AZURA Preparative HPLC System was used for the mass-directed purification of rosmarinic acid. The system consisted of an AZURA P 2.1L pump equipped with a 250 mL pump head and a three channel low pressure gradient (LPG) ternary module, a manual injection valve (1/8", 6 port 2 kanal) equipped with a 5 mL sample loop, an AZURA DAD 6.1L diode array detector equipped with a 3 mm PressureProof flow cell, a 4000 MiD mass spectrometer with the MiDas sampling unit, a AZURA V 2.1S equipped with a 6 port multi position valve for fractionation and an Eurospher II 100-10 C18 250 x 30 mm column. The gradient method run for 67 min at a flow rate of 21.3 mL/min under the same conditions as the analytical method described above. The wavelength of the DAD was set to 280 nm at a data rate of 10 Hz, while the mass selective detector was set to SIM mode monitoring the relevant mass of  $m/z$  359.2. The data trace of the mass selective detector was used for fractionation via the multi-position valve. 5 mL of the crude extract obtained under sonication from dried leaf material with 30 % isopropanol was injected.

## CONCLUSION


Rosmarinic acid was the main metabolite of the extracted lemon balm material. This target molecule was isolated in a short time with an AZURA Preparative HPLC system using the technique of mass-directed purification. By this, the number of fractions was reduced to one leading to a significant decrease of past analysis time.



Additional information

# Easy and fast isolation of rosmarinic acid from lemon balm with mass-directed purification

## ADDITIONAL RESULTS



**Fig. A1** Mass spectrum of rosmarinic acid ([M-H]).

## ADDITIONAL MATERIALS AND METHODS

**Tab. A1** Method parameters (preparative)

Eluent A	Water + 0.1 % formic acid		
Eluent B	Acetonitrile + 0.1 % formic acid		
Gradient	Time [min]	% A	% B
	0	80	20
	67	0	100
Flow rate	21.3 mL/min	System pressure	100 bar
Run temperature	RT	Run time	67 min
Injection volume	5 mL	Injection mode	Full loop
Detection wavelength	280 nm	Data rate	10 Hz
		Time constant	0.1 s

**Tab. A2** Method parameters (mass spectrometer analysis)

Scan mode	Interleave (Scan/SIM)
Mass range	100-400 m/z
Scan rate	1 Hz
Step	0.2
SIM	359.2 m/z
Ion mode	Negative
Gas flow	2.5 L/min

**Tab. A3** System configuration & data (analytical system)

Instrument	Description	Article No.
Pump	AZURA P 6.1L	<a href="#">APH34GA</a>
Autosampler	AZURA AS 6.1L	<a href="#">AAA00AA</a>
Detector	AZURA DAD 2.1L	<a href="#">ADC01</a>
Flow cell	PressureProof Flow cell 10 mm 10 µl	<a href="#">AMC38</a>
Column	Eurospher II 100-5 C18 with precolumn, Vertex Plus Column 150 x 4.6 mm	<a href="#">15VE181E2J</a>
Software	ClarityChrom	<a href="#">A1670</a>
	ClarityChrom 7.4.2 - PDA extension	<a href="#">A1676</a>

**Tab. A4** System configuration & data (preparative system)

Instrument	Description	Article No.
Pump	AZURA P 2.1L 250 mL AZURA LPG module for Pump P 2.1L	<a href="#">APE20LA</a> <a href="#">AZZ00AB</a>
Injection	AZURA V 2.1S Valve 6 Port 2 Position	<a href="#">A1359</a>
Sample loop	5 mL Sample loop	<a href="#">A0586-2</a>
Detector	AZURA DAD 6.1L	<a href="#">ADC11</a>
Flow cell	PressureProof flow cell 3 mm 2 µL	<a href="#">AMB18</a>
Mass spectrometer	4000 MiD with MiDas	<a href="#">A66900</a>
Fractionation	AZURA V 2.1S Valve 6 Port Multiposition	<a href="#">AWA10BC</a>
	ClarityChrom	<a href="#">A1670</a>
Software	ClarityChom PDA Extension	<a href="#">A1676</a>
	ClarityChom MS Extension	<a href="#">A1679</a>
	ClarityChrom FRC control module	<a href="#">A1682</a>
Column	VertexPlus AX Column 250x30 mm Eurospher II 100-10 C18	<a href="#">25QE181E2N</a>