

Alternaria alternata - determination of main metabolites

Juliane Böttcher, Marina Müller¹, Kate Monks; applications@knauer.net
KNAUER Wissenschaftliche Geräte GmbH, Hegauer Weg 38, 14163 Berlin; www.knauer.net
¹Leibniz-Zentrum für Agrarlandschaftsforschung e.V., Müncheberg, Germany



SUMMARY

Alternaria toxins represent a possible health-endangering group of mycotoxins produced mainly by the Alternaria species. These are a widespread group of fungi contaminating mainly fruits and vegetables, but also other crop plants, during growth as well as storage. The most important mycotoxin-producing species is Alternaria alternata which occurs mainly on cereals and seeds [1].

INTRODUCTION

Even though Alternaria toxins are normally associated with fruits and vegetables that are visibly infected by Alternaria species, they have also been found in cereals, such as wheat, rye, sorghum, rice, and even tobacco. Alternaria toxins have been shown to exhibit both acute and chronic effects and therefore represent a threat to animal and human health. The most studied mycotoxin in the group of toxins produced by the species Alternaria is tenuazonic acid. Its main function is the

inhibition of protein synthesis and results in antitumor, antiviral and antibacterial activity. Most of the other Alternaria toxins show cytotoxic activity in mammals, some of them are mutagenic like the altertoxins, while others are toxic to the unborn [1]. This application focusses on the determination of alternariol (AOH), alternariol monomethyl ether (AME), altenuene (ALT) and tenuazonic acid (TeA).

Alternaria alternata – determination of main metabolites

RESULTS

All samples and standards were provided from the Leibniz-Zentrum für Agrarlandschaftsforschung (ZALF) e.V. [2]. First a calibration was made using the standard ZALF 5 with five different injection volumes. AOH was calibrated in a range from 0.5 ng up to 5.0 ng (absolute), ALT and AME from 1 ng up to 10 ng (absolute) and TeA from 2ng up to 20 ng (absolute). AOH, ALT and AME were detected with a fluorescence detector. TeA was determined using a UV detector. **Tab 1** gives a short summary of the retention times of the substances and achieved correlation coefficients of calibration. As sample an extracted nutrient solution of

an *Alternaria* strain was used. The extract was divided into two fractions. One for the TeA determination (ZALF 7, dilution 1:20) and one for the other metabolites (ZALF 6, without dilution). **Fig 1** shows an overlay of the fluorescence traces of the standard ZALF 5 and sample ZALF 6. For ALT, a concentration of 1.18 ng/ μ L was calculated, for AOH 2.74 ng/ μ L and 0.36 ng/ μ L for AME. **Fig 2** Shows the UV traces of standard ZALF 5 and sample ZALF 7 for determination of TeA. In the second fraction, a value of 3.44 ng/ μ L TeA was calculated.

Tab.1 Retention times and correlation coefficients of altenuene, alternariol, alternariol monomethyl ether and tenuazonic acid calibration

Peak	Substance	Retention time	Correlation coefficient
1	ALT	4.063 min	0.99872
2	AOH	5.698 min	0.99853
3	TeA	8.424 min	0.99917
4	AME	10.664 min	0.99889




Fig. 1 Overlay of standard ZALF 5 (blue) and sample ZALF 6 (red), fluorescence detection, 1) ALT, 2) AOH, 3) TeA, 4) AME




Fig. 2 Overlay of standard ZALF 5 (blue) and sample ZALF 7 (red), UV-detection, 3) TeA

MATERIALS AND METHODS


An AZURA® Analytical HPLC Plus system was used for this application. The system consisted of an isocratic AZURA P 6.1L pump, an AZURA AS 6.1L autosampler, an AZURA DAD 6.1L, an AZURA CT 2.1 column thermostat and a RF 20 Axs fluorescence detector in combination with CBM 20 A under the Chromeleon™ software. The isocratic method [2] was applied for 30 minutes at a flow rate of 1 mL/min with a mixture of methanol and water in a ratio 70:30 (v/v). Furthermore 300 mg/mL zinc sulfate were added to the mobile phase. The column temperature was set to 30 °C. The substances were measured with an excitation at 253 nm and emission at 415 nm. The UV detector was set 280 nm. The used column, in a dimension 250 x 4.6 mm ID with precolumn, was filled with Prontosil Hypersorb 120-5 ODS silica.

CONCLUSION

It was possible to identify and quantify all *Alternaria alternata* metabolites with the described isocratic method. Using a fluorescence detector for enhanced sensitivity allows measurements of small amounts of mycotoxins even in a complex sample matrix. There are currently no statutory or guideline limits set for *Alternaria* mycotoxins because surveys to date have shown that their natural occurrence in foods is low and the possibility for human exposure is limited. The need for regulation is kept under review as new information becomes available [3].

REFERENCES

- [1] <https://www.romerlabs.com/en/knowledge-center/knowledge-library/articles/news/alternaria-toxins/>
- [2] Dr. Marina Müller, Leibniz-Zentrum für Agrarlandschaftsforschung (ZALF) e.V.
- [3] <http://www.micotoxinas.com.br/alertoxins.htm>



ADDITIONAL MATERIALS AND METHODS

Tab.A1 Method parameters

Eluent A	Methanol:Water 70:30 (v/v) with 300 mg/L ZnSO ₄ x 7 H ₂ O		
Gradient	isocratic		
Flow rate	1 mL/min	System pressure	ca. 160 bar
Column temperature	30 °C	Run time	30 min
Injection volume	1-10 µL	Injection mode	Partial loop
Detection UV	280 nm	Data rate	50 Hz
		Time constant	0.02 s
Detection FLD	Ex 253 nm / Em 415 nm	Data rate	100 Hz
		Time constant	0.01 s

Tab.A2 System configuration

Instrument	Description	Article No.
Pump	AZURA P6.1L, isocratic	APH30EA
Autosampler	AZURA AS 6.1L	AAA00AA
Detector	RF 20Axs with CBM-20A	A59201
Detector	AZURA DAD 6.1L	ADC11
Flow cell	LightGuide 50 mm, 6 µL	AMD59XA
Thermostat	AZURA CT 2.1	A05852
Column	Prontosil Hypersorb 120-5 ODS, VertexPlus Column 250 x 4.6 mm ID with precolumn	25VF180PYJ
Software	Chromeleon 7.2	

RELATED KNAUER APPLICATIONS

VFD0158 - Zearalenone and its major metabolites - a simple isocratic method

VFD0152 - Determination of Aflatoxin in milk