

Information for Installation of an Offline Version of PurityChrom®

1. PurityChrom operates under Microsoft Windows XP, Windows Vista, Windows 7, Windows 8 and Windows 10.
2. The user account settings have to be set to “Never notify” before installation. After the installation, the settings can be reset. (Start → Control Panel → User Accounts → User Accounts → Change User Account Control Settings → Never notify)
3. Restart your computer.
4. Insert the CD which includes PurityChrom and the configuration files of your system.
5. Start the installation as an administrator. The installer will ask you to choose the language of the software, your integrated UV-detector and the destination location for the files that will be installed. Do not change the given destination. Choosing “Demo-Version” in the installer is not necessary because when PurityChrom will be opened without a license file and dongle, the software switches automatically into demo mode.
6. If the installation is finished, copy and replace the PurityChrom.ini-file located on the CD (in the folder: *INI*) to C:\Windows
7. Copy and replace the PurityChrom.cfg-file located on the CD (in folder: *CFG*) to C:\PurityChrom.
8. Copy the background (*.bmp) located on the CD (in folder: *Visualisation\Background*) to C:\PurityChrom\Visualisation\Visualisation Backgrounds
9. Copy the visualization (*.vis) located on the CD (in the folder: *Visualisation\Files*) to C:\PurityChrom\Visualisation\Visualisation Files
10. **Installation under Windows versions older than 8 and 10:** Install PurityChrom. Go to C:\PurityChrom and right-click the PurityChrom.exe to open the properties. Go to compatibility and check “Run this program in compatibility mode for” and choose “Windows Vista (Service Pack 2). After this, check “Run this program as an administrator“ click “Apply“ and “Ok“, then close the window.
11. **Installation under Windows 8 or 10:** Install PurityChrom. Go to C:\PurityChrom and right-click the PurityChrom.exe to run the program as administrator. Open the visualization and the communication log once.
12. To open result or method files from your existing system, simply transfer them to C:\PurityChrom\Data.