

Important Information for the Installation of PurityChrom 5.9.36 Upgrade

Please refer to the Installation Guide and the KNAUER software manual for installation.

Don't install several versions of PurityChrom on one computer!

This Upgrade is a full version and will not update your existing PurityChrom software installation. If you install the **upgrade** on the same computer as the previous version, the previous version must be **completely removed** from the computer before the new version can be installed. The **existing folder must be deleted** before.

Before you delete the folder, make sure that your data have been copied / moved to a safe place.

Please copy or move minimum the following data to a safe place

PurityChrom.ini (C:\Windows),

PurityChrom.cfg, PurityChrom.lic (C:\PurityChrom)

Visualizations Backgrounds (C:\PurityChrom\Visualisation\Visualisation Backgrounds).

Visualizations C:\PurityChrom\Visualisation\Visualisation Files).

If the PurityChrom folder include data, please copy or move this data to a safe place.

Uninstall the PurityChrom completely and delete the existing folder.

The dongle license works independently from the installed PurityChrom version number.

The installation procedure, described below, is for installing the upgrade on the same computer as the previous version. **Instead of this, the software can also be installed on a new computer.**

A new configuration format was introduced with PurityChrom 5.9.x. This format is not backwards compatible: PurityChrom 5.9.28 can open the PurityChrom.cfg, acquired with older versions, but will **convert** them automatically into the new format. After the conversion of the PurityChrom.cfg, **older PurityChrom versions cannot open these file again.**

1. PurityChrom 5.9.36 operates under Microsoft Windows 10, 8.1, 8, 7, Windows Vista and Windows XP SP3.
2. Store PurityChrom.lic, PurityChrom.cfg and visualizations and all PurityChrom data you have already acquired, on a safe place outside of the PurityChrom installation folder.
3. Uninstall the "PurityChrom ..." from *Control Panel – Uninstall programs* (Windows 10 / 8.1 / 8 / 7) or *Control Panel - Add/Remove programs* (Windows XP).
4. Reboot your computer and delete the existing installation folder (C:\PurityChrom). Please check before, that you have already copied or moved your existing data.
5. Start the PurityChrom-Installer. Don't interrupt the installation.

6. If the installation is finished, move the
PurityChrom.ini to C:\Windows
PurityChrom.cfg to C:\PurityChrom
Visualisation.bmp to C:\PurityChrom\Visualisation\Visualisation Backgrounds
Visualisation.vis to C:\PurityChrom\Visualisation\Visualisation Files
7. Restart the PC.