



Read the corresponding technical documentation for handling and safety reasons.

AZURA tubing and fitting package for Foxy R1/R2 (A59113/A1407)



The standard version of Foxy R1 fraction collector can be used for flow rates up to 25 ml/min.

Use the 1/8" tubing and fitting package A59113 for flow rates up to 125 ml/min.

Use the adapter A1407 for 1/16" systems.



KNAUER

KNAUER

Wissenschaftliche Geräte GmbH

Hegauer Weg 38

14163 Berlin, Germany

Phone: +49 30 809727-0

Fax: +49 30 8015010

E-Mail: info@knauer.net

Internet: www.knauer.net